English

Instructions for Use: Leica 818/819 Blades & Leica DB80 LS/LX/HS Blades

Intended Use

The Leica 818, 819 and Leica DB80 series are highest quality blades specifically designed for cutting thin sections of human tissue specimens for use with Leica microtomes or Leica cryostats classified and labeled as in vitro diagnostic medical devices for histological medical diagnosis, e.g. cancer diagnosis.

Hazard warnings:

⚠	<ul style="list-style-type: none">Risk of injury. Blades are very sharp. Only trained laboratory personnel may use them. Never touch the blade edge! Only handle the blades by holding their sides between thumb and index finger crosswise to the blade edge. Wear cut-resistant safety gloves and protective eye wear when handling the blades. Do not use dispenser for disposal of the blades after cutting infectious or radioactive material.
----------------	---

1. Handling instructions:

- Inspect the dispenser upon receipt. If there is any damage, do not use and contact the manufacturer.
- If the dispenser does not eject blade, do not force the blade out of the dispenser. Do not use and contact the manufacturer.
- If the blade is not inserted directly into the blade holder, the blade dispenser must be held in one hand to remove a blade.
- To push out the blade, slowly push the blade slide all the way to the left until it stops.
- The blade must then be held by the other hand and removed from the blade dispenser.
- Insert the blade into the blade holder.
- If, after checking the additional instrument-specific, cut-relevant parameters, the cut result is not satisfactory, the blade must be moved or replaced.
- When cutting infectious or radioactive material, change the blade with every new specimen. This prevents unwanted contamination/transmission.
- Remove the blade from the blade holder.
- Observe the following disposal instructions.

2. Disposal:

- After being used to cut infectious or radioactive material, blades belong in the cut-resistant waste container provided for this purpose (not in the waste blade box on the bottom of the blade dispenser).
- Blades that have not come in contact with infectious or radioactive material can be disposed of into the waste blade box through the hole on the bottom of the blade dispenser.
- The holding capacity of the waste blade box is sufficient if all blades, without residues of section waste, are fed into the hole flat.

3. Technical data:

Dimensions:

- Leica 819 (low profile), DB 80 LS/LX: 80,00 x 8,00 x 0,254 mm
- Leica 818 (high profile), DB 80 HS: 80,00 x 14,00 x 0,317 mm

Material:

- High-alloyed stainless steel

Transport temperature/Storage humidity range:

- Transport temperature: -10 °C to +50 °C
- Storage humidity: < 80 % non-condensing

Operating temperature range:

- Blade dispenser: +15 °C to +40 °C
- Blade: -40 °C to +40 °C

4. Warranty:

Leica Biosystems Nussloch GmbH guarantees that the contractual product delivered has been subjected to a comprehensive quality control procedure based on the Leica in-house testing standards, and that the product is faultless and complies with all technical specifications and/or agreed characteristics warranted. The scope of the warranty is based on the scope of the concluded agreement. The warranty terms of your Leica sales organization or the organization from which you have purchased the contractual product shall apply exclusively.

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