Special Stain Kit

Modified Perls' Prussian Blue Iron Stain Kit

English	2
العربية (Arabic).....	5
简体中文 (Chinese Simplified).....	8
中國傳統的 (Chinese Traditional)	11
Dansk (Danish)	14
Nederlands (Dutch)	17
Français (French – Canada).....	21
Français (French – France)	25
Deutsch (German).....	29
Italiano (Italian).....	33
日本語 (Japanese)	36
한국어 (Korean)	39
Norsk (Norwegian)	42
Polski (Polish)	45
Português (Portuguese – Brazil).....	49
Português (Portuguese – Portugal)	52
Română (Romanian)	55
Русский (Russian)	58
Slovenski (Slovenian).....	62
Español (Spanish – Central America).....	65
Español (Spanish – Spain).....	68
Svensk (Swedish)	73
ไทย (Thai)	76
Türk (Turkish).....	79
Tiếng Việt (Vietnamese).....	82

Special Stain Kit

Modified Perls' Prussian Blue Iron Stain Kit

REF 38016SS7

Product Name

Modified Perls' Prussian Blue Iron Stain Kit

Intended Use

Detection/Measurement

The Leica Biosystems Modified Perls' Prussian Blue Iron Stain Kit does not detect or measure an analyte or marker. The Modified Perls' Prussian Blue Iron Stain Kit, when used with appropriate histological procedures, may be useful for the demonstration of hemosiderin in tissues sections and smears.

Product Function

The results obtained through use of Modified Perls' Prussian Blue Iron Stain Kit do not provide objective medical evidence. The coloration and contrast the Leica Biosystems Stain Kit provide to histologic specimens allows visualization of microscopic anatomy. This visualization, when interpreted by a trained professional, is utilized alongside other information such as the patient's medical history, physical condition, as well as results from other medical testing to render a medical diagnosis.

Specific Information Provided

Modified Perls' Prussian Blue Iron Stain Kit is not intended for the detection, definition or differentiation of a specific disorder, condition, or risk factor. The staining demonstrated with use of these products, when used as intended, provides trained professionals information which may define the physiological or pathological state of the tissue specimen.

Automation

Modified Perls' Prussian Blue Iron Stain Kit is not automated but can be used on automated staining platforms. Use on an automated platform should be validated at the point of use.

Qualitative/Quantitative

The Leica Biosystems Modified Perls' Prussian Blue Iron Stain Kit is qualitative stain.

Specimen Type

Modified Perls' Prussian Blue Iron Stain Kit can be used with any paraffin embedded human or animal specimen.

Testing Population

The Leica Biosystems Modified Perls' Prussian Blue Iron Stain Kit is intended for use with any patient requiring evaluation of biopsy or resection tissue for the assessment of a suspected pathology or disease.

Intended User

Modified Perls' Prussian Blue Iron Stain Kit is intended for use by qualified laboratory personnel and/or designee of the laboratory.

In Vitro Diagnostic

Modified Perls' Prussian Blue Iron Stain Kit is intended for *in vitro* diagnostics use only.

Test Principle

Hemosiderin is a large iron storage complex that is composed of ferric ion, the protein ferritin and denatured ferritin. It is insoluble in alkalis and water but soluble in strong acids even after fixation. Dilute hydrochloric acid releases loosely bound ferric iron (in the form of ferric chloride) from the protein complexes. The ferric ions can be easily demonstrated by Perls' method for ferric iron. The "Prussian Blue" reaction is the combination of potassium ferrocyanide with ferric (+3) salts in an acidic solution to form the bright blue ferric ferrocyanide complex known as Prussian Blue.

Calibrators & Controls

Modified Perls' Prussian Blue Iron Stain Kit does not require the use of any calibrators or controls.

Reagent Limitations

No reagent limitations are applicable to this product.

Applicable Products

Product Code	Material Description
38016SS7	Modified Perls' Prussian Blue Iron Stain Kit
38016SS7A	Potassium Ferrocyanide Solution, 500 mL
38016SS7B	Hydrochloric Acid Solution, 500 mL
38016SS3B	Nuclear Fast Red (Kernechtrot) Solution, 500 mL

Materials Not Included

Modified Perls' Prussian Blue Iron Stain Kit protocol requires the use of graded alcohols, xylene, or xylene substitutes, deionized or distilled water. Positive control slide(s), not included in this kit, should be included in each run.

Any iron rich tissue or spleen that normally contains small amounts of ferric iron might be a good positive control.

Devices Required

Modified Perls' Prussian Blue Iron Stain Kit can be used on any automated staining platform or with a manual staining method.

Special Stain Kit

Modified Perls' Prussian Blue Iron Stain Kit

REF 38016SS7

Storage and Stability

Store at room temperature (15-30 °C).

CAUTION: Do not use after the expiration date.

In Use Stability

User discretion should be utilized when determining in-use stability.

Sterility

Modified Perls' Prussian Blue Iron Stain Kit components are not sterile products.

Warnings/Precautions

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations. Refer to Material Safety Data Sheet and product labeling for any updated risk, hazard, or safety information.

Infectious Material Status

Modified Perls' Prussian Blue Iron Stain Kit does not include any infectious material. However, specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

Special Facilities

Modified Perls' Prussian Blue Iron Stain Kit should be used per facility guidelines.

Specimen Handling

- Fixation
Any general fixative including but not limited to 10% neutral buffered formalin, alcoholic formalin and alcohol may be used. Avoid iron containing materials and specimen containers while fixing since these may result in the contamination of tissue samples. Acidic fixatives and potassium dichromate should be avoided as they may remove some of the iron deposits.
- Paraffin Sections
Following processing and paraffin embedding section specimens at 4-6 microns. Smears and touch preparations may also be used.

Preparation for Use

Mix equal parts of the hydrochloric acid solution and the potassium ferrocyanide solution immediately before use.

Direction for Use

Conventional Staining Protocol

1. Deparaffinize with xylene and rehydrate through graded alcohols to deionized or distilled water.
2. Rinse in deionized or distilled water (6 changes).
Rinse well with deionized water before placing slides in the Potassium Ferrocyanide-Hydrochloric Acid working solution, as iron in the water or rust in the tap fixtures may cause false positive staining.
3. Place in working Potassium Ferrocyanide-Hydrochloric Acid working solution for 30 minutes.
4. Rinse in deionized water (5 changes).
Wash well as traces of iron will form granular red deposits with nuclear fast red.
5. Counterstain in Nuclear Fast Red (Kernechtrot) for 2-5 minutes.
6. Rinse in deionized or distilled water (2 rinses).
7. Dehydrate through graded alcohols, clear and mount.

Note: To demonstrate minute amounts of iron, a 5-minute pretreatment in 1% Potassium Ferrocyanide solution before treatment with the Potassium Ferrocyanide-Hydrochloric Acid solution can aid in its demonstration.

Table 1. Example of Conventional Modified Perls' Prussian Blue Iron Stain Protocol.

Steps	Action	Chemical	Time (mm: ss)
1-3	Deparaffinize	Xylene	3:00
4-5	Hydration	100% Alcohol	2:00
6	Hydration	80% or 95% Alcohol	1:00
7-12	Rinse	Deionized Water	0:30
13	Stain	Potassium Ferrocyanide-Hydrochloric Acid working solution	30:00
14-18	Rinse	Deionized Water	0:30
19	Stain	Nuclear Fast Red (Kernechtrot)	2:00 to 5:00
20-21	Rinse	Deionized Water	0:30
22	Dehydration	95% Alcohol	1:00
23-24	Dehydration	100% Alcohol	1:00
25-26	Clearing	Xylene	2:00

Special Stain Kit

Modified Perls' Prussian Blue Iron Stain Kit

REF 38016SS7

Microwave Staining Protocol

Exercise caution when using the microwave to heat any solution or reagent. The microwave must be properly ventilated to prevent the accumulation of fumes in the laboratory. Microwave transparent Coplin jars and caps should be used during the staining process. The caps should be loosely applied to prevent spills. Caps with ventilation holes also may be used. All microwaves should be used in accordance with the manufacturer's instructions. All microwave steps were conducted at a power setting of 800 watts unless otherwise noted.

1. Deparaffinize with xylene and rehydrate through graded alcohols to water.
2. Rinse well in deionized or distilled water. Rinse well with deionized water before placing slides in the Potassium Ferrocyanide-Hydrochloric Acid working solution, as iron in the water or rust in the tap fixtures may cause false positive staining.
3. Place slides in a plastic Coplin jar containing 40 mL of the Potassium Ferrocyanide-Hydrochloric Acid working solution and microwave at 400 watts for 30 seconds.
4. Gently mix the solution by gentle swirling and allow to stand for 10 seconds.
5. Rinse well in several changes of deionized water. Wash well as traces of iron will form granular red deposits with nuclear fast red.
6. Counter stain with Nuclear Fast Red for 5 minutes at room temperature.
7. Dehydrate through graded alcohols, clear and coverslip.

Readiness for Use

Once appropriate staining protocol is chosen and bath-layout is created, pour all the reagent into the reagent vessel. Place the reagent vessel back into the respective station.

Quality Control

Any iron rich tissue (hemochromatosis and hemosiderosis). Small amounts of ferric iron are found normally in the spleen.

Expected Results

- Ferric Iron deposits — Blue
- Nuclei — Red
- Cytoplasm — Pink

Analytical Performance

The Leica Biosystems Perl's Prussian Blue Iron Stain Kit is not used to detect a specific analyte or marker. This product is used in the demonstration of hemosiderin in tissues sections and smears. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

Clinical Performance

The Leica Biosystems Perl's Prussian Blue Iron Stain Kit is not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the Leica Biosystems Bluing Agents in a clinical setting.


Disposal

Spent or expired components of the Modified Perls' Prussian Blue Iron Stain Kit should be discarded in accordance with organizational, local, state, and federal regulations.

Leica Biosystems Richmond, Inc.
5205 Route 12
Richmond, IL 60071
USA
(1-844-534-2262)

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