# ST Infinity 2.0 H&E Staining System

<b>English</b>	<b>2</b>
<b>Arabic</b>	<b>8</b>
<b>Chinese</b>	<b>14</b>
<b>Croatian</b>	<b>20</b>
<b>Czech</b>	<b>26</b>
<b>Danish</b>	<b>32</b>
<b>French</b>	<b>38</b>
<b>German</b>	<b>44</b>
<b>Greek</b>	<b>50</b>
<b>Hungarian</b>	<b>56</b>
<b>Italian</b>	<b>62</b>
<b>Japanese</b>	<b>68</b>
<b>Korean</b>	<b>74</b>
<b>Lithuanian</b>	<b>80</b>
<b>Macedonian</b>	<b>86</b>
<b>Norwegian</b>	<b>92</b>
<b>Polish</b>	<b>98</b>
<b>Portuguese</b>	<b>104</b>
<b>Russian</b>	<b>110</b>
<b>Serbian</b>	<b>116</b>
<b>Slovakian</b>	<b>122</b>
<b>Spanish</b>	<b>128</b>
<b>Swedish</b>	<b>134</b>
<b>Turkish</b>	<b>140</b>

## ST Infinity 2.0 H&E Staining System Instructions for Use

Content: 5 bottles

Optimized for Use on the Leica ST5010 and ST5020 Staining Instruments.

For Laboratory Use

Item#: 3802098

The ST Infinity 2.0 H&E Staining System contains five proprietary reagent components and staining protocols that when used together provide quality, consistency and ease of use. This system has been optimized for the Leica ST5010 AutoStainer XL and ST5020 Multistainer. The ST Infinity 2.0 H&E Staining System reagents are intended for *in vitro* diagnostic use only.

This system delivers on a promise of being easy to use by allowing the technologist to initiate and terminate use of all of the components of the kit simultaneously while achieving a minimum of two thousand consistent high quality stained slides.

### Intended Use

The reagents in this kit are intended for *in vitro* use only. The reagents in this kit are intended for use on the Leica ST5010 AutoStainer XL and ST5020 Multistainer. The kit is used to stain histological samples with Hematoxylin and Eosin for subsequent microscopic evaluation.

### Reagents Provided:

- ST Hemalast 2.0 470 ml ..... Item #: 3802098B
- ST Hematoxylin 2.0 470 ml ..... Item #: 3802098A
- ST Differentiator 2.0 470 ml ..... Item #: 3802098C
- ST Bluing Agent 2.0 470 ml ..... Item #: 3802098E
- ST Eosin 2.0 470 ml ..... Item #: 3802098D

### Safety and Precautions:

**CAUTION:** Please follow local storage requirements for flammable products. For further information refer to Safety Data Sheet (SDS)

### Storage and Stability:

Kit components should be stored at room temperature (15 - 30° C). Do not use after the expiration date.

### Specimen Preparation:

#### Fixation

Any general fixative including but not limited to neutral buffered formalin may be used. **Note:** Underfixation or overfixation may result in compromised staining quality.

#### Paraffin Sections

Following processing and paraffin embedding, section tissues at standard thickness (2 - 5 µm).

#### Slides and Specimens

Apex Superior Adhesive or equivalent slides should be used for best results. Prior to staining, slides should be dried at 58° - 60° C for 20 - 30 minutes.

### Directions for Use of the ST Infinity 2.0 H & E Staining System on the ST5010 XL AutoStainer

Settings for operational parameters for the ST5010 XL AutoStainer as well as a recommended bath lay-out and protocol schedules are provided in order to optimize the performance and consistency of the ST Infinity 2.0 H&E Staining System.

**Failure to follow these instructions may result in compromised staining quality as well as staining capacity.**

#### ST5010 AutoStainer XL - Bath Lay-Out

The suggested Reagent Bath Lay-out for the ST5010 AutoStainer XL is shown below (Figure 1). Empty the entire contents of each of the bottles containing ST Hemalast 2.0, ST Hematoxylin 2.0, ST Differentiator 2.0, ST Bluing Agent 2.0 and ST Eosin 2.0 into the designated reagent vessels and positions.

**Note:** All of the reagents in each ST Infinity 2.0 Staining System kit must be used together to achieve the expected performance. Substitution of non-ST Infinity 2.0 H & E Staining System reagents or reagents from another ST Infinity 2.0 Staining System kit will potentially compromise performance.

**Note:** Reagent vessels should be covered when not in use to reduce evaporation. Discard the staining reagents (ST Hemalast 2.0, ST Hematoxylin 2.0, ST Differentiator 2.0, ST Bluing Agent 2.0, ST Eosin 2.0), if the 2000 slides staining capacity is not met within 8 calendar days after pouring into the reagent vessels.

#### ST5010 XL AutoStainer - Staining Protocol

The suggested immersion times and instrument settings for the ST5010 XL AutoStainer are detailed in Table 1. Please follow instructions for programming the ST5010 XL AutoStainer in the ST5010 XL AutoStainer Automated slide stainer instruction manual (manuals available online at [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)). A protocol designed to produce a moderate Hematoxylin intensity is shown in Table 1. The immersion times for the ST Hematoxylin 2.0 and ST Eosin 2.0 may be further varied within the parameters outlined in Table 2 to meet staining preferences. All other time requirements outlined in Table 1 should be adhered to in order to ensure quality and consistency of results.

Figure 1. Leica ST 5010 Bath Lay-Out												
1 Xylene	2 Xylene	3 Xylene	4 100% Reagent Alcohol	5 100% Reagent Alcohol	6 95% Reagent Alcohol	7	8 Tap Water Wash 1	9 Tap Water Wash 2	10 Tap Water Wash 3	11 Tap Water Wash 4	12 Tap Water Wash 5	Oven
Xylene Exit	18 Xylene	17 Xylene	16 100% Reagent Alcohol	15 100% Reagent Alcohol	14 100% Reagent Alcohol	13 ST Eosin 2.0	12 80% Reagent Alcohol	11 ST Bluing Agent 2.0	10 ST Differ- entiator 2.0	9 ST Hema- toxylin 2.0	8 ST Hema- last 2.0	Load

**Table 1: ST 5010- Staining Protocol**

ST 5010 AutoStainer				
Settings: Dips = 2, Withdrawal Speed = 9				
Step	Station	Reagent	Time	Exact
1	Load	N/A	0:00	N/A
2	1	Xylene	2:00	No
3	2	Xylene	2:00	No
4	3	Xylene	2:00	No
5	4	100% Reagent Alcohol	1:00	No
6	5	100% Reagent Alcohol	1:00	No
7	6	95% Reagent Alcohol	1:00	No
8	Water Rinse 5	Water Rinse	1:00	No
9	8	ST - HemaLast 2.0	1:00	No
10	9	ST-Hematoxylin 2.0	5:00	Yes
11	Water Rinse 4	Water Rinse	2:00	No
12	10	ST - Differentiator 2.0	0:30	Yes
13	Water Rinse 3	Water Rinse	1:00	No
14	11	ST - Bluing Agent 2.0	1:00	Yes
15	Water Rinse 2	Water Rinse	1:00	No
16	12	80% Reagent Alcohol	1:00	No
17	13	Eosin 2.0	0:45	Yes
18	14	100% Reagent Alcohol	1:00	Yes
19	15	100% Reagent Alcohol	1:00	Yes
20	16	100% Reagent Alcohol	1:00	Yes
21	17	Xylene	1:00	No
22	18	Xylene	1:00	No
23	Exit	Xylene	1:00	No

**Table 2. ST Hematoxylin 2.0 and ST Eosin 2.0 Variations**

Reagent	Staining Time Limits
Hematoxylin	2 – 10 min
Eosin	30 – 60 sec

**Alcohol and Xylene Rotation**

The alcohols and xylenes used for deparaffinization and rehydration should be rotated after every 600 slides. (Reagent rotation is defined as follows: Empty the most contaminated reagent of a series and refill the container. Move the remaining container(s) up one space in the sequence. Place the newly filled container into the last position of the series.) The alcohols and xylenes used for dehydration and clearing before coverslipping also should be rotated after every 600 slides as described above.

The 95% and 80% alcohols (steps 7 and 16) should be discarded and replaced with fresh 95% and 80% alcohol solutions after 600 slides.

**Failure to rotate and change reagents as advised may result in inadequate deparaffinization as well as excessive carry-over and dilution of reagents or solvents.**

**Expected Results**

By following the instructions for use, a single ST Infinity 2.0 Staining System Kit should provide at least 2,000 consistent high quality stained slides.

## Directions for Use of the ST Infinity 2.0 H & E Staining System on the ST5020 Multistainer

Settings for operational parameters for the ST5020 Multistainer as well as bath lay-out and protocol schedules are provided in order to optimize the performance and consistency of the ST Infinity 2.0 H&E Staining System.

**Failure to follow these instructions may result in compromised staining quality as well as staining capacity.**

### ST 5020 Multistainer - Bath Lay-Out

The suggested Reagent Bath Lay-out for the ST5020 Multistainer is shown below in Figure 2. Empty the entire contents of each of the bottles containing ST Hemalast 2.0, ST Hematoxylin 2.0, ST Differentiator 2.0, ST Bluing Agent 2.0 and ST Eosin 2.0 into the designated reagent vessels and positions.

**Note:** All of the reagents in each ST Infinity 2.0 Staining System Kit must be used together to achieve the expected performance. Substitution of non-ST Infinity 2.0 Staining System reagents or reagents from another ST Infinity 2.0 Staining System kit will potentially compromise performance.

**Note:** Reagent vessels should be covered when not in use to reduce evaporation. Discard the staining reagents (ST Hemalast 2.0, ST Hematoxylin 2.0, ST Differentiator 2.0, ST Bluing Agent 2.0, ST Eosin 2.0), if the 2000 slides staining capacity is not met within 8 calendar days after pouring into the reagent vessels.

Figure 2: Leica ST 5020 Bath Lay-Out											
1 Xylene	2 Xylene	3 Xylene	4 100% Reagent Alcohol	5 100% Reagent Alcohol	6 95% Reagent Alcohol	7	8	9 Tap Water Wash	10 Tap Water Wash	11 Tap Water Wash	12 Tap Water Wash
13 Xylene	14 Xylene	15 Xylene	16 100% Reagent Alcohol	17 100% Reagent Alcohol	18 100% Reagent Alcohol	19 ST Eosin 2.0	20 80% Reagent Alcohol	21 ST Bluing Agent 2.0	22 ST Differ- entia- tor 2.0	23 ST Hema- toxylin 2.0	24 ST Hema- last 2.0
25 Xylene Exit	26 Xylene Exit	27 Xylene Exit	28 Xylene Exit	29	30	31	32	L33 Load	L34 Load	L35 Load	L36 Load

The immersion times and instrument settings for the ST5020 Multistainer are included in Table 3. Please follow instructions for programming the ST5020 Multistainer in the ST5020 Multistainer Operating Instructions manual (manuals available online at [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)). A protocol designed to produce a moderate Hematoxylin intensity is shown in Table 3. The immersion times for the ST Hematoxylin and ST Eosin may be further varied within the parameters outlined in Table 2 to meet staining preferences. All other time requirements outlined in Table 3 should be adhered to in order to ensure quality and consistency of the staining results.

**Table 3. ST5020 Multistainer – Staining Protocol**

ST 5020 MultiStainer				
Settings: Dips = 2, Lift Speed = 1				
Step	Station	Reagent	Time	Exact
1	Load	N/A	0:00	N/A
2	1	Xylene	2:00	No
3	2	Xylene	2:00	No
4	3	Xylene	2:00	No
5	4	100% Reagent Alcohol	1:00	No
6	5	100% Reagent Alcohol	1:00	No
7	6	95% Reagent Alcohol	1:00	No
8	Water Rinse 12	Water Rinse	1:00	No
9	24	ST – HemaLast 2.0	1:00	No
10	23	ST- Hematoxylin 2.0	5:00	Yes
11	Water Rinse 11	Water Rinse	2:00	No
12	22	ST – Differentiator 2.0	0:30	Yes
13	Water Rinse 10	Water Rinse	1:00	No
14	21	ST - Bluing Agent 2.0	1:00	No
15	Water Rinse 9	Water Rinse	1:00	No
16	20	80% Reagent Alcohol	1:00	No
17	19	ST- Eosin 2.0	0:45	Yes
18	18	100% Reagent Alcohol	1:00	Yes
19	17	100% Reagent Alcohol	1:00	Yes
20	16	100% Reagent Alcohol	1:00	Yes
21	15	Xylene	1:00	No
22	14	Xylene	1:00	No
23	Exit	Xylene	1:00	No

**Alcohols and Xylene Rotation**  
 The alcohols and xylenes used for deparaffinization and rehydration should be rotated after every 600 slides. (Reagent rotation is defined as follows: Empty the most contaminated reagent of a series and refill the container. Move the remaining container(s) up one space in the sequence. Place the newly filled container into the last position of the series.) The alcohols and xylenes used for dehydration and clearing before coverslipping also should be rotated after 600 slides as described above.

The 95% and 80% alcohols (steps 7 and 16) should be discarded and replaced with fresh 95% and 80% alcohol solutions after 600 slides.

**Failure to rotate and change reagents as advised may result in inadequate deparaffinization as well as excessive carry-over and dilution of reagents or solvents.**

**Expected Results**  
 By following the instructions for use, a single ST Infinity 2.0 Staining System Kit should provide at least 2,000 consistent high quality stained slides.

Technical Notes:

1. Routine control slides containing tissue fixed and processed in a similar manner to the test specimens should be included to ensure that the reagents are performing and functioning properly.
2. Reagent vessels should be inspected at periodic intervals to insure the integrity and cleanliness of the vessels. If necessary, clean the vessels in accordance with the instructions provided in the operating instructions for the ST5010 XL AutoStainer and ST5020 Multistainer. If detergents are used, rinse vessels with adequate water to remove any residual detergent within the vessel prior to use.

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Leica Biosystems Richmond, Inc.
Sponsor: Leica Biosystems Melbourne Pty Ltd