Overview

Manual Input

Barcode Input

Managing Slide Levels

Job Queue and Pending Queue

Settings and Preferences

Multiple Printers

Caution!

The user needs to visually check each printed slide for legibility after the printing process.

Caution!

All laboratories must perform their own tests to ensure that the ink is resistant against the various reagents the slides will subsequently be exposed to.

A wide range of factors beyond Leica's control can have negative effects on the results. The laboratory operating the unit shall bear full responsibility for the legibility of the imprint after processing with reagents.

Manufactured by:
Primera Technology, Inc.
Two Carlson Parkway North, Suite 375
Plymouth, MN 55447-4446
U.S.A.

Distributed by:
Leica Biosystems Nussloch GmbH
Heidelberger Straße 17-19
DE-69226 Nussloch
Germany

Notices: The information in this document is subject to change without notice. NO WARRANTY OF ANY KIND IS MADE WITH REGARD TO THIS MATERIAL, INCLUDING, BUT NOT LIMITED TO, THE IMPLIED WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A PARTICULAR PURPOSE. No liability is assumed for errors contained herein or for incidental or consequential damages in connection with the furnishing, performance, or use of this material. This document contains proprietary information that is protected by copyright. All rights are reserved. No part of this document may be photocopied, reproduced, or translated into another language without prior written consent.

Trademark Acknowledgments: Windows is a registered trademark of Microsoft Corporation. All other trademarks are the property of their respective owners.

Printing History

Edition 1.3, #041714, Copyright 2012, All rights reserved.

FCC Compliance Statement: This device complies with part 15 of the FCC rules. Operation is subject to the following two conditions: (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

For Users in the United States: This product is intended to be supplied by a UL listed Direct Plug-In Power Supply marked "Class 2" or a UL listed ITE Power Supply marked "LPS" with output rated 12VDC, 4.5A or higher. This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to Part 15 of the FCC Rules. In a domestic environment this product may cause radio interference, in which case the user may be required to take adequate measures. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications.

However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Re-orient or relocate the receiving antenna.
- Increase the separation between the equipment and receiver.
- Connect the equipment into an outlet on a circuit different from that to which the receiver is connected.
- Consult the dealer or an experienced radio/TV technician for help. Use of shielded cables is required to comply with the Class A limits of Part 15 of the FCC Rules. You are cautioned that any changes or modifications not expressly approved in this manual could void your authority to operate and/or obtain warranty service for this equipment.

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CAUTION!

TO PREVENT FIRE OR SHOCK HAZARD, DO NOT EXPOSE THE UNIT TO RAIN OR MOISTURE. TO REDUCE THE RISK OF ELECTRIC SHOCK, DO NOT REMOVE EXTERIOR PANELS. NO USER-SERVICEABLE PARTS INSIDE. REFER SERVICING TO QUALIFIED SERVICE PERSONNEL. OPERATE THE UNIT WITH ONLY THE PROPER ELECTRICAL SPECIFICATIONS AS LABELED ON THE PRINTER AND AC ADAPTER.

CAUTION!

USE OF CONTROLS OR ADJUSTMENTS OR PERFORMANCE OF PROCEDURES OTHER THAN THOSE SPECIFIED HEREIN MAY RESULT IN HAZARDOUS RADIATION.

CAUTION!

THE USER NEEDS TO VISUALLY CHECK EACH PRINTED SLIDE FOR LEGIBILITY AFTER THE PRINTING PROCESS.

CAUTION!

ALL LABORATORIES MUST PERFORM THEIR OWN TESTS TO ENSURE THAT THE INK IS RESISTANT AGAINST THE VARIOUS REAGENTS THE SLIDES WILL SUBSEQUENTLY BE EXPOSED TO A WIDE RANGE OF FACTORS BEYOND LEICA'S CONTROL CAN HAVE NEGATIVE EFFECTS ON THE RESULTS. THE LABORATORY OPERATING THE UNIT SHALL BEAR FULL RESPONSIBILITY FOR THE LEGIBILITY OF THE IMPRINT AFTER PROCESSING WITH REAGENTS.

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
[7.2](#) Printing to Multiple Printers

Section 1: PTLab Overview

PTLab with a Slide Printer Connected.

Some hints and functions of the software might not be available with the current software version.


The printer can be connected to a network via the shared printer function of your computer. This shared printer function can e.g. be used to connect the printer e.g. to an LIS systems.

Section 2: Manual Input


This is the least difficult data entry method to implement but still presents the possibility for data errors to be introduced by the user. Enter the data by typing or by choosing from menus. The type of menu, the order, the labels and even the slide type are all customizable using the PTTemplate software.

The software comes preloaded with several template options. These options are meant to provide you with a sampling of the printer’s capability. They are not necessarily suited for practical use. Before using the printer in a lab you will probably want to edit an existing template or create an entirely new template using the PTTemplate software (See the PTTemplate software manual).

To quickly choose a preloaded template click the drop down arrow to view recent templates. Choose one.


You may also click the browse button to the right of the currently selected template to view the same list of available templates with a preview.


Once you have selected the template you must enter the information you would like to print.

1. For this example SampleTemplate1.svg was selected. Choose a color swatch to be printed.

2. Type the necessary information in the text fields that require manual input. The preview will update as you enter the information.

3. This template also includes two other automatic fields/features:


Barcode

A barcode encoded with the patient ID and the Slide Level will be printed.

Date

The current date will automatically be printed.

4. Finally, this template requires Slide Type 1. Slide type could be set to any number between 1 and 100 or removed entirely. Slide types can be used as internal designations in a lab for different types of slides such as charged slides (See section 7.6)

There are many other template fields/features available to preview in the other preloaded slide templates. Any of the fields found in the other preloaded templates can be added to a custom template that you can create in the PTTemplate software.


Section 3: Barcode Input

3.1 Overview (Sample / Demo Scan)

If the template supports it you can scan an existing linear or 2D barcode using the optional scanner from another label. This will allow you to directly transfer some or all the information on the label to the slide printer template. If all information to complete the template is present, you can either choose to automatically print the slide or send the job to the Pending Jobs queue.

If the scanned barcode contains the proper fields, the software can automatically choose the template needed, automatically adding the information to be printed, set a slide type (section 7.6), and even set the number of copies. Follow the instructions below to see an example of how this can work.

Connect your scanner and install the driver. Open PTLab. To verify that PTLab is communicating with your scanner, click the scanner icon at the top right part of your screen.

Scanner Status should report Connected and Running. If it cannot communicate you will see an error message and a “Reconnect” button will be present. Click “Reconnect”.


On this screen you may also display the character string which is actually encoded into any scanned barcode. Simply scan the barcode and the data string will appear.

Close this window before attempting a scan/print or scan/pending jobs operation.

2. At any point while the software is open, simply scan a barcode to print its corresponding template or send the template to the pending jobs queue for review. The following sample label is preconfigured to automatically print or send it to pending jobs.


3. If “Move all scanned jobs to pending queue” is **unchecked**, the associated template will automatically print the slide.


Tick box "Move all LIS jobs to the pending queue" is not supported with the current software version.

4. If “Move all scanned jobs to pending queue” is **checked**, the job will appear in the pending jobs queue.


5. Click on the job to display the template and review the data before printing.

6. Click Print to print the slide.

3.2 Scanned Barcode Requirements

In the sample case described in the previous section, the printed barcode on the specimen label required that the barcode be encoded with the following data in the order shown. This is not a requirement for all templates. You can edit or create a template to accept data in the format of an existing 2D barcode that is currently being printed. This can be accomplished using the PTemplate software.

3.3 Incomplete Templates

If the template has more input fields than are provided by the data string scanned in from the barcode, the job will automatically go into the Pending Jobs queue. This will be the case regardless of whether you have “Move all scanned jobs to pending queue” checked or unchecked.


Annotations:

- Scanned data will be shown here if fields are missing. (points to Pending Jobs list)
- Extra Field not found in Barcode Scan Data (points to the empty Level field)

You may now choose to fill in the missing information or print the slide without the extra field information.


Section 4: Managing Slide Levels

The slide level is tracked by the software and displayed on the status bar at the bottom.

Note: The printer does not have a physical slide level sensor. The software tracks when a slide has been printed and subtracts one from the total. If the original total was incorrect for some reason, the slide level reported by the software would also be incorrect.

To set the slide level, click on the slide icon. The following screen will appear.

This window will also appear if the slide cartridge is taken out of the printer with 2 or less slides remaining. When you replace the slide cartridge you will be prompted with the window above.

An explanation of the function of each button is below.

- **Reset Slide Level** resets the slide level to 100
- **Set Slide Level** sets to a specific number that you can set by typing the number to the left. This number will be remembered the next time this window appears. This is useful if your slides do not come in packs of 100. You can easily reset the slide counter to equal whatever number of slides your standard box includes.


- **Hide Slide Level** will stop this window from appearing after a cartridge change and remove the slide level icon from the status bar. You can always track the slide level by looking at the actual slides remaining through the clear blue slide cartridge. You can re-enable the slide icon in the settings area.
- **Ignore** dismisses the window and does not change the slide level. Use this if you removed and inserted the same slide cartridge but did not alter the slide count.

Section 5: Job Queue and Pending Jobs

5.1 Job Queue

The Job Queue displays all current jobs which are printing or waiting to be printed.

The currently printing job will show a status of “printing”. The waiting job will appear below the printing job.


When a job has been successfully printed or cancelled it will be moved to the Completed Jobs window. The results of the job will be displayed below the status.

You may also cancel a currently printing job, pause jobs while they are printing, clear all jobs from the list, or reprint jobs using the buttons at the bottom of the window.

Click “Reprint Multiples” to reprint specific slides in a job. Enter the slide range you would like to reprint including the starting slide, the total slides in the set and the number of slides you want to print. The slides the software will be printing are shown in the white box.

Note: Use this feature ONLY if your slide has a number range component.

Number Range

5.2 Pending Jobs

The Pending Jobs area is only used if data for your templates is received from a connected barcode scanner. You can set PTLab to automatically print any slide as soon as the necessary data is received by the template -or- you can set PTLab to automatically place all jobs from a scanner in the pending queue.

Note: If there is insufficient data to completely fill all fields in your template, the job will always go to the pending jobs queue no matter the setting. Click on the job to open the template and finish filling fields or just click “Print”.

If “Move all scanned jobs to pending queue” is **unchecked**, the associated template will automatically print the slide.

Tick box "Move all LIS jobs to the pending queue" is not supported with the current software version.

Click on the job to display the template and review the data before printing. You may enter additional data in empty fields, leave them blank or change data in any field. Click “Print” when you are ready to print.

Section 6: Settings and Preferences


The settings window contains the following information:

- Various software settings
- Printer firmware version and printing statistics
- Template management
- Logging options

Open the settings area by clicking the tool icon in the upper right corner of the software.


To open or close each settings area click on the ribbon header above each section.

In the following sections you will find descriptions of each of the functions within each settings area.


6.1 Settings


Tick box "Move all LIS jobs to the pending queue" is not supported with the current software version.

- **Full Screen Mode.** Check or uncheck to toggle between automatic full screen or standard window mode. You must restart the software for this setting to take effect.

- Require Login.** Check this box to require a login username when the software starts. This is not necessarily meant as a security feature. If this is checked, a user will be prompted to scan their 2D barcode on their security card or enter a username. This username can be used by a template to print on the slide, encode on the slide barcode or simply be logged by the software. The following screen will appear if this option is checked:

Once this setting is enabled you may add new users or manage current users by clicking the admin icon. Here you can switch the current user, delete a user, or setup a new user. Check the Administrator box to allow this user to manage users. You may set “Printed on Slide” to users initials or some other code that will discreetly identify the user on the slide.

- Show Slide Level.** Uncheck this box to remove the slide level icon from the start bar and stop tracking the number of slides in a cartridge.


- **Language.** Choose one of the many available languages. The text will immediately be changed to that language.
- **Move all scanned jobs to pending queue.** If unchecked, the associated template will automatically print the slide. If checked, the job will appear in the pending jobs queue. See section 6 for more information.
- **Color Print Heat / Black Print Heat.** Increase the heat settings to improve print quality. The minimum heat setting that achieves acceptable print quality should be used. Increasing print heat shortens print head life and increases the likelihood of a ribbon break. Separate heat settings for color and black ribbons can be stored. If using a color ribbon, the color heat setting will automatically be used. If using a black ribbon, the black heat setting will automatically be used.
- **Current Slide Type.** (Visible with Slide Printer Connected) See section 7.6
- **Enable auto-feed mode.** (Visible with Slide Printer Connected) See section 7.7

6.2 Printer Information

This window displays the current firmware version on the printer, total slides printed and other information that may be useful to tech support.

6.3 Template Management

Buttons "Clear LIS Default" and "Set LIS Default" are not supported with this software version

All current templates are displayed. Click "Remove Template" to remove a template from the list of templates available to the user via the "Choose Template" drop down menu or the browse button.

Any template that has been removed can be recovered using the "Recover Removed Templates" button. You will be brought to the deleted Templates folder to recover your template.


If you have exported a template (for another computer) from the PPTemplate editing software you can add that template by clicking the "Add Template" button.

Templates exported (for this computer) from the PPTemplate software installed on the same computer as the PTLab software do not need to be added manually in this way. They will be automatically added to the list.

You may also Set your Scanner Default in this location. See section 3 for more information.


6.4 Software Info

The Software Info area will tell you the current version of the software.

Buttons and License information referring to PTLab PE are not supported with this software version.

6.5 Logging


Several fields can be logged to the log file located where you choose.

- **Date.** Logs the date the slide was printed.
- **Time.** Logs the time the slide was printed.
- **Template Info.** Logs the Template fields printed for each slide.
- **Output Barcode.** Logs the information encoded in the barcode printed on the slide.
- **User Name.** Logs the Username associated with a template, if required.
- **User ID.** Logs the ID associated with a Username associated with a template.

6.6 Set the Slide Type on a Slide Cartridge

The slide type setting can be used to encode the blue slide cartridge with a unique number. Slide type can be set to any number between 1 and 100. Slides types can be used as internal designations in a lab for different types of slides such as charged slides. An encoded slide cartridge can be used in combination with a template that requires a specific slide type. In this way you can ensure certain procedures are done on the appropriate slides.

Note: To require a specific slide type your template must have either an image list or text list which has the “Set default copies and slide type” option checked in the “Edit Text/Image List” area within PPTemplate.

PPTemplate


When the cartridge is inserted, the slide type will be displayed below the slide icon on the status bar.

If a slide cartridge encoded with a slide type other than what the template requires is inserted into the printer the following error will appear.

Slide Cartridge


The default designation for any new slide cartridge is “1”.

If you would like to use the slide type setting to track different slide types, follow this procedure.

1. Purchase additional slide cartridges for each of the slide types you would like to track. Extra slide cartridges may be purchased from your printer supplier. You can track many different slide types. Below are a few examples:

- Charged slides
- Colored slides
- Various manufacturers of slides
- Slides with different coatings and treatments for specific procedures

2. Label the slide cartridges. Once you have received the extra slide cartridges, label them by placing a label or writing on the front or side of the blue cartridge. The label should be permanent so it will withstand daily handling. The label should specify the slide type number and a description of the slide type.
3. To encode the slide type number to the IC chip imbedded in the Slide Cartridge, insert the cartridge into the printer.
4. Open the settings area by clicking the tool icon in the upper right corner of the software.
5. The settings window will appear. Click the settings ribbon to view all settings.
6. Choose the slide type number you would like to encode (1-100). Click Set. The cartridge will be encoded.

Important Note!

If you switch between different slide cartridges encoded with different slide type numbers and you have the auto-feed mode enabled, the slide of a different type left inside the printer will be automatically ejected.

6.7 Preload slides in Auto-feed Mode

The auto-feed mode automatically feeds a slide from the cartridge into printing position before a print is sent to the printer. This greatly reduces the time it takes for the first slide in a batch of slides to print. This feature, if enabled, is especially useful if you only print one slide at a time.

Important Note!


1. If you switch between different slide types and you have the auto-feed mode enabled, a slide of a different type could be left inside the printer after you remove a slide cartridge and insert another one of a different type.

Press the Slide button  on the printer to eject the slide!

2. If the cartridges containing different slide types have been encoded with different slide type numbers, the slide left inside the printer will be automatically ejected (see section 7.6).

Enable or Disable Auto-Feed Mode.

1. Open the settings area by clicking the tool icon in the upper right corner of the software.
2. The settings window will appear. Click the settings ribbon to view all settings.

Tick box "Move all LIS jobs to the pending queue" is not supported with the current software version.

3. Check the Enable auto-feed mode box.

6.8 Advanced Settings

Press Shift +A on your keyboard to enable the advanced menu.

In this area you can set your print head resistance value, adjust print offsets, restore factory defaults or adjust heat offsets. Contact technical support before adjusting any of these settings.

Advanced Settings

Print Offsets

Horizontal Offset

Vertical Offset

Print Head Settings

Head Resistance

Type 2 Heat Offsets

Color

Black

Section 7: Multiple Slide Printers


7.1 Setup Multiple Printers

You can connect multiple Slide Printers to the same computer and control them all through PTLab. This is useful if you have two different slide types that you regularly print to and you do not want to continuously change slide cartridges.


To set this up click the gear icon above the picture of the currently connected slide printer.

The slide printer currently connected will be displayed only with the USB Serial. This is the number that Windows uses to differentiate between printers. If two different printers have the same USB serial number windows will recognize them as the same printer. If they both connected at the same time with the same serial number you may experience a Windows Blue Screen crash.

1. To get started click “Add Printer”.
2. You will be prompted to Unplug all slide printers from the system. Click OK when that has been completed.
3. Now you will be prompted to Plugin the new slide printer. Plug it in and wait 5 seconds. Now click “OK”.
4. The new printer will appear by default as Slide Printer (Copy 1). If successful it should have a USB Serial of 002.

Note: If the slide printer says “Not Connected” click the refresh button.

5. Connect the original Slide Printer. Click “Refresh”. You will now see both connected.

6. If you choose, you can now rename the slide printers. Select the one you want to rename and click “Rename Printer”. You will be presented with a window that allows you to do so. Consider naming the printers after the slide type which will be installed. You may want to further label the cartridge and the printer itself to match this new name.

7. Close this window. You will see both printers connected. Click the Refresh button if ever you have need to disconnect one printer and reconnect it while not restarting PTLab.

7.2 Printing to Multiple Slide Printers

To print to multiple slide printers simply click on the one you wish to print BEFORE clicking the “Print” button.

You will see the current number of prints remaining, slides remaining, printer name, current status and possibly the slide type change when you click on each different slide printer.

Note: Because both printers share the same PTLab printer queue, only one printer can be printing slides at any given time.

If you are in need of technical support, please contact your Leica representative or the Leica dealer where you purchased the instrument.

Please provide the following information:

- Model name and serial number of the instrument.
- Location of the instrument and name of a contact person.
- Reason for the service call.
- The delivery date of the instrument.