Leica Biosystems Bluing Agent products.

English.....	2
العربية (Arabic).....	5
简体中文 (Chinese Simplified).....	8
繁體中文 (Chinese Traditional).....	11
Dansk (Danish).....	14
Nederlands (Dutch).....	17
Français (French – Canada).....	20
Français (French – France).....	23
Deutsch (German).....	26
Italiano (Italian).....	29
日本語 (Japanese).....	32
한국어 (Korean).....	35
Norsk (Norwegian).....	38
Polski (Polish).....	41
Português (Portuguese – Brazil).....	44
Português (Portuguese – Portugal).....	47
Română (Romanian).....	50
Русский (Russian).....	53
Slovenščina (Slovenian).....	56
Español (Spanish – Central America).....	59
Español (Spanish – Spain).....	62
Svenska (Swedish).....	65
ภาษาไทย (Thai).....	68
Türkçe (Turkish).....	71
Tiếng Anh (Vietnamese).....	74

Bluing Agents

REF 3802900, 3802901, 3802901E, 3802915, 3802916, 3802918

Product Name

Leica Biosystems Bluing Agent products.

Intended Use

Detection/Measurement

The Leica Biosystems Bluing Agents do not detect or measure an analyte or marker. The Leica Biosystems Bluing Agents are used as a component of a hematoxylin and eosin (H&E) staining protocol. When used as recommended, the Bluing Agents convert the coloration of chromatin bound hematein-aluminum complexes from a dull red/purple hue to a bright purple/blue coloration.

Product Function

The results obtained through use of the Leica Biosystems Bluing Agents do not provide objective medical evidence. The Leica Biosystems Bluing Agents are used in a H&E staining protocol to enhance the contrast of hematoxylin stained nuclei with other components of cells and tissues. The H&E stained specimen, when interpreted by a trained professional, is utilized alongside other information such as the patient's medical history, physical condition, as well as results from other medical testing to render a medical diagnosis.

Specific Information Provided

Leica Biosystems Bluing Agents are not intended for the detection, definition or differentiation of a specific disorder, condition or risk factor. The staining demonstrated with use of these products, when used as intended, provides trained professionals information which may define the physiological or pathological state of the tissue specimen.

Automation

The Leica Biosystems Bluing Agents are not automated but can be used on automated staining platforms. Use on an automated platform should be validated at the point of use.

Qualitative/Quantitative

The Leica Biosystems Bluing Agents are used with qualitative stains.

Specimen Type

The Leica Biosystems Bluing Agents may be used with fixed or fresh histologic and cytologic specimens.

Testing Population

The Leica Biosystems Bluing Agents are intended for use with any patient requiring evaluation of biopsy or resection tissue as for the assessment of a suspected pathology or disease.

Intended User

The Leica Biosystems Bluing Agents are intended for use by qualified laboratory personnel and/or designee of the laboratory.

In Vitro Diagnostic

The Leica Biosystems Bluing Agents are intended for *in vitro* diagnostics use only.

Intended User

The Leica Biosystems Bluing Agent are intended for use by qualified laboratory personnel and/or designee.

Test Principle

The Leica Biosystems Bluing Agents work by converting the coloration of hematoxylin from a reddish to bluish hue in hematoxylin stains specimens. The bluing agents shall have an alkaline pH to promote a coloration change of hematoxylin in stained tissue specimens.

Calibrators & Controls

The Leica Biosystems Bluing Agents do not require the use of any calibrators or controls.

Reagent Limitation

No reagent limitations are applicable to these products.

Applicable Products

Product Code	Material Description
3802900	Surgipath Scott's Tap Water Concentrate (1pint)
3802901	Surgipath Scott's Tap Water Working Solution (1gal)
3802901E	Surgipath Scott's Tap Water Working Solution (5L)
3802915	SelecTech Blue Buffer 8 (500mL)
3802916	SelecTech Blue Buffer 8 (4-500mL)
3802918	SelecTech Blue Buffer 8 Ready-To-Use (1gal)

NOTE: Products listed here may not be available in all regions.

Materials Not Included

The Leica Biosystems Bluing Agents are designed to be used as part of a Hematoxylin & Eosin (H&E) stain protocol which require the use of graded alcohols, xylene or xylene substitutes, hematoxylin, differentiator, and eosin.

Bluing Agents

REF 3802900, 3802901, 3802901E, 3802915, 3802916, 3802918

Devices Required

Leica Biosystems Bluing Agents may be used on any open automated staining platform or with a manual staining method and should be validated at the point of use.

Storage and Stability

The product shall be stable for 24 months postproduction when stored at ambient temperature.

Store reagents at room temperature (15-30°C) in a well-ventilated place.

CAUTION: Do not use after the expiration date.

In Use Stability

User discretion should be utilized when determining in-use stability.

Sterility

The Leica Biosystems bluing agents are not sterile products.

Warnings/Precautions

This product and protocol(s) associated with the product, whether provided by Leica Biosystems in this instruction for use or developed by the user, shall be validated at the point of use by the user.

Infectious Material Status

The Leica Biosystems bluing agents do not include any infectious material. However, specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

Special Facilities

The Leica Biosystems bluing agents should be used per facility guidelines.

Specimen Handling

Specimens intended to be stained with H&E that includes Bluing Agent should be well fixed with neutral buffered formalin.

Following processing and paraffin embedding, section tissue at a standard thickness of (2 – 5µm).

Preparation for Use

- **Blue Buffer 8 RTU and Scott's Tap Water Working Solution** are ready-to-use formulas and therefore do not need to be altered or diluted in any way.
- **Blue Buffer Concentrate (20X)** must be diluted prior to use. Use one (1) part of concentrate, add nineteen (19) parts of deionized or distilled water to achieve working solution.
- **Scott's Tap Water Substitute Concentrate** must be diluted prior to use. Use one (1) part of concentrate to nine (9) parts of water.

Example: add 378.5 mL of concentrate to an empty gallon container and add tap water until fluid volume reaches one (1) gallon (3785 mL). Caution: Open Scott's Tap Water container slowly to release built up pressure.

Note - Inadequate washing after immersion in working solution of Blue Buffer 8 or Scott's Tap Water may affect desired staining intensity.

Protocol Set-up:

A typical H&E protocol (Table 1) should include water washes before and after immersion in a bluing agent.

Recommended immersion time in a bluing agent is no more than one minute. Water wash that follows bluing step removes any residual bluing solution and limits carry-over of alkaline solution into eosin.

Table 1. Example of H&E Staining Protocol

Steps	Action	Chemical	Time (mm:ss)
1	Deparaffinize	Xylene	3:00
2	Deparaffinize	Xylene	3:00
3	Deparaffinize	Xylene	3:00
4	Hydration	100% Alcohol	2:00
5	Hydration	100% Alcohol	1:00
6	Hydration	100% Alcohol	1:00
7	Hydration	80% or 95% Alcohol	1:00
8	Hydration	Water Wash	1:00
9	Stain	Progressive Hematoxylin	1:00 to 5:00
10	Wash	Water Wash	3:00
11	Differentiation	Differentiator	0:30 to 1:30

Bluing Agents

REF 3802900, 3802901, 3802901E, 3802915, 3802916, 3802918

12	Wash	Water Wash	1:00
13	Bluing	Bluing Buffer	0:30 to 1:00
14	Wash	Water Wash	2:00
15	Dehydration	80% to 95% Alcohol	1:00
16	Counterstaining	Alcoholic Eosin	0:30 to 1:30
17	Wash	Water Wash	2:00
18	Dehydration	95% to 100% Alcohol	1:00
19	Dehydration	100% Alcohol	1:00
20	Dehydration	100% Alcohol	1:00
21	Clearing	Xylene	2:00
22	Clearing	Xylene	2:00
23	Clearing	Xylene	2:00

Readiness for Use

Once concentrated version of a bluing agent is properly diluted, or if ready-to-use formula is used, pour all the reagent into the reagent vessel. Place the reagent vessel back into the respective station.

Quality Control

A routine quality control slide(s) containing tissue fixed and processed in a similar manner to the test specimens should be performed prior to routine use to ensure reagents are performing as intended.

Expected Results

By following the instruction for use, cell nuclei will show bright purple/blue coloration.

Analytical Performance

The Leica Biosystems Bluing Agents are not used to detect a specific analyte or marker. These products are used in conjunction with other products in a Hematoxylin & Eosin staining protocol system to stain cell nuclei blue and connective tissue, cytoplasm, muscle and erythrocytes various shades of orange, pink and red. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

Clinical Performance

The Leica Biosystems Bluing Agents are not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the Leica Biosystems Bluing Agents in a clinical setting.


Disposal

Bluing agents should be disposed in accordance with local governing regulations.

Leica Biosystems Richmond, Inc.
5205 Route 12
Richmond, IL 60071
USA
(1-844-534-2262)

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