Microtome Classic Blades

English.....	2
العربية (Arabic).....	5
简体中文 (Chinese Simplified).....	8
中國傳統的 (Chinese Traditional).....	11
Dansk (Danish).....	14
Nederlands (Dutch).....	17
Français (French – Canada).....	21
Français (French – France).....	25
Deutsch (German).....	29
Italiano (Italian).....	33
日本語 (Japanese).....	36
한국어 (Korean).....	39
Norsk (Norwegian).....	42
Polski (Polish).....	45
Português (Portuguese – Brazil).....	49
Português (Portuguese – Portugal).....	53
Română (Romanian).....	57
Русский (Russian).....	60
Slovenski (Slovenian).....	64
Español (Spanish – Central America).....	67
Español (Spanish – Spain).....	71
Svensk (Swedish).....	74
ไทย (Thai).....	77
Türk (Turkish).....	81
Tiếng Việt (Vietnamese).....	84

Microtome Classic Blades

REF 3802105, 3802106, 3802112, 3802118, 3802122, 3802110, 3802115, 3802123

Product Name

Leica Biosystems Microtome Classic Blades.

Intended Use

Detection/Measurement

The Leica Biosystems Microtome Classic Blades do not detect or measure an analyte or marker. These products are used for the sectioning of embedded tissues or cellular material in the microtomy process.

Product Function

The Leica Classic Blades function as aids to diagnosis. The Leica Classic Blades are designed for use in a laboratory environment for the sectioning of embedded tissue or cellular specimens in the microtomy. Following sectioning, the specimens may be evaluated by a range of diagnostic techniques including but not limited to general histology stains, immunohistochemistry and *in situ* hybridization. The specimens when assessed by a trained pathologist may provide valuable data that is useful in the assessment of various physiological or pathological states. This data alongside other information such as the patients' medical history, physical condition, as well as results from other medical testing is all considered when rendering a medical diagnosis.

Specific Information Provided

The Leica Classic Blades are not intended for the detection, definition or differentiation of a specific disorder, condition or risk factor. The Leica Classic Blades are intended for the sectioning of the embedded tissue and cellular specimens for the subsequent application of specific diagnostic techniques that when evaluated by a trained pathologist may provide valuable data that is useful in the assessment of various physiological or pathological states. The pathological states may include but are not limited to a multitude of disease processes such as cancer, dysplasia, metaplasia, hypertrophy, infection, inflammation, ischemia and scarring. This data alongside other information such as the patients' medical history, physical condition, as well as results from other medical testing is all considered when rendering a medical diagnosis.

Automation

The Leica Classic Blades are not automated but can be used in conjunction with automated equipment and procedures.

Quality/Quantity

The Leica Classic Blades are not used for the detection of a specific analyte or marker.

Specimen Type

The Leica Classic Blades are intended use for sectioning embedded tissue or cellular specimens that have been fixed and processed using the appropriate histologic techniques.

Testing Population

The Leica Classic Blades do not target a specific patient group. The Leica Classic Blades are intended for use with any patient requiring histopathologic evaluation of biopsy or resection tissue for the assessment of a suspected pathology or disease.

Intended User

The Leica Classic Blades are intended for use by qualified laboratory personnel and/or designee of the laboratory.

In Vitro Diagnostic

The Leica Biosystems Microtome Classic Blades are intended for *in vitro* diagnostics use only.

Test Principle

The Leica Biosystems Microtome Classic Blades function through using sharply ground, multi facet edge of metal to facilitate microns thick sectioning of material.

Calibrators & Controls

The Leica Classic Blades do not require the use of any calibrators or controls.

Reagent Limitations

No reagent limitations are applicable to these products.

Applicable Products

Product Code	Material Description
3802105	Classic Series Low-Profile, Triple Facet Microtome Blades, 1 pack
3802106	Classic Series Low-Profile, Triple Facet Microtome Blades, 1 case (10 packs)
3802112	Classic Series Low-Profile, Coated Microtome Blades, 1 case (10 packs)
3802118	Classic Series Low-Profile, Coated Microtome Blades, 1 pack
3802122	Classic Series Low-Profile, Thin, Triple Facet Microtome Blades, 1 case (10 packs)
3802110	Classic Series High-Profile, Double-Facet, Coated Microtome Blades, 1 case (10 packs)
3802115	Classic Series High-Profile, Double-Facet, Coated Microtome Blades, 1 pack
3802123	Classic Series High-Profile, Double-Facet, Thin Microtome Blades, 1 pack

Microtome Classic Blades

REF 3802105, 3802106, 3802112, 3802118, 3802122, 3802110, 3802115, 3802123

NOTE: Products listed here may not be available in all geographic regions.

Product Code	3802105 3802106	3802112 3802118	3802122	3802110 3802115	3802123
Triple Facet	x	x	x		
Double Facet				x	x
Blade angle	35 degrees			32 degrees	
Chromium coating	x	x	x	x	x
PTFE (Teflon)	x	x	x	x	
Dimensions	77.00 x 8.00 x 0.30 mm			77.00 x 14.00 x 0.30 mm	
Material	High-alloyed stainless steel				

Materials Not Included

The Leica Biosystems Microtome Classic Blades have been designed for use on microtomes.

Devices Required

Leica Biosystems Classic Blades can be used on any automated sectioning instrument.

Storage and Stability

Transport temperature: - 10 to + 50 °C.

Storage humidity: < 80 % non-condensing.

In Use Stability

User discretion should be utilized when determining in-use stability.

Sterility

The Leica Biosystems Microtome Classic Blades are not sterile products.

Warnings/Precautions

- Risk of injury. Blades are very sharp. Only trained laboratory personnel may use them.
- Never touch the blade edge!
- Only handle the blades by holding their sides between thumb and index finger crosswise to the blade edge.
- Wear cut-resistant safety gloves and protective eye wear when handling the blades.
- Do not use dispenser for disposal of the blades after cutting infectious or radioactive material.

Infectious Material Status

No infectious materials are applicable to these products.

Special Facilities

The Leica Biosystems Microtome Classic Blades should be used per facility guidelines.

Specimen Handling

Suggested fixatives include 10% neutral buffered formalin. Poor fixation, processing and rehydration will adversely affect the microtomy sectioning of paraffin embedded tissue. Tissue sections of 2 to 5-micron thickness are recommended.

Preparation for Use

- Inspect the dispenser upon receipt. If there is any damage, do not use and contact the manufacturer.
- If the dispenser does not eject blade, do not force the blade out of the dispenser. Do not use and contact the manufacturer.
- If the blade is not inserted directly into the blade holder, the blade dispenser must be held in the right hand to remove a blade.
- To push out the blade, slowly push the blade slide all the way to the left until it stops.
- The blade must then be held by the left hand and removed from the blade dispenser.

Direction for Use

- Insert the blade into the knife holder.
- If, after checking the additional instrument-specific, cut-relevant parameters, the cut result is not satisfactory, the blade must be moved or replaced.
- When cutting infectious or radioactive material, change the blade with every new specimen. This prevents unwanted contamination/transmission.
- Remove the blade from the knife holder.

Microtome Classic Blades

REF 3802105, 3802106, 3802112, 3802118, 3802122, 3802110, 3802115, 3802123

- Observe the following disposal instructions.

Readiness for Use

Once the blade is clamped securely in the blade holder of the microtome, sectioning may begin.

Quality Control

The need for quality controls may be determined at the point of use by the end user.

Expected Results

All blades are dispensed from the cartridge safely without damage to the cutting surface of the blade. Generated sections are free of defects such as knife lines.

Analytical Performance

The Leica Classic Blades are not used to detect a specific analyte or marker. The Leica Classic Blades are designed for use in a laboratory environment for the sectioning of embedded tissue or cellular specimens. Following sectioning, the specimens may be evaluated by a range of diagnostic techniques including but not limited to general histology stains, immunohistochemistry and *in situ* hybridization. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias) precision (repeatability and reproducibility) accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of these products.

A total of 22 dispensers were tested for the Leica Classic Blades. Of the 22 dispensers tested, 14 were low profile microtome blades and 8 high profile microtome blades. All 50 blades in each dispenser was ejected from the cartridge. Results demonstrated that all blades dispensed from the cartridge on the first try. After ejection, a total of 15 randomly selected blades from each cartridge were used for microtomy on blank paraffin blocks. Blank blocks were utilized as to not introduce any artificial artifacts or defects caused by the tissue and not reflective of the blade quality itself. Paraffin sections at a setting of 4µm were generated and evaluated for defects such as knife lines. All sections generated were free of defects.

These results support that the cartridge holds the blades securely and dispenses the blades in such a manner as to not compromise the quality of the ground honed surface necessary to achieve microtomy with the blades.

Clinical Performance

To evaluate clinical performance, testing occurred in a simulated laboratory environment. Testing utilized tissue specimens that were appropriately fixed, adequately hydrated and infiltrated and embedded in paraffin. A total of 12 blades, randomly selected from two lots were used for testing. The tissue blocks were faced and chilled for a minimum of 15 minutes prior to sectioning. Once chilled appropriately, the blocks were sectioned at 4µm and a ribbon of 4-8 sections was generated. The ribbon was inspected to ensure no defects were present and the results were documented. The testing supports the achievement of the intended use of the product and identifies that the product can be safely used in a clinical environment.



Disposal

After being used to cut infectious or radioactive material, the microtome blade should be disposed of in accordance with local governing regulations. Blades that have not come in contact with infectious or radioactive material should be disposed of in accordance with local governing regulations.

Leica Biosystems Richmond, Inc.
5205 Route 12
Richmond, IL 60071
USA
(1-844-534-2262)

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