## Gomori's Trichrome Stain Kit, Green Collagen Stain

**REF** 38016SS1

## Gomori's Trichrome Stain Kit, Blue Collagen Stain

**REF** 38016SS2

English.....	2
العربية (Arabic).....	6
简体中文 (Chinese Simplified).....	10
中國傳統的 (Chinese Traditional).....	14
Dansk (Danish).....	18
Nederlands (Dutch).....	22
Français (French – Canada).....	26
Français (French – France).....	30
Deutsch (German).....	34
Italiano (Italian).....	38
日本語 (Japanese).....	42
한국어 (Korean).....	46
Norsk (Norwegian).....	50
Polski (Polish).....	54
Português (Portuguese – Brazil).....	58
Português (Portuguese – Portugal).....	62
Română (Romanian).....	66
Русский (Russian).....	70
Slovenski (Slovenian).....	75
Español (Spanish – Central America).....	79
Español (Spanish – Spain).....	83
Svensk (Swedish).....	87
ไทย (Thai).....	91
Türk (Turkish).....	95
Tiếng Việt (Vietnamese).....	99

## Gomori's Trichrome Stain Kit, Green Collagen Stain

**REF** 38016SS1

## Gomori's Trichrome Stain Kit, Blue Collagen Stain

**REF** 38016SS2

### Product Name

Gomori's Trichrome Stain Kit, Green Collagen Stain.  
Gomori's Trichrome Stain Kit, Blue Collagen Stain.

### Intended Use

#### Detection/Measurement

Gomori's Trichrome Stain Kits do not detect or measure an analyte or marker.

The reagents of this kit, when used in histological protocols, are valuable for the differentiation of muscle and cytoplasm from collagen containing connective tissue.

#### Product Function

The results obtained through use of Gomori's Trichrome Stain Kits do not provide objective medical evidence. The coloration and contrast the Leica Biosystems Stain Kit provide to histologic specimens allows visualization of microscopic anatomy. This visualization, when interpreted by a trained professional, is utilized alongside other information such as the patient's medical history, physical condition, as well as results from other medical testing to render a medical diagnosis.

#### Specific Information Provided

Gomori's Trichrome Stain Kits are not intended for the detection, definition or differentiation of a specific disorder, condition, or risk factor. The staining demonstrated with use of these products, when used as intended, provides trained professionals information which may define the physiological or pathological state of the tissue specimen.

#### Automation

Gomori's Trichrome Stain Kits are not automated but can be used on automated staining platforms. Use on an automated platform should be validated at the point of use.

#### Qualitative/Quantitative

The Leica Biosystems Gomori's Trichrome Stain Kits are qualitative stains.

#### Specimen Type

Gomori's Trichrome Stain Kits can be used with any paraffin embedded human or animal specimen.

#### Testing Population

The Leica Biosystems Gomori's Trichrome Stain Kits are intended for use with any patient requiring evaluation of biopsy or resection tissue for the assessment of a suspected pathology or disease.

#### Intended User

Gomori's Trichrome Stain Kits are intended for use by qualified laboratory personnel and/or designee of the laboratory.

### In Vitro Diagnostic

Gomori's Trichrome Stain Kits are intended for *in vitro* diagnostics use only.

### Test Principle

The exact mechanism(s) by which trichrome stains differentially react with and stain various tissue components are not fully understood. Sulfated azo dyes such as chromotrope 2R, selectively stain the cytoplasm of most cell types including muscle tissue. This selectivity is due in part to phosphotungstic acid. Phosphotungstic acid binds preferentially to collagen and by doing so likely prevents or disrupts binding of the sulfated azo dyes such as chromotrope 2R. The acid triarylmethane dyes such as aniline blue or fast green also bind collagen but are not displaced by the phosphotungstic acid.

### Calibrators & Controls

Gomori's Trichrome Stain Kits do not require the use of any calibrators or controls.

### Reagent Limitations

No reagent limitations are applicable to this product.

### Applicable Products

Product Code	Material Description
38016SS1	Gomori's Trichrome Stain Kit, Green Collagen Stain
38016SS1A	Bouin's Solution, 500 mL
38016SS1B	1% Acetic Acid, 500 mL
38016SS1C	Weigert's Hematoxylin Solution A, 250 mL
38016SS1D	Weigert's Hematoxylin Solution B, 250 mL
38016SS1E	Gomori's Trichrome Green Collagen Stain, 500 mL

## Gomori's Trichrome Stain Kit, Green Collagen Stain

**REF** 38016SS1

## Gomori's Trichrome Stain Kit, Blue Collagen Stain

**REF** 38016SS2

Product Code	Material Description
38016SS2	Gomori's Trichrome Stain Kit, Blue Collagen Stain
38016SS1A	Bouin's Solution, 500 mL
38016SS1B	1% Acetic Acid, 500 mL
38016SS1C	Weigert's Hematoxylin Solution A, 250 mL
38016SS1D	Weigert's Hematoxylin Solution B, 250 mL
38016SS2E	Gomori's Trichrome Blue Collagen Stain, 500 mL

### Materials Not Included

Gomori's Trichrome Blue and Green Stain Kits protocols require the use of graded alcohols, xylene, or xylene substitutes, deionized or distilled water. Positive control slide(s) with tissue containing muscle and collagen (not included in this kit), should be included in each run.

### Devices Required

Gomori's Trichrome Stain Kits can be used on any automated staining platform or with a manual staining method.

### Storage and Stability

Store at room temperature (15-30 °C).

**CAUTION:** Do not use after the expiration date.

### In Use Stability

User discretion should be utilized when determining in-use stability.

### Sterility

Gomori's Trichrome Stain Kits components are not sterile products.

### Warnings/Precautions

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations. Refer to Material Safety Data Sheet and product labeling for any updated risk, hazard, or safety information.

Information.

### Infectious Material Status

Gomori's Trichrome Stain Kits do not include any infectious material. However, specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

### Special Facilities

Gomori's Trichrome Stain Kits should be used per facility guidelines.

### Specimen Handling

- **Fixation**  
Any general fixative including but not limited to neutral buffered formalin, alcoholic formalin or Bouin's solution may be used.
- **Paraffin Sections**  
Following processing and paraffin embedding, section tissue specimens at 4–6 microns.

### Preparation for Use

- Mix equal parts of Weigert's Hematoxylin Solutions A and B.
- Weigert's working solution is useable for 10 days.

### Direction for Use

#### Conventional Staining Protocol

1. Deparaffinize with xylene and rehydrate through graded alcohols to deionized or distilled water.
2. Place in Bouin's Solution at 56 °C for 1 hour (A).
3. Wash the sections in running tap water for 5 minutes or until the yellow coloration is removed.
4. Stain in the Weigert's Hematoxylin working solution for 5–10 minutes.
5. Rinse in running tap water for 5–7 minutes.
6. Stain in Gomori's Trichrome Blue Solution OR Gomori's Trichrome Green Solution for 15–20 minutes.
7. Place sections in 1% Acetic Acid solution for 1 minute (B).

## Gomori's Trichrome Stain Kit, Green Collagen Stain

**REF** 38016SS1

## Gomori's Trichrome Stain Kit, Blue Collagen Stain

**REF** 38016SS2

8. Rinse briefly in deionized water.
9. Dehydrate through two changes of alcohol for one minute each.
10. Clear sections in two changes of xylene for one minute each.
11. Coverslip using a mounting medium miscible with xylene.

Notes: (A). Bouin's Solution treatment at step 2 is unnecessary if the specimens initially are fixed in Bouin's solution.  
(B). The 1% Acetic Acid rinse solution should be discarded after each use.

Table 1. Example of Gomori's Trichrome Green Collagen Stain and Blue Collagen Stain.

Steps	Action	Chemical	Time (mm: ss)
1-3	Deparaffinize	Xylene	3:00
4-5	Hydration	100% Alcohol	2:00
6	Hydration	80% or 95% Alcohol	1:00
7-9	Rinse	Deionized Water	0:30
10	Pretreatment	Bouin's Solution	60:00 @ 56° C
11	Wash	Water Wash	5:00
12	Stain	Weigert's Hematoxylin	5:00 – 10:00
17	Wash	Water Wash	5:00 – 7:00
18	Stain	Gomori's Trichrome Blue or Green Solution	15:00 – 20:00
19	Rinse	1% Acetic Acid	1:00
20	Rinse	Deionized Water	0:30
21	Dehydration	95% Alcohol	1:00
22-23	Dehydration	100% Alcohol	1:00
24-25	Clearing	Xylene	1:00

### Microwave Staining Protocol

Exercise caution when using the microwave to heat any solution or reagent. The microwave must be properly ventilated to prevent the accumulation of fumes in the laboratory. Microwave transparent coplin jars and caps should be used during the staining process. The caps should be loosely applied to prevent spills. Caps with ventilation holes also may be used. All microwaves should be used in accordance with the manufacturer's instructions.

All microwave steps were conducted at a power setting of 800 watts unless otherwise noted. Because of differences in microwave power and frequencies among various models, it may be necessary to adjust power levels or times to achieve optimal results.

1. Deparaffinize with xylene and rehydrate through graded alcohols to deionized or distilled water.
2. Place sections in a loosely capped plastic Coplin jar containing Bouin's Solution (A). Microwave for 30 seconds.
3. Gently mix the Bouin's Solution by swirling and allow to stand for 10 minutes.
4. Rinse in running tap water for 5 minutes.
5. Place sections in a plastic Coplin jar with Weigert's Hematoxylin working solution and microwave for 20 seconds.
6. Gently mix the Weigert's Hematoxylin working solution by swirling and allow to stand for 2 minutes.
7. Rinse in running tap water for 2 minutes.
8. Place sections in a plastic Coplin jar with Gomori's Trichrome Green OR Blue Solution and microwave for 20 seconds at 600 watts.
9. Gently mix the Gomori's Trichrome Green OR Blue Solution by swirling and allow to stand for 6 minutes.
10. Place sections in Acetic Acid solution for 1 minute
11. Rinse sections in deionized or distilled water for 30 seconds.
12. Dehydrate sections in two changes of alcohol for one minute each.
13. Clear sections in two changes of xylene for one minute each.
14. Coverslip using a mounting medium that is miscible with xylene.

## Gomori's Trichrome Stain Kit, Green Collagen Stain

**REF** 38016SS1

## Gomori's Trichrome Stain Kit, Blue Collagen Stain

**REF** 38016SS2

**Note:** (A). Bouin's Solution treatment at step 2 is unnecessary if the specimens initially are fixed in Bouin's solution.  
(B). The 1% Acetic Acid rinse solution should be discarded after each use

### Readiness for Use

Once appropriate staining protocol is chosen and bath-layout is created, pour all the reagent into the reagent vessel. Place the reagent vessel back into the respective station.

### Quality Control

Tissue sections containing muscle fibers and collagen, fixed and processed in a similar manner to the test specimens, should be included in each staining assay to ensure Special Staining Kit is performing as intended.

A loss of staining intensity has been reported to occur in tissue sections stored for an extended period of time.

#### Expected Results for Collagen Blue Stain

- Muscle fiber and cytoplasm — Red
- Collagen — Blue
- Nuclei — Black

#### Expected Results for Collagen Green Stain

- Muscle fiber and cytoplasm — Red
- Collagen — Green
- Nuclei — Black

### Analytical Performance

The Leica Biosystems Gomori's Trichrome Stain Kits are not used to detect a specific analyte or marker. These products are used to provide coloration for the differentiation of muscle and cytoplasm from collagen containing connective tissue. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

### Clinical Performance

The Leica Biosystems Gomori's Trichrome Stain kits are not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the Leica Biosystems Bluing Agents in a clinical setting.


### Disposal

Spent or expired components of the Gomori's Trichrome Stain Kits should be discarded in accordance with organizational, local, state, and federal regulations.

Leica Biosystems Richmond, Inc.  
5205 Route 12  
Richmond, IL 60071  
USA  
(1-844-534-2262)

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