Alcian Blue/Periodic Acid Schiff Special Stain Kit

English	2
العربية (Arabic).....	6
简体中文 (Chinese Simplified).....	9
中國傳統的 (Chinese Traditional)	13
Dansk (Danish)	17
Nederlands (Dutch)	21
Français (French – Canada).....	25
Français (French – France)	29
Deutsch (German).....	33
Italiano (Italian).....	37
日本語 (Japanese)	41
한국어 (Korean)	44
Norsk (Norwegian)	48
Polski (Polish)	52
Português (Portuguese – Brazil).....	56
Português (Portuguese – Portugal)	60
Română (Romanian)	64
Русский (Russian).....	68
Slovenski (Slovenian).....	72
Español (Spanish – Central America).....	76
Español (Spanish – Spain).....	80
Svensk (Swedish)	84
ไทย (Thai).....	88
Türk (Turkish).....	92
Tiếng Việt (Vietnamese).....	96

Alcian Blue/Periodic Acid Schiff Special Stain Kit

REF 38016SS5

Product Name

Alcian Blue (AB)/Periodic Acid Schiff (PAS) Special Stain Kit

Intended Use

Detection/Measurement

The Leica Biosystems AB/PAS Special Stain Kit does not detect or measure an analyte or marker.

AB/PAS Special Stain Kit is used to identify acid mucin, neutral mucin, and glycogen in routinely prepared paraffin embedded tissue specimens.

Product Function

The results obtained through use of AB/PAS Special Stain Kit do not provide objective medical evidence. The coloration and contrast the Leica Biosystems AB/PAS Stain Kit provide to histologic specimens allows visualization of microscopic anatomy. This visualization, when interpreted by a trained professional, is utilized alongside other information such as the patient's medical history, physical condition, as well as results from other medical testing to render a medical diagnosis.

Specific Information Provided

The Leica Biosystems AB/PAS Special Stain Kit is not intended for the detection, definition or differentiation of a specific disorder, condition, or risk factor. The staining demonstrated with use of these products, when used as intended, provides trained professionals information which may define the physiological or pathological state of the tissue specimen.

Automation

AB/PAS Special Stain Kit is not automated but can be used on automated staining platforms. Use on an automated platform should be validated at the point of use.

Qualitative/Quantitative

The Leica Biosystems AB/PAS Special Stain kit is qualitative stain.

Specimen Type

AB/PAS Special Stain Kit can be used with any paraffin embedded human or animal specimen.

Testing Population

The Leica Biosystems AB/PAS Special Stain kit is intended for use with any patient requiring evaluation of biopsy or resection tissue for the assessment of a suspected pathology or disease.

Intended User

AB/PAS Special Stain Kit is intended for use by qualified laboratory personnel and/or designee of the laboratory.

In Vitro Diagnostic

AB/PAS Special Stain Kit is intended for *in vitro* diagnostics use only.

Test Principle

The Alcian blue molecule is a large conjugated dye molecule that contains a centrally located copper atom bound to a phthalocyanine ring and four isothiuronium groups. The isothiuronium groups impart an overall positive charge to the molecule. The cationic isothiuronium groups likely bind via electrostatic interactions to the anionic sulfate and carboxylate groups located within the carbohydrate moieties of acidic mucins. Neutral mucins that contain carbons with adjacent hydroxyl (1,2-glycol) groups are oxidized by periodic acid. The resultant dialdehyde groups react with the Schiff reagent.

Calibrators & Controls

AB/PAS Special Stain Kit does not require the use of any calibrators or controls.

Reagent Limitations

No reagent limitations are applicable to this product.

Applicable Products

Product Code	Material Description
38016SS5	Alcian Blue (AB)/Periodic Acid Schiff (PAS) Special Stain Kit
38016SS3A	Alcian Blue pH 2.5, 500 ml
38016SS4A	Periodic Acid, 0.5%, 500 ml
38016SS4B	Schiff Reagent, 500 ml

Materials Not Included

Alcian Blue/Periodic Acid Schiff Special Stain Kit protocol requires the use of graded alcohols, xylene, or xylene substitutes, deionized or distilled water. Positive PAS control slide(s), not included in this kit, should be included in each run.

Alcian Blue/Periodic Acid Schiff Special Stain Kit

REF 38016SS5

Devices Required

Leica Biosystems Alcian Blue/Periodic Acid Schiff Special Stain Kit can be used on any automated staining platform or with a manual staining method.

Storage and Stability

Store all components except Schiff reagent at room temperature (15-30 °C). The Schiff reagent should be stored at 2-8 °C.
CAUTION: Do not use after the expiration date.

In Use Stability

User discretion should be utilized when determining in-use stability.

Sterility

AB/PAS Special Stain Kit components are not sterile products.

Warnings/Precautions

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations. Refer to Material Safety Data Sheet and product labeling for any updated risk, hazard, or safety information.

Infectious Material Status

AB/PAS Special Stain Kit does not include any infectious material. However, specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

Special Facilities

AB/PAS Special Stain Kit should be used per facility guidelines.

Specimen Handling

Suggested fixatives include 10% neutral buffered formalin. Routine dehydration, clearing, and paraffin infiltration and embedding, and routine preparation of microtome sections. Poor fixation, processing, rehydration, and sectioning will adversely affect the staining quality. Tissue sections of 2 to 5-micron thickness are recommended.

Preparation for Use

Any general fixative including but not limited to neutral buffered formalin, alcoholic formalin and Bouin's may be used. Fixatives containing dialdehyde (glutaraldehyde) should be avoided, as free aldehyde groups may bind the Schiff reagent and produce non-specific staining. Following processing and paraffin embedding, cut sections at 4-6 microns.

Direction for Use

Conventional Staining Protocol

1. Deparaffinize and hydrate sections to deionized or distilled water.
2. Stain in Alcian Blue solution for 30 minutes.
3. Rinse in running tap water for 5 minutes.
4. Place in Periodic Acid solution for 5 minutes at room temperature (18-26 °C).
5. Rinse in several changes of deionized water.
6. Place slides in Schiff Reagent for 15 minutes at room temperature (18-26 °C).
7. Wash in gently running lukewarm tap water for 10 minutes.
8. Dehydrate through two changes of 95% alcohol and absolute alcohol, two minutes each.
9. Clear in two changes of xylene and mount in a mounting media that is miscible with xylene.

Table 1. Example of Conventional AB/PAS Staining Protocol.

Steps	Action	Chemical	Time (mm: ss)
1-3	Deparaffinize	Xylene	3:00
4-5	Hydration	100% Alcohol	2:00
6	Hydration	80% or 95% Alcohol	1:00
7	Hydration	Deionized Water	1:00
8	Stain	Alcian Blue	30:00
9	Wash	Water	5:00
10	Oxidize	Periodic Acid	5:00
11	Wash	Deionized Water	several changes, 10 sec each

Alcian Blue/Periodic Acid Schiff Special Stain Kit

REF 38016SS5

12	Stain	Schiff Reagent	15:00
13	Wash	Lukewarm Water	10:00
14-15	Dehydration	95% Alcohol	2:00
16-17	Dehydration	100% Alcohol	2:00
18-19	Clearing	Xylene	2:00

Note: When using xylene substitute, increase immersion times by approximately 50%.

Microwave Staining Protocol

Exercise caution when using the microwave to heat any solution or reagent. The microwave must be properly ventilated to prevent the accumulation of fumes in the laboratory. Microwave transparent Coplin jars and caps should be used during the staining process. The caps should be loosely applied to prevent spills. Caps with ventilation holes also may be used. All microwaves should be used in accordance with the manufacturer's instructions. All microwaved steps were conducted at a power setting of 800 watts.

1. Deparaffinize with xylene or a xylene substitute and rehydrate through graded alcohols to deionized water.
2. Place slides in a plastic Coplin jar containing the Alcian Blue solution (40-50 ml).
3. Microwave for 30 seconds.
4. Gently mix the Alcian Blue solution by swirling and allow to stand for 1 minute.
5. Wash in running tap water for 10 minutes.
6. Rinse briefly in deionized water.
7. Add Periodic Acid (40-50 ml) to a plastic Coplin jar and microwave for 10 seconds.
8. Gently mix the solution by swirling and allow to stand for 1 minute.
9. Rinse in several changes of deionized water.
10. Add the Schiff Reagent (40-50 ml) to a plastic Coplin jar and microwave at for 15 seconds.
11. Gently mix the Schiff Reagent by swirling and allow to stand for 1 minute.
12. Wash in gently running lukewarm tap water for 5 minutes.
13. Dehydrate through two changes of 95% alcohol and 100% alcohol, two minutes each.
14. Clear in two changes of xylene (two minutes each) and mount with a mounting media that is miscible with xylene.

Readiness for Use

Once appropriate staining protocol is chosen and bath-layout is created, pour all the reagent into the reagent vessel. Place the reagent vessel back into the respective station.

Quality Control

A quality control slide(s) containing liver or gastrointestinal epithelium specimen (small intestine, appendix, colon), fixed and processed in a similar manner to the test specimens should be included in each staining assay to ensure AB/PAS Special Staining Kit is performing as intended.

Expected Results

- Acidic mucins (sulfated and carboxylated), proteoglycans, hyaluronic acid – Blue
- Glycogen, neutral mucins, and various glycoproteins – Magenta/Red
- Mixed mucins (neutral and acidic) – Purple to Blue

Analytical Performance

The AB/ PAS Stain Kit is not used to detect a specific analyte or marker. These products used to identify acid mucin, neutral mucin, and glycogen in routinely prepared paraffin embedded tissue specimens. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

Clinical Performance

The AB/PAS Stain Kit is not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the Leica Biosystems Bluing Agents in a clinical setting.

Alcian Blue/Periodic Acid Schiff Special Stain Kit

REF 38016SS5

Disposal

Spent or expired components of the AB/PAS kit should be discarded in accordance with organizational, local, state, and federal regulations.

Leica Biosystems Richmond, Inc.
5205 Route 12
Richmond, IL 60071
USA
(1-844-534-2262)

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