Periodic Acid Schiff (PAS) Special Stain Kit

English	2
العربية (Arabic).....	6
简体中文 (Chinese Simplified).....	9
中國傳統的 (Chinese Traditional)	13
Dansk (Danish).....	17
Nederlands (Dutch)	21
Français (French – Canada).....	25
Français (French – France)	29
Deutsch (German)	33
Italiano (Italian).....	37
日本語 (Japanese).....	41
한국어 (Korean).....	44
Norsk (Norwegian)	48
Polski (Polish)	52
Português (Portuguese – Brazil)	56
Português (Portuguese – Portugal).....	60
Română (Romanian)	64
Русский (Russian).....	68
Slovenski (Slovenian).....	72
Español (Spanish – Central America).....	76
Español (Spanish – Spain)	80
Svensk (Swedish)	84
ไทย (Thai).....	88
Türk (Turkish).....	91
Tiếng Việt (Vietnamese)	95

Periodic Acid Schiff (PAS) Special Stain Kit

REF 38016SS4

Product Name

Periodic Acid Schiff (PAS) Special Stain Kit

Intended Use

Detection/Measurement

The Leica Biosystems PAS Special Stain Kit does not detect or measure an analyte or marker.

PAS Special Stain Kit is used for staining structures containing a high proportion of carbohydrates such as glycogen, glycoproteins, proteoglycans typically found in connective tissues, mucus, and basement membranes.

Product Function

The results obtained through use of PAS Special Stain Kit do not provide objective medical evidence. The coloration and contrast the Leica Biosystems PAS Stain Kit provide to histologic specimens allows visualization of microscopic anatomy. This visualization, when interpreted by a trained professional, is utilized alongside other information such as the patient's medical history, physical condition, as well as results from other medical testing to render a medical diagnosis.

Specific Information Provided

The Leica Biosystems PAS Special Stain Kit is not intended for the detection, definition or differentiation of a specific disorder, condition, or risk factor. The staining demonstrated with use of these products, when used as intended, provides trained professionals information which may define the physiological or pathological state of the tissue specimen.

Automation

PAS Special Stain Kit is not automated but can be used on automated staining platforms. Use on an automated platform should be validated at the point of use.

Qualitative/Quantitative

The Leica Biosystems PAS Special Stain kit is qualitative stain.

Specimen Type

PAS Special Stain Kit can be used with any paraffin embedded human or animal specimen.

Testing Population

The Leica Biosystems PAS Special Stain kit is intended for use with any patient requiring evaluation of biopsy or resection tissue for the assessment of a suspected pathology or disease.

Intended User

PAS Special Stain Kit is intended for use by qualified laboratory personnel and/or designee of the laboratory.

In Vitro Diagnostic

Periodic Acid Schiff (PAS) Special Stain Kit is intended for *in vitro* diagnostics use only.

Test Principle

This stain is used for the demonstration of glycogen. Tissue sections are first oxidized by periodic acid. The oxidative process results in the formation of aldehyde groupings through carbon-to-carbon bond cleavage. Free hydroxyl groups should be present for oxidation to take place. Oxidation is completed when it reaches the aldehyde stage. The aldehyde groups are detected by the Schiff reagent. A colorless, unstable dialdehyde compound is formed and then transformed to the colored final product by restoration of the quinoid chromophoric grouping which gives the magenta color.

Calibrators & Controls

PAS Special Stain Kit does not require the use of any calibrators or controls.

Reagent Limitations

No reagent limitations are applicable to this product.

Applicable Products

Product Code	Material Description
38016SS4	Periodic Acid Schiff (PAS) Special Stain Kit
38016SS4A	Periodic Acid, 0.5%, 500 ml
38016SS4B	Schiff Reagent, 500 ml
38016SS4C	Gill II Hematoxylin, 500 ml

Materials Not Included

Periodic Acid Schiff (PAS) Special Stain Kit protocol requires the use of graded alcohols, xylene, or xylene substitutes, deionized or distilled water. Positive PAS control slide(s), not included in this kit, should be included in each run.

Periodic Acid Schiff (PAS) Special Stain Kit

REF 38016SS4

Devices Required

Leica Biosystems Periodic Acid Schiff (PAS) Special Stain Kit can be used on any automated staining platform or with a manual staining method.

Storage and Stability

Periodic acid and Gill II Hematoxylin may be stored at room temperature. The Schiff Reagent should be stored at 2-8 °C. CAUTION: Do not use after the expiration date.

In Use Stability

User discretion should be utilized when determining in-use stability.

Sterility

Periodic Acid Schiff (PAS) Special Stain Kit components are not sterile products.

Warnings/Precautions

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations. Refer to Material Safety Data Sheet and product labeling for any updated risk, hazard, or safety information.

Infectious Material Status

Periodic Acid Schiff (PAS) Special Stain Kit does not include any infectious material. However, specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

Special Facilities

Periodic Acid Schiff (PAS) Special Stain Kit should be used per facility guidelines.

Specimen Handling

Suggested fixatives include 10% neutral buffered formalin. Routine dehydration, clearing, and paraffin infiltration and embedding, and routine preparation of microtome sections. Poor fixation, processing, rehydration, and sectioning will adversely affect the staining quality. Tissue sections of 2 to 5-micron thickness are recommended.

Preparation for Use

Any general fixative including but not limited to neutral buffered formalin, alcoholic formalin and Bouin's solution may be used. Fixatives containing dialdehyde (glutaraldehyde) should be avoided, as free aldehyde groups may bind the Schiff reagent and produce non-specific staining. Following processing and paraffin embedding, cut sections at 4-6 microns.

Direction for Use

Conventional Staining Protocol

1. Deparaffinize and hydrate sections to deionized water.
2. Place in Periodic Acid solution for 5 minutes at room temperature (18-26° C).
3. Rinse in several changes of deionized water.
4. Place slides in Schiff Reagent for 15 minutes at room temperature (18-26° C).
5. Wash in gently running lukewarm tap water for 10 minutes.
6. Counterstain with Gill II Hematoxylin for 3-4 minutes.
7. Rinse in running tap water for 5 minutes.
8. Dehydrate through two changes of 95% alcohol and absolute alcohol two minutes each.
9. Clear in two changes of xylene and mount in a mounting media that is miscible with xylene.

Table 1. Example of Conventional PAS Staining Protocol.

Steps	Action	Chemical	Time (mm: ss)
1-3	Deparaffinize	Xylene	3:00
4-5	Hydration	100% Alcohol	2:00
6	Hydration	80% or 95% Alcohol	1:00
7	Hydration	Deionized Water	1:00
8	Stain	Periodic Acid	5:00
9	Wash	Deionized Water	several changes, 10 sec each
10	Stain	Schiff Reagent	15:00

Periodic Acid Schiff (PAS) Special Stain Kit

REF 38016SS4

11	Wash	Lukewarm Water	10:00
12	Counterstain	Gill II Hematoxylin	3:00 to 4:00
13	Wash	Water Wash	5:00
14-15	Dehydration	95% Alcohol	2:00
16-17	Dehydration	100% Alcohol	2:00
18-19	Clearing	Xylene	2:00

Note: When using xylene substitute, increase immersion times by approximately 50%.

Microwave Staining Protocol

Exercise caution when using the microwave to heat any solution or reagent. The microwave must be properly ventilated to prevent the accumulation of fumes in the laboratory. Microwave transparent Coplin jars and caps should be used during the staining process. The caps should be loosely applied to prevent spills. Caps with ventilation holes also may be used. All microwaves should be used in accordance with the manufacturer's instructions.

1. Deparaffinize with xylene or a xylene substitute and rehydrate through graded alcohols to deionized water.
2. Place sections in a plastic Coplin jar containing Periodic Acid (40-50 ml) and microwave at 800 watts for 10 seconds.
3. Gently mix the solution by swirling and allow to stand for 1 minute.
4. Rinse in several changes of deionized water.
5. Add the Schiff Reagent (40-50 ml) to a plastic Coplin jar and microwave at 800 watts for 15 seconds.
6. Gently mix the Schiff Reagent by swirling and allow to stand for 1 minute.
7. Wash in gently running lukewarm tap water for 5 minutes.
8. Counterstain with Gill II Hematoxylin for 3-4 minutes.
9. Rinse in running tap water for 5 minutes.
10. Dehydrate through two changes of 95% alcohol and 100% alcohol, two minutes each.
11. Clear in two changes of xylene (two minutes each) and mount with a mounting media that is miscible with xylene.

Readiness for Use

Once appropriate staining protocol is chosen and bath-layout is created, pour all the reagent into the reagent vessel. Place the reagent vessel back into the respective station.

Quality Control

A quality control slide(s) containing liver or gastrointestinal epithelium specimen (small intestine, appendix, colon), fixed and processed in a similar manner to the test specimens should be included in each staining assay to ensure PAS Special Staining Kit is performing as intended.

Expected Results

Many tissues and cell types will demonstrate a bright magenta coloration indicative of a positive PAS reaction. This includes, but is not limited to: glycogen containing hepatocytes, mucin in the gastrointestinal tract, as well as basement membranes in various sites. Nuclei should stain purple/blue with hematoxylin.

Analytical Performance

The PAS Stain kit is not used to detect a specific analyte or marker. This product is used for staining structures containing a high proportion of carbohydrates such as glycogen, glycoproteins, proteoglycans typically found in connective tissues, mucus, and basement membranes. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

Clinical Performance

The PAS Stain Kit is not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the Leica Biosystems Bluing Agents in a clinical setting.

Disposal


Spent or expired components of the PAS kit should be discarded in accordance with organizational, local, state, and federal regulations.

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Leica Biosystems Richmond, Inc.