Biopsy Handling Devices

English.....	2
العربية (Arabic).....	5
简体中文 (Chinese Simplified).....	8
中國傳統的 (Chinese Traditional).....	11
Dansk (Danish).....	14
Nederlands (Dutch).....	17
Français (French – Canada).....	20
Français (French – France).....	23
Deutsch (German).....	26
Italiano (Italian).....	29
日本語 (Japanese).....	32
한국어 (Korean).....	35
Norsk (Norwegian).....	38
Polski (Polish).....	41
Português (Portuguese – Brazil).....	44
Português (Portugal).....	47
Română (Romanian).....	50
Русский (Russian).....	53
Slovenski (Slovenian).....	56
Español (Spanish – Central America).....	59
Español (Spanish – Spain).....	62
Svensk (Swedish).....	65
ไทย (Thai).....	68
Türk (Turkish).....	71
Tiếng Việt (Vietnamese).....	74

Biopsy Handling Devices

REF 3801090, 3801091, 3801092
 3800003, 3801000, 3801020, 3801040, 3801060, 3801080
 3801085, 3801087

Product Name

Leica Biosystems Biopsy Handling Devices:

- Bio-Wraps
- Biopsy Pads
- Mesh Biopsy Bags

Intended Use

Detection/Measurement

The Leica Biosystems Biopsy Handling Devices do not detect or measure an analyte or marker. These medical devices are intended to securely house tissue specimens within the cassette while undergoing tissue processing.

Product Function

Leica Biosystems Biopsy Handling Devices are used to secure small tissue specimens inside a histology cassette during processing. After processing, handling device can be extracted and opened, allowing the specimen to be recovered and embedded in paraffin. Following microtomy, the specimens may be evaluated by a range of diagnostic techniques including but not limited to general histology stains, immunohistochemistry and *in situ* hybridization. The specimens when assessed by a trained pathologist may provide valuable data that is useful in the assessment of various physiological or pathological states. This data alongside other information such as the patients' medical history, physical condition, as well as results from other medical testing is all considered when rendering a medical diagnosis.

Specific Information Provided

Leica Biosystems Biopsy Handling Devices are not intended for the detection, definition or differentiation of a specific disorder, condition, or risk factor. Leica Biosystems Biopsy Handling Devices are intended for the preparation of tissue specimens for the subsequent application of specific diagnostic techniques that when evaluated by a trained pathologist may provide valuable data that is useful in the assessment of various physiological or pathological states. The pathological states may include but are not limited to a multitude of disease processes such as cancer, dysplasia, metaplasia, hypertrophy, infection, inflammation, ischemia, and scarring. This data alongside other information such as the patients' medical history, physical condition, as well as results from other medical testing is all considered when rendering a medical diagnosis.

Automation

Leica Biosystems Biopsy Handling Devices are not automated.

Quality/Quantity

Leica Biosystems Biopsy Handling Devices are not used for the detection of a specific analyte or marker.

Specimen Type

Leica Biosystems Biopsy Handling Devices are intended for use with human tissue that are to be examined using the appropriate histologic techniques.

Testing Population

Leica Biosystems Biopsy Handling Devices do not target a specific patient group. Leica Biosystems Biopsy Handling Devices are intended for use with any human tissue requiring histopathologic evaluation of biopsy specimens and resection tissue for the assessment of a suspected pathology or disease.

Intended User

Leica Biosystems Biopsy Handling Devices are intended for use by qualified laboratory personnel or designee.

In Vitro Diagnostic

The Leica Biosystems Biopsy Handling Devices are intended for *in vitro* diagnostics use only.

Test Principle

The Leica Biosystems Biopsy Handling Devices hold and protect the tissue while it undergoes processing.

Calibrators & Controls

The Leica Biosystems Biopsy Handling Devices do not require the use of any calibrators or controls.

Reagent Limitations

No reagent limitations are applicable to these products.

Applicable Products

Product Code	Material Description
3801090	Bio Wraps, 2" x 3" (50 x 75 mm), Box/500
3801091	Bio Wraps, 5" x 3½" (127 x 89 mm), Box/500
3801092	Bio Wraps, 2½" x 4" (63 x 100 mm), Box/500
3801085	Mesh Biopsy Bags, 1¼" x 1¼" (50 x 30 mm), Bag/500
3801087	Mesh Biopsy Bags, 2¼" x 1¾" (60 x 45 mm), Bag/500
3800003	Biopsy Pad, Blue, 1"x1 1/4" (25 x 30 mm), 10,000/case bulk

Biopsy Handling Devices

REF 3801090, 3801091, 3801092

3800003, 3801000, 3801020, 3801040, 3801060, 3801080

3801085, 3801087

3801000	Biopsy Pad, Blue, 1"x1 1/4" (25 x 30 mm), 500/bag
3801020	Biopsy Pad, White, 1"x1 1/4" (25 x 30 mm), 500/bag
3801040	Biopsy Pad, Blue, 1" (35 mm) circle, 500/bag
3801060	Biopsy Pad, Blue, 1" (25 mm) circle, 500/bag
3801080	Biopsy Pad, Blue, 3/4" (19 mm) circle, 500/bag

NOTE: Products listed here may not be available in all geographic regions.

Materials Not Included

Products should be placed inside a standard cassette, which is not included.

Devices Required

Not Applicable.

Storage and Stability

The Leica Biosystems Biopsy Handling Devices are to be stored at room temperature (15-30°C).

In Use Stability

Stable under normal conditions. User discretion should be utilized when determining in-use stability.

Sterility

The Leica Biosystems Biopsy Handling Devices are not sterile products.

Warnings/Precautions

Not Applicable.

Infectious Material Status

The Leica Biosystems Biopsy Handling Devices do not include any infectious material. However, specimens, fresh or before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

Special Facilities

The Leica Biosystems Biopsy Handling Devices should be used per facility guidelines.

Specimen Handling

Fresh tissue grossing - all specimens should be handled as if capable of transmitting infection.

Preparation for Use

- Bio Wraps

The Bio Wraps shall fit within a standard cassette after folding such that the cassette can be securely latched. A colored pigment may help visualize the specimens on the bio wraps. Colored pigment applied to the specimens shall be validated at the point of use by the end user.

- Mesh Biopsy Bags

Remove a single bag from the dispenser. Use only with fixatives and processing reagents.

- Biopsy Pads

The Pads shall not extend beyond the internal side walls such that a standard size cassette can be securely latched.

Biopsy Pads, together with tissue specimens, should not exceed internal space of a standard cassette (4 mm in height).

To prevent floating and avoid artifacts, Biopsy Pads should be pre-saturated in fixative. Biopsy pads, when soaked in formalin, shall cover all pores of a cassette such that when moved from side to side, all pores remain covered.

Direction for Use

- Bio-Wraps

Remove single Wrap from a dispenser pack, place biopsy specimen in the center of it and fold the Wrap several times.

Place the Wrap inside the cassette making sure it fits entirely inside of it. Wrap can be also pre-wet to prevent drying artifact.

- Mesh Biopsy Bags

The fixative and specimen can be poured into the bag which filters fixative and leaves biopsies inside. The bag is then folded over at least once, placed inside a cassette, and processed as normal. Specimens can also be transferred from a vessel containing fixative into a Biopsy Bag using forceps. After processing, the bag is extracted and opened, allowing the specimen to be recovered and embedded in paraffin. Warning: Biopsy Bags cannot be used with acid based decalcifying solutions.

- Biopsy Pads

Use Pads to hold small biopsies in place during fixation and processing to prevent them from being lost. Samples can also be sandwiched between two foam pads and placed in cassettes with metal or plastic lids. Blue colored Pads assist with tissue visibility.

Readiness for Use

Biopsy Handling Devices are ready to use products.

Quality Control

Quality control should be assessed at the point of use.

Biopsy Handling Devices

REF

3801090, 3801091, 3801092

3800003, 3801000, 3801020, 3801040, 3801060, 3801080

3801085, 3801087

Expected Results

By following the instructions for use, Biopsy Handling Devices shall allow for flow of fixative and/or processing reagents while securing specimens inside a cassette.

Analytical Performance

The Leica Biosystems Biopsy Handling Devices do not detect or measure an analyte or marker. These medical devices are intended to securely house tissue specimens within the cassette while undergoing tissue processing.

Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

Clinical Performance

The Leica Biosystems Biopsy Handling Devices are not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the Leica Biosystems Bluing Agents in a clinical setting.


Disposal

Dispose of biopsy handling device in accordance with licensed collector's sorting instructions or per facility guidelines.

Leica Biosystems Richmond, Inc.
5205 Route 12
Richmond, IL 60071
USA
(1-844-534-2262)

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