# Leica Biosystems Special Stain Kit

English.....	2
العربية (Arabic).....	6
简体中文 (Chinese Simplified).....	10
中國傳統的 (Chinese Traditional).....	14
Dansk (Danish).....	18
Nederlands (Dutch).....	22
Français (French – Canada).....	27
Français (French – France).....	32
Deutsch (German).....	37
Italiano (Italian).....	42
日本語 (Japanese).....	47
한국어 (Korean).....	51
Norsk (Norwegian).....	56
Polski (Polish).....	60
Português (Portuguese – Brazil).....	65
Português (Portuguese – Portugal).....	70
Română (Romanian).....	75
Русский (Russian).....	80
Slovenski (Slovenian).....	85
Español (Spanish – Central America).....	90
Español (Spanish – Spain).....	95
Svensk (Swedish).....	100
ไทย (Thai).....	105
Türk (Turkish).....	110
Tiếng Việt (Vietnamese).....	115

# Special Stain Kit

## Modified Grocott's Methenamine Silver Stain

**REF** 38016SS12

### Product Name

Modified Grocott's Methenamine Silver Stain (GMS) Special Stain Kit.

### Intended Use

#### Detection/Measurement

Modified Grocott's Methenamine Silver Stain (GMS) Kit does not detect or measure an analyte or marker.

The Modified Grocott's Methenamine Silver Stain when used with appropriate histological protocols may be used for the demonstration of defined fungi and infectious agents such as *Aspergillus sp.*, *Pneumocystis carinii*, and *Cryptococcus neoformans* in formalin fixed, paraffin embedded tissue sections.

#### Product Function

The results obtained through use of GMS Special Stain Kit do not provide objective medical evidence. The coloration and contrast the Leica Biosystems GMS Special Stain Kit provide to histologic specimens allows visualization of microscopic anatomy. This visualization, when interpreted by a trained professional, is utilized alongside other information such as the patient's medical history, physical condition, as well as results from other medical testing to render a medical diagnosis.

#### Specific Information Provided

The Leica Biosystems GMS Kit is not intended for the detection, definition or differentiation of a specific disorder, condition, or risk factor. The staining demonstrated with use of these products, when used as intended, provides trained professionals information which may define the physiological or pathological state of the tissue specimen.

#### Automation

GMS Special Stain Kit is not automated but can be used on automated staining platforms. Use on an automated platform should be validated at the point of use.

#### Qualitative/Quantitative

The Leica Biosystems GMS Special Stain Kit is qualitative stain.

#### Specimen Type

GMS Special Stain Kit can be used with any paraffin embedded human or animal specimen.

#### Testing Population

The Leica Biosystems GMS Kit is intended for use with any patient requiring evaluation of biopsy or resection tissue for the assessment of a suspected pathology or disease.

#### Intended User

GMS Special Stain Kit is intended for use by qualified laboratory personnel and/or designee of the laboratory.

### In Vitro Diagnostic

GMS Kit is intended for *in vitro* diagnostics use only.

### Test Principle

The mechanism of action of Modified Grocott's Methenamine Silver Stain is based upon the capacity of aldehyde groups to reduce cationic silver (Ag<sup>+</sup>) to metallic silver. Chromic acid is used to generate aldehyde groups by the oxidation of 1-2 glycol groups within polysaccharide rich tissue components, e.g. glycogen, mucin, reticulin and fungal cell walls. When cationic silver is added to the section in the form of a methenamine-silver ion complex, the aldehyde groups reduce the silver ions to metallic silver. Sections are subsequently toned with Gold Chloride Solution to produce metallic gold which is more stable than metallic silver and produces superior contrast and clarity.

Because of the strong oxidizing potential of the Modified Chromic Acid Solution, many of the resultant aldehyde groups are further oxidized to carboxylic acid groups which are incapable of reducing silver. This capacity of the Modified Chromic Acid Solution has the advantage of reducing background reactions of collagen and basement membranes, and produces a strong impregnation with silver in only those structures that possess high levels of the reactive polysaccharide groups e.g. glycogen, mucin and fungal cell walls.

### Calibrators & Controls

GMS Kit does not require the use of any calibrators or controls.

### Reagent Limitations

No reagent limitations are applicable to this product.

### Applicable Products

Product Code	Material Description
38016SS12	Modified Grocott's Methenamine Silver Stain (GMS) Special Stain Kit
38016SS12A	Modified Chromic Acid Solution, 500 mL
38016SS12B	Sodium Bisulfite Solution, 500 mL
38016SS12C	Silver Nitrate Solution, 250 mL

# Special Stain Kit

## Modified Grocott's Methenamine Silver Stain

**REF** 38016SS12

38016SS12D	Methenamine/Borax Solution, 250 mL (A)
38016SS12E	Gold Chloride Solution, 500 mL
38016SS12F	Sodium Thiosulfate Solution, 500 mL
38016SS12G	Light Green SF, 500 mL

**Note: (A).** The refrigerated Methenamine/Borax Solution (Item No. 38016SS12D, 250 mL) is not included in the GMS Kit. Item must be ordered separately and will be shipped separately.

### Materials Not Included

GMS Special Stain Kit protocol requires the use of graded alcohols, xylene, or xylene substitutes, deionized or distilled water. Positive control slide(s) with tissue known to contain fungal elements (not included in this kit), should be included in each run.

### Devices Required

Leica Biosystems GMS Special Stain Kit can be used on any automated staining platform or with a manual staining method.

### Storage and Stability

The Methenamine/Borax Solution should be stored at 2–8 °C (36–46 °F). Store other components at room temperature 15–30 °C (59–86 °F).

**CAUTION:** Do not use after the expiration date.

### In Use Stability

User discretion should be utilized when determining in-use stability.

### Sterility

GMS Special Stain Kit components are not sterile products.

### Warnings/Precautions

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations. Refer to Material Safety Data Sheet and product labeling for any updated risk, hazard, or safety information.

### Infectious Material Status

GMS Special Stain Kit does not include any infectious material. However, specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

### Special Facilities

GMS Special Stain Kit should be used per facility guidelines.

### Specimen Handling

Suggested fixatives include 10% neutral buffered formalin. Routine dehydration, clearing, and paraffin infiltration and embedding,

and routine preparation of microtome sections. Poor fixation, processing, rehydration, and sectioning will adversely affect the staining quality. Following processing and paraffin embedding, cut sections at 4–6 microns.

### Preparation for Use

Prior to staining, place 20 mL of the Methenamine/Borax Solution and 20 mL of the Silver Nitrate Solution into separate acid cleaned beakers. Heat the solutions to 50–55 °C and just prior to silver impregnation combine the two solutions in an acid washed Coplin jar.

### Notes:

- Solutions may be heated using a water bath or laboratory oven. Avoid over-heating as break-down of Methenamine is accelerated at temperatures above 50 °C.
- Preheat a second Coplin jar with 40 mL of deionized or distilled water to 50–55 °C to be used as a rinse.

### Direction for Use

#### Conventional Staining Protocol

- Deparaffinize tissue sections with xylene and rehydrate through graded alcohols to deionized or distilled water.
- Place slides in Modified Chromic Acid Solution for 5 –10 minutes at room temperature. Please note: Exercise caution when handling the Modified Chromic Acid Solution.
- Rinse slides in two changes of tap water.
- Rinse slides in two changes of deionized or distilled water.
- Place slides in Sodium Bisulfite Solution for 1 minute.
- Rinse slides in running tap water for 30 seconds.

# Special Stain Kit

## Modified Grocott's Methenamine Silver Stain

**REF** 38016SS12

7. Rinse slides in two changes of deionized or distilled water.
8. Combine the prewarmed Silver Nitrate Solution and Methenamine/Borax Solution into a prewarmed acid cleaned Coplin jar.
9. Place slides in the Methenamine/Borax-Silver Nitrate Solution and incubate for 20–45 minutes at 50–55 °C. After 15–20 minutes, using non-metallic forceps, remove a control slide, dip in the prewarmed deionized water to rinse, and check microscopically for the completeness of silver deposition.
10. Rinse slides in 6 changes of deionized or distilled water.
11. Place slides in Gold Chloride Solution for 5 minutes.
12. Rinse slides in 3 changes of deionized or distilled or distilled water.
13. Place slides in Sodium Thiosulfate Solution for 2 minutes.
14. Rinse slides thoroughly in running tap water for 2 minutes.
15. Place slides in Light Green SF for 40 seconds.
16. Rinse slides briefly in deionized or distilled water.
17. Dehydrate slides in three changes of absolute alcohol.
18. Clear slides in two changes of xylene and mount in a xylene miscible medium.

**Table 1. Example of Conventional GMS Special Stain Kit Staining Protocol.**

Steps	Action	Chemical	Time (mm:ss)
1-3	Deparaffinize	Xylene	3:00
4-5	Hydration	100% Alcohol	2:00
6	Hydration	80% or 95% Alcohol	1:00
7	Hydration	Deionized Water	1:00
8	Stain	Modified Chromic Acid Solution	5:00-10:00
9-10	Wash	Tap Water	0:30
11-12	Rinse	Deionized Water	0:10
13	Neutralize	Sodium Bisulfite Solution	1:00
14	Rinse	Running Tap Water	0:30
15-16	Rinse	Deionized Water	0:10
17	Stain/Impregnate	Methenamine/Borax-Silver Nitrate Solution	20:00-45:00 @ 50–55 °C
18-23	Rinse	Deionized Water	0:10
24	Dehydration	Gold Chloride Solution	5:00
25-27	Rinse	Deionized Water	0:10
28	Stain	Sodium Thiosulfate Solution	2:00
29	Rinse	Running Tap Water	2:00
30	Counterstain	Light Green SF	0:40
31	Rinse	Deionized Water	0:30
32-34	Dehydration	100% Alcohol	2:00
35-36	Clearing	Xylene	2:00

Note: When using xylene substitute, increase immersion times by approximately 50%.

### Microwave Staining Protocol

Exercise caution when using the microwave to heat any solution or reagent. The microwave must be properly ventilated to prevent the accumulation of fumes in the laboratory. Microwave transparent coplin jars and caps should be used during the staining process. The caps should be loosely applied to prevent spills. Caps with ventilation holes also may be used. All microwaves should be used in accordance with the manufacturer's instructions.

1. Deparaffinize tissue sections with xylene and rehydrate through graded alcohols to deionized or distilled water.
2. Place slides in the Modified Chromic Acid Solution and microwave at 800 watts for 10 seconds.  
Please note: Exercise caution when handling the Modified Chromic Acid Solution.
3. Gently mix by swirling the Coplin jar and allow to stand for 1 minute.
4. Rinse slides in 6 changes of deionized or distilled water.
5. Place slides in Sodium Bisulfite Solution for 1 minute.

# Special Stain Kit

## Modified Grocott's Methenamine Silver Stain

REF 38016SS12

- Rinse slides in 6 changes of deionized or distilled water.
- Combine 20 mL of the Silver Nitrate Solution and 20 mL of the Methenamine/Borax Solution in an acid cleaned plastic Coplin jar.
- Place slides in the solution and microwave at 600 watts for 45–50 seconds.
- Gently mix by swirling the Coplin jar and allow to stand for 2 minutes.
- Microwave solution for 10 seconds at 600 watts.
- Gently agitate solution and using non-metallic forceps remove a control slide, dip in prewarmed deionized water to rinse and check microscopically for the completeness of silver impregnation. If necessary replace slides in the warm Methenamine/Borax Silver Nitrate Solution and swirl for 10–30 seconds.
- If necessary, repeat steps 10 and 11 until a satisfactory level of silver deposition is observed.
- Rinse slides in 6 changes of deionized or distilled water.
- Place slides in Gold Chloride Solution for 5 minutes.
- Rinse slides in 3 changes of deionized or distilled water.
- Place slides in Sodium Thiosulfate Solution for 2 minutes.
- Rinse slides thoroughly in running tap water for 2 minutes.
- Place slides in Light Green SF for 40 seconds.
- Rinse slides briefly in deionized or distilled water.
- Dehydrate slides in 3 changes of absolute alcohol.
- Clear slides in 2 changes of xylene and mount in a xylene miscible medium.

### Readiness for Use

Once appropriate staining protocol is chosen and bath-layout is created, pour all the reagent into the reagent vessel. Place the reagent vessel back into the respective station.

### Quality Control

A quality control slide(s) with tissue known to contain fungal elements, fixed and processed in a similar manner to the test specimens should be included in each staining assay to ensure GMS Special Stain Kit is performing as intended.

### Expected Results

- Fungi – sharply delineated black
- Inner parts of mycelia and hypha – taupe to old rose
- Mucin – taupe to dark grey
- Background – green

### Analytical Performance

The Leica Biosystems Modified Grocott's Methenamine Silver Stain Kits are not used to detect a specific analyte or marker. These products are used for demonstration of defined fungi and infectious agents such as *Aspergillus sp.*, *Pneumocystis carinii* and *Cryptococcus neoformans* in formalin fixed, paraffin embedded tissue sections. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

### Clinical Performance

The Leica Biosystems Modified Grocott's Methenamine Silver Stain Kits are not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the Leica Biosystems Bluing Agents in a clinical setting.


### Disposal

Spent or expired components of the GMS Special Stain Kit should be discarded in accordance with organizational, local, state, and federal regulations.

Leica Biosystems Richmond, Inc.  
5205 Route 12  
Richmond, IL 60071  
USA  
(1-844-534-2262)

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