## BOND Hybridization Solution

Catalog No: AR9037

Leica Biosystems Newcastle Ltd  
Balliol Business Park  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
☎ +44 191 215 4242

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## Catalog No: AR9037

### Intended Use

This reagent is for *in vitro* diagnostic use.

BOND Hybridization Solution is intended to be used for the dilution of *in situ* hybridization (ISH) probes for use on the automated BOND system (includes BOND-MAX system and BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

### Summary and Explanation

ISH allows visualization of specific nucleic acids in tissue sections (see "Using BOND Reagents" in your BOND user documentation).

The BOND system (includes BOND-MAX system and BOND-III system) provides an automated method of *in situ* hybridization on formalin-fixed, paraffin-embedded tissue sections. BOND Hybridization Solution has been developed for the dilution of *in situ* hybridization probes for use on the automated BOND system (includes BOND-MAX system and BOND-III system).

### Reagents Provided

BOND Hybridization Solution.

Total volume = 100 mL.

### Dilution and Mixing

BOND Hybridization Solution is ready to use. Reconstitution, mixing, dilution or titration of this reagent is not required.

### Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and ISH staining using the BOND system (includes BOND-MAX system and BOND-III system).

### Storage and Stability

Store at 2–8 °C. The product is stable under these conditions up to the expiry date indicated on the container label. Do not use after the expiration date indicated on the bottle label.

There are no obvious signs indicating contamination and/or instability. Appropriate positive and negative tissue controls should be run at the same time as test tissue.

BOND Hybridization Solution contains a small quantity of sodium dodecyl sulfate (SDS). This can appear as a white precipitate when the solution is first removed from 2–8 °C. The precipitate disappears when the solution reaches room temperature. The presence of the precipitate does not affect the performance of the solution.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

### Precautions

- This product is intended for *in vitro* diagnostic use.

#### BOND HYBRIDIZATION SOLUTION

Contains Formamide (<50%) and Dextran Sulphate (<30%).

GHS07: Exclamation mark.

GHS08: Health hazard.

Signal words: Danger.

H360D: May damage the unborn child.  
H315: Causes skin irritation.  
H319: Causes serious eye irritation.

P201: Obtain special instructions before use.  
P202: Do not handle until all safety precautions have been read and understood.  
P260: Do not breathe dust/fumes/gas/mist/vapours/spray.  
P308+313: If exposed or concerned: Get medical advice/attention.  
P314: Get medical advice/attention if you feel unwell.  
P264: Wash hands thoroughly after handling.  
P280: Wear protective gloves/protective clothing/eye protection/face protection.  
P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.  
P337+313: If eye irritation persists: Get medical advice/attention.  
P302+352: IF ON SKIN: Wash with plenty of soap and water.  
P332+313: If skin irritation occurs: Get medical advice/attention.  
P362+P364: Take off contaminated clothing and wash before reuse.

- To obtain a copy of the Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

## Instructions for Use

Dilute the ISH probe concentrate to the user-determined concentration with BOND Hybridization Solution. Transfer the diluted probe to a BOND Open Container and register as an ISH probe on the BOND system (includes BOND-MAX system and BOND-III system).

Probe dilutions must be tested using the appropriate positive control tissue for the particular ISH probe (see "Quality Control" in the "Using BOND Reagents" section of your user documentation).

When preparing dilutions for titration, ensure that allowance is made for the dead volume of the appropriate BOND Open Container (for values refer to your BOND user documentation).

## Product Specific Limitations

BOND Hybridization Solution was developed for use on the automated BOND system (includes BOND-MAX system and BOND-III system).

The recommended pre-treatment and staining protocol are determined by the specific probe (and detection) selected. Please refer to the relevant Instructions for Use for further information.

Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. Protocol times may vary due to variation in tissue type, fixation and processing. In addition BOND Enzyme concentration and incubation time may require optimization depending on tissue type, processing and fixation conditions.

Negative reagent controls should be used for CISH when optimizing pre-treatment conditions and protocol times.

The appropriate concentration of user's own ISH probe may vary, due to variation in tissue fixation and the effectiveness of enzymatic digestion, and must be determined empirically.

## Troubleshooting

Reference 3 may aid in remedial action.

Test samples should be complemented by the appropriate tissue and reagent controls.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

## Further Information

Further information on ISH with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

## Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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### Bibliographie

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In Situ Hybridization. A practical approach*. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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### Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In Situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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Leica Biosystems Inc
Pty Ltd