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Catalog No: PA0488

Intended Use

This reagent is for *in vitro* diagnostic use.

Inhibin Alpha (R1) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human inhibin alpha in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Inhibin Alpha (R1) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human inhibin alpha is achieved by first allowing the binding of Inhibin Alpha (R1) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

Inhibin Alpha (R1) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

R1

Immunogen

Synthetic peptide corresponding to amino acids 1-32 of the mature C-terminal region of the human inhibin alpha subunit. This corresponds to position 233-264 of the human inhibin alpha sequence starting at the signal peptide.

Specificity

The antibody will react with all size forms of the free alpha subunit which contain the mature sequence of the alpha subunit containing the epitope. It will also detect inhibin A or B forms which contain this region of the alpha subunit.

Ig Class

IgG2a

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 0.8 mg/L as determined by ELISA.

Dilution and Mixing

Inhibin Alpha (R1) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Inhibin Alpha (R1) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Inhibin Alpha (R1) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection. The recommended staining protocol for Inhibin Alpha (R1) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 2 for 40 minutes.

Results Expected

Normal Tissues

Clone R1 detected the inhibin alpha protein in the cytoplasm of epithelial cells in testis and epididymis, granulosa cells in ovary, and secretory cells in adrenal gland tissue. (Total number of normal cases = 126).

Tumor Tissues

Clone R1 stained 9/50 theca cell tumors (including 6/44 benign theca cell tumors and 3/6 malignant theca cell tumors), 2/10 cysts (including 1/1 follicular cyst, 1/1 corpus luteum cyst, 0/4 chocolate cysts, 0/3 simple cysts and 0/1 serous cyst), 2/2 tumors of the adrenal gland, 1/6 tumors of the head and neck (including 1/1 nasopharyngeal tumor, 0/1 adenocarcinoma of the oral cavity, 0/1 squamous cell carcinoma of the tongue, 0/1 pleomorphic adenoma of the salivary gland, 0/1 adenoid cystic carcinoma of the salivary gland and 0/1 melanoma of the nasal cavity), 1/4 hepatocellular carcinomas and 1/3 ovarian tumors (including 1/1 granulosa cell tumor, 0/1 ovarian adenocarcinoma and 0/1 endometrioid adenocarcinoma). No staining was observed in a variety of additional tumors evaluated, including intestinal tumors (0/9), breast tumors (0/5), metastatic tumors (0/5), thyroid tumors (0/5), brain tumors (0/4), lung tumors (0/4), esophageal tumors (0/3), stomach tumors (0/3), lymphomas (0/3), prostate tumors (0/3), squamous cell carcinomas of the cervix (0/2), endometrial tumors (0/2), bone tumors (0/2), kidney tumors (0/2), bladder carcinomas (0/2), a tumor of the testis (0/1), a pancreatic tumor (0/1), a pheochromocytoma (0/1), or a skin tumor (0/1). (Total number of abnormal cases = 133).

Inhibin Alpha (R1) is recommended for the detection of human inhibin (alpha) protein in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.

Product Specific Limitations

Inhibin Alpha (R1) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Robertson DM, Burger HG and Fuller PJ. Inhibin/activin and ovarian cancer. Endocrine-Related Cancer, 2004;11:35-49.
5. Zheng W, Senturk BZ and Parkash V. Immunohistochemical staining: A practical approach for the surgical pathologist in the diagnoses of ovarian sex cord stromal tumours. Advances in Anatomic Pathology, 2003;10(1):27-38.
6. Bernard DJ, Chapman SC and Woodruff TK. Mechanisms of inhibin signal transduction. Recent Progress in Hormone Research, 2001;56:417-450.
7. Groome N and Lawrence M. Preparation of monoclonal antibodies to the beta A subunit of ovarian inhibin using a synthetic peptide immunogen. Hybridoma, 1991;10(2):309-316.
8. Groome N, Hancock J, Betteridge A, et al. Monoclonal and polyclonal antibodies reactive with the 1-32 amino terminal sequence of the alpha subunit of human 32K inhibin. Hybridoma, 1990;9(1):31-42.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Robertson DM, Burger HG and Fuller PJ. Inhibin/activin and ovarian cancer. Endocrine-Related Cancer, 2004;11:35-49.
5. Zheng W, Senturk BZ and Parkash V. Immunohistochemical staining: A practical approach for the surgical pathologist in the diagnosis of ovarian sex cord stromal tumours. Advances in Anatomic Pathology, 2003;10(1):27-38.
6. Bernard DJ, Chapman SC and Woodruff TK. Mechanisms of inhibin signal transduction. Recent Progress in Hormone Research, 2001;56:417-450.
7. Groome N and Lawrence M. Preparation of monoclonal antibodies to the beta A subunit of ovarian inhibin using a synthetic peptide immunogen. Hybridoma, 1991;10(2):309-316.
8. Groome N, Hancock J, Betteridge A, et al. Monoclonal and polyclonal antibodies reactive with the 1-32 amino terminal sequence of the alpha subunit of human 32K inhibin. Hybridoma, 1990;9(1):31-42.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Robertson DM, Burger HG and Fuller P.J. Inhibin/activin and ovarian cancer. Endocrine-Related Cancer, 2004;11:35-49.
5. Zheng W, Senturk BZ and Parkash V. Immunohistochemical staining: A practical approach for the surgical pathologist in the diagnoses of ovarian sex cord stromal tumours. Advances in Anatomic Pathology, 2003;10(1):27-38.
6. Bernard DJ, Chapman SC and Woodruff TK. Mechanisms of inhibin signal transduction. Recent Progress in Hormone Research, 2001;56:417-450.
7. Groome N and Lawrence M. Preparation of monoclonal antibodies to the beta A subunit of ovarian inhibin using a synthetic peptide immunogen. Hybridoma, 1991;10(2):309-316.
8. Groome N, Hancock J, Betteridge A, et al. Monoclonal and polyclonal antibodies reactive with the 1-32 amino terminal sequence of the alpha subunit of human 32K inhibin. Hybridoma, 1990;9(1):31-42.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Robertson DM, Burger HG and Fuller PJ. Inhibin/activin and ovarian cancer. Endocrine-Related Cancer, 2004;11:35-49.
5. Zheng W, Senturk BZ and Parkash V. Immunohistochemical staining: A practical approach for the surgical pathologist in the diagnoses of ovarian sex cord stromal tumours. Advances in Anatomic Pathology, 2003;10(1):27-38.
6. Bernard DJ, Chapman SC and Woodruff TK. Mechanisms of inhibin signal transduction. Recent Progress in Hormone Research, 2001;56:417-450.

7. Groome N and Lawrence M. Preparation of monoclonal antibodies to the beta A subunit of ovarian inhibin using a synthetic peptide immunogen. *Hybridoma*, 1991;10(2):309-316.
8. Groome N, Hancock J, Betteridge A, et al. Monoclonal and polyclonal antibodies reactive with the 1-32 amino terminal sequence of the alpha subunit of human 32K inhibin. *Hybridoma*, 1990;9(1):31-42.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Robertson DM, Burger HG and Fuller PJ. Inhibin/activin and ovarian cancer. Endocrine-Related Cancer, 2004;11:35-49.
5. Zheng W, Senturk BZ and Parkash V. Immunohistochemical staining: A practical approach for the surgical pathologist in the diagnoses of ovarian sex cord tumours. Advances in Anatomic Pathology, 2003;10(1):27-38.
6. Bernard DJ, Chapman SC and Woodruff TK. Mechanisms of inhibin signal transduction. Recent Progress in Hormone Research, 2001;56:417-450.
7. Groome N and Lawrence M. Preparation of monoclonal antibodies to the beta A subunit of ovarian inhibin using a synthetic peptide immunogen. Hybridoma, 1991;10(2):309-316.
8. Groome N, Hancock J, Betteridge A, et al. Monoclonal and polyclonal antibodies reactive with the 1-32 amino terminal sequence of the alpha subunit of human 32K inhibin. Hybridoma, 1990;9(1):31-42.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd