## HistoCore SPECTRA H&E Staining System S1

English	2
العربية (Arabic)	5
中文 (简体) (Chinese Simplified)	8
中文(Chinese Traditional)	11
Dansk (Danish)	14
Français (French) Canada	17
Français (French) European	20
Deutsch (German)	23
Italiano (Italian)	26
日本語 (Japanese)	29
한국어 (Korean)	32
Norsk (Norwegian)	35
Polski (Polish)	38
Português (Portuguese) Brazil	41
Português (Portuguese) European	44
Română (Romanian)	47
Русский (Russian)	50
Slovenski (Slovenian)	54
Español (Spanish) Central America	57
Español (Spanish) European	61
Svensk (Swedish)	64
ไทย (Thai)	67
Türk (Turkish)	70
Tiếng Việt (Vietnamese)	73

# HistoCore SPECTRA H&E Staining System S1

**REF** 3801654

## HistoCore SPECTRA H&E Staining System S1

The HistoCore SPECTRA H&E Staining System S1 contains 5 proprietary reagent components and staining protocols that when used together provide consistent diagnostic quality and ease of use. This system has been optimized for use on HistoCore SPECTRA ST stainers.

### Intended Use

The reagents in this HistoCore SPECTRA H&E Staining System S1 are intended for *in vitro diagnostics* use only. The reagents in the HistoCore SPECTRA H&E Staining System S1 are intended for use on HistoCore SPECTRA ST Stainers. The HistoCore SPECTRA H&E Staining System S1 is used to stain histological samples with Hematoxylin and Eosin for subsequent microscopic evaluation.

### Reagents Provided

SPECTRA Hemalast S1 (REF 3801580, 380mL)  
SPECTRA Hematoxylin S1 (REF 3801581, 380mL)  
SPECTRA Differentiator S1 (REF 3801582, 380mL)  
SPECTRA Bluing Agent S1 (REF 3801583, 380mL)  
SPECTRA Eosin S1 (REF 3801584, 380mL)

### Storage and Stability

Store reagents at room temperature (15-30°C). Please follow local storage requirements for flammable products.

**CAUTION:** Do not use after the expiration date.

### Specimen Preparation:

#### Fixation

Any general fixative including but not limited to neutral buffered formalin may be used.

**NOTE:** Under fixation or over fixation may result in compromised staining quality.

#### Paraffin Sections

Following processing and paraffin embedding, section tissues at standard thickness (2 – 5µm).

### Directions for use of the HistoCore SPECTRA H&E Staining System S1:

Settings for operational parameters of the HistoCore SPECTRA ST stainer and protocol schedules are provided to optimize the performance and consistency of the HistoCore SPECTRA H&E Staining System S1.

**NOTE:** Staining protocols may be downloaded at [LeicaBiosystems.com/SPECTRAStainingProtocols](http://LeicaBiosystems.com/SPECTRAStainingProtocols)

### Slides and Specimens

If the oven step is not being utilized on the instrument prior to staining, slides should be dried at 66 - 74°C for 15-22 minutes.

### Oven Specifications:

To either activate the oven step, or deactivate the oven step for the SPECTRA H&E protocol, refer to the HistoCore SPECTRA ST Instructions for Use.

### Placing the HistoCore SPECTRA Staining System S1 on the Instrument:

Scan the HistoCore SPECTRA H&E staining system box and each individual bottle at the Radio Frequency Identification (RFID) scanner located on the left side of the HistoCore SPECTRA ST stainer. Follow the prompts as indicated by the instrument to place the reagents into the designated reagent vessels.

### HistoCore SPECTRA ST Bathlayout

When prompted by the user, the bathlayout is automatically generated by the instrument. The bathlayout includes specific protocols chosen by the user. The instrument will assign each reagent to a specific reagent station. The 95% alcohol in step 16 of the protocol must be duplicated manually as described below:

- In the station assignment view, press the "Adapt Bathlayout" button.
- To duplicate a reagent station, select it in the bathlayout by touching the screen.
- Finally, press the "Duplicate" button. Then, the reagent station will appear in the clipboard with an equivalence symbol (∞).
- After this, select the duplicate reagent station in the clipboard by touching the screen and place it in the desired (available) station, near the original station in the bathlayout.
- To save your changes, press the "Finish" button and confirm the subsequent message with a "Yes" or "No" to return to the clipboard.

**NOTE:** For additional instructions on adapting bathlayouts, refer to the SPECTRA ST Instructions for Use.

### Agitation Settings:

Agitation settings are predefined. Agitation will occur at a setting of 4. Refer to the HistoCore SPECTRA ST Instructions for Use.

### Immersion Times:

Immersion times and instrument settings for HistoCore SPECTRA ST stainer are detailed in Table 1 below. Please follow instructions for operating HistoCore SPECTRA ST as detailed in the HistoCore SPECTRA ST Instructions for Use.

**CAUTION:** Modifications to the validated staining protocol may result in compromised staining quality as well as staining capacity.

**CAUTION:** All of the reagents in each HistoCore SPECTRA H&E Staining System S1 must be used together to achieve the expected performance. Substitution of non-HistoCore SPECTRA H&E reagents or reagents from another HistoCore SPECTRA H&E Staining System will compromise performance.

**CAUTION:** Opening the lid of the HistoCore SPECTRA ST stainer during protocol execution may increase time in staining system reagent and impact the consistency in staining quality.

**CAUTION:** Reagent vessels should be covered when not in use to reduce evaporation.

# HistoCore SPECTRA H&E Staining System S1

**REF** 3801654

**CAUTION:** If the HistoCore SPECTRA ST is not connected to the HistoCore SPECTRA CV, the slides will be transferred to the unload drawers. In this occurrence, ensure there is the appropriate reagent in the unload vessel, otherwise quality of staining may be compromised.  
**CAUTION:** Any deviations from these instructions may alter expected results and are neither validated nor supported by Leica Biosystems.

## Technical Notes

- 100% alcohol and xylene reagents used for deparaffinization and rehydration should be rotated after every 800 slides. Reagent rotation is defined as follows:
  - Empty most contaminated reagent of a series and refill container.
  - Move remaining container(s) up one space in sequence.
  - Place newly filled container into the last position of the series.
- The 100% alcohol and xylene reagents used for dehydration and clearing before coverslipping also should be rotated after every 800 slides as described above.
- Two of the four 95% alcohols (steps 5, and 14) should be discarded and replaced with fresh 95% alcohol solution every 800 slides.
- The two 95% dehydrating alcohols (step 16) should be discarded and replaced with fresh 95% alcohol solution every 400 slides per reagent vessel.
- If using the Reagent Management System (RMS), after change and rotation, update the RMS with changes. Instruction for use of the RMS can be found in the HistoCore SPECTRA ST Instructions for Use.
- CAUTION:** Failure to rotate and change reagents as required may result in inadequate deparaffinization as well as excessive carry-over and dilution of reagents or solvent.
- NOTE:** Reagent vessels (including water rinse vessels) should be inspected at periodic intervals to ensure the integrity and cleanliness of the vessels. If cleaning solutions are used, rinse vessels with adequate water to remove any residual cleaning solutions prior to use. Clean the vessels in accordance with the instructions provided in the HistoCore SPECTRA ST Instructions for Use.
- NOTE:** The HistoCore SPECTRA ST Reagent Management System (RMS) tracks the number of slides that have utilized a specific station. However, for the ancillary reagents that have been rotated, the RMS will no longer capture the total number of slides that have utilized that specific ancillary reagent. For more information please refer to the SPECTRA ST Instructions for Use.

Table 1	Protocol Name: H&E Staining System S1		Abbreviation: HE1		Agitation Speed: 4		Oven Temperature: 70°C	
Step	Reagent	Abbreviation	Time (mm:ss)	Tolerance	Process Class	Heated Vessel	Days Max	Slide Max
Default: Include Oven Step	Oven	Oven	15:00	50%	NA	NA	NA	NA
1	Xylene Dewax 1	XylDwx 1	2:00	100%	Dewaxing	No	5	800
2	Xylene Dewax 2	XylDwx 2	2:00	100%	Dewaxing	No	5	800
3	100% Alcohol Dewax 1	100Dwx 1	2:00	100%	Dewaxing	No	5	800
4	100% Alcohol Dewax 2	100Dwx 2	2:00	100%	Dewaxing	No	5	800
5	95% Alcohol Dewax 1	95Dwx 1	2:00	100%	Dewaxing	No	5	800
6	Tap Water	Tap Water	2:00	100%	NA	No	NA	NA
7	SPECTRA Hemalast S1	L HLst S1	2:00	100%	Neutralizing	No	7	1600
8	SPECTRA Hematoxylin S1	L Hmtx S1	Refer to Hematoxylin Dial Settings	0%	Staining	No	7	1600
9	Tap Water	Tap Water	2:00	100%	NA	No	NA	NA
10	SPECTRA Differentiator S1	L Diff S1	1:30	0%	Differentiating	No	7	1600
11	Tap Water	Tap Water	2:00	100%	NA	No	NA	NA
12	SPECTRA Bluing Agent S1	L Blue S1	2:00	100%	Staining	No	7	1600
13	Tap Water	Tap Water	2:00	100%	NA	No	NA	NA
14	95% Alcohol Neutralizing H	95Neutr H	2:00	100%	Neutralizing	No	5	800
15	SPECTRA Eosin S1	L Eos S1	Refer to Eosin Dial Settings	0%	Staining	No	7	1600
16*	95% Alcohol Dehyd 1	95Dhy 1	2:00	0%	Dehydrating	No	5	400
17	100% Alcohol Dhyd 1 H	100Dhy 1H	2:00	0%	Dehydrating	No	5	800
18	100% Alcohol Dhyd 2 H	100Dhy 2H	2:00	0%	Dehydrating	No	5	800
19	Xylene Dehyd 1	XylDhy 1	2:00	100%	Dehydrating	No	5	800
20	Xylene Dehyd 2	XylDhy 2	2:00	100%	Dehydrating	No	5	800
21	Exit	NA	NA	NA	NA	NA	NA	NA

\* Indicates a manually duplicated step.

# HistoCore SPECTRA H&E Staining System S1

**REF** 3801654

## Protocol Adjustments

The validated HistoCore SPECTRA H&E Staining System S1 protocol comes pre-loaded on the HistoCore SPECTRA ST and will default to a dial setting of 4 for both the Hematoxylin and Eosin. Only the Hematoxylin and Eosin times may be adjusted when using this staining system. All other protocol steps and times outlined in Table 1 are fixed. To adjust the staining intensity of the Hematoxylin and/or Eosin, use the Dial Settings in Table 2 below.

Table 2: Dial Settings			
Hematoxylin Dial Settings		Eosin Dial Settings	
Setting	Timing (mm:ss)	Setting	Timing (mm:ss)
Dial 1	1:00	Dial 1	0:30
Dial 2	2:00	Dial 2	1:00
Dial 3	3:00	Dial 3	1:30
Dial 4	4:00	Dial 4	2:00
Dial 5	5:00	Dial 5	2:30
Dial 6	6:00	Dial 6	3:00
Dial 7	7:00	Dial 7	3:30
Dial 8	8:00	Dial 8	4:00
Dial 9	9:00	Dial 9	4:30

## Expected Results

By following the instruction for use, a single HistoCore SPECTRA H&E Staining System S1 should provide at least 1,600 stained slides of diagnostically consistent quality.

## Recommended Controls

Routine control slides containing tissue fixed and processed in a similar manner to the test specimens should be included to ensure reagents are performing and functioning properly.

# HistoCore SPECTRA H&E Staining System S1

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Leica Biosystems Richmond, Inc.