KREATECH™ FISH PROBES

- EN** - Instructions for use
- FR** - Mode d'emploi
- DE** - Gebrauchsanleitung
- IT** - Istruzioni per l'uso
- ES** - Instrucciones de uso
- PT** - Instruções de Utilização
- TR** - Kullanım talimatları

Table of Contents

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Instructions for use

Formalin-fixed paraffin-embedded tissue

Using Kreatech™ FISH probes

Fluorescent *in situ* hybridization (FISH) is used to identify, or label, target genomic sequences so that their location can be studied. DNA sequences from appropriate, chromosome specific probes are first labeled with reporter molecules. The labeled DNA probe is then hybridized to the target DNA in the nuclei of the tissue. After washing, the specimen is screened for the reporter molecules by fluorescence microscopy.

KREATECH REPEAT-FREE FISH probes do not contain Cot-1 DNA. Hybridization efficiency is therefore increased and background, due to unspecific binding, is greatly reduced.

For optimal results on paraffin embedded tissue sections it is advised to use Kreatech's complete Pre-treatment kits (KBI-60004, KBI-60007) which include a specially optimized "Instruction for Use" protocol.

For more info consult our website: <http://www.leicabiosystems.com/ihc-ish/kreatech-fish-probes/>

As an alternative the following protocols can also be used:

Slide pretreatment:

1. Mount 4 - 6 µm formalin-fixed paraffin-embedded (FFPE) tissue sections on positively charged slides.
2. Bake mounted FFPE tissue sections for 2 hours at 80 °C or 16 hours at 56 °C.

3. De-paraffinize slides in xylene or xylene substitute, incubate for 2 x 10 minutes (min) at room temperature (RT).
4. Re-hydrate slides in 100%, 85% and 70% ethanol, incubate for 3 min each at RT.
5. Place slides in deionized H₂O (dH₂O), incubate for 3 min at RT. Proceed with Protocol I or Protocol II.

Protocol I (standard protocol, if not satisfactory use protocol II)

1. Pre-warm 0.01 M sodium citrate pH 6.0 to 96 - 98 °C.
2. Pre-warm 0.01 M HCl to 37 °C (if not using Ready-to-Use (RtU) pepsin (LK-101A)).
3. Place slides in 0.01 M sodium citrate pH 6.0 at 96 - 98 °C, incubate for 15 min.
4. Place slides in dH₂O, incubate for 2 min at RT.
5. Add pepsin to pre-warmed 0.01 M HCl to reach a final concentration of 0.025%. (if not using RtU pepsin (LK-101A)).
6. Digest slides in 0.025% pepsin in 0.01 M HCl at 37 °C or cover tissue in RtU pepsin (LK-101A) at RT, incubate for 5 - 45 min, (time depending on tissue fixation and tissue type).
7. Place slides in dH₂O, incubate for 1 min at RT.
8. Place slides in 1 x PBS, pH 7.4 or 2 x SSC, pH 7.0, incubate for 5 min at RT.
9. Dehydrate slides in 70%, 85%, and 100% ethanol, incubate for 1 min each at RT. Air-dry at RT.

Proceed with probe preparation.

Protocol II (for heavily cross-linked samples)

1. Pre-warm 8% sodium thiocyanate to 80 °C.
2. Pre-treat slides with 0.2 M HCl, incubate for 20 min at RT.
3. Place slides in dH₂O, incubate for 3 min at RT.
4. Pre-warm 0.01 M HCl to 37 °C (if not using RtU pepsin (LK-101A)).
5. Place slides in pre-warmed 8% sodium thiocyanate in dH₂O at 80 °C, incubate for 30 min.

6. Place slides in 2 x SSC, pH 7.0, incubate for 3 min at RT.
7. Add pepsin to pre-warmed 0.01 M HCl to a final concentration of 0.025% (if not using RtU pepsin (LK-101A)).
8. Digest slides in 0.025% pepsin in 0.01 M HCl at 37 °C or cover tissue in RtU pepsin (LK-101A) at RT, incubate for 5 - 45 min, (time depending on tissue fixation and tissue type).
9. Place slides in dH₂O, incubate for 1 min at RT.
10. Place slides in 1 x PBS, pH 7.4 or 2 x SSC, pH 7.0, incubate for 5 min at RT.
11. Dehydrate slides in 70%, 85%, and 100% ethanol, incubate for 1 min each at RT. Air-dry at RT.

Proceed with probe preparation.

Probe preparation:

Briefly spin down probe vial, vortex probe vial and briefly spin down again before use. Allow probe to reach RT before use.

Kreatech FISH probes are supplied RtU unless specified otherwise in the product documentation. Consult label on vial and specific probe pack insert for dilution specifics.

Co-denaturation:

1. Apply 10 µl of probe or probe-mix per 22 x 22 mm field.
2. Cover with glass coverslip and seal with Fixogum (LK-071A) or rubber cement.
3. Denature sample and probe on a ThermoBrite (TS-01/02) for 5 min at 80 °C.

Hybridization:

Incubate overnight at 37 °C in a ThermoBrite (TS-01/02) or in a humidified chamber.

Post-hybridization wash:

1. Pre-warm Wash Buffer I (0.4 x SSC / 0.3% Igepal) (LK-102A) to 72 °C.
2. Remove rubber cement.

3. Place up to 14 slides in 200 ml of Wash Buffer II (2 x SSC / 0.1% Igepal) (LK-103A), incubate for 2 min at RT. Slide off coverslips. Re-use only once for a total of 28 slides.
4. Place up to 14 slides in 200 ml of pre-warmed Wash Buffer I (0.4 x SSC / 0.3% Igepal) (LK-102A), incubate for 2 min at 72 °C (± 1 °C) without agitation. Re-use only once for a total of 28 slides.
5. Place up to 14 slides in 200 ml of fresh Wash Buffer II (2 x SSC / 0.1% Igepal) (LK-103A), incubate for 1 min at RT without agitation. Re-use only once for a total of 28 slides.
6. Dehydrate in fresh 70%, 85% and 100% ethanol, incubate for 1 min each at RT. Air-dry at RT and proceed to Counterstaining.

Counterstaining:

Apply 15 μ l DAPI counterstain (LK-095A (0.1 μ g/ml) or LK-096A (1 μ g/ml)) and apply glass coverslip. DAPI can be diluted in counterstain diluent (LK-097A) to obtain desired concentration. Place slides in the dark and allow 10 - 15 min for counterstain to develop. Proceed with microscopy.

Recommendations for fluorescence microscopy:

For optimal visualization use a well maintained and regularly calibrated microscope equipped with a 100W mercury lamp or other appropriate light source and a 63x or 100x fluorescence objective. Triple band-pass filters (DAPI/FITC/Cy3 or DAPI/FITC/TRITC) are used to view multiple colors, single band-pass filters are used for individual color visualization.

Suitable excitation and emission range for Kreatech fluorophores:

Fluorophore	Excitation	Emission
DAPI	360 ±20 nm	460 ±30 nm
PlatinumBright™ 415	429 ±15 nm	470 ±15 nm
PlatinumBright™ 495	495 ±10 nm	517 ±10 nm
PlatinumBright™ 550	550 ±10 nm	580 ±10 nm

Troubleshooting

Check protein digestion and pre-treatment by applying 15 µl DAPI counterstain and evaluate slides using a fluorescence microscope equipped with a DAPI filter. A 15 min protein digestion is normally sufficient for a wide range of breast tumors. Remove coverslip and soak tissue in 2 x SSC, pH 7.4 for 2 min at RT. Prolong protein digestion by 2 - 20 min if sample is not sufficiently digested. Use a fresh sample and reduce protein digestion time if the sample is over-digested.

Alternative protocol for probe denaturation (separate probe and sample denaturation):

1. Denature slide in 70% formamide / 2 x SSC, pH 7.0 at 72 °C (±1 °C) for 2 min.
2. Dehydrate in ice cold (-20 °C) 70%, 85%, and 100% ethanol for 2 min each. Air-dry.
3. Denature probe mix at 90 °C for 10 min and then place on ice.
4. Briefly spin down probe vial, vortex probe vial and then spin down probe vial again.
5. Apply probe to denatured slide, cover with glass coverslip, seal with rubber cement and continue with Hybridization.

Frequently asked questions

I have weak or no signals

- Re-hybridize making sure that the probe has been mixed correctly

and that it is at room temperature before use.

- Re-hybridize making sure that the stringency wash (Wash Buffer I) is at the right temperature (72 °C (± 1 °C)).
- Account for temperature drop when adding slides to pre-warmed reagents.
- Tissue digestion not optimal (under-digested) – soak off cover slip in 2 x SSC, pH 7.0 for 2 min at RT. Prolong protein digestion by 2 - 20 min.
- Tissue digestion not optimal (over-digested) – Discard slide and start with fresh section. Reduce protein digestion by at least 10 min.
- Use a minimum of 10 μ l of probe per 22 x 22 mm coverslip.
- Check microscope filters and light source are correct and in full working order.

I have high background.

- Tissue digestion not optimal. Under-digested tissue increase background. Soak off cover slip in 2 x SSC, pH 7.0 for 2 min at RT. Prolong protein digestion by 2 - 20 min.
- Increase the temperature of the stringency wash (Wash Buffer I) when washing slides.

I have cross-hybridization when using Centromeric probes.

- Increase the temperature of the stringency wash (Wash Buffer I) when washing slides.

Minimum probe application:

Glass cover slip size	Minimum probe application
22 x 22 mm	10 μ l
22 x 32 mm	15 μ l
22 x 50 mm	23 μ l

Procedural recommendations:

Temperature and buffer concentration (stringency) of hybridization and washing are important, as lower stringency can result in non-specific binding of the probe to other sequences, and higher stringency can result in a lack of signal. Incomplete denaturation of target DNA and/or probe DNA can result in lack of signal.

Material required, but not supplied: Reagents:

1. Xylene
2. Formamide
3. Ethanol 100%, 85% and 70%
4. 0.01 M sodium citrate pH 6.0 or 8% sodium thiocyanate
5. 0.01 M HCl and 0.2M HCl
6. Pepsin Solution RtU (LK-101A)
7. 1 x PBS, pH 7.4
8. 2 x SSC, pH 7.0
9. Wash Buffer I (0.4 x SSC / 0.3% Igepal) (LK-102A)
10. Wash Buffer II (2 x SSC / 0.1% Igepal) (LK-103A)
11. DAPI counterstain (LK-095A (0.1 µg/ml) or LK-096A (1 µg/ml))
12. Counterstain diluent (LK-097A)
13. Fixogum (LK-071A) or rubber cement

Material required, but not supplied: equipment:

1. ThermoBrite (TS-01/02) or ThermoBrite Elite
(See www.Leicabiosystems.com for more information)
2. Incubator at 37 °C
3. Water bath with accurate temperature from 37 °C to 90 °C
4. Plastic or glass Coplin jar
5. Variable micropipettes (1 µl - 200 µl)
6. Fluorescence microscope equipped with suitable filters (see recommendations for fluorescence microscopy)

Warnings and precautions:

1. For *in vitro* use only. For professional use only. In case of emergencies check SDS sheets for safety information.
2. DNA probes and hybridization buffers contain formamide which is a teratogen; do not inhale or allow skin contact. Wear gloves and a lab coat when handling DNA probes and DAPI counterstain.
3. All materials should be disposed of according to your institution's guidelines for hospital waste disposal.
4. Labelling according Regulation (EC) No 1272/2008

Pictogram

Signal word: Danger

Hazard statement(s)

H351 Suspected of causing cancer.

H360D May damage the unborn child.

H373 May cause damage to organs (Blood) through prolonged or repeated exposure if swallowed.

Precautionary statement(s)

P201 Obtain special instructions before use.

P281 Use personal protective equipment as required.

P308 + P313 IF exposed or concerned: Get medical advice/ attention.

Supplemental Hazard Statements: none

Restricted to professional users

Patents: *These products or the use of these products is subject to proprietary rights. The probes in these products are produced using the Janssen Diagnostics, LLC REPEAT-FREE™ technology and labeled with*

the Universal Linkage System (ULS™).

The fluorophore used in the PlatinumBright - 415 labeling compound is subject to patents, owned or controlled, and manufactured by DYOMICS GmbH. US and International patents pending for Janssen Diagnostics, LLC REPEAT-FREE™ technology.

The ULS™ technology and products are covered by one or more of the following US patents 5,580,990; 5,714,327; 5,985,566; 6,133,038; US RE 40,557E, 6,797,818; 7,217,813. KREATECH is a trade name of KREATECH Biotechnology BV. ULS™, PlatinumBright™, are trademarks of KREATECH Biotechnology BV. REPEAT-FREE™ is a trademark of Janssen Diagnostics, LLC.

ThermoBrite™ is a registered trademark of Leica Biosystems.

Cell Preparations

Using Kreatech™ FISH probes

Fluorescent *in situ* hybridization (FISH) is used to identify, or label, target genomic sequences so that their location can be studied. DNA sequences from appropriate, chromosome specific probes are first labeled with reporter molecules. The labeled DNA probe is then hybridized to the metaphase chromosomes or interphase nuclei on a slide. After washing, the specimen is screened for the reporter molecules by fluorescence microscopy.

KREATECH REPEAT-FREE FISH probes do not contain Cot-1 DNA. Hybridization efficiency is therefore increased and background, due to unspecific binding, is greatly reduced.

For use on metaphase and interphase cells from peripheral blood and bone marrow cultures or direct preparations prepared by standard cytogenetic methods, see: The AGT genetics laboratory manual, 3rd ed. New York: Raven Press; 1996.

For optimal results it is advised to use Kreatech's complete Pre-treatment kits (KBI-60005, KBI-60006) which include a specially optimized "Instruction for Use" protocol.

For more info consult our website: <http://www.leicabiosystems.com/ihc-ish/kreatech-fish-probes/>

Pretreatment:

Protocol I (standard protocol, if not satisfactory use protocol II)

1. To age the slides, bake the slides for 10 minutes (min) at 80 °C, then allow slides to cool down to room temperature (RT).
2. Place slides in freshly prepared 70% acetic acid (in deionized H₂O)

(dH₂O)), incubate for 1 min at RT.

3. Place slides in 1 x PBS, pH 7.4, incubate for 30 seconds (s) at RT.
4. Place slides in 1 x PBS, pH7.4, incubate for 5 min at RT.
5. Dehydrate slide in 70%, 85% and 100% ethanol, incubate for 1 min each at RT. Air-dry at RT.

Proceed with Probe preparation.

Protocol II (for slides with high cytoplasmic background or direct preps (amniotic fluid or chronic villi))

1. Pre-warm 2 x SSC, pH 7.0 and 0.01 M HCl at 37 °C
2. Place dry sample slides in pre-warmed 2 x SSC, pH 7.0, incubate at 37 °C for 2 min.
3. Add pepsin to the pre-warmed 0.01 M HCl to a final concentration of 0.005%.
4. Place slides in 0.005% pepsin solution in 0.01 M HCl and incubate at 37 °C for 5 - 15 min (depending on the amount of cytoplasmic background).
5. Place slides in 1 x PBS, pH 7.4, incubate for 3 min at RT.
6. Post-fix slides by incubating in 1% buffered formaldehyde in 1 x PBS / 20 mM MgCl₂ for 10 min at RT.
7. Place slides in 1 x PBS, pH 7.4, incubate for 3 min at RT.
8. Dehydrate slide in 70%, 85% and 100% ethanol, incubate for 1 min each at RT. Air-dry at RT.

Proceed with Probe preparation.

Probe preparation:

Briefly spin down probe vial, vortex probe vial and briefly spin down again before use. Allow probe to reach RT before use.

Most Kreatech FISH probes are supplied Ready-to-Use (RtU). Centromere, Subtelomere and Whole chromosome FISH probes are provided in 5 x concentrated formats and must be diluted with the provided hybridization buffer and/or other concentrated probes. Consult label on vial and specific probe package insert for dilution specifics.

Co-denaturation:

1. Apply 10 μ l of probe or probe-mix per 22 x 22 mm field.
2. Cover with glass coverslip and seal with Fixogum (LK-071A) or rubber cement.
3. Denature sample and probe on a ThermoBrite (TS-01/02) for 5 min at 75 °C.

Hybridization:

Incubate overnight at 37 °C in a ThermoBrite (TS-01/02) or in a humidified chamber.

Post-hybridization wash:

1. Pre-warm Wash Buffer I (0.4 x SSC / 0.3% Igepal) (LK-102A) to 72 °C.
2. Remove rubber cement.
3. Place up to 14 slides in 200 ml of Wash Buffer II (2 x SSC / 0.1% Igepal) (LK-103A), incubate for 2 min at RT. Slide off coverslips. Re-use only once for a total of 28 slides.
4. Place up to 14 slides in 200 ml of pre-warmed Wash Buffer I (0.4 x SSC / 0.3% Igepal) (LK-102A), incubate for 2 min at 72 °C (± 1 °C) without agitation. Re-use only once for a total of 28 slides.
5. Place up to 14 slides in 200 ml of fresh Wash Buffer II (2 x SSC / 0.1% Igepal) (LK-103A), incubate for 1 min at RT without agitation. Re-use only once for a total of 28 slides.
6. Dehydrate slides in 70%, 85% and 100% ethanol, incubate for 1 min each at RT. Air-dry at RT.

Proceed to Counterstaining.

Counterstaining:

Apply 15 μ l DAPI counterstain (LK-095A (0.1 μ g/ml) or LK-096A (1 μ g/ml)) and apply glass coverslip. DAPI can be diluted in counterstain diluent (LK-097A) to obtain desired concentration. Place slides in the dark and allow 10 - 15 min for counterstain to develop. Proceed with microscopy.

Recommendations for fluorescence microscopy:

For optimal visualization use a well maintained and regularly calibrated microscope equipped with a 100 W mercury lamp or other appropriate light source and a 63x or 100x fluorescence objective. Triple band-pass filters (DAPI/FITC/Cy3 or DAPI/FITC/TRITC) are used to view multiple colors, single band-pass filters are used for individual color visualization.

Suitable excitation and emission range for Kreatech FISH fluorophores:

Fluorophore	Excitation	Emission
DAPI	360 ±20 nm	460 ±30 nm
PlatinumBright™ 415	429 ±15 nm	470 ±15 nm
PlatinumBright™ 495	495 ±10 nm	517 ±10 nm
PlatinumBright™ 550	550 ±10 nm	580 ±10 nm

Troubleshooting:

Alternative protocol for cytology specimens:

1. Pre-warm 50 ml of 2 x SSC / 0.5% Igepal (LK-105B) in a Coplin jar to 37 °C in a water bath. Place prepared slides in the Coplin jar and incubate for 15 minutes.
2. Dehydrate slides in 70%, 85%, and 100% ethanol, incubate for 1 min each at RT. Air-dry at RT.

Proceed with Probe preparation.

Alternative protocol for probe denaturation (separate probe and sample denaturation):

1. Prewarm 70% formamide / 2 x SSC, pH 7.0 to 72 °C.
2. Denature slides in 70% formamide / 2 x SSC, pH 7.0, incubate at 72 °C (±1 °C) for 2 min.
3. Dehydrate slides in ice cold (-20 °C) 70%, 85%, and 100% ethanol, incubate for 2 min each. Air-dry at RT.
4. Denature probe mix at 90 °C for 10 min and directly place on ice.
5. Briefly spin down probe vial, vortex probe vial and then spin down

probe vial again.

6. Apply probe to denatured slide, cover with glass coverslip, seal with rubber cement.

Proceed with Hybridization.

Frequently asked questions

I have weak or no signals

- Re-hybridize making sure that the probe has been mixed correctly and that is at room temperature before use.
- Re-hybridize making sure that the stringency wash (Wash Buffer I) is at the right temperature (72 °C (± 1 °C)).
- Account for temperature drop when adding slides to pre-warmed reagents.
- Use a minimum of 10 μ l of probe per 22 x 22 mm coverslip.
- Check microscope filters and light source are correct and in full working order.

I have high background or a-specific cross-hybridization when using Centromeric or Subtelomeric probes.

- Increase the temperature of the stringency wash (Wash Buffer I) when washing slides.

Minimum probe application:

Glass cover slip size	Minimum probe application
22 x 22 mm	10 μ l
22 x 32 mm	15 μ l
22 x 50 mm	23 μ l

Procedural recommendations:

Temperature and buffer concentration (stringency) of hybridization and washing are important, as lower stringency can result in non-specific binding of the probe to other sequences, and higher stringency can result in a lack of signal. Incomplete denaturation of target DNA

and/or probe DNA can result in lack of signal.

Material required, but not supplied: reagents:

1. FISH reagent Kit (KBI-60005)
2. FISH digestion Kit (KBI-60006)
3. 1% buffered formaldehyde / 1 x PBS / 20 mM MgCl₂
4. 1 x PBS, pH 7.4
5. 2 x SSC, pH 7.0
6. Pepsin solution (LK-101A)
7. Wash buffer I (0.4 x SSC / 0.3% Igepal) (LK-102A)
8. Wash buffer II (2 x SSC / 0.1% Igepal) (LK-103A)
9. Igepal
10. 0.01 M HCl
11. Carnoy's fixative (methanol : acetic acid = 3 : 1)
12. 70% acetic acid (in deionized H₂O)
13. DAPI counterstain (LK-095A (0.1 µg/ml) or LK-096A (1 µg/ml))
14. Counterstain diluent (LK-097A)
15. Formamide
16. Ethanol 100%, 85% and 70%
17. Fixogum (LK-071A) or rubber cement

Material required, but not supplied: equipment:

1. ThermoBrite (TS-01/02).
2. Incubator at 37 °C
3. Water bath with accurate temperature from 37 °C to 90 °C
4. Plastic or glass Coplin jar
5. Variable micropipettes (1 µl - 200 µl)
6. Fluorescence microscope equipped with suitable filters (see Recommendations for fluorescence microscopy)

Warnings and precautions:

1. For *in vitro* use only. For professional use only. In case of emergencies check SDS sheets for safety information.
2. DNA probes and hybridization buffers contain formamide which is

- a teratogen; do not inhale or allow skin contact. Wear gloves and a lab coat when handling DNA probes and DAPI counterstain.
3. All materials should be disposed of according to your institution's guidelines for hospital waste disposal.
 4. Labelling according Regulation (EC) No 1272/2008

Pictogram

Signal word: Danger

Hazard statement(s)

H351 Suspected of causing cancer.

H360D May damage the unborn child.

H373 May cause damage to organs (Blood) through prolonged or repeated exposure if swallowed.

Precautionary statement(s)

P201 Obtain special instructions before use.

P281 Use personal protective equipment as required.

P308 + P313 IF exposed or concerned: Get medical advice/ attention.

Supplemental Hazard Statements: none

Restricted to professional users

Patents: *These products or the use of these products is subject to proprietary rights. The probes in these products are produced using the Janssen Diagnostics, LLC REPEAT-FREE™ technology and labeled with the Universal Linkage System (ULS™).*

The fluorophore used in the PlatinumBright – 415 labeling compound is subject to patents, owned or controlled, and manufactured by DYOMICS GmbH. US and International patents pending for Janssen Diagnostics,

LLC REPEAT-FREE™ technology.

The ULS™ technology and products are covered by one or more of the following US patents 5,580,990; 5,714,327; 5,985,566; 6,133,038; US RE 40,557E, 6,797,818; 7,217,813. KREATECH is a trade name of KREATECH Biotechnology BV. ULS™, PlatinumBright™, are trademarks of KREATECH Biotechnology BV. REPEAT-FREE™ is a trademark of Janssen Diagnostics, LLC.

ThermoBrite™ is a registered trademark of Leica Biosystems

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Kreatech Biotechnology B.V.