# Leica Biosystems Eosins.

English.....	2
العربية (Arabic).....	6
简体中文 (Chinese Simplified).....	10
繁體中文 (Chinese Traditional).....	14
Dansk (Danish).....	18
Nederlands (Dutch).....	22
Français (French – Canada).....	26
Français (French – France).....	30
Deutsch (German).....	34
Italiano (Italian).....	38
日本語 (Japanese).....	41
한국어 (Korean).....	46
Norsk (Norwegian).....	50
Polski (Polish).....	54
Português (Portuguese – Brazil).....	58
Português (Portuguese – Portugal).....	62
Română (Romanian).....	66
Русский (Russian).....	70
Slovenščina (Slovenian).....	74
Español (Spanish – Central America).....	78
Español (Spanish – Spain).....	82
Svenska (Swedish).....	86
ภาษาไทย (Thai).....	90
Türkçe (Turkish).....	94
Tiếng Việt (Vietnamese).....	98

# Eosins

**REF 3801590BBE, 3801590E, 3801591E, 3801592E, 3801600, 3801600E, 3801600BBE, 3801601E, 3801602, 3801602E, 3801606, 3801610, 3801611, 3801615, 3801616, 3801619**

**Product Name**

Leica Biosystems Eosins.

**Intended Use**

**Detection/Measurement**

The Leica Biosystems Eosins do not detect or measure an analyte or marker. The Leica Biosystems Eosins are used to provide coloration and contrast to histologic, cytologic and fresh frozen specimens mounted on microscope slides. Eosins primarily are used to provide contrast of the cell cytoplasm, muscle tissue and connective tissue elements with cell nuclei when using nuclear stain such as hematoxylin. The Leica Biosystems Eosins stain cell cytoplasm, muscle tissue and connective tissue elements varying shades of pink, red and orange.

**Product Function**

The results obtained through use of the Leica Biosystems Eosins do not provide objective medical evidence. The coloration and contrast the Leica Biosystems Eosins provide to histologic and cytologic specimens allows visualization of microscopic anatomy. This visualization, when interpreted by a trained professional, is utilized alongside other information such as the patient’s medical history, physical condition, as well as results from other medical testing to render a medical diagnosis.

**Specific Information Provided**

Leica Biosystems Eosins are not intended for the detection, definition or differentiation of a specific disorder, condition or risk factor. The staining demonstrated with use of these products, when used as intended, provides trained professionals information which may define the physiological or pathological state of the tissue specimen.

**Automation**

The Leica Biosystems Eosins are not automated but can be used on automated staining platforms. Use on an automated platform should be validated at the point of use.

**Qualitative/Quantitative**

The Leica Biosystems Eosins are qualitative stains.

**Specimen Type**

The Leica Biosystems Eosins may be used with fixed or fresh histologic and cytologic specimens.

**Testing Population**

The Leica Biosystems Eosins are intended for use with any patient requiring evaluation of biopsy or resection tissue as well as cytologic preparations for the assessment of a suspected pathology or disease.

**Intended User**

The Leica Biosystems Eosins are intended for use by qualified laboratory personnel and/or designee of the laboratory.

**In Vitro Diagnostic**

The Leica Biosystems Eosins are intended for *in vitro* diagnostics use only.

**Intended User**

The Leica Biosystems Eosins are intended for use by qualified laboratory personnel and/or designee.

**Test Principle**

The Leica Biosystems Eosins shall stain cytoplasm, muscle and connective tissue various shades red, pink and orange. The Eosin shall provide contrast to nuclear stains such as hematoxylin.

**Calibrators & Controls**

The Leica Eosins do not require the use of any calibrators or controls.

**Reagent Limitations**

No reagent limitations are applicable to these products.

**Applicable Products**

Product Code	Material Description
3801590BBE	Surgipath Aqueous Eosin 1% (5L bag in a box)
3801590E	Surgipath Aqueous Eosin 1% (5L)
3801591E	Surgipath Aqueous Eosin 1% (2.5L)
3801592E	Surgipath Aqueous Eosin 1% (1L)
3801600	Surgipath Eosin (1gal)
3801600E	Surgipath Eosin (5L)
3801600BBE	Surgipath Eosin (5L bag in a box)
3801601E	Surgipath Eosin (2.5L)
3801602	Surgipath Eosin (1qt)
3801602E	Surgipath Eosin (1L)
3801606	SelecTech Eosin 515 Phloxine (500mL)
3801610	SelecTech Eosin Trichrome (500mL)

# Eosins

**REF 3801590BBE, 3801590E, 3801591E, 3801592E, 3801600, 3801600E, 3801600BBE, 3801601E  
3801602, 3801602E, 3801606, 3801610, 3801611, 3801615, 3801616, 3801619**

3801611	SelecTech Eosin Trichrome (4-500mL)
3801615	SelecTech Alcoholic Eosin Y 515 (500mL)
3801616	SelecTech Alcoholic Eosin Y 515 Y (4-500mL)
3801619	SelecTech Eosin 515 LT (4-500mL)

NOTE: Products listed here may not be available in all regions.

## Materials Not Included

The Leica Biosystems Eosins are designed to be used as part of a Hematoxylin & Eosin (H&E) stain protocol which require the use of graded alcohols, xylene or xylene substitutes, hematoxylin, differentiator, and bluing agent.

## Devices Required

Leica Biosystems Eosin may be used on an open automated staining platform or with a manual staining method and should be validated at the point of use by the user.

## Storage and Stability

The product shall be stable for 24 months postproduction when stored at ambient temperature.

Store reagents at room temperature (15-30°C) in a well-ventilated place.

**CAUTION:** Do not use after the expiration date.

## In Use Stability

User discretion should be utilized when determining in-use stability.

## Sterility

The Leica Biosystems Eosins are not sterile products.

## Warnings/Precautions

This product and protocol(s) associated with the product, whether provided by Leica Biosystems in this instruction for use or developed by the user, shall be validated at the point of use by the user.

## Infectious Material Status

The Leica Biosystems Eosins do not include any infectious material. However, specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

## Special Facilities

The Leica Biosystems Eosins should be used per facility guidelines.

## Specimen Handling

Suggested fixatives include 10% neutral buffered formalin. Routine dehydration, clearing, and paraffin infiltration and embedding, and routine preparation of microtome sections. Poor fixation, processing, rehydration, and sectioning will adversely affect the staining quality. Sections of 2-5micron thickness are recommended.

## Preparation for Use

- **SelecTech Alcoholic Eosin Y 515** - most commonly used cytoplasmic stain for hematoxylin and eosin staining. Stains cytoplasmic components: erythrocytes, collagen, muscle and epithelial cells with three varying shades of red/pink. Prominent nucleoli may be defined in an orange hue when tissue is properly fixed. Recommended staining time is from 30-90 seconds.
- **SelecTech Eosin 515 LT** - produces a relatively subtle stain of the cytoplasm. Designed to provide contrast between cytoplasm and nuclear chromatin. Designed to work specifically with Hematoxylin 560 MX, Eosin 515 LT enhances darker nuclear staining. Recommended staining time is from 30-90 seconds.
- **SelecTech Eosin Phloxine 515** - Alcoholic Eosin Y with addition of Phloxine, formulated to produce contrast with 560 and 560MX Hematoxylin. Phloxine provides a high vibrancy, more reddish staining and greater range of pink hues of connective tissue and cytoplasm. Recommended staining time is from 30-90 seconds.
- **SelecTech Eosin Trichrome 515** - Generates three distinct colors: stains muscle a vibrant red-pink, connective tissue orange-pink and erythrocytes red. Designed to provide contrast between cytoplasm, nuclear chromatin and connective tissue components. Recommended exposure time is from 30-90 seconds. Post rinsing in 70% or 80% alcohol may increase staining intensity of muscle fibers.
- **Surgipath Eosin** – Eosin Y is an alcohol-based, ready-to-use secondary stain used primarily in histology for routine hematoxylin and eosin staining.
- **Surgipath Aqueous Eosin 1%** - produces lighter shades of pink, red and orange. Convenient, water-based Eosin can be used in histology for routine hematoxylin and eosin staining.

# Eosins

**REF 3801590BBE, 3801590E, 3801591E, 3801592E, 3801600, 3801600E, 3801600BBE, 3801601E, 3801602, 3801602E, 3801606, 3801610, 3801611, 3801615, 3801616, 3801619**

**Protocol Set-up when using Alcoholic Eosin (Table 1):**

1. De-paraffinize and rehydrate section as appropriate using routine procedure.
2. Stain in hematoxylin following manufacturer's protocol.
3. Wash in running tap water.
4. Differentiate in a proper acid alcohol, wash in running tap water.
5. Immerse in bluing agent, wash in running tap water.
6. Dehydrate in reagent alcohol
7. Counterstain with selected alcoholic eosin stain until desired intensity is reached.
8. Dehydrate, clear and mount in mounting medium.

**Results:** The nuclei and their components should be stained in varying shades of blue/purple. Muscle should be stained dark pink, erythrocytes are orange/pink, and connective tissue light pink.

**Table 1. Example of H&E Staining Protocol with Alcoholic Eosin**

Steps	Action	Chemical	Time (mm:ss)
1	Deparaffinize	Xylene	3:00
2	Deparaffinize	Xylene	3:00
3	Deparaffinize	Xylene	3:00
4	Hydration	100% Alcohol	2:00
5	Hydration	100% Alcohol	1:00
6	Hydration	100% Alcohol	1:00
7	Hydration	80% or 95% Alcohol	1:00
8	Hydration	Water Wash	1:00
9	Stain	<b>Progressive Hematoxylin</b>	1:00 to 5:00
10	Wash	Water Wash	3:00
11	Differentiation	<b>Differentiator</b>	0:30 to 1:30
12	Wash	Water Wash	1:00
13	Bluing	<b>Bluing Buffer</b>	0:30 to 1:00
14	Wash	Water Wash	2:00
15	Dehydration	80% to 95% Alcohol	1:00
16	Counterstaining	<b>Alcoholic Eosin</b>	0:30 to 1:30
17	Wash	Water Wash	2:00
18	Dehydration	95% to 100% Alcohol	1:00
19	Dehydration	100% Alcohol	1:00
20	Dehydration	100% Alcohol	1:00
21	Clearing	Xylene	2:00
22	Clearing	Xylene	2:00
23	Clearing	Xylene	2:00

**Protocol Set-up when using Aqueous Eosin 1% (Table 2):**

1. De-paraffinize and rehydrate section as appropriate using routine procedure.
2. Stain in hematoxylin following manufacturer's protocol.
3. Wash in running tap water.
4. Differentiate in a proper acid alcohol, wash in running tap water.
5. Immerse in bluing agent, wash in running tap water.
6. Counterstain in Aqueous Eosin 1% for 30 seconds or until desired intensity is reached.
7. Wash in running tap water.
8. Dehydrate, clear and mount in mounting medium.

**Results:** The nuclei and their components should now be stained in varying shades of blue/purple. Erythrocytes are orange/pink, muscle is pink, and connective tissue is light pink.

# Eosins

**REF 3801590BBE, 3801590E, 3801591E, 3801592E, 3801600, 3801600E, 3801600BBE, 3801601E, 3801602, 3801602E, 3801606, 3801610, 3801611, 3801615, 3801616, 3801619**

**Table 2. Example of H&E Staining Protocol with Aqueous Eosin.**

Steps	Action	Chemical	Time (mm:ss)
1	Deparaffinize	Xylene	3:00
2	Deparaffinize	Xylene	3:00
3	Deparaffinize	Xylene	3:00
4	Hydration	100% Alcohol	2:00
5	Hydration	100% Alcohol	1:00
6	Hydration	100% Alcohol	1:00
7	Hydration	80% or 95% Alcohol	1:00
8	Hydration	Water Wash	1:00
9	Stain	Progressive Hematoxylin	1:00 to 5:00
10	Wash	Water Wash	3:00
11	Differentiation	Differentiator	0:30 to 1:30
12	Wash	Water Wash	1:00
13	Bluing	Bluing Buffer	0:30 to 1:00
14	Wash	Water Wash	2:00
15	Counterstaining	Aqueous Eosin	0:30
16	Wash	Water Wash	2:00
17	Dehydration	70% to 100% Alcohol	1:00
18	Dehydration	100% Alcohol	1:00
19	Dehydration	100% Alcohol	1:00
20	Clearing	Xylene	2:00
21	Clearing	Xylene	2:00
22	Clearing	Xylene	2:00

## Readiness for Use

Once appropriate staining protocol is chosen and bath layout is created, pour all the reagent into the reagent vessel. Place the reagent vessel back into the respective station.

## Quality Control

A routine quality control slide(s) containing tissue fixed and processed in a similar manner to the test specimens should be performed prior to routine use to ensure Eosin is performing as intended.

## Expected Results

By following the instruction for use, Eosin shall stain cytoplasm, muscle and connective tissue various shades red, pink and orange.

## Analytical Performance

The Leica Biosystems Eosins are not used to detect a specific analyte or marker. These products are used in conjunction with other products in a Hematoxylin & Eosin staining protocol system to stain cell nuclei blue and connective tissue, cytoplasm, muscle and erythrocytes various shades of orange, pink and red. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

## Clinical Performance

The Leica Biosystems Eosins are not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the Leica Biosystems Eosins in a clinical setting.

## Disposal

Eosin should be disposed in accordance with local governing regulations.

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