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BOND™ Ready-To-Use Primary Antibody Cytokeratin 7 (RN7)

Catalog No: PA0942

Intended Use

This reagent is for *in vitro* diagnostic use.

Cytokeratin 7 (RN7) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human cytokeratin 7 intermediate filament protein in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Cytokeratin 7 (RN7) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human cytokeratin 7 intermediate filament protein is achieved by first, allowing the binding of Cytokeratin 7 (RN7) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

Cytokeratin 7 (RN7) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

RN7.

Immunogen

Prokaryotic recombinant protein corresponding to part of the C-terminal region of the cytokeratin 7 intermediate filament molecule.

Specificity

Human cytokeratin 7 intermediate filament protein.

Ig Class

IgG1.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 0.4 mg/L as determined by ELISA.

Dilution and Mixing

Cytokeratin 7 (RN7) primary antibody is optimally diluted for use on the BOND system. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system.

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Cytokeratin 7 (RN7) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Cytokeratin 7 (RN7) primary antibody was developed for use on the automated BOND system in combination with BOND Polymer Refine Detection. The recommended staining protocol for Cytokeratin 7 (RN7) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 1 for 20 minutes.

Results Expected

Normal Tissues

Clone RN7 detects cytokeratin 7 in the cytoplasm and membrane of a large number of epithelial cell types; simple, complex and transitional, including ductal and glandular epithelia. It consistently stains endometrium, lung alveoli, bile ducts in liver, bronchiolar epithelium, thyroid, distal convoluted tubules and collecting ducts in kidney. Staining was also seen in the glands of the endocervix, bronchus, breast, salivary, sweat glands of the skin, myoepithelium (variably), ovarian mesothelium, trophoblasts of the placenta, all cell layers of the urothelium and the surface epithelium of the stomach. Non-epithelial tissues such as smooth muscle, connective tissue, lymphoid tissue and blood vessels are consistently negative. (Total number of normal cases evaluated = 85).

Tumor Tissues

Clone RN7 stained 58/118 tumors evaluated, including urothelial carcinomas (20/22), colon adenocarcinomas (1/17), ovarian tumors (7/8, including 4/4 ovarian adenocarcinomas, 2/2 ovarian clear cell carcinomas, 1/1 ovarian endometrioid carcinomas and 0/1 ovarian germ cell tumors), lung carcinomas (5/8, including 4/4 adenocarcinomas of the lung, 1/2 non-small cell carcinomas of the lung and 0/2 squamous cell carcinomas of the lung), thyroid carcinomas (6/6, including 4/4 thyroid papillary carcinomas, 1/1 thyroid follicular carcinomas and 1/1 thyroid medullary carcinomas), breast carcinomas (4/5, including 3/4 ductal carcinomas of the breast and 1/1 lobular carcinomas of the breast), liver carcinomas (2/4), pancreatic adenocarcinomas (3/3), squamous cell carcinomas of the cervix (2/2), metastatic tumors of unknown origin (2/2), adenocarcinomas of the stomach (1/2), adenocarcinomas of the rectum (1/2), skin tumors (1/2), endocervical carcinomas (1/1), endometrioid carcinomas (1/1), intestinal carcinoid tumors (1/1), lymphomas (0/8), renal cell carcinomas (0/4), testicular seminomas (0/4), prostate adenocarcinomas (0/2), brain tumors (0/2), squamous cell carcinomas of the esophagus (0/2), soft tissue tumors (0/2), squamous cell carcinomas of the tongue (0/2), sarcomas (0/2), melanomas (0/2), squamous cell carcinomas of the larynx (0/1) and atypical carcinoid tumors of the thymus (0/1). (Total number of tumor cases evaluated = 118).

Cytokeratin 7 (RN7) is recommended for the detection of human CK7 protein in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.

Product Specific Limitations

Cytokeratin 7 (RN7) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. van de Molengraff FJMM, van Niekerk CC, Jap PHK, et al. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. Histopathology. 1993; 22:35–38.
5. van Niekerk CC, Jap PH, Ramaekers FC, et al. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. Journal of Pathology. 1991; 165(2):145–152.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991;7(9). Référence à commander : M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4ème édition. Churchill Livingstone, New York. 1996.
4. van de Molengraft FJJM, van Niekerk CC, Jap PHK, et al. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. Histopathology 1993;22:35-38.
5. van Niekerk CC, Jap PH, Ramaekers FC, et al. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. Journal of Pathology 1991;165(2):145-152.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. van de Molengraft FJJM, van Niekerk CC, Jap PHK, et al. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. Histopathology. 1993; 22:35–38.
5. van Niekerk CC, Jap PH, Ramaekers FC, et al. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. Journal of Pathology. 1991; 165(2):145–152.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 28. February 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD und Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. van de Molengraff FJJM, van Niekerk CC, Jap PHK, et al. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. Histopathology. 1993; 22:35–38.
5. van Niekerk CC, Jap PH, Ramaekers FC, et al. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. Journal of Pathology. 1991; 165(2):145–152.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. van de Molengraff FJJM, van Niekerk CC, Jap PHK, et al. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. *Histopathology*. 1993; 22:35–38.
5. van Niekerk CC, Jap PH, Ramaekers FC, et al. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. *Journal of Pathology*. 1991; 165(2):145–152.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. van de Molengraaf FJJM, van Niekerk CC, Jap PHK, et al. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. Histopathology. 1993; 22:35–38.
5. van Niekerk CC, Jap PH, Ramaekers FC, et al. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. Journal of Pathology. 1991; 165(2):145–152.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. van de Molengraff FJJM, van Niekerk CC, Jap PHK, et al. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. Histopathology. 1993; 22:35–38.
5. van Niekerk CC, Jap PH, Ramaekers FC, et al. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. Journal of Pathology. 1991; 165(2):145–152.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. van de Molengraff FJMM, van Niekerk CC, Jap PHK, et al. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. *Histopathology*. 1993; 22:35–38.
5. van Niekerk CC, Jap PH, Ramaekers FC, et al. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. *Journal of Pathology*. 1991; 165(2):145–152.

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Kaynakça

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. van de Molengraft FJJM, van Niekerk CC, Jap PHK, et al. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. Histopathology. 1993; 22:35-38.
5. van Niekerk CC, Jap PH, Ramaekers FC, et al. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. Journal of Pathology. 1991; 165(2):145-152.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. van de Molengraaf FJMM, van Niekerk CC, Jap PHK, et al. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. Histopathology. 1993; 22:35–38.
5. van Niekerk CC, Jap PH, Ramaekers FC, et al. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. Journal of Pathology. 1991; 165(2):145–152.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York, 1996.
4. van de Molengraaf FJJM, van Niekerk CC, Jap PHK, et al. OV-TL 12/30 (keratin 7 antibody) is a marker of glandular differentiation in lung cancer. Histopathology. 1993; 22:35–38.
5. van Niekerk CC, Jap PH, Ramaekers FC, et al. Immunohistochemical demonstration of keratin 7 in routinely fixed paraffin-embedded human tissues. Journal of Pathology. 1991; 165(2):145–152.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd