Bond™ Polymer Refine Red Detection

Catalog No: DS9390

Leica Biosystems Newcastle Ltd
Balliol Business Park
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

EN FR IT DE ES PT
SV EL DA NL NO TR

Instructions for Use

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Catalog No: DS9390

Intended Use

This detection system is for in vitro diagnostic use

Bond Polymer Refine Red Detection is a biotin-free, polymeric alkaline phosphatase (AP)-linker antibody conjugate system for the detection of tissue-bound mouse and rabbit IgG and some mouse IgM primary antibodies. It is intended for staining sections of formalin-fixed, paraffin-embedded tissue on the BOND automated system.

The clinical interpretation of any staining or its absence should be complemented by morphological studies. Proper controls should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

The Bond Polymer Refine Red Detection must be used with laboratory best practice in the use of tissue controls. For assurance, laboratories should stain each patient sample in conjunction with positive, negative and other tissue specific controls as needed.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation).

Bond Polymer Refine Red Detection utilizes a novel controlled polymerization technology to prepare polymeric AP-linker antibody conjugates. The detection system avoids the use of streptavidin and biotin, and therefore eliminates nonspecific staining as a result of endogenous biotin.

Bond Polymer Refine Red Detection works as follows:

- A user-supplied specific primary antibody is applied.
- Post Primary IgG linker reagent localizes mouse antibodies.
- Poly-AP IgG reagent localizes rabbit antibodies.
- The substrate chromogen, Fast Red, visualizes the complex via a red precipitate.
- Hematoxylin (blue) counterstaining allows the visualization of cell nuclei.

Using the Bond Polymer Refine Red Detection in combination with the BOND automated system reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

The reagents provided are sufficient for 20 individual BOND staining runs, a maximum of 100 slides.

To achieve a maximum of 100 slides from this detection system, slides must be batched in quantities of 5 or greater per Slide Staining Assembly. Batching in quantities of less than 5 will result in fewer stained slides.

1. Post Primary AP (15 mL) Rabbit anti-mouse IgG (<10 µg/mL) in Tris-buffered saline containing 10% (v/v) animal serum and 0.09% ProClin™ 950.
2. Polymer AP (15 mL) Poly-AP anti-rabbit IgG (<25 µg/mL) in Tris-buffered saline containing 10% (v/v) animal serum and 0.09% ProClin™ 950.
3. Red Part A (4.5 mL) Activator containing 0.5% ProClin™ 950.
4. Red Part B (1.0 mL) Substrate.
5. Red Part C (1.0 mL) Substrate.
6. Red Part D (32 mL) Buffer solution containing 0.5% ProClin™ 950.
7. Hematoxylin (15 mL) <0.1% Hematoxylin.

Dilution and Mixing

Bond Polymer Refine Red Detection is optimized for use on the BOND system.

Reconstitution, mixing, dilution, or titration of these reagents is not required.

Materials Required But Not Provided

Refer to "Using Bond Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system.

Storage and Stability

Store at 2–8 °C. Do not freeze. Do not use after the expiration date indicated on the tray handle label. Return to 2–8 °C immediately after use.

There are no obvious signs to indicate instability of this product; therefore positive and negative controls should be run simultaneously with unknown specimens (refer to "Quality Control" in the "Using BOND Reagents" section of your BOND user documentation).

If unexpected staining is observed that cannot be explained by variations in laboratory procedures, and a problem with the detection system is suspected, contact your local distributor or the regional office of Leica Biosystems immediately.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- WARNING - some mounting media may cause the red chromogen to crystallize or fade
- DO NOT USE DPX, Entellan®, Leica CV Mount, Eukitt® quick-hardening mounting medium, Surgipath MM 24®, Surgipath Sub-X® and Pertex mounting media as they are not suitable for use with this product.
- Fading or crystallization has not been observed with: Leica CV Ultra Mounting Media (Cat 14070937891), Poly-Mount® (Polysciences Inc. Cat 08381), VectaMount™ (Vector Laboratories Cat H-5000) and Limonene Mount (Electron Microscopy Sciences Cat 17987-01).

- This detection system is intended for in vitro diagnostic use.
- The activator and substrate solutions may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents. Dispose of according to local regulations.

Red Part C

Contains 1.0 M (<10%)

Hydrochloric Acid.

GHS05: Corrosion.

Signal words: Warning.

H290: May be corrosive to metals.

P234: Keep only in original container.

P390: Absorb spillage to prevent material damage.

- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems. Alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user¹.
- Do not mix reagents from different detection systems.

Instructions for Use

The Bond Polymer Refine Red Detection was developed for use on the BOND automated system using *IHC Protocol J. Operating parameters for application of the detection system reagents on the BOND Processing Module have been optimized by Leica Biosystems. These can be displayed by following the instructions in your BOND user documentation.

Product Specific Limitations

Bond Polymer Refine Red Detection has been optimized at Leica Biosystems for use with BOND ancillary reagents. Laboratories may use their own primary antibodies provided they have been diluted to an appropriate concentration with Bond Primary Antibody Diluent (Catalog No. AR9352). Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances.

The appropriate concentration of user's own primary antibodies may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and primary antibody concentrations.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls.

They should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

The Bond Polymer Refine Red Detection must be used with laboratory best practice in the use of tissue controls. For assurance, laboratories should stain each patient sample in conjunction with positive, negative and other tissue specific controls as needed.

Some mounting media are not compatible with Bond Polymer Refine Red Detection (see Precautions).

Rapid alcohol/xylene dehydration is recommended.

Contact your local distributor or regional office of Leica Biosystems for further information.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 est une marque commerciale de Supelco, membre du groupe Sigma-Aldrich Corporation.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 è un marchio di fabbrica di Supelco, società del gruppo Sigma-Aldrich.

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Bibliografie

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 ist eine Marke von Supelco, einem Unternehmen der Sigma-Aldrich Corporation.

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- Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 - Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 - Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
- ProClin® 950 es una marca comercial de Supelco, parte de Sigma-Aldrich Corporation.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin: 950 é uma marca comercial de Supelco, parte da Sigma-Aldrich Corporation.

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1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
Το ProClin[®] 950 είναι εμπορικό σήμα της Supelco, τμήμα της Sigma-Aldrich Corporation.

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Bibliografi

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 er et varemærke tilhørende Supelco, en del af Sigma-Aldrich Corporation.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin™ 950 er et varmerke fra Supelco, et selskap i Sigma-Aldrich Corporation.

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Kaynakça

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163, 28 Şubat 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Sipariş kodu M29-P.
3. Bancroft JD ve Stevens A. Theory and Practice of Histological Techniques. 4. Baskı. Churchill Livingstone, New York. 1996. ProClin™ 950, bir Sigma-Aldrich Corporation şirketi olan Supelco'nın markasıdır.

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09 Şubat 2015

Leica Biosystems Newcastle Ltd
Balliol Business Park
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

Leica Biosystems Canada
71 Four Valley Drive
Concord, Ontario L4K 4V8
Canada
☎ +1 800 248 0123

Leica Biosystems Inc
1700 Leider Lane
Buffalo Grove IL 60089
USA
☎ +1 800 248 0123

Leica Biosystems Melbourne Pty Ltd
495 Blackburn Road
Mt Waverley VIC 3149
Australia
☎ +61 2 8870 3500