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BOND™ Ready-to-Use ISH RNA Negative Control Probe

Catalog No: PB0809

Intended Use

This reagent is for *in vitro* diagnostic use.

RNA Negative Control Probe is intended for use in the identification of background staining resulting from non-specific interactions in formalin-fixed, paraffin-embedded tissue by *in situ* hybridization (ISH) using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

RNA Negative Control Probe is a single oligonucleotide, designed from zebra fish DNA and analysed using Basic Local Alignment Search Tool (BLAST) analysis to confirm that the sequence bears no homology with any human sequences. The RNA Negative Control Probe is generated with a fluorescein label using the same procedures as applied to other oligonucleotide probes used in the detection of RNA on BOND. Therefore, RNA Negative Control Probe is ideal as a negative control probe for RNA ISH on BOND.

Reagents Provided

RNA Negative Control Probe is a fluorescein-conjugated oligonucleotide probe (200 ng/mL) supplied in hybridization solution.
Total volume = 5.5 mL

Dilution and Mixing

RNA Negative Control Probe is ready to use. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Needed but not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and *in situ* hybridization staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8 °C. The product is stable under these conditions up to the expiry date indicated on the container label.

There are no obvious signs indicating contamination and/or instability. Appropriate positive and negative tissue controls should be run at the same time as test tissue.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.

RNA NEGATIVE CONTROL PROBE

Contains Formamide (<50%) and Dextran Sulphate (<30%).

GHS07: Exclamation mark.

GHS08: Health hazard.

Signal words: Danger.

H360D: May damage

the unborn child.

H315: Causes skin irritation.

H319: Causes serious eye irritation.

P201: Obtain special instructions before use.

P202: Do not handle until all safety precautions have been read and understood.

P260: Do not breathe dust/fumes/gas/mist/vapours/spray.

P308+313: IF exposed or concerned: Get medical advice/attention.

P314: Get medical advice/attention if you feel unwell.

P264: Wash hands thoroughly after handling.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+313: If eye irritation persists: Get medical advice/attention.

P302+352: IF ON SKIN: Wash with plenty of soap and water.

P332+313: If skin irritation occurs: Get medical advice/attention.

P362+364: Take off contaminated clothing and wash before reuse. Restricted to professional users.

- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

RNA Negative Control Probe was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with Anti-Fluorescein Antibody, BOND Polymer Refine Detection and BOND ancillary reagents. The recommended staining protocol for RNA Negative Control Probe is ISH Protocol A. Enzyme retrieval is recommended using the BOND Enzyme Pre-treatment Kit, Enzyme 1, for 15 minutes.

Appropriate tissue and reagent controls should always be used. The protocol for the tissue and reagent controls should correspond to that of the RNA test probe.

Quality Control

Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls in addition to the following procedures.

Positive Tissue Control

Used to indicate correctly prepared tissues and proper staining techniques. One positive tissue control should be included for each set of test conditions in each staining run.

A tissue with weak positive staining is more suitable than a tissue with strong positive staining for optimal quality control and to detect minor levels of reagent degradation.

Negative Tissue Control

Should be examined after the positive tissue control to verify the specificity of the labeling of probe to the target. Alternatively, the variety of different cell types present in most tissue sections frequently offers negative control sites, but this should be verified by the user.

Negative Reagent Control

Use RNA Negative Control Probe PB0809 in place of the RNA test probe with a section of each patient specimen to evaluate non-specific staining and allow better interpretation of specific staining at the target.

Positive Reagent Control

Use RNA Positive Control Probe PB0785 in place of the RNA test probe with a section of each patient specimen to provide information on the preservation of nucleic acids in the tissue as well as accessibility of nucleic acids to the probe. If the RNA Positive Control Probe fails to demonstrate positive staining, results with the test specimens should be considered invalid.

Patient Tissue

Examine patient specimens stained with RNA test probe last. Positive staining intensity should be assessed within the context of any non-specific background staining of the RNA Negative Control Probe PB0809.

Results Expected

Normal Tissues

Using PB0809 no staining was observed in a wide range of tissues. (Total number of normal cases evaluated = 96).

Abnormal Tissues

Using PB0809 no staining was observed in ovarian tumors (0/3), thyroid tumors (0/3), lung tumors (0/3), liver tumors (0/3), brain tumors (0/2), esophageal tumors (0/2), breast tumors (0/2), skin tumors (0/2), metastatic tumors of unknown origin (0/2), stomach tumors (0/1), colon tumors (0/1), rectal tumors (0/1), a tumor of the larynx (0/1), and a tumor of the thymus (0/1) (Total number of abnormal cases evaluated = 29).

PB0809 is recommended to aid as a screening tool to evaluate staining of an RNA test probe on a patient specimen.

Product Specific Limitations

RNA Negative Control Probe has been optimized at Leica Biosystems for use with Anti-Fluorescein antibody, BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of enzymatic digestion, and must be determined empirically. RNA Negative Control Probe should be used when optimising retrieval conditions and protocol times.

Troubleshooting

Reference 3 may aid in remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on *in situ* hybridization with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

References

1. Pringle JH, Primrose L, Kind CN, et al. *In situ* hybridization demonstration of poly-adenylated RNA sequences in formalin-fixed paraffin sections using a biotinylated oligonucleotide poly d(T) probe. *Journal of Pathology*. 1989;158(4):279—286.
2. Lewin B. Units of transcription and translation: the relationship between heterogeneous nuclear RNA and messenger RNA. *Cell*. 1975;4:11—20.
3. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
4. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
5. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18—20.

Date of Issue

25 September 2018

PB0809

Page 3

BOND™ Ready-to-Use ISH RNA Negative Control Probe

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Referanslar

1. Pringle JH, Primrose L, Kind CN, et al. *In situ* hybridization demonstration of poly-adenylated RNA sequences in formalin-fixed paraffin sections using a biotinylated oligonucleotide poly d(T) probe. *Journal of Pathology*. 1989;158(4):279—286.
2. Lewin B. Units of transcription and translation: the relationship between heterogeneous nuclear RNA and messenger RNA. *Cell*. 1975;4:11—20.

3. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
4. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
5. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18—20.

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1. Pringle JH, Primrose L, Kind CN, et al. *In situ* hybridization demonstration of poly-adenylated RNA sequences in formalin-fixed paraffin sections using a biotinylated oligonucleotide poly d(T) probe. *Journal of Pathology*. 1989;158(4):279—286.
2. Lewin B. Units of transcription and translation: the relationship between heterogeneous nuclear RNA and messenger RNA. *Cell*. 1975;4:11—20.
3. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
4. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
5. Wilkinon DG. The theory and practice of *in situ* hybridization. In: Wilkinon DG. (ed.) *In situ* Hybridization A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18—20.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd