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Peroxidase Detection System (250 tests)

Product No: RE7110-K

Peroxidase Detection System (500 tests)

Product No: RE7120-K

Streptavidin-HRP

Product No: RE7104

Intended Use

For in vitro diagnostic use.

Novocastra™ Peroxidase Detection Systems (250 tests) RE7110-K and (500 tests) RE7120-K are for the visualization of mouse IgG, mouse IgM and rabbit IgG primary antibodies. Novocastra™ Streptavidin-HRP RE7104 are component reagents of these systems that are intended for use in the procedure described below. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Principle of Procedure

The first immunoperoxidase technique was reported by Nakane and Pierce¹, since then many developments have occurred leading to the immunoperoxidase techniques commonly used today. The Novocastra™ Peroxidase Detection Systems are based upon the technique described by Guesdon.² These products are used in an immunohistochemical (IHC) procedure, which allows the qualitative identification by light microscopy of Epitopes in sections of formalin-fixed, paraffin-embedded tissue, via sequential steps with interposed washing steps.

If required by the primary antibody, sections are subjected to epitope retrieval prior to staining. Endogenous peroxidase activity is neutralized using the Novocastra™ Peroxidase Block RE7101. This is followed by application of the Novocastra™ Protein Block RE7102 to reduce non-specific binding of primary and secondary antibodies. The section is subsequently incubated with optimally diluted primary antibody. Novocastra™ Biotinylated Secondary Antibody RE7103 a biotin-conjugated secondary antibody formulation that recognizes mouse and rabbit immunoglobulins, is used to detect any tissue-bound primary antibody. Novocastra™ Streptavidin-HRP RE7104, a streptavidin-peroxidase conjugate is then applied and binds to the biotin present on the secondary antibody. Sections are further incubated with the substrate/chromogen, 3,3' - diaminobenzidine (DAB), prepared from Novocastra™ DAB Chromogen RE7105 and Novocastra™ DAB Substrate Buffer RE7106 as described below. Reaction catalyzed by the peroxidase produces a visible brown precipitate at the Epitope site. Sections are counterstained with Novocastra™ Hematoxylin RE7107 and coverslipped. Results are interpreted using a light microscope and aid in the differential diagnosis of pathophysiological processes, which may or may not be associated with a particular Epitope.

Reagents Provided

Details of those reagents from the following list that are supplied in each product are given in the table below.

1. Peroxidase Block RE7101. 3-4% Hydrogen peroxide.
2. Protein Block RE7102. 0.4% Casein in phosphate-buffered saline, with stabilizers, surfactant, and 0.2% Bronidox L as a preservative.
3. Biotinylated Secondary Antibody RE7103. Affinity purified goat anti-mouse IgM (5µg/ml), horse anti-mouse IgG (10.5µg/mL) and horse anti-rabbit IgG (10.5µg/ml) in Tris-buffered saline with protein stabilizer and 0.35% ProClin™ 950.
4. Streptavidin-HRP RE7104. 1 to 1.4µg/ml in tris-buffered saline with protein stabilizer and 0.35% ProClin™ 950.
5. DAB Chromogen RE7105. 1.74% w/v 3,3' - diaminobenzidine, in a stabilizer solution.
6. DAB Substrate Buffer RE7106. Buffered-solution containing ≤ 0.1% hydrogen peroxide and preservative.
7. Hematoxylin RE7107 <0.1% hematoxylin.

Reagent	Product Number	Peroxidase Detection System (250 tests) RE7110-K	Peroxidase Detection System (500 tests) RE7120-K	Streptavidin-HRP RE7104
Peroxidase Block	RE7101	1 x 25ml	2 x 25ml	
Protein Block	RE7102	1 x 25ml	2 x 25ml	
Biotinylated Secondary Antibody	RE7103	1 x 25ml	2 x 25ml	
Streptavidin-HRP	RE7104	1 x 25ml	2 x 25ml	1 x 25ml
DAB Chromogen	RE7105	1 x 3ml	1 x 3ml	
DAB Substrate Buffer	RE7106	1 x 30ml	2 x 30ml	
Hematoxylin	RE7107	1 x 25ml	2 x 25ml	

Reconstitution, Mixing, Dilution, Titration

The Peroxidase Block RE7101, Protein Block RE7102, Biotinylated Secondary Antibody RE7103 Streptavidin-HRP RE7104 and Hematoxylin RE7107 solutions are prediluted. Reconstitution, mixing, dilution, or titration of these reagents is not recommended. Further dilution may result in loss of Epitope staining. The user must validate any such change.

The DAB Chromogen RE7105 requires dilution to 1/20 in DAB Substrate Buffer RE7106 prior to use. Further dilution may result in loss of Epitope staining. The user must validate any such change.

Storage and Stability

Store at 2–8 °C. Do not freeze. Return to 2–8 °C immediately after use. Do not use after expiration date indicated on the product label. Storage conditions other than those specified must be verified by the user. There are no obvious signs to indicate instability of this product therefore positive and negative controls should be run simultaneously with patient samples.

Specimen Preparation

The recommended fixative is 10% neutral-buffered formalin for paraffin-embedded tissue sections.

Warnings and Precautions

DAB CHROMOGEN

Contains <10% Biphenyl-3,3',4,4'-Tetrahydroamine.

GHS08: Health hazard.

Signal words: Danger.

H341: Suspected of causing genetic defects.

H350: May cause cancer.

P201: Obtain special instructions before use.

P202: Do not handle until all safety precautions have been read and understood.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P308+313: IF exposed or concerned: Get medical advice/attention.

P501: Special waste of contents/container to hazardous or special waste collection point.

A Material Safety Data Sheet is available upon request or available from www.LeicaBiosystems.com

For professional users.

Do not mix reagents from different detection systems.

Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions.³

Never pipette reagents by mouth and avoid contacting the skin and mucus membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water.

Consult Federal, State or local regulations for disposal of any potentially toxic components.

Minimize microbial contamination of reagents or an increase in non-specific staining may occur. Incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Procedure

A. Reagents required but not supplied

1. Standard solvents used in immunohistochemistry.
2. 50mM tris-buffered saline (TBS) pH7.6.
3. Epitope retrieval solution(s).
4. Enzyme retrieval solution(s).
5. Antibody diluent.
6. Primary antibody.
7. Mounting medium.

B. Equipment required but not supplied

1. Equipment required for Epitope retrieval, if recommended for the primary antibody.
2. General immunohistochemistry laboratory equipment.

C. Methodology

Prior to undertaking this methodology, users must be trained in immunohistochemical techniques.

The combination of the primary antibody, its dilution, together with the detection system should be validated by the user on a series of known positive and negative controls.

Unless indicated, all steps are performed at room temperature (25 °C).

1. Cut and mount sections on slides coated with a suitable tissue adhesive.
2. De-paraffinize sections in xylene or xylene substitutes.
3. Re-hydrate through graded alcohols.
4. Wash slides in running tap water.
5. Perform Epitope retrieval as required (see **Recommendations for Use** for primary antibody).
6. Wash slides in de-ionized water.
7. Neutralize endogenous peroxidase using Peroxidase Block RE7101 for 5 minutes.
8. Wash in TBS for 2 x 5 minutes.
9. Incubate with Protein Block RE7102 for 5 minutes.
10. Wash in TBS for 2 x 5 minutes.
11. Incubate with optimally diluted primary antibody (see **Recommendations for Use** for primary antibody).
12. Wash in TBS for 2 x 5 minutes.

13. Incubate with Biotinylated Secondary Antibody RE7103 for 30 minutes.
14. Wash in TBS for 2 x 5 minutes.
15. Incubate with Streptavidin-HRP RE7104 for 30 mins.
16. Wash in TBS for 2 x 5 minutes with gentle rocking.
17. Develop peroxidase activity with DAB working solution (see **DAB Working Solution**) for 5 minutes.
18. Rinse slides in water.
19. Counterstain with Hematoxylin RE7107.
20. Rinse slides in water for 5 minutes.
21. Dehydrate, clear and mount sections.

DAB Working Solution

Add 50µl of DAB Chromogen RE7105 to 1ml of DAB Substrate Buffer RE7106. Use within six hours of preparation.

Quality Control

Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls in addition to the following procedures. Controls should be fresh autopsy/biopsy/ surgical specimens, formalin-fixed, processed and paraffin wax-embedded as soon as possible in the same manner as the patient sample(s).

Positive Tissue Control

Used to indicate correctly prepared tissues and proper staining techniques. One positive tissue control should be included for each set of test conditions/primary antibody in each staining run. A tissue with weak positive staining is more suitable than a tissue with strong positive staining for optimal quality control and to detect minor levels of reagent degradation.⁴ For recommended positive control tissue see primary antibody Instructions for Use. If the positive tissue control fails to demonstrate positive staining, results with the test specimens should be considered invalid.

Negative Tissue Control

Should be examined after the positive tissue control to verify the specificity of the labeling of the target Epitope by the primary antibody. For recommended negative control tissue see primary antibody Instructions for Use. Alternatively, the variety of different cell types present in most tissue sections frequently offers negative control sites, but this should be verified by the user. Non-specific staining, if present, usually has a diffuse appearance. Sporadic staining of connective tissue may also be observed in sections from excessively formalin-fixed tissues. Use intact cells for interpretation of staining results. Necrotic or degenerated cells often stain non-specifically.⁵ False-positive results may be seen due to non-immunological binding of proteins or substrate reaction products. They may also be caused by endogenous enzymes such as pseudoperoxidase (erythrocytes), endogenous peroxidase (cytochrome C), or endogenous biotin⁶ (eg. liver, breast, brain, kidney). To differentiate endogenous enzyme activity or non-specific binding of enzymes from specific immunoreactivity, additional patient tissues may be stained exclusively with substrate-chromogen or Streptavidin-HRP, and substrate-chromogen, respectively. If specific staining occurs in the negative tissue control, results with the patient specimens should be considered invalid.

Negative Reagent Control

Use a non-specific negative reagent control in place of the primary antibody with a section of each patient specimen to evaluate non-specific staining and allow better interpretation of specific staining at the Epitope site.

Patient Tissue

Examine stained patient specimens last. Positive staining intensity should be assessed within the context of any non-specific background staining of the negative reagent control. As with any immunohistochemical test, a negative result means that the Epitope was not detected, not that the Epitope was absent in the cells/tissue assayed. If necessary, use a panel of antibodies to identify false-negative reactions.

Limitations

Immunohistochemistry is a multistep diagnostic process that consists of specialized training in the selection of the appropriate reagents; tissue selection, fixation, and processing; preparation of the IHC slide; and interpretation of the staining results.

Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning or contamination with other tissues or fluids may produce artifacts, antibody trapping, or false negative results. Inconsistent results may be due to variations in fixation and embedding methods, or to inherent irregularities within the tissue.⁷

Excessive or incomplete counterstaining may compromise proper interpretation of results.

The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Novocastra™ Peroxidase Detection Systems RE7110-K or RE7120-K are for use on paraffin-embedded sections with specific fixation requirements. Unexpected Epitope expression may occur, especially in neoplasms. The clinical interpretation of any stained tissue section must include morphological analysis and the evaluation of appropriate controls.

Performance Characteristics

The performance of Novocastra™ Peroxidase Detection Systems (250 tests)/ (500 tests) RE7110-K/RE7120-K, Novocastra™ Streptavidin-HRP RE7104 have been validated using a range of Novocastra™ mouse IgG, mouse IgM and rabbit IgG primary antibodies. These products are stable up to the expiry date(s) printed on the product label.

Bibliography

1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. Journal of Histochemistry and Cytochemistry. 1967; 14 :929-931.
2. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. Journal of Histochemistry and Cytochemistry.1979; 27:1131.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.

4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL and Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Non-immunologic binding of horseradish peroxidase to hepatitis B surface Epitope: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

Amendments to Previous Issue

Not applicable.

Date of Issue

15 July 2019

Peroxidase Detection System (250 tests)

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Peroxidase Detection System (500 tests)

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Bibliographie

1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. *Journal of Histochemistry and Cytochemistry*. 1967; 14 :929-931.
2. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. *Journal of Histochemistry and Cytochemistry*. 1979; 27:1131.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL and Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Non-immunologic binding of horseradish peroxidase to hepatitis B surface Epitope: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Reagente	Codice prodotto	Peroxidase Detection System (250 tests) RE7110-K	Peroxidase Detection System (500 tests) RE7120-K	Streptavidin-HRP RE7104
Peroxidase Block	RE7101	1 x 25 ml	2 x 25 ml	
Protein Block	RE7102	1 x 25 ml	2 x 25 ml	
Biotinylated Secondary Antibody	RE7103	1 x 25 ml	2 x 25 ml	
Streptavidin-HRP	RE7104	1 x 25 ml	2 x 25 ml	1 x 25 ml
DAB Chromogen	RE7105	1 x 3 ml	1 x 3 ml	
DAB Substrate Buffer	RE7106	1 x 30 ml	2 x 30 ml	
Hematoxylin	RE7107	1 x 25 ml	2 x 25 ml	

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1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. Journal of Histochemistry and Cytochemistry. 1967; 14 :929-931.

2. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. *Journal of Histochemistry and Cytochemistry*. 1979; 27:1131.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL and Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Non-immunologic binding of horseradish peroxidase to hepatitis B surface Epitope: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Reagenz	Produktnummer	Peroxidase Detection System (250 tests) RE7110-K	Peroxidase Detection System (500 tests) RE7120-K	Streptavidin-HRP RE7104
Peroxidase Block	RE7101	1 x 25 ml	2 x 25 ml	
Protein Block	RE7102	1 x 25 ml	2 x 25 ml	
Biotinylated Secondary Antibody	RE7103	1 x 25 ml	2 x 25 ml	
Streptavidin-HRP	RE7104	1 x 25 ml	2 x 25 ml	1 x 25 ml
DAB Chromogen	RE7105	1 x 3 ml	1 x 3 ml	

DAB Substrate Buffer	RE7106	1 x 30 ml	2 x 30 ml	
Hematoxylin	RE7107	1 x 25 ml	2 x 25 ml	

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Bibliografia

1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. *Journal of Histochemistry and Cytochemistry*. 1967; 14 :929-931.
2. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. *Journal of Histochemistry and Cytochemistry*.1979; 27:1131.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL and Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Non-immunologic binding of horseradish peroxidase to hepatitis B surface Epitope: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Reagens	Produkt-numme	Peroxidase Detection System (250 tests) RE7110-K	Peroxidase Detection System (500 tests) RE7120-K	Streptavidin-HRP RE7104
Peroxidase Block	RE7101	1 x 25ml	2 x 25 ml	
Protein Block	RE7102	1 x 25ml	2 x 25 ml	
Biotinylated Secondary Antibody	RE7103	1 x 25ml	2 x 25 ml	
Streptavidin-HRP	RE7104	1 x 25ml	2 x 25ml	1 x 25ml
DAB Chromogen	RE7105	1 x 3ml	1 x 3ml	
DAB Substrate Buffer	RE7106	1 x 30ml	2 x 30ml	
Hematoxylin	RE7107	1 x 25ml	2 x 25ml	

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Bibliografi

1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. Journal of Histochemistry and Cytochemistry. 1967; 14 :929-931.
2. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. Journal of Histochemistry and Cytochemistry.1979; 27:1131.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
6. Mount SL and Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. Current Diagnostic Pathology. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Non-immunologic binding of horseradish peroxidase to hepatitis B surface Epitope: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.

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Bibliografi

1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. Journal of Histochemistry and Cytochemistry. 1967; 14 :929-931.
2. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. Journal of Histochemistry and Cytochemistry. 1979; 27:1131.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
6. Mount SL and Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. Current Diagnostic Pathology. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Non-immunologic binding of horseradish peroxidase to hepatitis B surface Epitope: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.

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Reagens	Productnummer	Peroxidase Detection System (250 tests) RE7110-K	Peroxidase Detection System (500 tests) RE7120-K	Streptavidin-HRP RE7104
Peroxidase Block	RE7101	1 x 25 ml	2 x 25 ml	
Protein Block	RE7102	1 x 25 ml	2 x 25 ml	
Biotinylated Secondary Antibody	RE7103	1 x 25 ml	2 x 25 ml	
Streptavidin-HRP	RE7104	1 x 25 ml	2 x 25 ml	1 x 25 ml
DAB Chromogen	RE7105	1 x 3 ml	1 x 3 ml	
DAB Substrate Buffer	RE7106	1 x 30 ml	2 x 30 ml	
Hematoxylin	RE7107	1 x 25 ml	2 x 25 ml	

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1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. Journal of Histochemistry and Cytochemistry. 1967; 14 :929-931.
2. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. Journal of Histochemistry and Cytochemistry. 1979; 27:1131.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.

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Reagens	Produkt Antall	Peroxidase Detection System (250 tests) RE7110-K	Peroxidase Detection System (500 tests) RE7120-K	Streptavidin-HRP RE7104
Peroxidase Block	RE7101	1 x 25ml	2 x 25ml	
Protein Block	RE7102	1 x 25ml	2 x 25ml	
Biotinylated Secondary Antibody	RE7103	1 x 25ml	2 x 25ml	
Streptavidin-HRP	RE7104	1 x 25ml	2 x 25ml	1 x 25ml
DAB Chromogen	RE7105	1 x 3ml	1 x 3ml	
DAB Substrate Buffer	RE7106	1 x 30ml	2 x 30ml	
Hematoxylin	RE7107	1 x 25ml	2 x 25ml	

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Bibliografi

8. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. Journal of Histochemistry and Cytochemistry. 1967; 14 :929-931.
9. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. Journal of Histochemistry and Cytochemistry. 1979; 27:1131.
10. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
11. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
12. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
13. Mount SL and Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. Current Diagnostic Pathology. 2001; 7:161-167.
14. Omata M, Liew CT, Ashcavai M, Peters RL. Non-immunologic binding of horseradish peroxidase to hepatitis B surface Epitope: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.

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1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. Journal of Histochemistry and Cytochemistry. 1967; 14 :929-931.
2. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. Journal of Histochemistry and Cytochemistry.1979; 27:1131.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
6. Mount SL and Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. Current Diagnostic Pathology. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavi M, Peters RL. Non-immunologic binding of horseradish peroxidase to hepatitis B surface Epitope: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.

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Bibliografia

1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. *Journal of Histochemistry and Cytochemistry*. 1967; 14 :929-931.
2. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. *Journal of Histochemistry and Cytochemistry*. 1979; 27:1131.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL and Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Non-immunologic binding of horseradish peroxidase to hepatitis B surface Epitope: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Reagent	Šifra izdelka	Peroxidase Detection System (250 tests) RE7110-K	Peroxidase Detection System (500 tests) RE7120-K	Streptavidin-HRP RE7104
Peroxidase Block	RE7101	1 x 25 ml	2 x 25 ml	
Protein Block	RE7102	1 x 25 ml	2 x 25 ml	
Biotinylated Secondary Antibody	RE7103	1 x 25 ml	2 x 25 ml	
Streptavidin-HRP	RE7104	1 x 25 ml	2 x 25 ml	1 x 25 ml
DAB Chromogen	RE7105	1 x 3ml	1 x 3ml	
DAB Substrate Buffer	RE7106	1 x 30ml	2 x 30ml	
Hematoxylin	RE7107	1 x 25 ml	2 x 25 ml	

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Literatura

1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. Journal of Histochemistry and Cytochemistry. 1967; 14 :929-931.
2. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. Journal of Histochemistry and Cytochemistry. 1979; 27:1131.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
6. Mount SL and Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. Current Diagnostic Pathology. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Non-immunologic binding of horseradish peroxidase to hepatitis B surface Epitope: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.

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Literatura

1. Nakane PK and Pierce GB. Enzyme labeled antibodies : Preparations and applications for the localization of Epitopes. Journal of Histochemistry and Cytochemistry. 1967; 14 :929-931.
2. Guesdon JL, Ternynck T and Avrameas S. The use of avidin-biotin interaction in immunoenzymatic techniques. Journal of Histochemistry and Cytochemistry.1979; 27:1131.

3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
6. Mount SL and Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. Current Diagnostic Pathology. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Non-immunologic binding of horseradish peroxidase to hepatitis B surface Epitope: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd