BIOSYSTEMS

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BOND™ Ready-to-Use Primary Antibody Programmed Death Ligand 1 (73-10)

Catalog No: PA0832

Intended Use

This reagent is for *in vitro* diagnostic use.

Programmed Death Ligand 1 (73-10) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of endogenous PD-L1 protein in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Programmed Death Ligand 1 (73-10) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of endogenous PD-L1 protein is achieved by first allowing the binding of Programmed Death Ligand 1 (73-10) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

Programmed Death Ligand 1 (73-10) is a recombinant rabbit anti-human monoclonal antibody supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

73-10

N.B. This PD-L1 antibody has been created by Epitomics Inc., using Epitomics' proprietary rabbit monoclonal antibody technology covered under Patent No.'s 5,675,063 and 7,402,409.

Immunogen

Peptide corresponding to a region of the C-terminal cytoplasmic domain of PD-L1.

Specificity

Human PD-L1 protein containing the immunized target sequence.

Ig Class

Rabbit IgG

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 0.2 mg/L as determined by ELISA.

Dilution and Mixing

Programmed Death Ligand 1 (73-10) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Programmed Death Ligand 1 (73-10) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Programmed Death Ligand 1 (73-10) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection. The recommended staining protocol for Programmed Death Ligand 1 (73-10) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 1 for 20 minutes.

Results Expected

Normal Tissues

Clone 73-10 detected Programmed Death Ligand 1 protein in the crypts and germinal centres of tonsil, occasional cells in lung including alveolar macrophages, macrophages in small intestine, Kupffer cells in liver, occasional megakaryocytes in bone marrow, some immune cells of the oesophagus, colon, rectum and spleen, weak membrane staining in a proportion of cells in the pituitary gland, parathyroid chief cells, Hassall's corpuscles and immune cells in thymus, occasional cells of the endometrium (possibly macrophages), cells of the glomerulus and a proportion of glandular epithelium in prostate. (Total number of normal cases evaluated = 120).

Tumor Tissues

Clone 73-10 stained 22/42 lung cancers (including 16/24 squamous cell carcinomas, 4/16 adenocarcinomas, 1/1 large cell carcinoma, 1/1 small cell carcinoma), 1/2 bladder transitional cell carcinomas, 1/3 oesophageal squamous cell carcinomas, 1/3 stomach adenocarcinomas, 1/2 cervical squamous cell carcinoma, 1/1 Hodgkin lymphoma, 1/1 anaplastic large cell lymphoma, 1/1 nasopharyngeal carcinoma, 1/1 pancreatic adenocarcinoma, 1/1 follicular carcinoma of the thyroid, 1/1 squamous cell carcinoma of the tongue. No staining was observed in a variety of additional abnormal tissues evaluated, including breast tumours (0/5), metastatic tumours (0/4), tumours of the colon (0/4), liver tumours (0/4), brain tumours (0/4), ovarian tumours (0/3), thyroid adenomas (0/3), tumours of the adrenal gland (0/2), tumours of the small intestine (0/2), tumours of the rectum (0/2), kidney tumours (0/2), tumours of the salivary gland (0/2), testis (0/2), and endometrium (0/2), a follicular papillary adenocarcinoma of the thyroid (0/1), a prostatic adenocarcinoma (0/1), a hyperplastic prostate, (0/1), a pulmonary oedema (0/1), a pancreatic adenocarcinoma (0/1), a hard palate adenocarcinoma (0/1), a non-Hodgkin B-cell lymphoma (0/1), a chondrosarcoma (0/1) and a melanoma (0/1). (Total number of abnormal tissues evaluated = 107).

Programmed Death Ligand 1 (73-10) is recommended for the detection of endogenous levels of total PD-L1 protein, in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.

Product Specific Limitations

Programmed Death Ligand 1 (73-10) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. International Association for the Study of Lung Cancer. IASLC Atlas of PD-L1 Immunohistochemistry testing in lung cancer. 2017
5. Patel SP and Kurzrock R. PD-L1 Expression as a Predictive Biomarker in Cancer Immunotherapy. Molecular Cancer Therapeutics. 2015; 14(4):847-56
6. Wang X, Teng F, Kong L et al. PD-L1 expression in human cancers and its association with clinical outcomes. OncoTargets and Therapy. 2016; 9:5023-5039

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. International Association for the Study of Lung Cancer. IASLC Atlas of PD-L1 Immunohistochemistry testing in lung cancer. 2017.
5. Patel SP and Kurzrock R. PD-L1 Expression as a Predictive Biomarker in Cancer Immunotherapy. Molecular Cancer Therapeutics. 2015; 14(4):847-56
6. Wang X, Teng F, Kong L et al. PD-L1 expression in human cancers and its association with clinical outcomes. OncoTargets and Therapy. 2016; 9:5023-5039

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. International Association for the Study of Lung Cancer. IASLC Atlas of PD-L1 Immunohistochemistry testing in lung cancer. 2017
5. Patel SP and Kuzrock R. PD-L1 Expression as a Predictive Biomarker in Cancer Immunotherapy. Molecular Cancer Therapeutics. 2015; 14(4):847-56
6. Wang X, Teng F, Kong L et al. PD-L1 expression in human cancers and its association with clinical outcomes. OncoTargets and Therapy. 2016; 9:5023-5039

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd