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BOND™ Ready-To-Use Anti-Fluorescein Antibody 15 mL Catalog No: AR0222

Intended Use

This reagent is for *in vitro* diagnostic use.

Anti-Fluorescein Antibody is intended for use with fluorescein-labeled nucleic acid probes bound to specific nucleotide sequences in formalin-fixed, paraffin-embedded tissue by *in situ* hybridization (ISH) using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

ISH allows the detection and visualization of specific nucleic acids in tissues sections (see "Using BOND Reagents" in your BOND user documentation).

The Anti-Fluorescein Antibody allows linking of the fluorescein-labeled nucleic acid probe with the detection reagents, and consequently visualization of a chromogenic product by light microscopy.

Reagents Provided

Anti-Fluorescein Antibody is a purified IgG fraction of a mouse monoclonal antibody (108 µg/mL).

Total volume = 15 mL

Dilution and Mixing

Anti-Fluorescein Antibody is ready to use. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and *in situ* hybridization staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8 °C. The product is stable under these conditions up to the expiry date indicated on the container label.

The signs indicating contamination and/or instability of Anti-Fluorescein Antibody are: turbidity of the solution, odor development, and presence of precipitate.

Appropriate positive and negative tissue controls should be run at the same time as test tissue.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Anti-Fluorescein Antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with fluorescein-labeled nucleic acid probes and BOND detection. The recommended pre-treatment and staining protocol are determined by the specific probe and/or detection selected. Please refer to the relevant Instructions for Use for further information.

Product Specific Limitations

Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. Protocol times may vary due to variation in tissue type, fixation and processing. In addition, BOND Enzyme concentration and incubation time may require optimization depending on tissue type, processing and fixation conditions. Negative reagent controls should be used when optimizing pre-treatment conditions and protocol times.

Troubleshooting

Reference 3 may aid in remedial action.

Test samples should be complemented by the appropriate tissue and reagent controls.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on *in situ* hybridization with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

References

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In Situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

U.S. Patent No. 5, 597, 692 licensed from Syngene Ltd.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In Situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp00,18–20.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Wilkinon DG. The theory and practice of *in situ* hybridization. In: Wilkinon DG. (ed.) *In Situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In Situ Hybridization. A practical approach*. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In Situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.
- U.S. Patent No. 5, 597, 692 onder licentie van Syngene Ltd.

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Irodalomjegyzék

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In Situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.
- 5., 597., 692. számú USA szabadalom, melynek licenctulajdonosa a Syngene Kft.

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Referințe

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In Situ* Hybridization. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.
- Nr. brevet SUA 5, 597, 692 licențiat de la Syngene Ltd.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In Situ Hybridization*. A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd