BOND™ Ready-to-Use Primary Antibody Carcinoembryonic Antigen (CD66e) (COL-1)

Catalog No: PA0848

Leica Biosystems Newcastle Ltd
Balliol Business Park
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

[EN](#)

[FR](#)

[IT](#)

[DE](#)

[ES](#)

[PT](#)

[SV](#)

[EL](#)

[DA](#)

[NL](#)

[NO](#)

[TR](#)

[BG](#)

[HU](#)

[RO](#)

[RU](#)

[PL](#)

[SL](#)

[CS](#)

[SK](#)

[AR](#)

Instructions for Use

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Catalog No: PA0848

Intended Use

This reagent is for *in vitro* diagnostic use.

Carcinoembryonic Antigen (CD66e) (COL-1) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of carcinoembryonic antigen in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Carcinoembryonic Antigen (CD66e) (COL-1) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of carcinoembryonic antigen is achieved by first allowing the binding of Carcinoembryonic Antigen (CD66e) (COL-1) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

Carcinoembryonic Antigen (CD66e) (COL-1) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

COL-1

Immunogen

BALB/C mice injected with the extract of colon carcinoma cells.

Specificity

Human carcinoembryonic antigen.

Ig Class

IgG2a.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 0.2 mg/L as determined by ELISA.

Dilution and Mixing

Carcinoembryonic Antigen (CD66e) (COL-1) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Carcinoembryonic Antigen (CD66e) (COL-1) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Carcinoembryonic Antigen (CD66e) (COL-1) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection. The recommended staining protocol for Carcinoembryonic Antigen (CD66e) (COL-1) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 2 for 20 minutes.

Results Expected

Normal Tissues

Clone COL-1 detected the carcinoembryonic antigen in epithelial cells of the colon, rectum, esophagus and larynx, and in the squamous mucosa of the cervix. (Total number of normal cases = 148).

Tumor Tissues

Clone COL-1 stained 59/65 colon tumors (including 21/23 metastatic adenocarcinomas, 18/22 adenocarcinomas, 6/6 adenomas, 5/5 polyps, 4/4 mucinous adenocarcinomas, 3/3 metastatic mucinous adenocarcinomas and 2/2 metastatic signet-ring cell carcinomas), 34/77 lung tumors (including 16/28 adenocarcinomas, 11/28 squamous cell carcinomas, 3/8 small cell carcinomas, 1/6 atypical carcinoids, 1/3 large cell carcinomas, 1/2 metastatic squamous cell carcinomas and 1/2 metastatic adenocarcinomas), 9/49 ovarian tumors (including 3/12 endometrioid adenocarcinomas, 2/3 mucinous adenocarcinomas, 2/6 serous papillary adenocarcinomas, 1/3 mucinous cystadenocarcinomas, 1/3 mucinous cystadenomas, 0/5 serous adenocarcinomas, 0/4 theca cell tumors, 0/3 granulosa cell tumors, 0/3 undifferentiated carcinomas, 0/2 endodermal sinus carcinomas, 0/1 dysgerminoma, 0/1 immature teratoma, 0/1 leiomyosarcoma, 0/1 borderline mucinous cystadenoma and 0/1 adenocarcinoma), 5/11 metastatic tumors (including 3/4 metastatic gastric carcinomas, 2/4 metastatic breast carcinomas, 0/1 metastatic nasopharyngeal carcinoma, 0/1 metastatic thyroid adenocarcinoma and 0/1 metastatic esophageal squamous cell carcinoma), 5/5 chronic colitis, 4/4 Crohn's disease, 3/3 stomach adenocarcinomas, 2/2 tumors of the small intestine (including 1/1 adenoma and 1/1 adenocarcinoma), 2/2 squamous cell carcinomas of the cervix, 1/3 esophageal squamous cell carcinomas, 1/2 endometrial adenocarcinomas, 1/1 squamous cell carcinoma of the skin, 1/1 pancreatic tumor and 1/1 squamous cell carcinoma of the tongue. No staining was observed in breast tumors (0/5), thyroid tumors (0/5), brain tumors (0/4), hepatocellular carcinomas (0/4), lymphomas (0/3), prostatic tumors (0/2), tumors of the head and neck (0/2), tumors of the adrenal gland (0/2), bladder tumors (0/2), bone tumors (0/2), tumors of the salivary gland (0/2), seminomas (0/2), kidney tumors (0/2), a prostatic hyperplasia (0/1), a melanoma (0/1), a case of tuberculosis of the colon (0/1) and a pheochromocytoma (0/1). (Total number of abnormal cases = 267).

Carcinoembryonic Antigen (CD66e) (COL-1) is recommended for the detection of carcinoembryonic antigen in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.

Product Specific Limitations

Carcinoembryonic Antigen (CD66e) (COL-1) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Chhieng DC, Rodriguez-Burford C, Telley LI, et al. Expression of CEA, Tag-72 and Lewis-Y antigen in primary and metastatic lesions of ovarian carcinoma. Human Pathology. 2003; 34(10): 1016-1021.
5. Kass ES, Greiner JW, Kantor JA, et al. Carcinoembryonic antigen as a target for specific antitumor immunotherapy of head and neck cancer. Cancer Research. 2002; 62(17): 5049-5057.
6. Tendler A, Kaufman HL, Kadish AS. Increased carcinoembryonic antigen expression in cervical intraepithelial neoplasia grade 3 and in cervical squamous cell carcinoma. Human Pathology. 2000; 31(11): 1357-1362.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code : M9-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Chhieng DC, Rodriguez-Burford C, Telley LI, et al. Expression of CEA, Tag-72 and Lewis-Y antigen in primary and metastatic lesions of ovarian carcinoma. Human Pathology. 2003; 34(10): 1016-1021.
5. Kass ES, Greiner JW, Kantor JA, et al. Carcinoembryonic antigen as a target for specific antitumor immunotherapy of head and neck cancer. Cancer Research. 2002; 62(17): 5049-5057.
6. Tendler A, Kaufman HL, Kadish AS. Increased carcinoembryonic antigen expression in cervical intraepithelial neoplasia grade 3 and in cervical squamous cell carcinoma. Human Pathology. 2000; 31(11): 1357-1362.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 28. February 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Chhieng DC, Rodriguez-Burford C, Telley LI, et al. Expression of CEA, Tag-72 and Lewis-Y antigen in primary and metastatic lesions of ovarian carcinoma. Human Pathology. 2003; 34(10): 1016-1021.
5. Kass ES, Greiner JW, Kantor JA, et al. Carcinoembryonic antigen as a target for specific antitumor immunotherapy of head and neck cancer. Cancer Research. 2002; 62(17): 5049-5057.
6. Tendler A, Kaufman HL, Kadish AS. Increased carcinoembryonic antigen expression in cervical intraepithelial neoplasia grade 3 and in cervical squamous cell carcinoma. Human Pathology. 2000; 31(11): 1357-1362.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code : M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Chhieng DC, Rodriguez-Burford C, Telley LI, et al. Expression of CEA, Tag-72 and Lewis-Y antigen in primary and metastatic lesions of ovarian carcinoma. Human Pathology. 2003; 34(10): 1016-1021.
5. Kass ES, Greiner JW, Kantor JA, et al. Carcinoembryonic antigen as a target for specific antitumor immunotherapy of head and neck cancer. Cancer Research. 2002; 62(17): 5049-5057.
6. Tendler A, Kaufman HL, Kadish AS. Increased carcinoembryonic antigen expression in cervical intraepithelial neoplasia grade 3 and in cervical squamous cell carcinoma. Human Pathology. 2000; 31(11): 1357-1362.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Chhieng DC, Rodriguez-Burford C, Tellez LI, et al. Expression of CEA, Tag-72 and Lewis-Y antigen in primary and metastatic lesions of ovarian carcinoma. Human Pathology. 2003; 34(10): 1016-1021.
5. Kass ES, Greiner JW, Kantor JA, et al. Carcinoembryonic antigen as a target for specific antitumor immunotherapy of head and neck cancer. Cancer Research. 2002; 62(17): 5049-5057.
6. Tendler A, Kaufman HL, Kadish AS. Increased carcinoembryonic antigen expression in cervical intraepithelial neoplasia grade 3 and in cervical squamous cell carcinoma. Human Pathology. 2000; 31(11): 1357-1362.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Chhieng DC, Rodriguez-Burford C, Tolley LI, et al. Expression of CEA, Tag-72 and Lewis-Y antigen in primary and metastatic lesions of ovarian carcinoma. Human Pathology. 2003; 34(10): 1016-1021.
5. Kass ES, Greiner JW, Kantor JA, et al. Carcinoembryonic antigen as a target for specific antitumor immunotherapy of head and neck cancer. Cancer Research. 2002; 62(17): 5049-5057.
6. Tendler A, Kaufman HL, Kadish AS. Increased carcinoembryonic antigen expression in cervical intraepithelial neoplasia grade 3 and in cervical squamous cell carcinoma. Human Pathology. 2000; 31(11): 1357-1362.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Chhieng DC, Rodriguez-Burford C, Telley LI, et al. Expression of CEA, Tag-72 and Lewis-Y antigen in primary and metastatic lesions of ovarian carcinoma. Human Pathology. 2003; 34(10): 1016-1021.

5. Kass ES, Greiner JW, Kantor JA, et al. Carcinoembryonic antigen as a target for specific antitumor immunotherapy of head and neck cancer. *Cancer Research*. 2002; 62(17): 5049-5057.
6. Tendler A, Kaufman HL, Kadish AS. Increased carcinoembryonic antigen expression in cervical intraepithelial neoplasia grade 3 and in cervical squamous cell carcinoma. *Human Pathology*. 2000; 31(11): 1357-1362.

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Literatúra

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova P. Rational Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Chhieng DC, Rodriguez-Burford C, Telley LI, et al. Expression of CEA, Tag-72 and Lewis-Y antigen in primary and metastatic lesions of ovarian carcinoma. Human Pathology. 2003; 34(10): 1016-1021.
5. Kass ES, Greiner JW, Kantor JA, et al. Carcinoembryonic antigen as a target for specific antitumor immunotherapy of head and neck cancer. Cancer Research. 2002; 62(17): 5049-5057.
6. Tendler A, Kaufman HL, Kadish AS. Increased carcinoembryonic antigen expression in cervical intraepithelial neoplasia grade 3 and in cervical squamous cell carcinoma. Human Pathology. 2000; 31(11): 1357-1362.

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Leica Biosystems Inc
Pty Ltd