BIO SYSTEMS

BOND Primary Antibody Diluent

Catalog No: AR9352

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

EN FR IT DE ES PT SV EL DA NL
NO TR BG HU RO RU PL SL CS SK AR

Instructions for Use

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Catalog No: AR9352

Intended Use

This reagent is for in vitro diagnostic use.

BOND Primary Antibody Diluent is a ready to use reagent for the dilution of individual primary antibody concentrates for use on the BOND system.

Not to be used for the reconstitution of lyophilized antibodies.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation).

The BOND system is an automated method of immunostaining formalin-fixed, paraffin embedded tissue sections. BOND Primary Antibody Diluent has been developed to ensure the optimal performance of primary antibody concentrates on the BOND system, together with the appropriate BOND detection system and ancillary reagents.

Reagents Provided

BOND Primary Antibody Diluent contains Tris-buffered saline, surfactant, protein stabilizer and 0.35% ProClin™ 950. Total volume = 500 mL.

Dilution and Mixing

BOND Primary Antibody Diluent is ready to use. Dilution of this reagent is not recommended.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system.

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the bottle label. The signs indicating contamination and/or instability of BOND Primary Antibody Diluent are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methylisothiazol-3(2H)-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Other buffers should not be used to replace BOND Primary Antibody Diluent for the dilution of primary antibodies on the BOND system.
- BOND Primary Antibody Diluent was developed to facilitate the use of primary antibodies from the user's laboratory on the BOND system. The user must conduct preliminary experiments to determine the appropriate concentration of primary antibody with the specific BOND detection system to be employed (see Instructions For Use following).

Instructions for Use

Dilute primary antibody concentrate to optimal concentration with BOND Primary Antibody Diluent. Transfer the diluted antibody to a BOND Open container and use as the primary antibody reagent. It must be registered by the BOND system using the hand held bar code scanner according to the instructions given in your BOND user documentation.

The user must conduct preliminary experiments to optimize the concentration of primary antibody to be used with the appropriate BOND detection system. Serial dilutions must be tested using the appropriate positive control tissue for the particular primary antibody (see "Quality Control" in the "Using BOND Reagents" section of your user documentation).

Titration of concentrated antibodies is usually accomplished by serial two-fold dilutions. When preparing dilutions for titration, ensure that allowance is made for the dead volume of the appropriate BOND Open container (for values refer to your BOND user documentation).

Product Specific Limitations

The appropriate concentration of user's own primary antibodies may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

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Leica Biosystems Inc
Pty Ltd