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BOND™ Ready-To-Use Primary Antibody Human Growth Hormone (Polyclonal)

Catalog No: PA0704

Intended Use

This reagent is for in vitro diagnostic use.

Human Growth Hormone (Polyclonal) antibody is intended to be used for the qualitative identification by light microscopy of human growth hormone in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Human Growth Hormone (Polyclonal) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human growth hormone is achieved by first, allowing the binding of Human Growth Hormone (Polyclonal) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

Human Growth Hormone (Polyclonal) is a rabbit anti-human polyclonal supplied as rabbit serum diluted in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

Polyclonal.

Specificity

Human growth hormone.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Not Applicable.

Dilution and Mixing

Human Growth Hormone (Polyclonal) primary antibody is optimally diluted for use on the BOND system. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system.

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Human Growth Hormone (Polyclonal) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Human Growth Hormone (Polyclonal) primary antibody was developed for use on the automated BOND system in combination with BOND Polymer Refine Detection. The recommended staining protocol for Human Growth Hormone (Polyclonal) primary antibody is IHC Protocol F. No pretreatment is recommended.

Results Expected

Normal Tissues

The Human Growth Hormone polyclonal antibody detected human growth hormone in somatotrophic cells of the anterior pituitary gland. No other normal tissues stained positively (Total number of cases = 64).

Tumor Tissues

The Human Growth Hormone polyclonal antibody gave no tumor specific staining in the cases evaluated, including squamous cell carcinomas, adenocarcinomas, sarcomas, melanomas and lymphomas (Total number of cases stained = 44).

Human Growth Hormone (Polyclonal) is recommended for the detection of human growth hormone by qualitative immunohistochemistry on normal and neoplastic formalin-fixed, paraffin-embedded tissue sections, to be viewed by light microscopy.

Product Specific Limitations

Human Growth Hormone (Polyclonal) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Dejneka NS, Auricchio A, Maguire AM, et al. Pharmacologically regulated gene expression in the retina following transduction with viral vectors. Gene Therapy 2001; 8(6):442-446.
5. Osamura RY, Tahara S, Kurotani R, et al. Contributions of immunohistochemistry and in situ hybridization to the functional analysis of pituitary adenomas. Journal of Histochemistry and Cytochemistry. 2000; 48(4):445-458.

ProClin 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Dejneka NS, Auricchio A, Maguire AM, et al. Pharmacologically regulated gene expression in the retina following transduction with viral vectors. Gene Therapy 2001; 8(6):442–446.
5. Osamura RY, Tahara S, Kurotani R, et al. Contributions of immunohistochemistry and in situ hybridization to the functional analysis of pituitary adenomas. Journal of Histochemistry and Cytochemistry. 2000; 48(4):445–458.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 1763 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Dejneka NS, Auricchio A, Maguire AM, et al. Pharmacologically regulated gene expression in the retina following transduction with viral vectors. Gene Therapy 2001; 8(6):442–446.
5. Osamura RY, Tahara S, Kurotani R, et al. Contributions of immunohistochemistry and in situ hybridization to the functional analysis of pituitary adenomas. Journal of Histochemistry and Cytochemistry. 2000; 48(4):445–458.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd