BOND Ready-to-Use ISH

Kappa Probe

Catalog No: PB0645

Leica Biosystems Newcastle Ltd
Balliol Business Park
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
+44 191 215 4242

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Catalog No: PB0645

Intended Use

This reagent is for *in vitro* diagnostic use.

Kappa Probe is intended for use in the qualitative identification of Kappa light chain messenger RNA (mRNA) in formalin-fixed, paraffin-embedded tissue by *in situ* hybridization (ISH) using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunoglobulins are glycoproteins produced in mature B-cells against a specific antigen. Each individual immunoglobulin molecule comprises two heavy and two light polypeptide chains. There are five classes of immunoglobulin, determined by the type of heavy chain. In contrast, there are only two types of light chain: Kappa or Lambda. Each individual immunoglobulin molecule is composed of one of five classes of heavy chains and either Kappa or Lambda light chains. In normal human lymphoid populations, the ratio of Kappa to Lambda light chains is approximately 2:1¹.

B-cell neoplasms are thought to arise from a single transformed cell (monoclonal)². In contrast, reactive states result in proliferation of a number of B-cells (polyclonal). Since immunoglobulins from the same B-cell contain either Kappa or Lambda light chains, light chain restriction or monoclonality can be used to make the distinction between reactive and neoplastic B cell proliferations³.

Kappa Probe is used in conjunction with Lambda Probe for the detection of antibody producing B-cells in formalin-fixed, paraffin-embedded tissue. The method is reproducible and should result in intense, brown cytoplasmic staining of cells containing Kappa mRNA.

Reagents Provided

Kappa Probe is a fluorescein-conjugated oligonucleotide probe supplied in hybridization solution.

Total volume = 5.5 mL

Dilution and Mixing

Kappa Probe is ready to use. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Needed but not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and *in situ* hybridization staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8 °C. The product is stable under these conditions up to the expiry date indicated on the vial label.

There are no obvious signs indicating contamination and/or instability. Appropriate positive and negative tissue controls should be run at the same time as test tissue.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user⁴.

Precautions

- This product is intended for *in vitro* diagnostic use.

KAPPA PROBE

Contains Formaldehyde (<50%) and Dextran Sulphate (<30%).

GHS07: Exclamation mark.

GHS08: Health hazard.

Signal words: Danger.

H360D: May damage the unborn child.

H315: Causes skin irritation.

H319: Causes serious eye irritation.

P201: Obtain special instructions before use.

P202: Do not handle until all safety precautions have been read and understood.

P260: Do not breathe dust/fumes/gas/mist/vapours/spray.

P308+313: IF exposed or concerned: Get medical advice/attention.

P314: Get medical advice/attention if you feel unwell.

P264: Wash hands thoroughly after handling.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P305+351+338: IF IN EYES: Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to do. Continue rinsing.

P337+313: If eye irritation persists: Get medical advice/attention.

P302+352: IF ON SKIN: Wash with plenty of soap and water.

P332+313: If skin irritation occurs: Get medical advice/attention.

P362+P364: Take off contaminated clothing and wash it before reuse. Restricted to professional users.

- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions⁵. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.

- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Kappa Probe was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), in combination with Anti-Fluorescein Antibody and BOND Polymer Refine Detection. The recommended staining protocol for Kappa Probe is ISH Protocol A. Enzyme retrieval is recommended using the BOND Enzyme Pretreatment Kit, Enzyme 1, for 15 minutes.

Appropriate tissue and reagent controls should always be used. The protocol for the tissue and reagent controls should correspond to that of the RNA test probe.

Quality Control

Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls in addition to the following procedures.

Positive Tissue Control

Used to indicate correctly prepared tissues and proper staining techniques.

One positive tissue control should be included for each set of test conditions in each staining run. A tissue with weak positive staining is more suitable than a tissue with strong positive staining for optimal quality control and to detect minor levels of reagent degradation.

Negative Tissue Control

Should be examined after the positive tissue control to verify the specificity of the labeling of probe to the target.

Alternatively, the variety of different cell types present in most tissue sections frequently offers negative control sites, but this should be verified by the user.

Negative Reagent Control

Use RNA Negative Control Probe PB0809 in place of the Kappa Probe with a section of each patient specimen to evaluate non-specific staining and allow better interpretation of specific staining at the target.

Positive Reagent Control

Use RNA Positive Control Probe PB0785 in place of the Kappa Probe with a section of each patient specimen to provide information on the preservation of nucleic acids in the tissue as well as accessibility of nucleic acids to the probe. If the RNA Positive Control Probe fails to demonstrate positive staining, results with the test specimens should be considered invalid.

Patient Tissue

Examine patient specimens stained with Kappa Probe last. Positive staining intensity should be assessed within the context of any non-specific background staining of the RNA Negative Control Probe PB0809.

Results Expected

Normal Tissues

PB0645 detected Kappa mRNA in B-cells in spleen and tonsil as well as the occasional B-cell in a wide range of tissues including breast, cervix, liver, salivary gland, parotid gland, small intestine and large intestine. Precursor B-cells were also positive for Kappa mRNA in bone marrow. Additional staining was also seen in epithelial cells of the thyroid, prostate, larynx, esophagus and endometrium, tubules of the kidney, secretory cells in the adrenal glands and Leydig cells. (Total number of normal cases evaluated = 110).

Abnormal Tissues

PB0645 stained 2/6 plasma cell myelomas. No staining was seen in diffuse large B-cell lymphomas (0/3), T-cell lymphomas (0/3), Hodgkin's lymphomas (0/2), a Burkitt's lymphoma (0/1), lung tumors (0/3), liver tumors (0/3), ovarian tumors (0/3), thyroid tumors (0/3), esophageal tumors (0/2), breast tumors (0/2), metastatic tumors of unknown origin (0/2), brain tumors (0/2), testicular tumors (0/2), skin tumors (0/2), a colon tumor (0/1), a stomach tumor (0/1), a rectal tumor (0/1), a tumor of the larynx (0/1) and a tumor of the thymus (0/1). (Total number of abnormal cases evaluated = 44).

PB0645 is recommended for the detection of Kappa light chain mRNA.

Product Specific Limitations

Kappa Probe has been optimised at Leica Biosystems for use with Anti-Fluorescein Antibody, BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of enzymatic digestion, and must be determined empirically. RNA Negative Control Probe should be used when optimising retrieval conditions and protocol times.

Troubleshooting

Reference 6 may aid in remedial action.

Test samples should be complemented by the appropriate tissue and reagent controls.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on *in situ* hybridization with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

References

1. Levy R, Warnke R, Dorfman RF, Haimovich J. The monoclonality of human B-cell lymphomas. *J Exp Med.* 1977;145:1014–1024.
2. Friedman JM, Fialow PJ. Cell marker studies of human tumorigenesis. *Transplant Rev.* 1976;28:17–33.
3. McNicol AM, Farquharson MA, Lee FD, Foulis AK. Comparison of *in situ* hybridization and polymerase chain reaction in the diagnosis of B cell lymphoma. *J Clin Pathol.* 1998;51:229–233.
4. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
5. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
6. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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Bibliografia

1. Levy R, Warnke R, Dorfman RF, Haimovich J. The monoclonality of human B-cell lymphomas. *J Exp Med.* 1977;145:1014–1024.
2. Friedman JM, Fialow PJ. Cell marker studies of human tumorigenesis. *Transplant Rev.* 1976;28:17–33.
3. McNicol AM, Farquharson MA, Lee FD, Foulis AK. Comparison of *in situ* hybridization and polymerase chain reaction in the diagnosis of B cell lymphoma. *J Clin Pathol.* 1998;51:229–233.
4. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
5. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
6. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ Hybridization A practical approach.* 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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1. Levy R, Warnke R, Dorfman RF, Haimovich J. The monoclonality of human B-cell lymphomas. *J Exp Med*. 1977;145:1014–1024.
2. Friedman JM, Fialow PJ. Cell marker studies of human tumorigenesis. *Transplant Rev*. 1976;28:17–33.
3. McNicol AM, Farquharson MA, Lee FD, Foulis AK. Comparison of *in situ* hybridization and polymerase chain reaction in the diagnosis of B cell lymphoma. *J Clin Pathol*. 1998;51:229–233.
4. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
5. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
6. Wilkinson DG. The theory and practice of *in situ* hybridization. In: *In situ Hybridization A practical approach*. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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1. Levy R, Warnke R, Dorfman RF, Haimovich J. The monoclonality of human B-cell lymphomas. J Exp Med. 1977;145:1014–1024.
2. Friedman JM, Fialow PJ. Cell marker studies of human tumorigenesis. Transplant Rev. 1976;28:17–33.
3. McNicol AM, Farquharson MA, Lee FD, Foulis AK. Comparison of *in situ* hybridization and polymerase chain reaction in the diagnosis of B cell lymphoma. J Clin Pathol. 1998;51:229–233.
4. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
5. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
6. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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1. Levy R, Warnke R, Dorfman RF, Haimovich J. The monoclonality of human B-cell lymphomas. J Exp Med. 1977;145:1014–1024.
2. Friedman JM, Fialow PJ. Cell marker studies of human tumorigenesis. Transplant Rev. 1976;28:17–33.
3. McNicol AM, Farquharson MA, Lee FD, Foulis AK. Comparison of *in situ* hybridization and polymerase chain reaction in the diagnosis of B cell lymphoma. J Clin Pathol. 1998;51:229–233.
4. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.

5. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
6. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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Leica Biosystems Inc
Pty Ltd