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www.LeicaBiosystems.com

Estrogen Receptor Clone 6F11

Ready-To-Use Primary Antibody For Bond™

Catalog No: PA0151 (7 mL) / PA0009 (30 mL)

Intended Use

This reagent is for *in vitro* diagnostic use.

Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human estrogen receptor (ER) in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond specifically binds to the ER antigen located in the nucleus of ER positive normal and neoplastic cells.

ER (6F11) is indicated as an aid in the management, prognosis and prediction of therapy outcome of breast cancer. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

ER content of breast cancer tissue is an important parameter in the prediction of prognosis and response to endocrine therapy. The introduction of monoclonal antibodies to ER has allowed the determination of receptor status of breast tumors to be carried out in routine histopathology laboratories. ER (6F11) is a mouse monoclonal antibody directed against the human estrogen receptor molecule. A prokaryotic recombinant protein, corresponding to the full length human ER molecule was used as the immunogen. ER (6F11) has been shown to react with a 66 kD protein from MCF-7 cell lysates via Western blot.⁴

Principle of Procedure

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond is a ready to use product that has been specifically optimized for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with Bond Polymer Refine Detection. The recommended staining protocol for Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond is IHC Protocol F. Heat induced epitope retrieval is recommended using Bond Epitope Retrieval Solution 1 for 20 minutes. Bond Polymer Refine Detection utilizes a novel controlled polymerization technology to prepare polymeric HRP-linker antibody conjugates. The detection system avoids the use of streptavidin and biotin, and therefore eliminates nonspecific staining as a result of endogenous biotin.

Bond Polymer Refine Detection works as follows:

- The specimen is incubated with hydrogen peroxide to quench endogenous peroxidase activity
- Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond is applied
- A post primary antibody solution enhances penetration of the subsequent polymer reagent
- A poly-HRP anti-mouse/rabbit IgG reagent localizes the primary antibody
- The substrate chromogen, 3,3'-diaminobenzidine (DAB), visualizes the complex via a brown precipitate
- Hematoxylin (blue) counterstaining allows the visualization of cell nuclei.

Using Bond Polymer Refine Detection in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagent Provided

Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL / 30 mL.

Clone

6F11.

Immunogen

Prokaryotic recombinant protein corresponding to the full length alpha form of the human estrogen receptor molecule.

Specificity

Human estrogen receptor.

Subclass

IgG1.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 0.88 mg/L as determined by ELISA.

Method

ER (6F11) was raised against recombinant ER protein derived from the mRNA of MCF-7 cells. Balb/c mice were immunized with the resulting (His)6-ER recombinant antigen. Screening was conducted by ELISA, with ELISA positive supernatants tested in IHC on formalin-fixed, paraffin-embedded sections of breast carcinoma of known receptor status. Colonies demonstrating positive immunohistochemical staining were cloned by limiting dilution.

Dilution and Mixing

Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond is optimally diluted for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with Bond Polymer Refine Detection. Further dilution may result in loss of antigen staining. The user must validate any such change. Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results necessitating regular performance of in-house controls. Refer to "Using BOND Reagents" in your BOND user documentation.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user.

Specimen Preparation

The recommended fixative is 10% neutral-buffered formalin for paraffin-embedded tissue sections.

Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with Bond Polymer Refine Detection. The recommended staining protocol for Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond is IHC Protocol F. Heat induced epitope retrieval is recommended using Bond Epitope Retrieval Solution 1 for 20 minutes.

Quality Control

Refer to "Using BOND Reagents" in your BOND user documentation.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Interpretation of Staining

Refer to "Using BOND Reagents" in your BOND user documentation.

General Limitations

- Immunohistochemistry is a multistep diagnostic process that consists of specialized training in the selection of the appropriate reagents; tissue selection, fixation, and processing; preparation of the IHC slide; and interpretation of the staining results.
- Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning or contamination with other tissues or fluids may produce artifacts, antibody trapping, or false negative results. Inconsistent results may be due to variations in fixation and embedding methods, or to inherent irregularities within the tissue.
- Excessive or incomplete counterstaining may compromise proper interpretation of results.

- The clinical interpretation of any positive or negative staining should be evaluated within the context of clinical presentation, morphology and other histopathological criteria. The clinical interpretation of any positive or negative staining should be complemented by morphological studies using proper positive and negative internal and external controls as well as other diagnostic tests. It is the responsibility of a qualified pathologist who is familiar with the proper use of IHC antibodies, reagents and methods to interpret all of the steps used to prepare and interpret the final IHC preparation.
- The manufacturer provides these antibodies/reagents at optimal dilution for use following the provided instructions for IHC on prepared tissue sections or cytologic preparation. Any deviation from recommended test procedures may invalidate declared expected results; appropriate controls must be employed and documented. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results.
- This product is not intended for use in flow cytometry. Performance characteristics have not been determined for flow cytometry. Tissues from persons infected with hepatitis B virus and containing hepatitis B surface antigen (HBsAg) may exhibit nonspecific staining with horseradish peroxidase.
- Reagents may demonstrate unexpected reactions in previously untested tissues. The possibility of unexpected reactions even in tested tissue groups cannot be completely eliminated due to biological variability of antigen expression in neoplasms, or other pathological tissues.
- Normal/nonimmune sera from the same animal source as secondary antisera used in blocking steps may cause false-negative or false-positive results due to autoantibodies or natural antibodies.
- False-positive results may be seen due to non-immunological binding of proteins or substrate reaction products. They may also be caused by pseudoperoxidase activity (erythrocytes), endogenous peroxidase activity (cytochrome C), or endogenous biotin (e.g. liver, breast, brain, kidney) depending on the type of immunostain used.

Product Specific Limitations

Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond has been optimized at Leica Biosystems for use with Bond Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Performance Characteristics

Reproducibility

Intra run reproducibility of staining with Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond was determined by staining 10 sections of the same tissue using Leica Biosystems Bond Polymer Refine Detection (DS9800). 10 of 10 slides stained positively. All slides stained with similar staining specificity and intensity (varied by <1).

Inter-run reproducibility of staining with Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond was determined by staining 10 sections of the same tissue, on 3 different staining runs using Leica Biosystems Bond Polymer Refine Detection (DS9800). 10 of 10 slides stained positively on each run. All slides stained with similar staining specificity and intensity (varied by <1).

Immunoreactivity

Table 1: Reactivity of Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond on Normal Tissues

Tissue	Number of Cases	Description of Staining	Staining Intensity (0-3+)
Adrenal	3	No staining of tissue elements	0
Brain, Cerebellum	3	No staining of tissue elements	0
Brain, Cerebrum	3	No staining of tissue elements	0
Breast	3	Duct nuclei in 2/3 tissues	3+
Cervix	3	Variable nuclei staining in the stratified squamous epithelium of the ectocervix, cervical stromal cells and glandular tissue	1–2+
Colon	3	No staining of tissue elements	0
Esophagus	3	No staining of tissue elements	0
Heart	3	No staining of tissue elements	0
Kidney	3	No staining of tissue elements	0
Liver	3	No staining of tissue elements	0
Lung	3	No staining of tissue elements	0
Mesothelial cells	1	No staining of tissue elements	0
Ovary	3	Faint variable nuclear staining of a small percentage of ovarian stromal cells. Positive nuclear staining of the <i>theca externa</i> branching into the ovarian stroma in a percentage of ovarian follicles.	1–2+

Tissue	Number of Cases	Description of Staining	Staining Intensity (0-3+)
Pancreas	3	No staining of tissue elements	0
Peripheral nerve	3	No staining of tissue elements	0
Pituitary	3	No staining of tissue elements	0
Prostate	3	No staining of tissue elements	0
Salivary/ Submandibular Gland	3	No staining of tissue elements	0
Skeletal Muscle	3	No staining of tissue elements	0
Skin	3	No staining of tissue elements	0
Small Intestine	3	No staining of tissue elements	0
Spleen	3	No staining of tissue elements	0
Stomach	3	No staining of tissue elements	0
Testis	3	No staining of tissue elements	0
Thymus	2	No staining of tissue elements	0
Thyroid	3	No staining of tissue elements	0
Tonsil	3	No staining of tissue elements	0
Uterus	3	Variable nuclei staining of endometrial/myometrial stromal cells	1+
Bone Marrow	3	No staining of tissue elements	0

Key to Staining Intensity:

- 0 - Negative
- 1+ - Weak
- 2+ - Moderate
- 3+ - Strong

Normal Tissues

Clone 6F11 detects the estrogen receptor alpha antigen in the nuclei of cells that express high levels of ER, a proportion of endometrial, ovarian and myometrial cells, and normal breast ductal cells.

Published Immunoreactivity

Characterization of ER (6F11) during antibody development included a comparative evaluation of a series of 55 sequential breast carcinomas. The tissues evaluated were routinely processed formalin-fixed, paraffin-embedded specimens stained using both ER (6F11) and ER (1D5). There was an observed concordance of staining for 50/55 cases⁴.

Estrogen receptor status was evaluated in 592 cases using routinely prepared paraffin-embedded tissue samples from primary breast carcinomas with ER (6F11) and anti-ER (1D5). Overall, ER (1D5) and ER (6F11) showed a 97.5% concordance rate¹⁶.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Bevitt DJ, Milton ID, Piggot N et al. New monoclonal antibodies to oestrogen and progesterone receptors effective for paraffin section immunohistochemistry. *Journal of Pathology* 1997; 183(2), 228–232.
5. Diaz L, Sahin A and Sneige N. Immunohistochemical detection of estrogen receptor in breast cancer: a laboratory quality improvement study. United States and Canadian Academy of Pathology (Annual Meeting Abstracts March 22–28), 2003; 27A.
6. Dabbs DJ, Landrenau RJ, Liu Y et al. Detection of estrogen receptor by immunohistochemistry in pulmonary adenocarcinoma. *Ann Thorac Surg* 2002; 73(2), 403–405.
7. Khan SA, Yee KA, Kaplan C et al. Estrogen receptor alpha expression in normal human breast epithelium is consistent over time. *International Journal of Cancer* 2002; 102(4), 334–337.
8. Radzikowska E, Langfort R and Giedronowicz D. Estrogen and progesterone receptors in non small cell lung cancer patients. *Ann Thorac Cardiovasc Surg* 2002; 8(2), 69–73.
9. Leav I, Lau KM, Adams JY, et al. Comparative studies of the estrogen receptors beta and alpha and the androgen receptor in normal human prostate glands, dysplasia, and in primary and metastatic carcinoma. *American Journal of Pathology* 2001; 159(1), 79–92.
10. Braidman IP, Baris C, Selby PL et al. Preliminary report of impaired oestrogen receptor-alpha expression in bone, but no involvement of androgen receptor, in male idiopathic osteoporosis. *Journal of Pathology* 2000; 192, 90–96.
11. de las Mulas JM, van Niel M, Millan Y et al. Immunohistochemical analysis of estrogen receptors in feline mammary gland benign and malignant lesions: comparison with biochemical assay. *Domest Anim Endocrinol* 2000; 18(1), 111–125.
12. Khan SA, Rogers MA, Khuruna KK et al. Oestrogen receptor expression in normal breast epithelium. *European Journal of Cancer* 2000; 36(Suppl 4), S27-S28.
13. Leake R, Barnes D, Pinder S et al. Immunohistochemical detection of steroid receptors in breast cancer: a working protocol. *Journal of Clinical Pathology* 2000; 53(8), 634–635.
14. Zafrani B, Aubriot MH, Mouret E et al. High sensitivity and specificity of immunohistochemistry for the detection of hormone receptors in breast carcinoma: comparison with biochemical determination in a prospective study of 793 cases. *Histopathology* 2000; 37(6), 536–545.
15. Harvey JM, Clark GM, Osborne CK et al. Estrogen receptor status by immunohistochemistry is superior to the ligand-binding assay for predicting response to adjuvant endocrine therapy in breast cancer. *Journal of Clinical Oncology* 1999; 17(5), 1474–1481.
16. Kaplan, P.A. et al. 1D5 and 6F11 an immunohistochemical comparison of two monoclonal antibodies for the evaluation of estrogen receptor status in primary breast carcinoma. *Am J Clin Pathol* 2005: 276–280.

Date of Issue

27 June 2016

Estrogen Receptor Clone 6F11 Ready-To-Use Primary Antibody For Bond™

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991;7(9). Référence à commander : M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4ème édition. Churchill Livingstone, New York. 1996.
4. Bevitt DJ, Milton ID, Piggot N, et al. New monoclonal antibodies to oestrogen and progesterone receptors effective for paraffin section immunohistochemistry. *Journal of Pathology* 1997;183(2):228–232.
5. Diaz L, Sahin A and Sneige N. Immunohistochemical detection of estrogen receptor in breast cancer: a laboratory quality improvement study. *United States and Canadian Academy of Pathology* 2003;27A (résumés de la conférence annuelle, 22–28 mars).
6. Dabbs DJ, Landrenau RJ, Liu Y, et al. Detection of estrogen receptor by immunohistochemistry in pulmonary adenocarcinoma. *Ann Thorac Surg* 2002;73(2):403–405.
7. Khan SA, Yee KA, Kaplan C, et al. Estrogen receptor alpha expression in normal human breast epithelium is consistent over time. *International Journal of Cancer* 2002;102(4):334–337.
8. Radzikowska E, Langfort R and Giedronowicz D. Estrogen and progesterone receptors in non small cell lung cancer patients. *Ann Thorac Cardiovasc Surg* 2002;8(2):69–73.
9. Leav I, Lau KM, Adams JY, et al. Comparative studies of the estrogen receptors beta and alpha and the androgen receptor in normal human prostate glands, dysplasia, and in primary and metastatic carcinoma. *American Journal of Pathology* 2001;159(1):79–92.
10. Braidman IP, Baris C, Selby PL, et al. Preliminary report of impaired oestrogen receptor-alpha expression in bone, but no involvement of androgen receptor, in male idiopathic osteoporosis. *Journal of Pathology* 2000;192:90–96.
11. de las Mulas JM, van Niel M, Millan Y, et al. Immunohistochemical analysis of estrogen receptors in feline mammary gland benign and malignant lesions: comparison with biochemical assay. *Domest Anim Endocrinol* 2000;18(1):111–125.
12. Khan SA, Rogers MA, Khuruna KK, et al. Oestrogen receptor expression in normal breast epithelium. *European Journal of Cancer* 2000;36 suppl 4:S27–S28.
13. Leake R, Barnes D, Pinder S, et al. Immunohistochemical detection of steroid receptors in breast cancer: a working protocol. *Journal of Clinical Pathology* 2000;53(8):634–635.
14. Zafrani B, Aubriot MH, Mouret E, et al. High sensitivity and specificity of immunohistochemistry for the detection of hormone receptors in breast carcinoma: comparison with biochemical determination in a prospective study of 793 cases. *Histopathology* 2000;37(6):536–545.
15. Harvey JM, Clark GM, Osborne CK, et al. Estrogen receptor status by immunohistochemistry is superior to the ligand-binding assay for predicting response to adjuvant endocrine therapy in breast cancer. *Journal of Clinical Oncology* 1999;17(5):1474–1481.
16. Kaplan PA, et al. 1D5 and 6F11 an immunohistochemical comparison of two monoclonal antibodies for the evaluation of estrogen receptor status in primary breast carcinoma. *Am J Clin Pathol* 2005:276–280.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 28. February 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Bestellnummer M29-P.
3. Bancroft JD and Stevens A. Stevens. Theory and Practice of Histological Techniques. 4. Auflage. Churchill Livingstone, New York. 1996.
4. Bevitt DJ, Milton ID, Piggot N et al. Piggot et al. New monoclonal antibodies to oestrogen and progesterone receptors effective for paraffin section immunohistochemistry. *Journal of Pathology* 1997; 183(2), 228–232.
5. Diaz L, Sahin A and Sneige N. Immunohistochemical detection of estrogen receptor in breast cancer: a laboratory quality improvement study. United States and Canadian Academy of Pathology (Annual Meeting Abstracts March 22–28), 2003; 27A.
6. Dabbs DJ, Landrenau RJ, Liu Y et al. Detection of estrogen receptor by immunohistochemistry in pulmonary adenocarcinoma. *Ann Thorac Surg* 2002; 73(2), 403–405.
7. Khan SA, Yee KA, Kaplan C et al. Estrogen receptor alpha expression in normal human breast epithelium is consistent over time. *International Journal of Cancer* 2002; 102(4), 334–337.
8. Radzikowska E, Langfort R and Giedronowicz D. Estrogen and progesterone receptors in non small cell lung cancer patients. *Ann Thorac Cardiovasc Surg* 2002; 8(2), 69–73.
9. Leav I, Lau KM, Adams JY, et al. Comparative studies of the estrogen receptors beta and alpha and the androgen receptor in normal human prostate glands, dysplasia, and in primary and metastatic carcinoma. *American Journal of Pathology* 2001; 159(1), 79–92.
10. Braidman IP, Baris C, Selby PL et al. Preliminary report of impaired oestrogen receptor-alpha expression in bone, but no involvement of androgen receptor, in male idiopathic osteoporosis. *Journal of Pathology* 2000; 192, 90–96.
11. de las Mulas JM, van Niel M, Millan Y et al. Immunohistochemical analysis of estrogen receptors in feline mammary gland benign and malignant lesions: comparison with biochemical assay. *Domest Anim Endocrinol* 2000; 18(1), 111–125.
12. Khan SA, Rogers MA, Khuruna KK et al. Oestrogen receptor expression in normal breast epithelium. *European Journal of Cancer* 2000; 36(Suppl 4), S27-S28.
13. Leake R, Barnes D, Pinder S et al. Immunohistochemical detection of steroid receptors in breast cancer: a working protocol. *Journal of Clinical Pathology* 2000; 53(8), 634–635.
14. Zafrani B, Aubriot MH, Mouret E et al. High sensitivity and specificity of immunohistochemistry for the detection of hormone receptors in breast carcinoma: comparison with biochemical determination in a prospective study of 793 cases. *Histopathology* 2000; 37(6), 536–545.
15. Harvey JM, Clark GM, Osborne CK et al. Estrogen receptor status by immunohistochemistry is superior to the ligand-binding assay for predicting response to adjuvant endocrine therapy in breast cancer. *Journal of Clinical Oncology* 1999; 17(5), 1474–1481.
16. Kaplan, P.A. et. al. 1D5 and 6F11 an immunohistochemical comparison of two monoclonal antibodies for the evaluation of estrogen receptor status in primary breast carcinoma. *Am J Clin Pathol* 2005: 276–280.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. Bevitt DJ, Milton ID, Piggot N et al. New monoclonal antibodies to oestrogen and progesterone receptors effective for paraffin section immunohistochemistry. Journal of Pathology 1997; 183(2), 228–232.
4. Diaz L, Sahin A and Sneige N. Immunohistochemical detection of estrogen receptor in breast cancer: a laboratory quality improvement study. United States and Canadian Academy of Pathology (Annual Meeting Abstracts March 22–28), 2003; 27A.

5. Dabbs DJ, Landrenau RJ, Liu Y et al. Detection of estrogen receptor by immunohistochemistry in pulmonary adenocarcinoma. *Ann Thorac Surg* 2002; 73(2), 403–405.
6. Khan SA, Yee KA, Kaplan C et al. Estrogen receptor alpha expression in normal human breast epithelium is consistent over time. *International Journal of Cancer* 2002; 102(4), 334–337.
7. Radzikowska E, Langfort R and Giedronowicz D. Estrogen and progesterone receptors in non small cell lung cancer patients. *Ann Thorac Cardiovasc Surg* 2002; 8(2), 69–73.
8. Leav I, Lau KM, Adams JY, et al. Comparative studies of the estrogen receptors beta and alpha and the androgen receptor in normal human prostate glands, dysplasia, and in primary and metastatic carcinoma. *American Journal of Pathology* 2001; 159(1), 79–92.
9. Braidman IP, Baris C, Selby PL et al. Preliminary report of impaired oestrogen receptor-alpha expression in bone, but no involvement of androgen receptor, in male idiopathic osteoporosis. *Journal of Pathology* 2000; 192, 90–96.
10. de las Mulas JM, van Niel M, Millan Y et al. Immunohistochemical analysis of estrogen receptors in feline mammary gland benign and malignant lesions: comparison with biochemical assay. *Domest Anim Endocrinol* 2000; 18(1), 111–125.
11. Khan SA, Rogers MA, Khurana KK et al. Oestrogen receptor expression in normal breast epithelium. *European Journal of Cancer* 2000; 36(Suppl 4), S27-S28.
12. Leake R, Barnes D, Pinder S et al. Immunohistochemical detection of steroid receptors in breast cancer: a working protocol. *Journal of Clinical Pathology* 2000; 53(8), 634–635.
13. Zafrani B, Aubriot MH, Mouret E et al. High sensitivity and specificity of immunohistochemistry for the detection of hormone receptors in breast carcinoma: comparison with biochemical determination in a prospective study of 793 cases. *Histopathology* 2000; 37(6), 536–545.
14. Harvey JM, Clark GM, Osborne CK et al. Estrogen receptor status by immunohistochemistry is superior to the ligand-binding assay for predicting response to adjuvant endocrine therapy in breast cancer. *Journal of Clinical Oncology* 1999; 17(5), 1474–1481.
15. Kaplan, P.A. et. al. 1D5 and 6F11 an immunohistochemical comparison of two monoclonal antibodies for the evaluation of estrogen receptor status in primary breast carcinoma. *Am J Clin Pathol* 2005: 276–280.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Bevitt DJ, Milton ID, Piggot N et al. New monoclonal antibodies to oestrogen and progesterone receptors effective for paraffin section immunohistochemistry. *Journal of Pathology* 1997; 183(2), 228–232.
5. Diaz L, Sahin A and Sneige N. Immunohistochemical detection of estrogen receptor in breast cancer: a laboratory quality improvement study. United States and Canadian Academy of Pathology (Annual Meeting Abstracts March 22–28), 2003; 27A.
6. Dabbs DJ, Landrenau RJ, Liu Y et al. Detection of estrogen receptor by immunohistochemistry in pulmonary adenocarcinoma. *Ann Thorac Surg* 2002; 73(2), 403–405.
7. Khan SA, Yee KA, Kaplan C et al. Estrogen receptor alpha expression in normal human breast epithelium is consistent over time. *International Journal of Cancer* 2002; 102(4), 334–337.
8. Radzikowska E, Langfort R and Giedronowicz D. Estrogen and progesterone receptors in non small cell lung cancer patients. *Ann Thorac Cardiovasc Surg* 2002; 8(2), 69–73.
9. Leav I, Lau KM, Adams JY, et al. Comparative studies of the estrogen receptors beta and alpha and the androgen receptor in normal human prostate glands, dysplasia, and in primary and metastatic carcinoma. *American Journal of Pathology* 2001; 159(1), 79–92.
10. Braidman IP, Baris C, Selby PL et al. Preliminary report of impaired oestrogen receptor-alpha expression in bone, but no involvement of androgen receptor, in male idiopathic osteoporosis. *Journal of Pathology* 2000; 192, 90–96.
11. de las Mulas JM, van Niel M, Millan Y et al. Immunohistochemical analysis of estrogen receptors in feline mammary gland benign and malignant lesions: comparison with biochemical assay. *Domest Anim Endocrinol* 2000; 18(1), 111–125.
12. Khan SA, Rogers MA, Khuruna KK et al. Oestrogen receptor expression in normal breast epithelium. *European Journal of Cancer* 2000; 36(Suppl 4), S27-S28.
13. Leake R, Barnes D, Pinder S et al. Immunohistochemical detection of steroid receptors in breast cancer: a working protocol. *Journal of Clinical Pathology* 2000; 53(8), 634–635.
14. Zafrani B, Aubriot MH, Mouret E et al. High sensitivity and specificity of immunohistochemistry for the detection of hormone receptors in breast carcinoma: comparison with biochemical determination in a prospective study of 793 cases. *Histopathology* 2000; 37(6), 536–545.
15. Harvey JM, Clark GM, Osborne CK et al. Estrogen receptor status by immunohistochemistry is superior to the ligand-binding assay for predicting response to adjuvant endocrine therapy in breast cancer. *Journal of Clinical Oncology* 1999; 17(5), 1474–1481.
16. Kaplan, P.A. et. al. 1D5 and 6F11 an immunohistochemical comparison of two monoclonal antibodies for the evaluation of estrogen receptor status in primary breast carcinoma. *Am J Clin Pathol* 2005: 276–280.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Bevitt DJ, Milton ID, Piggot N et al. New monoclonal antibodies to oestrogen and progesterone receptors effective for paraffin section immunohistochemistry. Journal of Pathology 1997; 183(2), 228–232.
5. Diaz L, Sahin A and Sneige N. Immunohistochemical detection of estrogen receptor in breast cancer: a laboratory quality improvement study. United States and Canadian Academy of Pathology (Annual Meeting Abstracts March 22–28), 2003; 27A.

6. Dabbs DJ, Landrenau RJ, Liu Y et al. Detection of estrogen receptor by immunohistochemistry in pulmonary adenocarcinoma. *Ann Thorac Surg* 2002; 73(2), 403–405.
7. Khan SA, Yee KA, Kaplan C et al. Estrogen receptor alpha expression in normal human breast epithelium is consistent over time. *International Journal of Cancer* 2002; 102(4), 334–337.
8. Radzikowska E, Langfort R and Giedronowicz D. Estrogen and progesterone receptors in non small cell lung cancer patients. *Ann Thorac Cardiovasc Surg* 2002; 8(2), 69–73.
9. Leav I, Lau KM, Adams JY, et al. Comparative studies of the estrogen receptors beta and alpha and the androgen receptor in normal human prostate glands, dysplasia, and in primary and metastatic carcinoma. *American Journal of Pathology* 2001; 159(1), 79–92.
10. Braidman IP Baris C, Selby PL et al. Preliminary report of impaired oestrogen receptor-alpha expression in bone, but no involvement of androgen receptor, in male idiopathic osteoporosis. *Journal of Pathology* 2000; 192, 90–96.
11. de las Mulas JM, van Niel M, Millan Y et al. Immunohistochemical analysis of estrogen receptors in feline mammary gland benign and malignant lesions: comparison with biochemical assay. *Domest Anim Endocrinol* 2000; 18(1), 111–125.
12. Khan SA, Rogers MA, Khuruna KK et al. Oestrogen receptor expression in normal breast epithelium. *European Journal of Cancer* 2000; 36(Suppl 4), S27-S28.
13. Leake R, Barnes D, Pinder S et al. Immunohistochemical detection of steroid receptors in breast cancer: a working protocol. *Journal of Clinical Pathology* 2000; 53(8), 634–635.
14. Zafrani B, Aubriot MH, Mouret E et al. High sensitivity and specificity of immunohistochemistry for the detection of hormone receptors in breast carcinoma: comparison with biochemical determination in a prospective study of 793 cases. *Histopathology* 2000; 37(6), 536–545.
15. Harvey JM, Clark GM, Osborne CK et al. Estrogen receptor status by immunohistochemistry is superior to the ligand-binding assay for predicting response to adjuvant endocrine therapy in breast cancer. *Journal of Clinical Oncology* 1999; 17(5), 1474–1481.
16. Kaplan, P.A. et. al. 1D5 and 6F11 an immunohistochemical comparison of two monoclonal antibodies for the evaluation of estrogen receptor status in primary breast carcinoma. *Am J Clin Pathol* 2005: 276–280.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Bevit DJ, Milton ID, Piggot N et al. New monoclonal antibodies to oestrogen and progesterone receptors effective for paraffin section immunohistochemistry. *Journal of Pathology* 1997; 183(2), 228–232.
5. Diaz L, Sahin A and Sneige N. Immunohistochemical detection of estrogen receptor in breast cancer: a laboratory quality improvement study. United States and Canadian Academy of Pathology (Annual Meeting Abstracts March 22–28), 2003; 27A.
6. Dabbs DJ, Landrenau RJ, Liu Y et al. Detection of estrogen receptor by immunohistochemistry in pulmonary adenocarcinoma. *Ann Thorac Surg* 2002; 73(2), 403–405.
7. Khan SA, Yee KA, Kaplan C et al. Estrogen receptor alpha expression in normal human breast epithelium is consistent over time. *International Journal of Cancer* 2002; 102(4), 334–337.
8. Radzikowska E, Langfort R and Giedronowicz D. Estrogen and progesterone receptors in non small cell lung cancer patients. *Ann Thorac Cardiovasc Surg* 2002; 8(2), 69–73.
9. Leav I, Lau KM, Adams JY, et al. Comparative studies of the estrogen receptors beta and alpha and the androgen receptor in normal human prostate glands, dysplasia, and in primary and metastatic carcinoma. *American Journal of Pathology* 2001; 159(1), 79–92.
10. Braidman IP, Baris C, Selby PL et al. Preliminary report of impaired oestrogen receptor-alpha expression in bone, but no involvement of androgen receptor, in male idiopathic osteoporosis. *Journal of Pathology* 2000; 192, 90–96.
11. de las Mulas JM, van Niel M, Millan Y et al. Immunohistochemical analysis of estrogen receptors in feline mammary gland benign and malignant lesions: comparison with biochemical assay. *Domest Anim Endocrinol* 2000; 18(1), 111–125.
12. Khan SA, Rogers MA, Khuruna KK et al. Oestrogen receptor expression in normal breast epithelium. *European Journal of Cancer* 2000; 36(Suppl 4), S27-S28.
13. Leake R, Barnes D, Pinder S et al. Immunohistochemical detection of steroid receptors in breast cancer: a working protocol. *Journal of Clinical Pathology* 2000; 53(8), 634–635.
14. Zafrani B, Aubriot MH, Mouret E et al. High sensitivity and specificity of immunohistochemistry for the detection of hormone receptors in breast carcinoma: comparison with biochemical determination in a prospective study of 793 cases. *Histopathology* 2000; 37(6), 536–545.
15. Harvey JM, Clark GM, Osborne CK et al. Estrogen receptor status by immunohistochemistry is superior to the ligand-binding assay for predicting response to adjuvant endocrine therapy in breast cancer. *Journal of Clinical Oncology* 1999; 17(5), 1474–1481.
16. Kaplan, P.A. et. al. 1D5 and 6F11 an immunohistochemical comparison of two monoclonal antibodies for the evaluation of estrogen receptor status in primary breast carcinoma. *Am J Clin Pathol* 2005: 276–280.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Bevtit DJ, Milton ID, Piggot N et al. New monoclonal antibodies to oestrogen and progesterone receptors effective for paraffin section immunohistochemistry. *Journal of Pathology* 1997; 183(2), 228–232.
5. Diaz L, Sahin A and Sneige N. Immunohistochemical detection of estrogen receptor in breast cancer: a laboratory quality improvement study. United States and Canadian Academy of Pathology (Annual Meeting Abstracts March 22–28), 2003; 27A.
6. Dabbs DJ, Landrenau RJ, Liu Y et al. Detection of estrogen receptor by immunohistochemistry in pulmonary adenocarcinoma. *Ann Thorac Surg* 2002; 73(2), 403–405.
7. Khan SA, Yee KA, Kaplan C et al. Estrogen receptor alpha expression in normal human breast epithelium is consistent over time. *International Journal of Cancer* 2002; 102(4), 334–337.
8. Radzikowska E, Langfort R and Giedronowicz D. Estrogen and progesterone receptors in non small cell lung cancer patients. *Ann Thorac Cardiovasc Surg* 2002; 8(2), 69–73.
9. Leav I, Lau KM, Adams JY, et al. Comparative studies of the estrogen receptors beta and alpha and the androgen receptor in normal human prostate glands, dysplasia, and in primary and metastatic carcinoma. *American Journal of Pathology* 2001; 159(1), 79–92.
10. Braidman IP, Baris C, Selby PL et al. Preliminary report of impaired oestrogen receptor-alpha expression in bone, but no involvement of androgen receptor, in male idiopathic osteoporosis. *Journal of Pathology* 2000; 192, 90–96.
11. de las Mulas JM, van Niel M, Millan Y et al. Immunohistochemical analysis of estrogen receptors in feline mammary gland benign and malignant lesions: comparison with biochemical assay. *Domest Anim Endocrinol* 2000; 18(1), 111–125.
12. Khan SA, Rogers MA, Khuruna KK et al. Oestrogen receptor expression in normal breast epithelium. *European Journal of Cancer* 2000; 36(Suppl 4), S27-S28.
13. Leake R, Barnes D, Pinder S et al. Immunohistochemical detection of steroid receptors in breast cancer: a working protocol. *Journal of Clinical Pathology* 2000; 53(8), 634–635.
14. Zafrani B, Aubriot MH, Mouret E et al. High sensitivity and specificity of immunohistochemistry for the detection of hormone receptors in breast carcinoma: comparison with biochemical determination in a prospective study of 793 cases. *Histopathology* 2000; 37(6), 536–545.
15. Harvey JM, Clark GM, Osborne CK et al. Estrogen receptor status by immunohistochemistry is superior to the ligand-binding assay for predicting response to adjuvant endocrine therapy in breast cancer. *Journal of Clinical Oncology* 1999; 17(5), 1474–1481.
16. Kaplan, P.A. et al. 1D5 and 6F11 an immunohistochemical comparison of two monoclonal antibodies for the evaluation of estrogen receptor status in primary breast carcinoma. *Am J Clin Pathol* 2005: 276–280.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Bevitt DJ, Milton ID, Piggot N et al. New monoclonal antibodies to oestrogen and progesterone receptors effective for paraffin section immunohistochemistry. *Journal of Pathology* 1997; 183(2), 228–232.
5. Diaz L, Sahin A and Sneige N. Immunohistochemical detection of estrogen receptor in breast cancer: a laboratory quality improvement study. United States and Canadian Academy of Pathology (Annual Meeting Abstracts March 22–28), 2003; 27A.
6. Dabbs DJ, Landrenau RJ, Liu Y et al. Detection of estrogen receptor by immunohistochemistry in pulmonary adenocarcinoma. *Ann Thorac Surg* 2002; 73(2), 403–405.
7. Khan SA, Yee KA, Kaplan C et al. Estrogen receptor alpha expression in normal human breast epithelium is consistent over time. *International Journal of Cancer* 2002; 102(4), 334–337.
8. Radzikowska E, Langfort R and Giedronowicz D. Estrogen and progesterone receptors in non small cell lung cancer patients. *Ann Thorac Cardiovasc Surg* 2002; 8(2), 69–73.
9. Leav I, Lau KM, Adams JY, et al. Comparative studies of the estrogen receptors beta and alpha and the androgen receptor in normal human prostate glands, dysplasia, and in primary and metastatic carcinoma. *American Journal of Pathology* 2001; 159(1), 79–92.
10. Braidman IP, Baris C, Selby PL et al. Preliminary report of impaired oestrogen receptor-alpha expression in bone, but no involvement of androgen receptor, in male idiopathic osteoporosis. *Journal of Pathology* 2000; 192, 90–96.
11. de las Mulas JM, van Niel M, Millan Y et al. Immunohistochemical analysis of estrogen receptors in feline mammary gland benign and malignant lesions: comparison with biochemical assay. *Domest Anim Endocrinol* 2000; 18(1), 111–125.
12. Khan SA, Rogers MA, Khurana KK et al. Oestrogen receptor expression in normal breast epithelium. *European Journal of Cancer* 2000; 36(Suppl 4), S27-S28.
13. Leake R, Barnes D, Pinder S et al. Immunohistochemical detection of steroid receptors in breast cancer: a working protocol. *Journal of Clinical Pathology* 2000; 53(8), 634–635.
14. Zafrani B, Aubriot MH, Mouret E et al. High sensitivity and specificity of immunohistochemistry for the detection of hormone receptors in breast carcinoma: comparison with biochemical determination in a prospective study of 793 cases. *Histopathology* 2000; 37(6), 536–545.
15. Harvey JM, Clark GM, Osborne CK et al. Estrogen receptor status by immunohistochemistry is superior to the ligand-binding assay for predicting response to adjuvant endocrine therapy in breast cancer. *Journal of Clinical Oncology* 1999; 17(5), 1474–1481.
16. Kaplan, P.A. et. al. 1D5 and 6F11 an immunohistochemical comparison of two monoclonal antibodies for the evaluation of estrogen receptor status in primary breast carcinoma. *Am J Clin Pathol* 2005: 276–280.

Yayın Tarihi

27 Haziran 2016

Leica Biosystems Newcastle Ltd 
Balliol Business Park
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

Leica Biosystems Canada
71 Four Valley Drive
Concord, Ontario L4K 4V8
Canada
☎ +1 800 248 0123

Leica Biosystems Inc
1700 Leider Lane
Buffalo Grove IL 60089
USA
☎ +1 800 248 0123

Leica Biosystems Melbourne
Pty Ltd
495 Blackburn Road
Mt Waverley VIC 3149
Australia
☎ +61 2 8870 3500