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### Catalog No: PA0364

#### Intended Use

*This reagent is for in vitro diagnostic use.*

Thyroid Transcription Factor-1 (SPT24) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human thyroid transcription factor-1 in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

#### Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Thyroid Transcription Factor-1 (SPT24) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human thyroid transcription factor-1 is achieved by first allowing the binding of Thyroid Transcription Factor-1 (SPT24) to the section and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

#### Reagents Provided

Thyroid Transcription Factor-1 (SPT24) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL.

#### Clone

SPT24.

#### Immunogen

Prokaryotic recombinant protein corresponding to a 123 amino acid fragment of the N-terminal region of the TTF-1 molecule.

#### Specificity

Human thyroid transcription factor-1 (TTF-1).

#### Subclass

IgG1, kappa.

#### Total Protein Concentration

Approx 10 mg/mL.

#### Antibody Concentration

Greater than or equal to 2.6 mg/L as determined by ELISA.

#### Dilution and Mixing

Thyroid Transcription Factor-1 (SPT24) primary antibody is optimally diluted for use on the BOND system. Reconstitution, mixing, dilution or titration of this reagent is not required.

#### Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system.

#### Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Thyroid Transcription Factor-1 (SPT24) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

#### Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

## Instructions for Use

Thyroid Transcription Factor-1 (SPT24) primary antibody was developed for use on the automated BOND system in combination with BOND Polymer Refine Detection. The recommended staining protocol for Thyroid Transcription Factor-1 (SPT24) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 1 for 20 minutes.

## Results Expected

### Normal Tissues

Clone SPT24 detected the thyroid transcription factor-1 protein in the nucleus of follicular epithelial cells of the thyroid (13/13), type II pneumocytes and Clara cells of the lung (9/9). Some reactivity was also noted in glial cells in cerebrum (3/3).

Clone SPT24 did not stain TTF-1 in a variety of other normal tissues (n=140) including colon (0/12), brain (0/8), placenta (0/4), kidney (0/6), prostate (0/4), testis (0/4), heart (0/4), spleen (0/4), cervix (0/5), lymph node (0/2), thymus (0/13), tongue (0/1), endometrium (0/2), skeletal muscle (0/4), tonsil (0/4), salivary gland (0/4), myometrium (0/1), umbilical cord (0/1), ureter (0/1), bronchus (0/1), ovary (0/4), fallopian tube (0/1), liver (0/4), breast (0/4), esophagus (0/4), stomach (0/4), spinal cord (0/1), eye (0/1), pancreas (0/4), ileum (0/4), cecum (0/1), rectum (0/1), adrenal (0/4), mesothelium (0/1), parathyroid (0/1), peripheral nerve (0/3), pituitary (0/3), bone marrow (0/3), uterus (0/3) and skin (0/4).

### Abnormal Tissues

Clone SPT24 detected the thyroid transcription factor-1 protein in the nucleus of 26/32 lung adenocarcinoma, 4/7 lung small cell carcinoma, 2/3 lung bronchioalveolar carcinoma, 6/6 thyroid papillary carcinoma, 4/4 thyroid medullary carcinoma, 3/3 thyroid follicular carcinoma, 1/1 thyroid follicular adenoma, 1/1 Hashimoto's thyroiditis, 19/56 thymoma, 1/1 desmoplastic small round cell tumor, 3/3 moderately differentiated colon adenocarcinoma and 1/1 rectal carcinoma.

Clone SPT24 did not stain lung squamous cell carcinoma (0/18), lung large cell carcinoma (0/5), poorly differentiated colon adenocarcinoma (0/4), well differentiated colon adenocarcinoma (0/3), ungraded colon adenocarcinoma (0/3), mesothelioma (0/5), small bowel carcinoid (0/4), thymus atypical carcinoid (1/5), metastatic thymic tumor (0/1), breast tumor (0/6), liver tumor (0/8), kidney renal cell carcinoma (0/5), kidney transitional cell carcinoma (0/1), ovary thecoma (0/1), ovary granulosa cell tumor (0/1), ovary juvenile granulosa (0/1), ovary serous carcinoma (0/3), ovary mucinous carcinoma (0/2), ovary germ cell tumor (0/1), ovary clear cell carcinoma (0/1), stomach adenocarcinoma (0/4), pancreas adenocarcinoma (0/2), pancreas papillary mucinous carcinoma (0/1), pancreas islet cell tumor (0/1), pancreas glucagonoma (0/1), testis seminoma (0/3), testis mixed germ cell tumor (0/1), testis embryonal carcinoma (0/1), brain astrocytoma (0/1), brain choroid plexus papilloma (0/1), melanoma (0/3), skin basal cell carcinoma (0/1), skin squamous cell carcinoma (0/2), penis squamous cell carcinoma (0/2), esophagus squamous cell carcinoma (0/2), larynx squamous cell carcinoma (0/1), tongue squamous cell carcinoma (0/2), cervix squamous cell carcinoma (0/1), small bowel carcinoma (0/1), GIST (0/1), synovial sarcoma (0/1), leiomyosarcoma (0/1), Ewing's sarcoma (0/1), spindle cell rhabdomyosarcoma (0/1), omental fibrous tumor (0/1), bladder transitional cell carcinoma (0/2), bladder small cell carcinoma (0/1), large B-cell lymphoma (0/1), adrenal oncocytoma (0/2), adrenal adenoma (0/1), ganglioneuroma (0/1), prostate adenocarcinoma (0/1), benign prostatic hyperplasia (0/1), endometrial stromal sarcoma (0/1), endometrial adenocarcinoma (0/1), endometrial clear cell carcinoma (0/1), pheochromocytoma (0/1), paraganglioma (0/1) and cervix tumor (0/2).

**Thyroid Transcription Factor-1 (SPT24) is recommended for use in the identification of TTF-1 in normal and neoplastic tissues and for use in diagnosis as part of a panel of antibodies.**

## Product Specific Limitations

Thyroid Transcription Factor-1 (SPT24) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

## Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

## Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

## **Bibliography**

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Berghmans T, Mascaux C, Martin B et al. Prognostic role of thyroid transcription factor-1 in stage III non-small cell lung cancer. Lung Cancer. 2006; 52(2): 219-224.
5. Penman D, Downie I, Roberts F. Positive immunostaining for thyroid transcription factor-1 in primary and metastatic colonic adenocarcinoma: a note of caution. Journal of Clinical Pathology. 2006; 59:663-664.
6. Comperat E, Zhang F, Perrotin C et al. Variable sensitivity and specificity of TTF-1 antibodies in lung metastatic adenocarcinoma of colorectal origin. Modern Pathology. 2005; 18(10):1371-1376.
7. Pan C-C, Chen PC-H, Tsay S-H et al. Cytoplasmic immunoreactivity for thyroid transcription factor-1 in hepatocellular carcinoma: a comparative immunohistochemical analysis of four commercial antibodies using a tissue array technique. American Journal of Clinical Pathology. 2004; 121(3):343-349.

ProClin™ 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Berghmans T, Mascaloux C, Martin B et al. Prognostic role of thyroid transcription factor-1 in stage III non-small cell lung cancer. Lung Cancer. 2006; 52(2): 219-224.
5. Penman D, Downie I, Roberts F. Positive immunostaining for thyroid transcription factor-1 in primary and metastatic colonic adenocarcinoma: a note of caution. Journal of Clinical Pathology. 2006; 59:663-664.
6. Comperat E, Zhang F, Perrotin C et al. Variable sensitivity and specificity of TTF-1 antibodies in lung metastatic adenocarcinoma of colorectal origin. Modern Pathology. 2005; 18(10):1371-1376.
7. Pan C-C, Chen PC-H, Tsay S-H et al. Cyttoplasmic immunoreactivity for thyroid transcription factor-1 in hepatocellular carcinoma: a comparative immunohistochemical analysis of four commercial antibodies using a tissue array technique. American Journal of Clinical Pathology. 2004; 121(3):343-349.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd