BIO SYSTEMS

# Novocastra™ IHC Diluent

## Product Code: RE7133

Leica Biosystems Newcastle Ltd  
Balliol Business Park West  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
☎ +44 191 215 4242

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### Intended Use

*For in vitro diagnostic use.*

Novocastra™ IHC Diluent RE7133 is intended for use as a diluent for Novocastra™ primary antibodies, Novocastra™ Concentrated Biotinylated Secondary Antibody RE7108 and Novocastra™ Concentrated Streptavidin-HRP RE7109 in immunohistochemical (IHC) procedures. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

### Principle of Procedure

The use of an antibody to detect an antigen in tissue cells was first described by Coons AH and Kaplan MH in 1950.<sup>1</sup> Since then, many developments have occurred leading to the immunoperoxidase procedures commonly used today (for IHC staining Principles of Procedure see appropriate detection system Instructions for Use). To obtain optimal staining in IHC procedures primary antibodies require dilution. The binding of antibodies to their epitopes can be readily dissociated by extreme pH and ionic strength. Therefore, an appropriate diluent must be used to obtain optimal staining. This product is used as a diluent for primary antibodies in IHC procedures, which allow the qualitative identification by light microscopy of antigens in sections of formalin-fixed, paraffin-embedded tissue, via sequential steps with interposed washing steps. Novocastra™ IHC Diluent RE7133 helps maintain the specificity and stability of antibodies.

### Reagents Provided

IHC Diluent RE7133 (500ml). Tris-buffered saline, surfactant and protein stabilizer with 0.35% ProClin™ 950.

### Reconstitution, Mixing, Dilution, Titration

IHC Diluent RE7133 is a ready to use reagent. Reconstitution, mixing, dilution, or titration of this reagent is not recommended. Further dilution may result in loss of antigen staining. The user must validate any such change.

### Storage and Stability

Store at 2–8 °C. Do not freeze. Return to 2–8 °C immediately after use. Do not use after expiration date indicated on the product label. Storage conditions other than those specified must be verified by the user. There are no obvious signs to indicate instability of this product, therefore positive and negative controls should be run simultaneously with patient samples. Indications of contamination of IHC Diluent RE7133 are turbidity, odor development or the presence of a precipitate.

### Specimen Preparation

The recommended fixative is 10% neutral-buffered formalin for paraffin-embedded tissue sections.

### Warnings and Precautions

For professional users.

**Do not use IHC Diluent RE7133 for the reconstitution of lyophilized antibodies.**

**Do not use IHC Diluent RE7133 if it becomes turbid, develops odor or contains precipitate.**

The concentration of ProClin™ 950, the preservative used in IHC Diluent RE7133 is 0.35%. It may cause irritation to the skin, eyes, mucus membranes and upper respiratory tract.

A Material Safety Data Sheet (MSDS) is available upon request.

Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions.<sup>2</sup>

Never pipette reagents by mouth and avoid contacting the skin and mucus membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water.

Consult Federal, State or local regulations for disposal of any potentially toxic components.

Minimize microbial contamination of reagents or an increase in non-specific staining may occur. Incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

### Procedure

#### A. Reagents required but not supplied

1. Standard solvents used in immunohistochemistry.
2. Antigen retrieval solution(s) (see **Recommendations on Use** for primary antibody).
3. Enzyme retrieval solution(s) (see **Recommendations on Use** for primary antibody).
4. Primary antibody.
5. Detection system.
6. Mounting medium.

#### B. Equipment required but not supplied

1. Equipment required for antigen retrieval, if recommended for the primary antibody.
2. General immunohistochemistry laboratory equipment.

## C. Methodology

Prior to undertaking this methodology, users must be trained in immunohistochemical techniques.

The combination of the primary antibody, its dilution and optimum conditions for Epitope Retrieval, together with the detection system should be validated by the user on a series of known positive and negative controls.

Mix appropriate volumes of primary antibody, Concentrated Biotinylated Secondary Antibody RE7108 or Concentrated Streptavidin-HRP RE7109 and IHC Diluent RE7133 to achieve the validated dilutions.

Store diluted antibody at 2-8 °C. Return all reagents to 2-8 °C immediately after use.

### Quality Control

Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls in addition to the following procedures. Controls should be fresh autopsy/biopsy/surgical specimens, formalin-fixed, processed and paraffin wax-embedded as soon as possible in the same manner as the patient sample(s).

### Positive Tissue Control

Used to indicate correctly prepared tissues and proper staining techniques. One positive tissue control should be included for each set of test conditions/primary antibody in each staining run. A tissue with weak positive staining is more suitable than a tissue with strong positive staining for optimal quality control and to detect minor levels of reagent degradation.<sup>3</sup> For recommended positive control tissue see primary antibody Instructions for Use. If the positive tissue control fails to demonstrate positive staining, results with the test specimens should be considered invalid.

### Negative Tissue Control

Should be examined after the positive tissue control to verify the specificity of the labeling of the target antigen by the primary antibody. For recommended negative control tissue see primary antibody Instructions for Use. Alternatively, the variety of different cell types present in most tissue sections frequently offers negative control sites, but this should be verified by the user. Non-specific staining, if present, usually has a diffuse appearance. Sporadic staining of connective tissue may also be observed in sections from excessively formalin-fixed tissues. Use intact cells for interpretation of staining results. Necrotic or degenerated cells often stain non-specifically.<sup>4</sup> False-positive results may be seen due to non-immunological binding of proteins or substrate reaction products. They may also be caused by endogenous enzymes such as pseudoperoxidase (erythrocytes), endogenous peroxidase (cytochrome C), or endogenous biotin<sup>5</sup> (eg. liver, breast, brain, kidney). To differentiate endogenous enzyme activity or non-specific binding of enzymes from specific immunoreactivity, additional patient tissues may be stained exclusively with substrate-chromogen, Streptavidin-HRP or labeled polymer and substrate-chromogen, respectively. If specific staining occurs in the negative tissue control, results with the patient specimens should be considered invalid.

### Negative Reagent Control

Use a non-specific negative reagent control in place of the primary antibody with a section of each patient specimen to evaluate non-specific staining and allow better interpretation of specific staining at the antigen site.

### Patient Tissue

Examine stained patient specimens last. Positive staining intensity should be assessed within the context of any non-specific background staining of the negative reagent control. As with any immunohistochemical test, a negative result means that the antigen was not detected, not that the antigen was absent in the cells/tissue assayed. If necessary, use a panel of antibodies to identify false-negative reactions.

### Limitations

Immunohistochemistry is a multistep diagnostic process that consists of specialized training in the selection of the appropriate reagents; tissue selection, fixation, and processing; preparation of the IHC slide; and interpretation of the staining results.

Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning or contamination with other tissues or fluids may produce artifacts, antibody trapping, or false negative results. Inconsistent results may be due to variations in fixation and embedding methods, or to inherent irregularities within the tissue.<sup>7</sup>

Excessive or incomplete counterstaining may compromise proper interpretation of results.

The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Novocastra™ Epitope Retrieval Solutions are for use on paraffin-embedded sections with specific fixation requirements. Unexpected antigen expression may occur, especially in neoplasms. The clinical interpretation of any stained tissue section must include morphological analysis and the evaluation of appropriate controls.

### Performance Characteristics

The performance of Novocastra™ IHC Diluent RE7133 has been validated using a range of Novocastra™ primary antibodies and detection systems.

These products are stable up to the expiry date(s) indicated on the product labels.

### Bibliography

1. Coons AH and Kaplan MH. Localization of antigen in tissue cells II. Improvements in method of detection of antigen by means of fluorescent antibody. *Journal of Experimental Medicine* 1950; 91: 1-13.
2. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
3. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
4. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
5. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.

6. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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## **Bibliographie**

1. Coons AH and Kaplan MH. Localization of antigen in tissue cells II. Improvements in method of detection of antigen by means of fluorescent antibody. *Journal of Experimental Medicine* 1950; 91: 1-13.
2. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
3. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
4. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
5. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
6. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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1. Coons AH and Kaplan MH. Localization of antigen in tissue cells II. Improvements in method of detection of antigen by means of fluorescent antibody. *Journal of Experimental Medicine* 1950; 91: 1-13.
2. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
3. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
4. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
5. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
6. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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#### Bibliografia

1. Coons AH and Kaplan MH. Localization of antigen in tissue cells II. Improvements in method of detection of antigen by means of fluorescent antibody. *Journal of Experimental Medicine* 1950; 91: 1-13.
2. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.

3. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
4. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
5. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
6. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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### Kaynakça

1. Coons AH and Kaplan MH. Localization of antigen in tissue cells II. Fluorescent in method of detection of antigen by means of fluorescent antibody. Journal of Experimental Medicine 1950; 91: 1-13.
2. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
3. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
4. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.

5. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
6. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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### Bibliografie

1. Coons AH and Kaplan MH. Localization of antigen in tissue cells II. Improvements in method of detection of antigen by means of fluorescent antibody. Journal of Experimental Medicine 1950; 91: 1-13.
2. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.

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Leica Biosystems Inc
Pty Ltd