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### Catalog No: PA0051

#### Intended Use

This reagent is for *in vitro* diagnostic use.

Muc-1 Glycoprotein (Ma695) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of Muc-1 Glycoprotein in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

#### Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Muc-1 Glycoprotein (Ma695) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of Muc-1 Glycoprotein is achieved by first allowing the binding of Muc-1 Glycoprotein (Ma695) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

#### Reagents Provided

Muc-1 Glycoprotein (Ma695) diluted in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

#### Clone

Ma695.

#### Immunogen

Human breast cancer cell line ZR-75-1.

#### Specificity

Carbohydrate epitope of the human Muc-1 glycoprotein.

#### Ig Class

IgG1.

#### Total Protein Concentration

Approx 10 mg/mL.

#### Antibody Concentration

Greater than or equal to 0.29 mg/L as determined by ELISA.

#### Dilution and Mixing

Muc-1 Glycoprotein (Ma695) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

#### Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

#### Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Muc-1 Glycoprotein (Ma695) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

#### Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

## Instructions for Use

Muc-1 Glycoprotein (Ma695) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection. The recommended staining protocol for Muc-1 Glycoprotein (Ma695) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 2 for 20 minutes.

## Results Expected

### Normal Tissues

Clone Ma695 detected Muc-1 glycoprotein on the luminal aspects of ducts in breast and thyroid tissue, pneumocytes, alveolar walls and macrophages in lung, glandular cells in endometrium, gastric glands in stomach, collecting ducts and convoluted tubules in kidney, umbrella cells and superficial keratinocytes in bladder, follicular dendritic cells in appendix and rectum, a subset of cells in parathyroid, a proportion of keratinocytes in esophagus, brush borders of acinar cells in salivary gland, luminal aspect of acinar cells in pancreas, sebaceous glands in skin, follicular dendritic cells and squamous mucosa in tonsil and fibroblasts in colon. (Total number of normal cases evaluated = 123).

### Tumor Tissues

Clone Ma695 stained 5/5 breast tumors (including 3/3 invasive ductal carcinomas and 2/2 fibroadenomas and), 4/5 metastatic tumors, 4/4 lung tumors (including 2/2 squamous cell carcinoma, 1/1 adenocarcinomas and 1/1 small cell carcinomas), 2/8 bowel tumors (including 2/6 adenocarcinomas and 0/2 adenomas), 2/5 thyroid tumors (including 1/1 follicular carcinoma, 1/1 follicular variant of papillary carcinoma and 0/3 adenomas), 2/3 adenocarcinomas of the stomach, 2/3 squamous cell carcinomas of the esophagus, 2/2 transitional cell carcinomas of the bladder, 2/2 clear cell carcinomas of the kidney, 2/2 adenocarcinomas of the prostate, 2/2 tumors of the salivary gland (including 1/1 pleomorphic adenoma and 1/1 adenoid cystic carcinoma), 2/2 squamous cell carcinomas of the cervix, 2/2 adenocarcinomas of the endometrium, 1/3 ovarian tumors (including 1/2 endometrioid carcinoma and 0/1 granulosa cell tumor), 1/1 squamous cell carcinoma of the tongue and 1/1 adenocarcinoma of the pancreas. No staining was observed in a variety of additional abnormal tissues evaluated, including brain tumors (0/4), liver tumors (0/4), lymphomas (0/3), tumors of the adrenal gland (0/2), bone tumors (0/2), tumors of the head and neck (0/2), seminomas (0/2), a prostatic hyperplasia (0/1) and a melanoma (0/1). (Total number of abnormal tissues evaluated = 71).

**Muc-1 glycoprotein (Ma695) is recommended for the detection of Muc-1 glycoprotein in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.**

## Product Specific Limitations

Muc-1 Glycoprotein (Ma695) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

## Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

## Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

## Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Khoury T, Tan D, Wang J, et al. Inclusion of MUC1 (Ma695) in a panel of immunohistochemical markers is useful for distinguishing between endocervical and endometrial mucinous adenocarcinoma. BMC Clinical Pathology. 2006; 12:6:1.
5. Lau SK, Weiss LM and Chu PG. Differential expression of MUC1, MUC2 and MUC5AC in carcinomas of various sites. Anatomic Pathology. 2004; 122:61-69.
6. Rakha EA, Boyce RWG, Abd El-Rehim D et al. Expression of mucins (MUC1, MUC2, MUC3, MUC4, MUC5AC and MUC6) and their prognostic significance in human breast cancer. Modern Pathology. 2005; 18:1295-1304.
7. Yuan S, Li K, Wang L et al. Expression of MUC1 and its significance in hepatocellular and cholangiocarcinoma tissue. World Journal of Gastroenterology. 2005; 11(30):4661-4666.
8. Langner C, Ratschek M, Rehak P, et al. Expression of MUC1 (EMA) and E-cadherin in renal cell carcinoma: a systematic immunohistochemical analysis of 188 cases. Modern Pathology. 2004; 17(2):180-188.

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## Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Khoury T, Tan D, Wang J, et al. Inclusion of MUC1 (Ma695) in a panel of immunohistochemical markers is useful for distinguishing between endocervical and endometrial mucinous adenocarcinoma. BMC Clinical Pathology. 2006; 12:6:1.
5. Lau SK, Weiss LM and Chu PG. Differential expression of MUC1, MUC2 and MUC5AC in carcinomas of various sites. Anatomic Pathology. 2004; 12:2:61-69.
6. Rakha EA, Boyce RWG, Abd El-Rehim D et al. Expression of mucins (MUC1, MUC2, MUC3, MUC4, MUC5AC and MUC6) and their prognostic significance in human breast cancer. Modern Pathology. 2005; 18:1295-1304.
7. Yuan S, Li K, Wang L et al. Expression of MUC1 and its significance in hepatocellular and cholangiocarcinoma tissue. World Journal of Gastroenterology. 2005; 11(30):4661-4666.
8. Langner C, Ratschek M, Rehak P, et al. Expression of MUC1 (EMA) and E-cadherin in renal cell carcinoma: a systematic immunohistochemical analysis of 188 cases. Modern Pathology. 2004; 17(2):180-188.

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## Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 1163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Khoury T, Tan D, Wang J, et al. Inclusion of MUC1 (Ma695) in a panel of immunohistochemical markers is useful for distinguishing between endocervical and endometrial mucinous adenocarcinoma. BMC Clinical Pathology. 2006; 12:6:1.
5. Lau SK, Weiss LM and Chu PG. Differential expression of MUC1, MUC2 and MUC5AC in carcinomas of various sites. Anatomic Pathology. 2004; 122:61-69.
6. Rakha EA, Boyce RWG, Abd El-Rehim D et al. Expression of mucins (MUC1, MUC2, MUC3, MUC4, MUC5AC and MUC6) and their prognostic significance in human breast cancer. Modern Pathology. 2005; 18:1295-1304.
7. Yuan S, Li K, Wang L et al. Expression of MUC1 and its significance in hepatocellular and cholangiocarcinoma tissue. World Journal of Gastroenterology. 2005; 11(30):4661-4666.
8. Langner C, Ratschek M, Rehak P et al. Expression of MUC1 (EMA) and E-cadherin in renal cell carcinoma: a systematic immunohistochemical analysis of 188 cases. Modern Pathology. 2004; 17(2):180-188.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd