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### Catalog No: PA0575

#### Intended Use

This reagent is for *in vitro* diagnostic use.

Enhancer of Zeste Homolog 2 (Drosophila) (EZHZ) (6A10) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human EZHZ antigen in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

#### Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Enhancer of Zeste Homolog 2 (Drosophila) (EZHZ) (6A10) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of EZHZ antigen is achieved by first allowing the binding of Enhancer of Zeste Homolog 2 (Drosophila) (EZHZ) (6A10) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

#### Reagents Provided

Enhancer of Zeste Homolog 2 (Drosophila) (EZHZ) (6A10) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

#### Clone

6A10

#### Immunogen

Prokaryotic recombinant protein corresponding to a portion of the N-terminus domain of the human EZHZ molecule.

#### Specificity

Human EZHZ (Enhancer of Zeste Homolog 2 (Drosophila)), (85 kD). Region common to isoforms A and B of the protein.

#### Ig Class

IgG1

#### Total Protein Concentration

Approx 10 mg/mL.

#### Antibody Concentration

Greater than or equal to 0.4 mg/L as determined by ELISA.

#### Dilution and Mixing

Enhancer of Zeste Homolog 2 (Drosophila) (EZHZ) (6A10) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

#### Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

#### Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Enhancer of Zeste Homolog 2 (Drosophila) (EZHZ) (6A10) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

#### Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

## Instructions for Use

Enhancer of Zeste Homolog 2 (*Drosophila*) (EZHZ) (6A10) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection. The recommended staining protocol for Enhancer of Zeste Homolog 2 (*Drosophila*) (EZHZ) (6A10) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 2 for 40 minutes.

## Results Expected

### Normal Tissues

Clone 6A10 detected the EZHZ antigen in the nucleus of a variety of cell types in normal tissues, including lymphocytes in tonsil, lymph node, spleen and thymus and epithelium in tonsil, skin, esophagus, ovary, small bowel, large bowel, breast, cervix and larynx. Staining was also noted in seminiferous tubules of testis, and hepatocytes. (Total number of normal cases = 140).

### Tumor Tissues

Clone 6A10 stained 97/130 lymphomas (including 46/56 diffuse B-cell lymphomas, 19/30 T-cell lymphomas, 19/20 diffuse large B-cell lymphomas, 4/8 follicular lymphomas, 2/4 lymphocytic predominance Hodgkin's lymphomas, 2/2 diffuse small cell lymphomas, 2/2 T-cell rich B-cell lymphomas, 1/5 mixed cellularity Hodgkin's lymphomas, 1/1 non-Hodgkin's B-cell lymphoma, 1/1 anaplastic large cell lymphoma and 0/1 Hodgkin's lymphoma), 14/42 prostate adenocarcinomas, 9/18 breast tumors (including 4/8 invasive ductal carcinomas, 2/2 medullary carcinomas, 1/3 fibroadenomas, 1/2 invasive lobular carcinomas, 1/1 Paget's disease and 0/2 adenomas), 7/9 bowel tumors (including 5/7 adenocarcinomas and 2/2 adenomas), 3/4 lung tumors (including 1/2 squamous cell carcinomas, 1/1 adenocarcinoma and 1/1 small cell carcinoma), 2/5 metastatic tumors (including 1/1 metastatic breast invasive ductal carcinoma, 1/1 metastatic esophageal squamous cell carcinoma, 0/1 metastatic colon adenocarcinoma, 0/1 metastatic tumor from gastrointestinal site and 0/1 metastatic colon signet ring carcinoma), 2/5 hepatocellular carcinomas, 2/3 ovarian tumors (including 1/1 adenocarcinoma, 1/1 endometrioid adenocarcinoma and 0/1 granulosa cell tumor), 2/3 squamous cell carcinomas of the esophagus, 2/3 stomach adenocarcinomas, 2/2 squamous cell carcinomas of the cervix, 2/2 endometrial adenocarcinomas, 1/4 brain tumors (including 1/3 meningiomas and 0/1 astrocytoma), 1/2 transitional cell carcinoma of the bladder, 1/2 seminomas, 1/2 bone tumors (including 1/1 chondrosarcoma and 0/1 osteosarcoma), 1/2 tumors of the salivary gland (including 1/1 adenoid cystic carcinoma and 0/1 pleomorphic adenoma), 1/1 melanoma, 1/1 hyalinosis of breast fibrous tissue and 1/1 pancreatic adenocarcinoma. No staining was detected in thyroid tumors (0/5), breast adenosis (0/3), tumors of the adrenal gland (0/2), renal clear cell carcinomas (0/2), a nasopharyngeal carcinoma (0/1), a squamous cell carcinoma of the tongue (0/1), a prostatic hyperplasia (0/1), a squamous cell carcinoma of the skin (0/1) and a pheochromocytoma (0/1). (Total number of abnormal cases = 258).

**Enhancer of Zeste Homolog 2 (*Drosophila*) (EZHZ) (6A10) is recommended for the detection of human EZHZ antigen in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.**

## Product Specific Limitations

Enhancer of Zeste Homolog 2 (*Drosophila*) (EZHZ) (6A10) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

## Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

## Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

## Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Leinonen K, Saramäki O, Furusato B, et al. Loss of PTEN is associated with aggressive behavior in ERG-positive prostate cancer. Cancer Epidemiology, Biomarkers and Prevention. 2013; 22 (12):2333-2344.
5. Gong Y, Huo L, Liu P, et al. Polycomb group protein EZHZ2 is frequently expressed in inflammatory breast cancer and is predictive of worse clinical outcome. Cancer. 2011; 117:5476-5484.
6. Hebbard L, Cecena G, Golas J, et al. Control of mammary tumor differentiation by SKI-606 (bosutinib). Oncogene. 2011; 30 (3):301-312.

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## Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Leinonen K, Saramäki O, Furusato B, et al. Loss of PTEN is associated with aggressive behavior in ERG-positive prostate cancer. Cancer Epidemiology, Biomarkers and Prevention. 2013; 22 (12):2333-2344.
5. Gong Y, Huo L, Liu P, et al. Polycomb group protein EZH2 is frequently expressed in inflammatory breast cancer and is predictive of worse clinical outcome. Cancer. 2011; 117:5476-5484.
6. Hebbard L, Cecena G, Golas J, et al. Control of mammary tumor differentiation by SKI-606 (bosutinib). Oncogene. 2011; 30 (3):301-312.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Leinonen K, Saramäki O, Furusato B, et al. Loss of PTEN is associated with aggressive behavior in ERG-positive prostate cancer. *Cancer Epidemiology, Biomarkers and Prevention*. 2013; 22 (12):2333-2344.
5. Gong Y, Huo L, Liu P, et al. Polycomb group protein EZH2 is frequently expressed in inflammatory breast cancer and is predictive of worse clinical outcome. *Cancer*. 2011; 117:5476-5484.
6. Hebbard L, Cecena G, Golas J, et al. Control of mammary tumor differentiation by SKI-606 (bosutinib). *Oncogene*. 2011; 30 (3):301-312.

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## Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

4. Leinonen K, Saramäki O, Furusato B, et al. Loss of PTEN is associated with aggressive behavior in ERG-positive prostate cancer. *Cancer Epidemiology, Biomarkers and Prevention*. 2013; 22 (12):2333-2344.
5. Gong Y, Huo L, Liu P, et al. Polycomb group protein EZH2 is frequently expressed in inflammatory breast cancer and is predictive of worse clinical outcome. *Cancer*. 2011; 117:5476-5484.
6. Hebbard L, Cecena G, Golas J, et al. Control of mammary tumor differentiation by SKI-606 (bosutinib). *Oncogene*. 2011; 30 (3):301-312.

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## Literatura

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Leinonen K, Saramäki O, Furusato B, et al. Loss of PTEN is associated with aggressive behavior in ERG-positive prostate cancer. Cancer Epidemiology, Biomarkers and Prevention. 2013; 22 (12):2333-2344.

5. Gong Y, Huo L, Liu P, et al. Polycomb group protein EZH2 is frequently expressed in inflammatory breast cancer and is predictive of worse clinical outcome. *Cancer*. 2011; 117:5476-5484.
6. Hebbard L, Cecena G, Golas J, et al. Control of mammary tumor differentiation by SKI-606 (bosutinib). *Oncogene*. 2011; 30 (3):301-312.

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## Literatúra

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Leinonen K, Saramäki O, Furusato B, et al. Loss of PTEN is associated with aggressive behavior in ERG-positive prostate cancer. Cancer Epidemiology, Biomarkers and Prevention. 2013; 22 (12):2333-2344.
5. Gong Y, Huo L, Liu P, et al. Polycarbonyl group protein EZHZ2 is frequently expressed in inflammatory breast cancer and is predictive of worse clinical outcome. Cancer. 2011; 117:5476-5484.
6. Hebbard L, Cecena G, Golas J, et al. Control of mammary tumor differentiation by SKI-606 (bosutinib). Oncogene. 2011; 30 (3):301-312.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd