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BOND™ Ready-To-Use Primary Antibody CD57 (NK-1)

Catalog No: PA0443

Intended Use

This reagent is for in vitro diagnostic use.

CD57 (NK-1) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human CD57 molecule in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). CD57 (NK-1) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human CD57 molecule is achieved by first, allowing the binding of CD57 (NK-1) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

CD57 (NK-1) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin® 950 as a preservative.

Total volume = 7 mL.

Clone

NK-1.

Immunogen

Human peripheral blood mononuclear cells.

Specificity

Human CD57 molecule.

Subclass

IgM.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 0.18 mg/L as determined by ELISA.

Dilution and Mixing

CD57 (NK-1) primary antibody is optimally diluted for use on the BOND system. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system.

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of CD57 (NK-1) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin® 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

CD57 (NK-1) primary antibody was developed for use on the automated BOND system in combination with BOND Polymer Refine Detection. The recommended staining protocol for CD57 (NK-1) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 1 for 20 minutes.

Results Expected

Normal Tissues

Clone NK-1 detected the CD57 antigen on the membrane of large granular lymphocytes in tonsil, spleen and in a variety of non-lymphoid tissues (n=90). In the tissue tested, staining was also observed in glandular epithelium of prostate and endometrium, epithelial cells of thymus, nerve elements of the peripheral and central nervous system, neuroendocrine cells in gastrointestinal tract, gastric glands and serous cells of the submandibular gland. Variable staining was seen in the colloid of thyroid, kidney tubules, adrenal medulla, hepatocytes and breast ducts. Additionally, some staining was also observed in, capillaries of parathyroid, gastric glands, a percentage of pancreatic islet cells and spermatogonia in testis.

Tumor Tissues

Clone NK-1 stained 5/12 schwannomas, 2/4 neurofibromas, 1/1 neuroendocrine tumor, 5/16 liposarcomas, 1/11 malignant fibrous histiocytomas, 2/8 leiomyosarcomas, 2/4 rhabdomyosarcomas, 1/4 synovial sarcomas, 1/1 glomus tumor, 0/2 alveolar soft part sarcomas, 0/1 osteosarcoma, 0/1 giant cell tumor, 1/1 hemangiopericytoma, 0/2 hemangioendothelial sarcomas, 0/1 breast carcinoma, 0/1 papillary carcinoma of thyroid, 0/1 leiomyoma and 0/1 adenocarcinoma.

CD57 (NK-1) is recommended for the identification of cells expressing CD57 in normal and neoplastic tissues.

Product Specific Limitations

CD57 (NK-1) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Villadsen LS, Skov L, Dam TN, et al. In situ depletion of CD4+ T cells in human skin by Zanolimumab. Archives of Dermatological Research. 2007; 298:449-455.
5. Villadsen LS, Schuurman J, Beurskens F, et al. Resolution of psoriasis upon blockade of IL-15 biological activity in a xenograft mouse model. Journal of Clinical Investigation. 2003; 112:1571-1580.
6. Moreira PR, Santos DFM, Martins RD, et al. CD57+ cells in radicular cyst. International Endodontic Journal. 2000; 33:99-102.
7. Wongweragiat S, Searle RF and Bulmer JN. Decidual T lymphocyte activation in hydatidiform mole. Journal of Clinical Pathology. 1999; 52:888-894.

ProClin® 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Villadsen LS, Skov L, Dam TN, et al. In situ depletion of CD4+ T cells in human skin by Zanolimumab. Archives of Dermatological Research. 2007; 298:449–455.
5. Villadsen LS, Schuurman J, Beurskens F, et al. Resolution of psoriasis upon blockade of IL-15 biological activity in a xenograft mouse model. Journal of Clinical Investigation. 2003; 112:1571–1580.
6. Moreira PR, Santos DFM, Martins RD, et al. CD57+ cells in radicular cyst. International Endodontic Journal. 2000; 33:99–102.
7. Wongweragiat S, Searle RF and Bulmer JN. Decidual T lymphocyte activation in hydatidiform mole. Journal of Clinical Pathology. 1999; 52:888–894.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Villadsen LS, Skov L, Dam TN, et al. In situ depletion of CD4+ T cells in human skin by Zanolimumab. Archives of Dermatological Research. 2007; 298:449–455.
5. Villadsen LS, Schuurman J, Beurskens F, et al. Resolution of psoriasis upon blockade of IL-15 biological activity in a xenograft mouse model. Journal of Clinical Investigation. 2003; 112:1571–1580.
6. Moreira PR, Santos DFM, Martins RD, et al. CD57+ cells in radicular cyst. International Endodontic Journal. 2000; 33:99–102.
7. Wongweraijat S, Searle RF and Bulmer JN. Decidual T lymphocyte activation in hydatidiform mole. Journal of Clinical Pathology. 1999; 52:888–894.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Villadsen LS, Skov L, Dam TN, et al. In situ depletion of CD4+ T cells in human skin by Zanolimumab. Archives of Dermatological Research. 2007; 298:449–455.
5. Villadsen LS, Schuurman J, Beurskens F, et al. Resolution of psoriasis upon blockade of IL-15 biological activity in a xenograft mouse model. Journal of Clinical Investigation. 2003; 112:1571–1580.
6. Moreira PR, Santos DFM, Martins RD, et al. CD57+ cells in radicular cyst. International Endodontic Journal. 2000; 33:99–102.
7. Wongweraqiatt S, Searle RF and Bulmer JN. Decidual T lymphocyte activation in hydatidiform mole. Journal of Clinical Pathology. 1999; 52:888–894.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Villadsen LS, Skov L, Dam TN, et al. In situ depletion of CD4+ T cells in human skin by Zanolimumab. Archives of Dermatological Research. 2007; 298:449–455.
5. Villadsen LS, Schuurman J, Beurskens F, et al. Resolution of psoriasis upon blockade of IL-15 biological activity in a xenograft mouse model. Journal of Clinical Investigation. 2003; 112:1571–1580.
6. Moreira PR, Santos DFM, Martins RD, et al. CD57+ cells in radicular cyst. International Endodontic Journal. 2000; 33:99–102.
7. Wongweraqiāt S, Searle RF and Bulmer JN. Decidual T lymphocyte activation in hydatidiform mole. Journal of Clinical Pathology. 1999; 52:888–894.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Villadsen LS, Skov L, Dam TN, et al. In situ depletion of CD4+ T cells in human skin by Zanilimumab. Archives of Dermatological Research. 2007; 298:449–455.
5. Villadsen LS, Schuurman J, Beurskens F, et al. Resolution of psoriasis upon blockade of IL-15 biological activity in a xenograft mouse model. Journal of Clinical Investigation. 2003; 112:1571–1580.
6. Moreira PR, Santos DFM, Martins RD, et al. CD57+ cells in radicular cyst. International Endodontic Journal. 2000; 33:99–102.
7. Wongweragiat S, Searle RF and Bulmer JN. Decidual T lymphocyte activation in hydatidiform mole. Journal of Clinical Pathology. 1999; 52:888–894.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Villadsen LS, Skov L, Dam TN, et al. In situ depletion of CD4+ T cells in human skin by Zanolimub. Archives of Dermatological Research. 2007; 298:449–455.
5. Villadsen LS, Schuurman J, Beurskens F, et al. Resolution of psoriasis upon blockade of IL-15 biological activity in a xenograft mouse model. Journal of Clinical Investigation. 2003; 112:1571–1580.
6. Moreira PR, Santos DFM, Martins RD, et al. CD57+ cells in radicular cyst. International Endodontic Journal. 2000; 33:99–102.
7. Wongweragiat S, Searle RF and Bulmer JN. Decidual T lymphocyte activation in hydatidiform mole. Journal of Clinical Pathology. 1999; 52:888–894.

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Leica Biosystems Canada
71 Four Valley Drive
Concord, Ontario L4K 4V8
Canada
J +1 800 248 0123

Leica Biosystems Inc
1700 Leider Lane
Buffalo Grove IL 60089
USA
J +1 800 248 0123

Leica Biosystems Melbourne
Pty Ltd
495 Blackburn Road
Mt Waverley VIC 3149
Australia
J +61 2 8870 3500

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Leica Biosystems Newcastle Ltd