Bond™ Ready-To-Use Primary Antibody Vimentin (SRL33)

Catalog No: PA0033

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Instructions for Use

Please read before using this product.

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Check the integrity of the packaging before use.

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Catalog No: PA0033

Intended Use

This reagent is for in vitro diagnostic use.

Vimentin (SRL33) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human vimentin intermediate filament in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated Bond™ system. The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using Bond Reagents" in your Bond user documentation). Vimentin (SRL33) primary antibody is a ready to use product that has been specifically optimized for use with Bond Polymer Refine Detection. The demonstration of human vimentin intermediate filament is achieved by first, allowing the binding of Vimentin (SRL33) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated Bond system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

Vimentin (SRL33) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

SRL33.

Immunogen

Prokaryotic recombinant protein corresponding to a portion of the C-terminus of the human vimentin molecule.

Specificity

Human vimentin.

Subclass

IgG2a.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 0.58 mg/L as determined by ELISA.

Dilution and Mixing

Vimentin (SRL33) primary antibody is optimally diluted for use on the Bond system. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using Bond Reagents" in your Bond user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the Bond system.

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Vimentin (SRL33) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Vimentin (SRL33) primary antibody was developed for use on the automated Bond system in combination with Bond Polymer Refine Detection. The recommended staining protocol for Vimentin (SRL33) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using Bond Epitope Retrieval Solution 1 for 20 minutes.

Results Expected

Normal Tissues

Clone SRL33 detected the intermediate filament protein, vimentin, in the cytoplasm of cells of mesenchymal origin. Staining was observed in a variety of cell types including vascular endothelial cells, adipose cells, smooth muscle cells, pancreatic acinar cells, peripheral nerve cells, macrophages (including Kupffer), Langerhans cells and lymphoid cells (n=76).

Tumor Tissues

Clone SRL33 stained 46/135 tumors evaluated including lymphomas, sarcomas, renal cell carcinomas, melanomas, rhabdomyosarcomas, soft tissue tumors some adenocarcinomas and squamous cell carcinomas, and some thyroid, neural and ovarian tumors.

Vimentin (SRL33) is recommended for use as part of an antibody panel for the classification of tumors of mesenchymal origin.

Product Specific Limitations

Vimentin (SRL33) has been optimized at Leica Biosystems for use with Bond Polymer Refine Detection and Bond ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with Bond reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using Bond Reagents" in your Bond user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Guorong L, Cuilleron M, Gentil-Perret A, et al. Rapid and sensitive detection of messenger RNA expression for molecular differential diagnosis of renal cell carcinoma. Clinical Cancer Research 2003; 9:6441–6446.
5. Langa F, Kress C, Colucci-Guyon E, et al. Teratocarcinomas induced by embryonic stem (ES) cells lacking vimentin: an approach to study the role of vimentin in tumorigenesis. Journal of Cell Science 2000; 113:3463–3472.
6. Paramio J, Jorcano J. Beyond structure: do intermediate filaments modulate cell signalling? BioEssays 2000; 24:836–844.
7. Evans R. Vimentin: the conundrum of the intermediate filament gene family. BioEssays 1998; 20:79–86.

ProClin 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Guorong L, Cuilleron M, Gentil-Perret A, et al. Rapid and sensitive detection of messenger RNA expression for molecular differential diagnosis of renal cell carcinoma. *Clinical Cancer Research* 2003; 9:6441–6446.
5. Langa F, Kress C, Colucci-Guyon E, et al. Teratocarcinomas induced by embryonic stem (ES) cells lacking vimentin: an approach to study the role of vimentin in tumorigenesis. *Journal of Cell Science* 2000; 113:3463–3472.
6. Paramio J, Jorcano J. Beyond structure: do intermediate filaments modulate cell signalling? *BioEssays* 2000; 24:836–844.
7. Evans R. Vimentin: the conundrum of the intermediate filament gene family. *BioEssays* 1998; 20:79–86.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Guorong L, Cuilleron M, Gentil-Perret A, et al. Rapid and sensitive detection of messenger RNA expression for molecular differential diagnosis of renal cell carcinoma. Clinical Cancer Research 2003; 9:6441-6446.
 5. Langa F, Kress C, Colucci-Guyon E, et al. Teratocarcinomas induced by embryonic stem (ES) cells lacking vimentin: an approach to study the role of vimentin in tumorigenesis. Journal of Cell Science 2000; 113:3463-3472.
 6. Paramio J, Jorcano J. Beyond structure: do intermediate filaments modulate cell signalling? BioEssays 2000; 24:836-844.
 7. Evans R. Vimentin: the conundrum of the intermediate filament gene family. BioEssays 1998; 20:79-86.
- ProClin® 950 è un marchio di fabbrica di Supelco, società del gruppo Sigma-Aldrich.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Guorong L, Cuilleron M, Gentil-Perret A, et al. Rapid and sensitive detection of messenger RNA expression for molecular differential diagnosis of renal cell carcinoma. Clinical Cancer Research 2003; 9:6441–6446.
 5. Langa F, Kress C, Colucci-Guyon E, et al. Teratocarcinomas induced by embryonic stem (ES) cells lacking vimentin: an approach to study the role of vimentin in tumorigenesis. Journal of Cell Science 2000; 113:3463–3472.
 6. Paramio J, Jorcano J. Beyond structure: do intermediate filaments modulate cell signalling? BioEssays 2000; 24:836–844.
 7. Evans R. Vimentin: the conundrum of the intermediate filament gene family. BioEssays 1998; 20:79–86.
- ProClin: 950 ist eine Marke von Supelco, einem Unternehmen der Sigma-Aldrich Corporation.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Guorong L, Cuilleron M, Gentil-Perret A, et al. Rapid and sensitive detection of messenger RNA expression for molecular differential diagnosis of renal cell carcinoma. Clinical Cancer Research 2003; 9:6441–6446.
 5. Langa F, Kress C, Colucci-Guyon E, et al. Teratocarcinomas induced by embryonic stem (ES) cells lacking vimentin: an approach to study the role of vimentin in tumorigenesis. Journal of Cell Science 2000; 113:3463–3472.
 6. Paramio J, Jorcano J. Beyond structure: do intermediate filaments modulate cell signalling? BioEssays 2000; 24:836–844.
 7. Evans R. Vimentin: the conundrum of the intermediate filament gene family. BioEssays 1998; 20:79–86.
- ProClin® 950 es una marca registrada de Supelco, parte de Sigma-Aldrich Corporation.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Guorong L, Cuilleron M, Gentil-Perret A, et al. Rapid and sensitive detection of messenger RNA expression for molecular differential diagnosis of renal cell carcinoma. Clinical Cancer Research 2003; 9:6441–6446.
 5. Langa F, Kress C, Colucci-Guyon E, et al. Teratocarcinomas induced by embryonic stem (ES) cells lacking vimentin: an approach to study the role of vimentin in tumorigenesis. Journal of Cell Science 2000; 113:3463–3472.
 6. Paramio J, Jorcano J. Beyond structure: do intermediate filaments modulate cell signalling? BioEssays 2000; 24:836–844.
 7. Evans R. Vimentin: the conundrum of the intermediate filament gene family. BioEssays 1998; 20:79–86.
- ProClin® 950 é uma marca registada de Supelco, parte da Sigma-Aldrich Corporation.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD och Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Guorong L, Cuilleron M, Gentil-Perret A, et al. Rapid and sensitive detection of messenger RNA expression for molecular differential diagnosis of renal cell carcinoma. Clinical Cancer Research 2003; 9:6441–6446.
5. Langa F, Kress C, Colucci-Guyon E, et al. Teratocarcinomas induced by embryonic stem (ES) cells lacking vimentin: an approach to study the role of vimentin in tumorigenesis. Journal of Cell Science 2000; 113:3463–3472.
6. Paramio J, Jorcano J. Beyond structure: do intermediate filaments modulate cell signalling? BioEssays 2000; 24:836–844.
7. Evans R. Vimentin: the conundrum of the intermediate filament gene family. BioEssays 1998; 20:79–86.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Guorong L, Cuilleron M, Gentil-Perret A, et al. Rapid and sensitive detection of messenger RNA expression for molecular differential diagnosis of renal cell carcinoma. Clinical Cancer Research 2003; 9:6441–6446.
 5. Langa F, Kress C, Colucci-Guyon E, et al. Teratocarcinomas induced by embryonic stem (ES) cells lacking vimentin: an approach to study the role of vimentin in tumorigenesis. Journal of Cell Science 2000; 113:3463–3472.
 6. Paramio J, Jorcano J. Beyond structure: do intermediate filaments modulate cell signalling? BioEssays 2000; 24:836–844.
 7. Evans R. Vimentin: the conundrum of the intermediate filament gene family. BioEssays 1998; 20:79–86.
- ProClin[®] 950 er et varemærke tilhørende Supelco, en del af Sigma-Aldrich Corporation.

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20 januar 2014

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
J +44 191 215 4242