


## Kreatech<sup>™</sup> FISH probes Product Information Sheet

KBI-40008 RCAN1 (21q22), SE X, SE Y






Kreatech Biotechnology B.V. Vlierweg 20 1032 LG Amsterdam The Netherlands www.LeicaBiosystems.com

Published December 2021

PI-KBI-40008\_D3.0







## Kreatech™ RCAN1 (21q22), SE X, SE Y FISH probe

| Intended purpose:   | using uncultured amniotic cells<br>assay to determine the following<br>The <b>RCAN1 (21q22)</b> specific FI<br>The <b>SE X (DXZ1)</b> (Satellite Enur<br>cells.   | The SE Y (DYZ3) (Satellite Enumeration) FISH probe) is optimized to detect copy numbers of chromosome Y at Yp11-Yq11 on uncultured amniotic   |  |              |  |
|---|---|---|--|--------------|--|
| Warnings and Limitation   | standalone test, but always in co   | This product is not intended for use on any type of cells/tissue other than uncultured amniotic cells. This test should never be used as a<br>standalone test, but always in conjunction with other results/follow-up testing.<br>This FISH assay will not detect the presence of structural chromosome abnormalities that can also result in birth defects.  |  |              |  |
|   | We do not recommend this proc   | duct for the detection of the iAMP21.   |  |              |  |
| Introduction:   | of phenotypic features that inclu-<br>appears to contain the gene(s<br>chromosomes) are slightly less<br>people with sex chromosome a<br>inherit only one X chromosome<br>more rarely XXXX or XXXXX.                            | Trisomy 21 is one of the most common chromosomal abnormalities in live born children and causes Down syndrome, a particular combination<br>of phenotypic features that includes mental retardation and characteristic facies. Molecular analysis has revealed that the 21q22.1-q22.3 region<br>appears to contain the gene(s) responsible for the congenital heart. Chromosomal abnormalities involving the X and Y chromosome (sex<br>chromosomes) are slightly less common than autosomal abnormalities and are usually much less severe in their effects. The high frequency of<br>people with sex chromosome aberrations is partly due to the fact that they are rarely lethal conditions. Turner syndrome occurs when females<br>inherit only one X chromosome; their genotype is X0. Metafemales or triple-X females, inherit three X chromosomes; their genotype is XXX or<br>more rarely XXXX or XXXXX. Kineflater syndrome males inherit one or more extra X chromosome; their genotype is XXY or more rarely<br>XXXY, XXXXY, or XY/XXY mosaic. XYY syndrome males inherit an extra Y chromosome; their genotype is XYY. |  |              |  |
| Critical region 1 (red):<br>Critical region 2 (green):<br>Critical region 3 (blue): | The <b>RCAN1 (21q22)</b> FISH probe is direct-labeled with PlatinumBright™550.<br>The <b>SE X</b> FISH probe is direct-labeled with PlatinumBright™495<br>The <b>SE Y</b> FISH probe is direct-labeled with PlatinumBright™415. |   |  |              |  |
| Reagent:  | Kreatech probes are direct-labeled DNA probes provided in a ready-to-use format.<br>Apply 10 µl of probe to a sample area of approximately 22 x 22 mm.  |   |  |              |  |
|   | Please refer to the Instruction<br>Kreatech FISH probes are REF<br>background, due to unspecifio  |   | rotocol.<br>ain Cot-1 DNA. Hybridization efficiency is i | ncreased and |  |
| Interpretation:   | hybridization assay. In females<br>two single colors red (R) and or   | The RCAN1 (21q22), SE X, SE Y FISH probe is designed as a triple-color assay to detect aneuploidies of chromosome 21, X, and Y in a single<br>hybridization assay. In females two single color red (R) and green (G) signals will identify the normal chromosomes 21 and X (2R2G). In males<br>two single colors red (R) and one green (G) and one blue (B) signals will identify the normal chromosomes 21, x, and Y (2R1GB). Deviations<br>in the number of sex chromosomes will be detected by more or less signals of the X- or Y- chromosome than normal (as shown below); see also<br>Indications for use.  |  |              |  |
| Γ   | Expected Signals  | Female  | Male   | 1            |  |
| F   | Normal  | 2R2G  | 2R1G1B   | 1            |  |
| F   | Trisomy 21  | 3R2G  | 3R1G1B   | 1            |  |
| F   | XÓ  | 2R1G  | -  | 1            |  |
|   | XXX   | 2R3-5G  | -  | 1            |  |
|   | XXY   |   | 2R2G1B<br>2R3-4G1B<br>2R1G1B/2R2G1B in mosaics           | ]            |  |
| References:   | Korenberg J. et al, 1994, Proc. N<br>Spathas D et al, 1994, Prenat D<br>Tepperberg et al, 2001, Prenat I  | iagn. 14(11); 1049-1054   |  |              |  |

Warning and precautions: In case of emergencies check SDS sheets for medical advice. SDS sheets may be obtained by either contacting Leica Technical Support or visiting <u>www.LeicaBiosystems.com</u>. DNA probes contain formamide which is a teratogen; do not inhale or allow skin contact. Wear gloves and a lab coat when handling DNA probes. All materials should be disposed of according to your institution's guidelines for hospital waste disposal.

| Reagent Storage and<br>Handling: | Store at 2-8 °C. Reagents should not be used after the expiration date on the vial label.   |
|----------------------------------|---|
| TECHNICAL SUPPORT                | Technical support is available at <u>www.LeicaBiosystems.com</u> or +31 20 6919181or via e-mail:<br><u>kreatech-support@leicabiosystems.com</u> . |
| CUSTOMER SERVICE                 | Kreatech probes may be ordered through Leica Customer Service +31 20 6919181 or order via e-mail:<br>purchase.orders@leica-microsystems.com.      |