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Catalog No: PA0016

Intended Use

This reagent is for *in vitro* diagnostic use.

p16 (6H12) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human p16 protein in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). p16 (6H12) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human p16 protein is achieved by first allowing the binding of p16 (6H12) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

p16 (6H12) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

6H12

Immunogen

Prokaryotic recombinant fusion protein corresponding to the entire human p16 molecule.

Specificity

Human p16 protein.

Ig Class

IgG2b

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 5 mg/L as determined by ELISA.

Dilution and Mixing

p16 (6H12) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using a BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of p16 (6H12) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

p16 (6H12) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection. The recommended staining protocol for p16 (6H12) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 2 for 20 minutes.

Results Expected

Normal Tissues

Clone 6H12 detects the p16 protein that is expressed in a highly restricted pattern in certain cells within specific tissues. The staining pattern is nuclear, but cytoplasmic staining is observed in most cell types that express p16. Immunoreactivity was observed in reticulated crypt epithelium and follicular dendritic cells in tonsil, islet cells in the pancreas, Leydig cells of the testis, scattered urothelial cells and in occasional glandular epithelial cells in endometrium. Occasional weak staining was seen in antral glands of the stomach, glandular and myo-epithelium in the breast, and in the adrenal gland. Staining of macrophages and some fibroblasts was also seen. No additional staining was observed in a variety of other normal tissues evaluated (Total number of normal cases = 149).

Tumor Tissues

Clone 6H12 stained 30/42 cervical cancers (including 25/31 squamous cell carcinomas of the cervix, 2/4 adenocarcinomas of the cervix, 2/2 endometrioid adenocarcinomas, 1/2 endocervical adenocarcinomas, 0/2 intestinal type adenocarcinomas and 0/1 serous papillary carcinomas), 21/39 oropharyngeal squamous cell carcinomas, 3/6 breast tumors (including 2/4 invasive ductal carcinomas and 1/2 fibroadenomas), 2/10 tumors of the bowel (including 2/4 adenocarcinomas of the colon, 0/3 adenocarcinomas of the rectum, 0/2 adenomas of the bowel and 0/1 adenocarcinomas of the small intestine), 2/5 hepatocellular carcinomas, 2/5 head and neck tumors (including 1/1 adenocarcinoma of the oral cavity, 1/1 pleomorphic adenoma of the salivary gland, 0/1 adenoid cystic carcinoma, 0/1 squamous cell carcinoma of the tongue and 0/1 nasopharyngeal carcinoma), 2/4 brain tumors (including 1/3 meningiomas and 1/1 astrocytoma), 2/4 adenocarcinomas of the stomach, 2/4 lung tumors (including 1/1 adenocarcinoma, 1/1 small cell carcinoma and 0/2 squamous cell carcinoma), 2/3 adenocarcinomas of the prostate, 2/2 adenocarcinomas of the endometrium, 1/5 thyroid tumors (including 1/1 follicular carcinoma, 0/1 follicular papillary adenocarcinoma and 0/3 adenomas), 1/5 metastatic tumors, 1/4 squamous cell carcinomas of the esophagus, 1/4 ovarian tumors (including 1/2 endometrioid adenocarcinomas, 0/1 granulosa cell tumor and 0/1 adenocarcinoma), 1/2 bone tumors (including 1/1 chondrosarcoma and 0/1 osteosarcoma), and 1/1 adenocarcinoma of the pancreas. No staining was detected in lymphomas (0/3), bladder tumors (0/3), renal tumors (0/3), skin tumors (0/3), adrenal tumors (0/2), seminomas (0/2) and a prostatic hyperplasia (0/1). (Total number of abnormal cases evaluated = 162).

p16 (6H12) is recommended for the detection of p16 protein in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.

Product Specific Limitations

p16 (6H12) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Schache AG, Powell NG, Cuschieri KS et al. HPV-Related Oropharynx cancer in the United Kingdom: An Evolution in the understanding of disease etiology. Cancer Research. 2016; 76(22): 6598-6606.
5. Klaes R, Friedrich T, Spitkovsky D et al. Overexpression of p16(INK4A) as a specific marker for dysplastic and neoplastic epithelial cells of the cervix uteri. Int J Cancer. 2001 Apr 15;92(2):276-84.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JG and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Schache AG, Powell NG, Cuschieri KS et al. HPV-Related Oropharynx cancer in the United Kingdom: An Evolution in the understanding of disease etiology. Cancer Research. 2016; 76(22): 6598-6606.
5. Klaes R, Friedrich T, Spitkovsky D et al. Overexpression of p16(INK4A) as a specific marker for dysplastic and neoplastic epithelial cells of the cervix uteri. Int J Cancer. 2001 Apr 15;92(2):276-84.

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