ST Infinity H&E Staining System

English.....	2
العربية (Arabic).....	8
简体中文 (Chinese Simplified).....	13
中國傳統的 (Chinese Traditional).....	18
Dansk (Danish).....	24
Nederlands (Dutch).....	30
Français (French – Canada).....	36
Français (French – France).....	42
Deutsch (German).....	48
Italiano (Italian).....	54
日本語 (Japanese).....	60
한국어 (Korean).....	66
Norsk (Norwegian).....	72
Polski (Polish).....	78
Português (Portuguese – Brazil).....	84
Português (Portuguese – Portugal).....	90
Română (Romanian).....	96
Русский (Russian).....	102
Slovenski (Slovenian).....	108
Español (Spanish – Central America).....	114
Español (Spanish – Spain).....	120
Svensk (Swedish).....	126
ไทย (Thai).....	132
Türk (Turkish).....	138
Tiếng Việt (Vietnamese).....	144

ST Infinity H&E Staining System

REF 3801698

Product Name

Leica Biosystems ST Infinity H&E Staining System.

ST Infinity H&E Staining System contains 5 proprietary reagent components and staining protocols that when used together provide consistent diagnostic quality and ease of use. This system has been optimized and validated specifically for the Leica ST5010 and ST5020 stainers.

Intended Use

Detection/Measurement

ST Infinity H&E Staining System does not detect or measure an analyte or marker. Staining System is used to provide coloration and contrast to tissue sections mounted on microscopic slides from paraffin embedded specimens. The components of the stain system stain cell nuclei blue and color the connective tissue, cytoplasm, muscle, and erythrocytes various shades of orange, pink and red.

Product Function

The results obtained through use of the ST Infinity H&E Staining System are utilized as an aid to diagnosis. The coloration and contrast the ST Infinity H&E Staining System provides to tissue sections, allows visualization to microscopic anatomy to aid in diagnosis. This visualization, when interpreted by a trained professional, is utilized alongside other information such as the patient's medical history, physical condition, as well as results from other medical testing to render a medical diagnosis.

Specific Information Provided

ST Infinity H&E Staining System is not intended for the detection, definition or differentiation of a specific disorder, condition, or risk factor. The staining demonstrated with use of this product, when used as intended, provides trained professionals information which may define the physiological or pathological state of the tissue specimen.

Automation

ST Infinity H&E Staining System is not automated but are optimized and validated for use on the Leica ST5010 and ST5020 automated stainers.

Qualitative/Quantitative

The ST Infinity H&E Staining System is quantitative.

Specimen Type

The ST Infinity H&E Staining System can be used with any paraffin embedded human or animal specimen.

Testing Population

The ST Infinity H&E Staining System is intended for use with any patient requiring histopathologic evaluation of biopsy or resection tissue for the assessment of a suspected pathology or disease.

Intended User

The ST Infinity H&E Staining System is intended for use by qualified laboratory personnel and/or designee of the laboratory.

In Vitro Diagnostic

The ST Infinity H&E Staining System is intended for *in vitro* diagnostics use only.

Test Principle

The ST Infinity H&E Staining System functions through the attraction of charged dyes to oppositely charged molecules within cells and tissues.

Calibrators & Controls

The ST Infinity H&E Staining System does not require the use of any calibrators or controls. It is recommended that control slides containing tissue fixed and processed in a similar manner to the patient specimens should be included to ensure reagents are performing and functioning properly prior to staining the test specimens.

Reagent Limitations

The ST Infinity H&E Staining System contains 5 proprietary reagent components. The reagents have been developed and optimized specifically on the Leica ST5010 and ST5020 stainers. The ST Infinity H&E Staining System is designed to be utilized with a tap water hook up to the Leica ST5010 and ST5020 stainers. The ST Infinity H&E Staining System is developed and optimized for use with reagent grade alcohols and xylene.

Applicable Products

Product Code	Material Description
3801698	ST Infinity Staining System
3801698A	ST Hematoxylin, 470 mL
3801698B	ST HemaLast, 450 mL
3801698C	ST Differentiator, 450 mL
3801698D	ST Eosin, 470 mL
3801698E	ST Bluing Agent, 450 mL

NOTE: Products listed here may not be available in all geographic regions.

Materials Not Included

The ST Infinity H&E Staining System has been designed for tissue specimens placed on non-adhesive microscope slides. The ST Infinity H&E System requires the use of reagent grade alcohols and xylene.

ST Infinity H&E Staining System

REF 3801698

Devices Required

The ST Infinity H&E Staining System was developed and optimized for use with the Leica ST5010 and ST5020 stainers.

Storage and Stability

Store reagents at room temperature (15-30°C). Please follow facility guidelines for storage of flammable products.

CAUTION: Do not use after the expiration date.

In Use Stability

The ST Infinity H&E Staining System has a 14 day in use stability and a slide capacity of 2,000 slides. After 14 days on the ST 5010 or ST5020 stainer or after staining 2,000 slides, whichever occurs first, the reagents of the ST Infinity H&E Staining System should be disposed in accordance with local governing regulations.

When the reagents are not in use, the reagent vessels should be covered to reduce evaporation to maintain performance characteristics.

Sterility

The ST Infinity H&E Staining System components are not sterile products.

Warnings/Precautions

A routine quality control slide(s) containing tissue fixed and processed in a similar manner to the test specimens should be performed prior to routine use to ensure reagents are performing as intended. Reagents should be disposed of per facility guidelines for hazardous material.

Infectious Material Status

The ST Infinity H&E Staining System does not include any infectious material. However, specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

Special Facilities

The ST Infinity H&E Staining System does not require any use of special facilities.

Specimen Handling

Specimens intended for use with the ST Infinity H&E Staining System should be well fixed with neutral buffered formalin. Following processing and paraffin embedding, section tissue at a standard thickness of (2 – 5µm).

Poor fixation, processing, and sectioning will adversely affect the staining quality.

Preparation for Use

Fixation

Any general fixative including but not limited to neutral buffered formalin may be used.

Note: Underfixation or overfixation may result in compromised staining quality.

Paraffin Sections

Following processing and paraffin embedding, section tissues at standard thickness (2 – 5µm).

Slides and Specimens

Apex Superior Adhesive or equivalent slides should be used for best results. Prior to staining, slides should be dried at 58° - 60° C for 20 - 30 minutes.

Directions for use of the ST Infinity H&E Staining System on the ST5010 XL AutoStainer.

Settings for operational parameters for the ST5010 XL AutoStainer as well as a recommended bath lay-out and protocol schedules are provided in order to optimize the performance and consistency of the ST Infinity H&E Staining System. Failure to follow these instructions may result in compromised staining quality as well as staining capacity.

ST5010 AutoStainer XL: Bath Lay-Out

The required Reagent Bath Lay-out for the ST5010 AutoStainer XL is shown below.

Empty the entire contents of each of the bottles into the designated reagent vessels and positions. All of the reagents in each ST Infinity Kit must be used together to achieve the expected performance.

Substitution of non-ST Infinity reagents or reagents from another ST Infinity kit will potentially compromise performance.

Note: Reagent vessels should be covered when not in use to reduce evaporation. Discard the staining reagents, if the minimal staining capacity (2000 slides) is not met within 14 calendar days after pouring into the reagent vessels.

ST Infinity H&E Staining System

REF 3801698

Figure 1. ST5010 Reagent Bath Lay-Out

1	2	3	4	5	6	7	Tap Water Wash 1	ST Hemalast	Tap Water Wash 3	Tap Water Wash 4	Tap Water Wash 5	Oven
Xylene	Xylene	Xylene	100% Reagent Alcohol	100% Reagent Alcohol	100% Reagent Alcohol	80% Reagent Alcohol						
Exit	18	17	16	15	14	13	12	11	10	9	8	Load
Xylene	Xylene	Xylene	Xylene	100% Reagent Alcohol	100% Reagent Alcohol	100% Reagent Alcohol	ST Eosin	80% Reagent Alcohol	ST Bluing Agent	ST Differentiator	ST Hematoxylin	

ST5010 AutoStainer XL: Staining Protocol

The immersion times and instrument settings for the ST5010 XL AutoStainer are detailed in Table 1. Three protocols for varying contrasts of the H&E stain are included in Table 1. The incubation times for the ST Hematoxylin and ST Eosin may be further varied within the parameters outlined in Table 2 to meet staining preferences. All other time requirements outlined in Table 1 must be adhered to in order to ensure quality and consistency of the staining results.

Table 1. ST5010 Staining Protocol

Settings: Dips = 2, Withdrawal Speed = 9					
Step	Station	Reagent	Time	Exact	
1	Load	N/A	0:00	N/A	
2	1	Xylene	2:00	No	
3	2	Xylene	2:00	No	
4	3	Xylene	2:00	No	
5	4	100% Reagent Alcohol	1:00	No	
6	5	100% Reagent Alcohol	1:00	No	
7	6	100% Reagent Alcohol	1:00	No	
8	7	80% Reagent Alcohol	1:00	No	
9	Water Rinse 1	Water Rinse	1:00	No	
10	Water Rinse 2*	ST HemaLast	0:30	Yes	
11	8	ST Hematoxylin	Protocol 1	2:00	Yes
			Protocol 2	3:00	Yes
			Protocol 3	5:00	Yes
12	Water Rinse 5	Water Rinse	2:00	Yes	
13	9	ST Differentiator	0:45	Yes	
14	Water Rinse 4	Water Rinse	1:00	No	
15	10	ST Bluing Agent	1:00	Yes	
16	Water Rinse 3	Water Rinse	1:00	No	
17	11	80% Reagent Alcohol	1:00	No	
18	12	ST Eosin	Protocol 1	0:30	Yes
			Protocol 2	0:50	Yes
			Protocol 3	1:00	Yes
19	13	100% Reagent Alcohol	1:00	Yes	
20	14	100% Reagent Alcohol	1:00	Yes	
21	15	100% Reagent Alcohol	1:00	Yes	
22	16	Xylene	1:00	No	
23	17	Xylene	1:00	No	
24	18	Xylene	1:00	No	

*Replace Water Rinse station 2 with a standard reagent container. The opening within the rinse manifold should be sealed with a (screw) blocking plug (part number 14045635393) to prevent a loss of water pressure to the remaining water rinse stations.

ST Infinity H&E Staining System

REF 3801698

Incubation times for ST Hematoxylin and ST Eosin may be modified within the limits outlined below (Table 2) in order to achieve staining preferences with respect to color contrast.

Table 2: ST Hematoxylin and ST Eosin Variations

Reagent	Staining time range
ST Hematoxylin	2 – 5 min
ST Eosin	30 – 90 sec

Alcohol and Xylene Rotation

The alcohols and xylenes used for deparaffinization and rehydration should be rotated after every 600 slides. Reagent rotation is defined as follows:

1. Empty the most contaminated reagent of a series and refill the container.
2. Move the remaining container(s) up one space in the sequence.
3. Place the newly filled container into the last position of the series.

The alcohols and xylenes used for dehydration and clearing before coverslipping also should be rotated after every 600 slides as described above.

The two 80% alcohols (steps 8 and 17) should be discarded and replaced with fresh 80% alcohol solutions after 600 slides.

Failure to rotate and change reagents as advised may result in inadequate deparaffinization as well as excessive carry-over and dilution of reagents or solvents.

Expected Results

By following the instructions for use, a single ST Infinity Kit should provide at least 2,000 consistent high-quality stained slides.

Directions for use of the ST Infinity H&E Staining System on the ST5020 Multistainer.

Settings for operational parameters for the ST5020 Multistainer as well as bath lay-out and protocol schedules are provided in order to optimize the performance and consistency of the ST Infinity H&E Staining System. Failure to follow these instructions may result in compromised staining quality as well as staining capacity.

ST 5020 Multistainer: Bath Lay-Out

The required Reagent Bath Lay-out for the ST5020 Multistainer is shown below in Figure 2. Empty the entire contents of each of the bottles into the designated reagent vessels and positions. All of the reagents in each ST Infinity Kit must be used together to achieve the expected performance. Substitution of non-ST Infinity reagents or reagents from another ST Infinity kit will potentially compromise performance.

Note: Reagent vessels should be covered when not in use to reduce evaporation. Discard the staining reagents if the minimal staining capacity (2000 slides) is not met within 14 calendar days after pouring into the reagent vessels.

Figure 2. ST5020 Reagent Bath Lay-Out

1	2	3	4	5	6	7	8	9	10	11	12
Xylene	Xylene	Xylene	100% Reagent Alcohol	100% Reagent Alcohol	100% Reagent Alcohol	80% Reagent Alcohol	Tap Water Wash	Tap Water Wash	Tap Water Wash	Tap Water Wash	Tap Water Wash
13	14	15	16	17	18	18	20	21	22	23	24
Xylene	Xylene	100% Reagent Alcohol	100% Reagent Alcohol	100% Reagent Alcohol	ST Eosin	80% Reagent Alcohol	ST Bluing Agent	ST Differen- tiator	ST Hema- toxylin	ST HemaLast	
U25	U26	U27	28	29	30	31	32	L33	L34	L35	L36
Xylene Exit	Xylene Exit	Xylene Exit	Xylene Exit					Load	Load	Load	Load

ST Infinity H&E Staining System

REF 3801698

ST5020 Multistainer: Staining Protocol

The immersion times and instrument settings for the ST5020 Multistainer are included in Table 3. Three protocols for varying contrasts of the H&E stain are included. The incubation times for the ST Hematoxylin and ST Eosin may be further varied within the parameters outlined in Table 4 to meet staining preferences. All other time requirements outlined in Table 3 must be adhered to in order to ensure quality and consistency of the staining results.

Table 3. ST5020 Staining Protocol

Settings: Dips = 2, Lift = 1					
Step	Station	Reagent	Time	Exact	
1	Load	N/A	0:00	N/A	
2	1	Xylene	2:00	No	
3	2	Xylene	2:00	No	
4	3	Xylene	2:00	No	
5	4	100% Reagent Alcohol	1:00	No	
6	5	100% Reagent Alcohol	1:00	No	
7	6	100% Reagent Alcohol	1:00	No	
8	7	80% Reagent Alcohol	1:00	No	
9	Water Rinse 11	Water Rinse	1:00	No	
10	23	ST HemaLast	0:30	Yes	
11	22	ST Hematoxylin	Protocol 1	2:00	Yes
			Protocol 2	3:00	Yes
			Protocol 3	5:00	Yes
12	Water Rinse 10	Water Rinse	2:00	Yes	
13	21	ST Differentiator	0:45	Yes	
14	Water Rinse 9	Water Rinse	1:00	No	
15	20	ST Bluing Agent	1:00	Yes	
16	Water Rinse 8	Water Rinse	1:00	No	
17	19	80% Reagent Alcohol	1:00	No	
18	18	ST Eosin	Protocol 1	0:30	Yes
			Protocol 2	0:50	Yes
			Protocol 3	1:00	Yes
19	17	100% Reagent Alcohol	1:00	Yes	
20	16	100% Reagent Alcohol	1:00	Yes	
21	15	100% Reagent Alcohol	1:00	Yes	
22	14	Xylene	1:00	No	
23	13	Xylene	1:00	No	
24	Exit	Xylene	1:00	No	

Incubation times for ST Hematoxylin and ST Eosin may be modified within the limits outlined below (Table 4) in order to achieve staining preferences with respect to color contrast.

Table 4: ST Hematoxylin and ST Eosin Variations

Reagent	Staining time range
ST Hematoxylin	2 – 5 min
ST Eosin	30 – 90 sec

Alcohol and Xylene Rotation

The alcohols and xylenes used for deparaffinization and rehydration should be rotated after every 600 slides. Reagent rotation is defined as follows:

1. Empty the most contaminated reagent of a series and refill the container.

ST Infinity H&E Staining System

REF 3801698

2. Move the remaining container(s) up one space in the sequence.
3. Place the newly filled container into the last position of the series.

The alcohols and xylenes used for dehydration and clearing before coverslipping also should be rotated after every 600 slides as described above.

The two 80% alcohols (steps 8 and 17) should be discarded and replaced with fresh 80% alcohol solutions after 600 slides. Failure to rotate and change reagents as advised may result in inadequate deparaffinization as well as excessive carry-over and dilution of reagents or solvents.

Expected Results

By following the instructions for use, a single ST Infinity Kit should provide at least 2,000 consistent high-quality stained slides.

Readiness for Use

Once the ST Infinity H&E Staining System is placed onto the ST5010 or ST5020 Stainer, the system is ready to be used. After the staining protocol is chosen and bath-layout is created, pour all the reagents into specified reagent vessels.

Quality Control

A routine control slides containing tissue fixed and processed in a similar manner to the test specimens should be included to ensure that the reagents are performing and functioning properly.

Expected Results

By following the instructions for use, a single ST Infinity Kit can provide at least 2,000 consistent high-quality stained slides as based on internal performance characteristics. The components of the stain system will stain cell nuclei blue and color the connective tissue, cytoplasm, muscle, and erythrocytes various shades of orange, pink and red.

Analytical Performance

The ST Infinity H&E Staining System is not used to detect a specific analyte or marker. The system is used to stain cell nuclei blue and connective tissue, cytoplasm, muscle and erythrocytes various shades of orange, pink and red. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

Clinical Performance

The ST Infinity H&E Staining System is not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the ST Infinity H&E Staining System in a clinical setting.


Disposal

All reagents of the ST Infinity H&E Staining System should be disposed in accordance with local governing regulations.

Leica Biosystems Richmond, Inc.
5205 Route 12
Richmond, IL 60071
USA
(1-844-534-2262)

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