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BOND Ready-to-Use Primary Antibody LAG-3 (12H6), 7 mL

Catalog No: PA0300

Intended Use

This reagent is for *in vitro* diagnostic use.

LAG-3 (12H6) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of LAG-3 protein in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). LAG-3 (12H6) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human LAG-3 protein is achieved by first allowing the binding of LAG-3 (12H6) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

LAG-3 (12H6) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

12H6

Immunogen

Prokaryotic recombinant protein corresponding to the 3 Immunoglobulin superfamily C2 external domains of the C-terminal region of the LAG-3 molecule.

Specificity

Human lymphocyte activation gene-3 protein.

Ig Class

IgG1, kappa

Total Protein Concentration

Approx 10 mg/mL

Antibody Concentration

Greater than or equal to 0.9 mg/L as determined by ELISA.

Dilution and Mixing

LAG-3 (12H6) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8°C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of LAG-3 (12H6) are: turbidity of the solution, odor development, and presence of precipitate. Return to 2–8°C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methylisothiazol-3(2H)-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

LAG-3 (12H6) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection. The recommended staining protocol for LAG-3 (12H6) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 2 (AR9640) for 20 minutes.

Results Expected

Normal Tissues

Clone 12H6 detects the LAG-3 protein in the cytoplasm/membrane of a proportion of immune cells. Staining was observed within immune cells in most tissue types, and particularly in lymph node, tonsil, spleen and thymus gland. Staining was also observed in alveolar macrophages. (Total number of normal cases evaluated = 122).

Tumor Tissues

Clone 12H6 stained 1/51 lymphomas (including 1/37 B-cell lymphomas, 0/8 T-cell lymphomas, 0/2 small lymphocytic lymphomas, 0/1 mantle cell lymphoma, 0/1 follicular lymphoma, 0/1 Hodgkin lymphoma and 0/1 anaplastic large cell lymphoma). No staining (other than infiltrating immune cells) was detected in bowel tumors (0/8) (including 0/6 adenocarcinomas and 0/2 adenomas), head and neck tumors (0/5) (including 0/1 adenocarcinoma, 0/1 nasopharyngeal carcinoma, 0/1 squamous cell carcinoma, 0/1 adenoid cystic carcinoma and 0/1 adenoma), breast tumors (0/5) (including 0/3 invasive ductal carcinomas and 0/2 fibroadenomas), metastatic tumors (0/5), thyroid tumors (0/5) (including 0/3 follicular adenomas, 0/1 follicular carcinoma, and 0/1 papillary thyroid carcinoma), brain tumors (0/4) (including 0/3 meningiomas and 0/1 astrocytoma), hepatocellular carcinomas (0/4), lung tumors (0/4) (including 0/2 squamous cell carcinomas, 0/1 adenocarcinoma and 0/1 small cell carcinoma), squamous cell carcinomas of the esophagus (0/3), adenocarcinomas of the stomach (0/3), ovarian tumors (0/3) (including 0/2 adenocarcinomas and 0/1 granulosa cell tumor), adrenal tumors (0/2) (including 0/1 adrenocortical carcinoma and 0/1 adenoma), transitional cell carcinomas of the bladder (0/2), bone tumours (0/2) (including 0/1 osteosarcoma and 0/1 chondrosarcoma), prostatic adenocarcinomas (0/2), renal clear cell carcinomas (0/2), seminomas (0/2), squamous cell carcinoma of the cervix (0/2), endometrial adenocarcinomas (0/2), a melanoma (0/1), an adenocarcinoma of the pancreas (0/1), a squamous cell carcinoma of the skin (0/1) and a prostatic hyperplasia (0/1). (Total number of abnormal cases evaluated = 120).

LAG-3 (12H6) is recommended for the detection of LAG-3 protein in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.

Product Specific Limitations

LAG-3 (12H6) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gandhi MK, Lambley E, Duraiswamy J et al. Expression of LAG-3 by tumor-infiltrating lymphocytes is coincident with the suppression of latent antigen-specific CD8+ T-cell function in Hodgkin lymphoma patients. Blood. 2006; 108(7):2280-2289.
5. He Y, Yu H, Rozeboom L et al. LAG-3 protein expression in non-small cell lung cancer and its relationship with PD-1/PD-L1 and tumor-infiltrating lymphocytes. Journal of Thoracic Oncology. 2017; 12(5):814-823.
6. Woo S, Li N, Bruno TC et al. Differential subcellular localization of the regulatory T-cell protein LAG-3 and the coreceptor CD4. European Journal of Immunology 2010; 40(6):1768-1777.
7. Hannier S and Triebel F. The MHC class II ligand lymphocyte activation gene-3 is co-distributed with CD8 and CD3-TCR molecules after their engagement by mAb or peptide-MHC class I complexes. International Immunology. 1999; 11(11):1745-1752.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gandhi MK, Lambley E, Duraiswamy J et al. Expression of LAG-3 by tumor-infiltrating lymphocytes is coincident with the suppression of latent antigen-specific CD8+ T-cell function in Hodgkin lymphoma patients. Blood. 2006; 108(7):2280-2289.
5. He Y, Yu H, Rozeboom L et al. LAG-3 protein expression in non-small cell lung cancer and its relationship with PD-1/PD-L1 and tumor-infiltrating lymphocytes. Journal of Thoracic Oncology. 2017; 12(5):814-823.
6. Woo S, Li N, Bruno TC et al. Differential subcellular localization of the regulatory T-cell protein LAG-3 and the coreceptor CD4. European Journal of Immunology 2010; 40(6):1768-1777.
7. Hanner S and Triebel F. The MHC class II ligand lymphocyte activation gene-3 is co-distributed with CD8 and CD3-TCR molecules after their engagement by mAb or peptide-MHC class I complexes. International Immunology. 1999; 11(11):1745-1752.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gandhi MK, Lambley E, Duraiswamy J et al. Expression of LAG-3 by tumor-infiltrating lymphocytes is coincident with the suppression of latent antigen-specific CD8+ T-cell function in Hodgkin lymphoma patients. *Blood*. 2006; 108(7):2280-2289.
5. He Y, Yu H, Rozeboom L et al. LAG-3 protein expression in non-small cell lung cancer and its relationship with PD-1/PD-L1 and tumor-infiltrating lymphocytes. *Journal of Thoracic Oncology*. 2017; 12(5):814-823.
6. Woo S, Li N, Bruno TC et al. Differential subcellular localization of the regulatory T-cell protein LAG-3 and the coreceptor CD4. *European Journal of Immunology* 2010; 40(6):1768-1777.
7. Hannier S and Triebel F. The MHC class II ligand lymphocyte activation gene-3 is co-distributed with CD8 and CD3-TCR molecules after their engagement by mAb or peptide-MHC class I complexes. *International Immunology*. 1999; 11(11):1745-1752.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gandhi MK, Lambley E, Duraiswamy J et al. Expression of LAG-3 by tumor-infiltrating lymphocytes is coincident with the suppression of latent antigen-specific CD8+ T-cell function in Hodgkin lymphoma patients. Blood. 2006; 108(7):2280-2289.

5. He Y, Yu H, Rozeboom L et al. LAG-3 protein expression in non-small cell lung cancer and its relationship with PD-1/PD-L1 and tumor-infiltrating lymphocytes. *Journal of Thoracic Oncology*. 2017; 12(5):814-823.
6. Woo S, Li N, Bruno TC et al. Differential subcellular localization of the regulatory T-cell protein LAG-3 and the coreceptor CD4. *European Journal of Immunology* 2010; 40(6):1768-1777.
7. Hannier S and Triebel F. The MHC class II ligand lymphocyte activation gene-3 is co-distributed with CD8 and CD3-TCR molecules after their engagement by mAb or peptide-MHC class I complexes. *International Immunology*. 1999; 11(11):1745-1752.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gandhi MK, Lambley E, Duraiswamy J et al. Expression of LAG-3 by tumor-infiltrating lymphocytes is coincident with the suppression of latent antigen-specific CD8+ T-cell function in Hodgkin lymphoma patients. Blood. 2006; 108(7):2280-2289.
5. He Y, Yu H, Rozeboom L et al. LAG-3 protein expression in non-small cell lung cancer and its relationship with PD-1/PD-L1 and tumor-infiltrating lymphocytes. Journal of Thoracic Oncology. 2017; 12(5):814-823.
6. Woo S, Li N, Bruno TC et al. Differential subcellular localization of the regulatory T-cell protein LAG-3 and the coreceptor CD4. European Journal of Immunology 2010; 40(6):1768-1777.
7. Hanner S and Triebel F. The MHC class II ligand lymphocyte activation gene-3 is co-distributed with CD8 and CD3-TCR antigens after their engagement by mAb or peptide-MHC class I complexes. International Immunology. 1999; 11(11):1745-1752.

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Literatura

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gandhi MK, Lambley E, Duraiswamy J et al. Expression of LAG-3 by tumor-infiltrating lymphocytes is coincident with the suppression of latent antigen-specific CD8+ T-cell function in Hodgkin lymphoma patients. Blood. 2006; 108(7):2280-2289.
5. He Y, Yu H, Rozeboom L et al. LAG-3 protein expression in non-small cell lung cancer and its relationship with PD-1/PD-L1 and tumor-infiltrating lymphocytes. Journal of Thoracic Oncology. 2017; 12(5):814-823.
6. Woo S, Li N, Bruno TC et al. Differential subcellular localization of the regulatory T-cell protein LAG-3 and the coreceptor CD4. European Journal of Immunology 2010; 40(6):1768-1777.
7. Hamier S and Triebel F. The MHC class II ligand lymphocyte activation gene-3 is co-distributed with CD8 and CD3-TCR molecules after their engagement by mAb or peptide-MHC class I complexes. International Immunology. 1999; 11(11):1745-1752.

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Literatúra

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. American Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gandhi MK, Lambley E, Duraiswamy J et al. Expression of LAG-3 by tumor-infiltrating lymphocytes is coincident with the suppression of latent antigen-specific CD8+ T-cell function in Hodgkin lymphoma patients. Blood. 2006; 108(7):2280-2289.
5. He Y, Yu H, Rozeboom L et al. LAG-3 protein expression in non-small cell lung cancer and its relationship with PD-1/PD-L1 and tumor-infiltrating lymphocytes. Journal of Thoracic Oncology. 2017; 12(5):814-823.
6. Woo S, Li N, Bruno TC et al. Differential subcellular localization of the regulatory T-cell protein LAG-3 and the coreceptor CD4. European Journal of Immunology 2010; 40(6):1768-1777.
7. Hannier S and Triebel F. The MHC class II ligand lymphocyte activation gene-3 is co-distributed with CD8 and CD3-TCR molecules after their engagement by mAb or peptide-MHC class I complexes. International Immunology. 1999; 11(11):1745-1752.

Dátum vydania

03 apríl 2020

Leica Biosystems Newcastle Ltd 
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

Leica Biosystems Canada
71 Four Valley Drive
Concord, Ontario L4K 4V8
Canada
☎ +1 800 248 0123

Leica Biosystems Inc
1700 Leider Lane
Buffalo Grove IL 60089
USA
☎ +1 800 248 0123

Leica Biosystems Melbourne
Pty Ltd
495 Blackburn Road
Mt Waverley VIC 3149
Australia
☎ +61 2 8870 3500