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BOND™ Ready-To-Use ISH DNA Positive Control Probe

Catalog No: PB0682

Intended Use

This reagent is for *in vitro* diagnostic use.

The DNA Positive Control Probe is intended for use as a positive control in formalin-fixed, paraffin-embedded tissue by DNA *in situ* hybridization (ISH) using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

The DNA Positive Control Probe is a DNA probe designed to specifically hybridize to the genomic ALU repeat sequences, which represent approximately 10% of the human genome¹. ALU repeats are a part of the short interspersed elements superfamily. Although their function is unclear, they are known to be responsible for genomic insertional mutagenesis via retrotransposition and are associated with human disease^{2,3}. The DNA Positive Control Probe is generated with a biotin label using the same procedures as applied to the BOND HPV Probes.

DNA ISH using DNA Positive Control Probe is used as a control; by the detection of human ALU repeat sequences in formalin-fixed, paraffin-embedded tissue. The method is reproducible and should result in brown nuclear staining of most human nucleated cells.

Reagents Provided

DNA Positive Control Probe is a DNA probe supplied in hybridization solution.

Total volume = 6.25 mL

Dilution and Mixing

DNA Positive Control Probe is ready to use. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Needed but Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and *in situ* hybridization staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8 °C. The product is stable under these conditions up to the expiry date indicated on the container label.

There are no obvious signs indicating contamination and/or instability. Appropriate positive and negative tissue controls should be run at the same time as test tissue.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user⁴.

Precautions

- This product is intended for *in vitro* diagnostic use.

DNA POSITIVE CONTROL PROBE

Contains Formamide (<70%).

GHS08: Health hazard.

Signal words: Danger.

H360D: May damage the unborn child.

P201: Obtain special instructions before use.

P202: Do not handle until all safety precautions have been read and understood.

P260: Do not breathe dust/fumes/gas/mist/vapours/spray.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P308+313: If exposed or concerned: Get medical advice/attention.

P314: Get medical advice/attention if you feel unwell. Restricted to professional users.

- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimise microbial contamination of reagents or an increase in non-specific staining may occur.
- Enzyme digestion incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

DNA Positive Control Probe was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with Anti-Biotin Antibody, Stringency Wash Solution and BOND Polymer Refine Detection. The recommended staining protocol for the DNA Positive Control Probe is ISH Protocol B. Enzyme pre-treatment is recommended using

the BOND Enzyme Pretreatment Kit, Enzyme 1, for 15 minutes. Appropriate tissue and reagent controls should always be used. The protocol for the tissue and reagent controls should correspond to that of the DNA test probe.

To gain optimal performance from the DNA Positive Control Probe tissue specimens should be placed at least 5 mm from the frosted end of the slide and 9 mm from the non-frosted end of the slide.

Quality Control

Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls in addition to the following procedures.

Positive Tissue Control

Used to indicate correctly prepared tissues and proper staining techniques. One positive tissue control should be included for each set of test conditions in each staining run. A tissue with weak positive staining is more suitable than a tissue with strong positive staining for optimal quality control and to detect minor levels of reagent degradation.

Negative Tissue Control

Should be examined after the positive tissue control to verify the specificity of the labeling of probe to the target. Alternatively, the variety of different cell types present in most tissue sections frequently offers negative control sites, but this should be verified by the user.

Negative Reagent Control

Use DNA Negative Control PB0731 in place of the DNA test probe with a section of each patient specimen to evaluate non-specific staining and allow better interpretation of specific staining at the target.

Positive Reagent Control

Use DNA Positive Control Probe PB0682 in place of the DNA test probe with a section of each patient specimen to provide information on the preservation of nucleic acids in the tissue as well as accessibility of nucleic acids to the probe. If the DNA Positive Control Probe fails to demonstrate positive staining, results with the test specimens should be considered invalid.

Patient Tissue

Examine patient specimens stained with DNA test probe last. Positive staining intensity should be assessed within the context of any non-specific background staining of the DNA Negative Control Probe PB0731.

Results Expected

Normal Tissues

Where measurable, PB0682 detected DNA in a wide range of tissues. (Total number of normal cases evaluated = 118).

Abnormal Tissues

PB0682 detected DNA in 4/4 lung tumors, 4/4 ovarian tumors, 3/3 thyroid tumors, 2/2 esophageal tumors, 2/2 breast tumors, 2/2 stomach tumors, 2/2 soft tissue tumors, 2/2 tongue tumors, 2/2 metastatic tumors of unknown origin, 2/2 liver tumors, 2/2 kidney tumors, 2/2 cervical tumors, 2/2 testicular tumors, 2/2 colon tumors, 2/2 rectal tumors, 2/2 skin tumors, 1/1 larynx tumor 1/1 thymus tumor and 1/1 brain tumor. (Total number of abnormal cases evaluated = 40).

PB0682 is recommended as a screening tool to detect the preservation of DNA in cells.

Product Specific Limitations

DNA Positive Control Probe has been optimized by Leica Biosystems for use with Anti-Biotin Antibody, Stringency Wash Solution, BOND Polymer Refine Detection and BOND ancillary reagents on formalin-fixed, paraffin-embedded cervical tissue. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue type, fixation and processing. In addition, BOND Enzyme concentration and incubation time may require optimization depending on tissue type, processing and fixation conditions. Negative reagent controls should be used when optimizing pre-treatment conditions and protocol times.

Troubleshooting

Reference 6 may aid in remedial action.

Test samples should be complemented by the appropriate tissue and reagent controls.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on *in situ* hybridization with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. International Human Genome Sequencing Consortium., 2001. Initial sequencing and analysis of the human genome. *Nature* 409 pp. 860–921.
2. Batzer MA and Deininger PL., 2002 Alu repeats and human genomic diversity. *Nat Rev Genet.* 372 pp. 370–9.
3. Schmid CW., 1998. Does SINE evolution preclude Alu function? *Nucleic Acids Res.* 26 pp.4541–4550.
4. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
5. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
6. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

Date of Issue

24 September 2018

PB0682
Page 3

BOND™ Ready-To-Use ISH DNA Positive Control Probe

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Bibliographie

1. International Human Genome Sequencing Consortium., 2001. Initial sequencing and analysis of the human genome. Nature 409 pp. 860–921.
2. Batzer MA and Deininger PL., 2002 Alu repeats and human genomic diversity. Nat Rev Genet. 372 pp. 370–9.
3. Schmid CW., 1998. Does SINE evolution preclude Alu function? Nucleic Acids Res. 26 pp.,4541–4550.
4. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
5. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
6. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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Bibliografia

1. International Human Genome Sequencing Consortium., 2001. Initial sequencing and analysis of the human genome. *Nature* 409 pp. 860–921.
2. Batzer MA and Deininger PL., 2002 Alu repeats and human genomic diversity. *Nat Rev Genet.* 372 pp. 370–9.
3. Schmid CW., 1998. Does SINE evolution preclude Alu function? *Nucleic Acids Res.* 26 pp. 4541–4550.
4. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
5. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
6. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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Bibliografie

1. International Human Genome Sequencing Consortium., 2001. Initial sequencing and analysis of the human genome. Nature 409 pp. 860–921.
2. Batzer MA and Deininger PL., 2002 Alu repeats and human genomic diversity. Nat Rev Genet. 372 pp. 370–9.
3. Schmid CW., 1998. Does SINE evolution preclude Alu function? Nucleic Acids Res. 26 pp.4541–4550.
4. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
5. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
6. Wilkinson DG. The theory and practice of *in situ* hybridization. In: Wilkinson DG. (ed.) *In situ* Hybridization A practical approach. 2nd Edition. New York: Oxford University Press, 1998, pp.18–20.

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1. International Human Genome Sequencing Consortium.,2001. Initial sequencing and analysis of the human genome.Nature 409 pp.860–921.
2. Batzer MA and Deininger PL.,2002 Alu repeats and human genomic diversity.Nat Rev Genet.372 pp.370–9.
3. Schmid CW.,1998. Does SINE evolution preclude Alu function?Nucleic Acids Res.26 pp.4541–4550.
4. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
5. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
6. Wilkinson DG.The theory and practice of *in situ* hybridization.In:Wilkinson DG.(ed.)*In situ* Hybridization A practical approach.2nd Edition.New York:Oxford University Press, 1998, pp.18–20.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd