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# BOND Ready-to-Use Primary Antibody S-100 (EP32)

## Catalog No: PA0031

### Intended Use

*This reagent is for in vitro diagnostic use.*

S-100 (EP32) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of S-100B protein in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

### Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). S-100 (EP32) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of S-100B is achieved by first allowing the binding of S-100 (EP32) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

### Reagents Provided

S-100 (EP32) is a rabbit anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

### Clone

EP32

### Immunogen

A synthetic peptide corresponding to residues of human S100 Beta protein.

### Specificity

Human.

### Ig Class

Rabbit IgG

### Total Protein Concentration

Total Protein

Approx 10 mg/mL

### Antibody Concentration

Greater than or equal to 0.15 mg/L as determined by ELISA.

### Dilution and Mixing

S-100 (EP32) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

### Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

### Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of S-100 (EP32) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

### Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

## Instructions for Use

S-100 (EP32) primary antibody is recommended for use on an automated BOND system in combination with either BOND Polymer Refine Detection (DS9800) or BOND Polymer Refine Red Detection (DS9390). The recommended staining protocols for S-100 (EP32) primary antibody are IHC Protocol F when using BOND Polymer Refine Detection and IHC Protocol J when using BOND Polymer Refine Red Detection. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 1 for 20 minutes.

## Results Expected

### Normal Tissues

Clone EP32 detects S-100 protein in the nucleus and cytoplasm of cells of neuroectodermal origin in a variety of tissues. Immunoreactivity was observed in nerve fibres, dendritic cells, adipocytes and a percentage of macrophages, lymphocytes and plasma cells.

Staining was also observed in parenchymal cells in the cerebrum and cerebellum, islets of Langerhans in pancreas, myoepithelial cells of salivary glands, retinal cells of the eye and myoepithelial and ductal cells of the breast. (Total number of normal cases = 117)

### Tumor Tissues

Clone EP32 stained 20/42 peripheral nerve tumours (including 10/12 malignant Schwannomas, 4/4 ganglioneuromas, 3/21 malignant peripheral nerve sheath tumours, 3/4 neurofibromas and 0/1 primitive neuroectodermal tumours), 9/11 malignant melanomas. No staining was detected in other various tumours including GI tract tumours (0/12), breast tumours (0/6), lung tumours (0/6), ovarian tumours (0/3), hepatocellular carcinomas (0/3), cervical tumours (0/3), endometrial tumours (0/3), bladder tumours (0/2), kidney clear cell carcinomas (0/2), prostate adenocarcinomas (0/1), prostate hyperplasia (0/1), and squamous cell carcinoma of the skin (0/1). (Total number of abnormal cases evaluated = 96)

**S100 (EP32) is recommended for the detection of S-100 B protein in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.**

## Product Specific Limitations

S-100 (EP32) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

## Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

## Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

## Bibliography - General

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Chen H, Xu C, Jin Q et al. S100 protein family in human cancer. American Journal of Cancer Research 2014;4(2):89-115
5. Torlakovic EE, Nielsen S, Francis G et al. Standardization of Positive Controls in Diagnostic Immunohistochemistry. Applied Immunohistochemistry & Molecular Morphology. 2015; 23(1):1-18

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Chen H, Xu C, Jin Q et al. S100 protein family in human cancer. American Journal of Cancer Research 2014;4(2):89-115
5. Tortlakovic EE, Nielsen S, Francis G et al. Standardization of Positive Controls in Diagnostic Immunohistochemistry. Applied Immunohistochemistry & Molecular Morphology. 2015; 23(1):1-18

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd