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BOND Ready-to-Use Primary Antibody TIM-3 (EPR22241), 7 mL

Catalog No: PA0360

Intended Use

This reagent is for *in vitro* diagnostic use.

TIM-3 (EPR22241) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of TIM-3 protein in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). TIM-3 (EPR22241) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human TIM-3 protein is achieved by first allowing the binding of TIM-3 (EPR22241) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

TIM-3 (EPR22241) is a rabbit anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

EPR22241

N.B. This TIM-3 (EPR22241) antibody has been created by Epitomics Inc., using Epitomics' proprietary rabbit monoclonal antibody technology covered under Patent No.'s 5,675,063 and 7,402,409.

Immunogen

Recombinant fragment within Human TIM 3 aa 200-300.

Specificity

Human TIM-3 Protein.

Ig Class

IgG

Total Protein Concentration

Approx 10 mg/mL

Antibody Concentration

Greater than or equal to 0.5 mg/L.

Dilution and Mixing

TIM-3 (EPR22241) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Storage and Stability

Store at 2–8°C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of TIM-3 (EPR22241) are: turbidity of the solution, odor development, and presence of precipitate. Return to 2–8°C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methylisothiazol-3(2H)-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

TIM-3 (EPR22241) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection. The recommended staining protocol for TIM-3 (EPR22241) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 2 (AR9640) for 20 minutes.

Results Expected

Normal Tissues

Clone EPR22241 detects the TIM-3 protein in the cytoplasm/membrane of a proportion of immune cells. Staining was observed within immune cells in most tissue types, particularly in lymph node, tonsil, spleen, bone marrow, and thymus gland. Staining was also observed in breast epithelial cells, Kupffer cells in liver, alveolar macrophages, Langerhans cells in skin and proximal tubules and brush borders in kidney. Weak staining was observed in cells within the medulla and cortex of the adrenal gland. (Total number of normal cases evaluated = 132).

Tumor Tissues

Clone EPR22214 stained 4/4 brain tumors (including 3/3 meningiomas and 1/1 astrocytoma), 3/78 bone tumors (including 3/50 osteosarcomas and 0/28 chondrosarcomas), 2/3 lymphomas (including 1/1 non-Hodgkin B-cell lymphoma, 1/1 anaplastic large cell lymphoma and 0/1 Hodgkin lymphoma), 2/2 renal clear cell carcinomas, 1/5 breast tumors (including 1/3 invasive ductal carcinomas and 0/2 fibroadenomas) and 1/4 hepatocellular carcinomas. No staining (other than infiltrating immune cells) was detected in pancreatic adenocarcinomas (0/11), prostate adenocarcinomas (0/8), bowel tumors (0/8) (including 0/6 adenocarcinomas and 0/2 adenomas), thyroid tumors (0/5) (including 0/3 follicular adenomas, 0/1 follicular carcinoma and 0/1 papillary thyroid carcinoma), metastatic tumors (0/5), lung tumors (0/4) (including 0/2 squamous cell carcinomas, 0/1 adenocarcinoma and 0/1 small cell carcinoma), squamous cell carcinomas of the esophagus (0/3), adenocarcinomas of the stomach (0/3), ovarian tumors (0/3) (including 0/1 granulosa cell tumor, 0/1 adenocarcinoma and 0/1 endometrioid adenocarcinoma), squamous cell carcinomas of the cervix (0/2), endometrial adenocarcinomas (0/2), adrenal gland tumors (0/2), (including 0/1 cortical adenoma and 0/1 adrenocortical carcinoma), transitional cell carcinomas of the bladder (0/2), salivary gland tumors (0/2) (including 0/1 pleomorphic adenoma, 0/1 adenoid cystic carcinoma), seminomas (0/2), melanoma (0/1), oral adenocarcinoma (0/1), squamous cell carcinoma of the tongue (0/1), squamous cell carcinoma of skin (0/1) and prostatic hyperplasia (0/1) (Total number of abnormal cases evaluated = 163).

TIM-3 (EPR22241) is recommended for the detection of TIM-3 protein in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.

Product Specific Limitations

TIM-3 (EPR22241) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gorris MAJ, Hallilovic A, Rabold K et al. Eight color multiplex immunohistochemistry for simultaneous detection of multiple immune checkpoint molecules within the tumor microenvironment. The Journal of Immunology. 2017; 200, 1701262.
5. Wu J, Lin G, Zhu Y et al. Low TIM-3 expression indicated poor prognosis of metastatic prostate cancer and acts as an independent predictor of castration resistant status. Scientific Reports 2017; 7:8869.
6. Anderson AC. TIM-3: an emerging target in the cancer immunotherapy landscape. Cancer Immunology Research. 2014; 2:393- 398
7. Gao X, Zhu Y, Gang L et al. TIM-3 expression characterises regulatory T cells in tumour tissues and is associated with lung cancer progression. PLoS ONE. 2012; 7(2): e30676.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gorris MAJ, Hallilovic A, Rabold K et al. Eight color multiplex immunohistochemistry for simultaneous detection of multiple immune checkpoint molecules within the tumor microenvironment. The Journal of Immunology. 2017; 200, 1701262.
5. Wu J, Lin G, Zhu Y et al. Low TIM-3 expression indicated poor prognosis of metastatic prostate cancer and acts as an independent predictor of castration resistant status. Scientific Reports 2017; 7:8869.

6. Anderson AC. TIM-3: an emerging target in the cancer immunotherapy landscape. *Cancer Immunology Research*. 2014; 2:393- 398
7. Gao X, Zhu Y, Gang L et al. TIM-3 expression characterises regulatory T cells in tumour tissues and is associated with lung cancer progression. *PLoS ONE*. 2012; 7(2): e30676.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gorris MAJ, Hallilovic A, Rabold K et al. Eight color multiplex immunohistochemistry for simultaneous detection of multiple immune checkpoint molecules within the tumor microenvironment. The Journal of Immunology. 2017; 200, 1701262.

5. Wu J, Lin G, Zhu Y et al. Low TIM-3 expression indicated poor prognosis of metastatic prostate cancer and acts as an independent predictor of castration resistant status. *Scientific Reports* 2017; 7:8869.
6. Anderson AC. TIM-3: an emerging target in the cancer immunotherapy landscape. *Cancer Immunology Research*. 2014; 2:393- 398
7. Gao X, Zhu Y, Gang L et al. TIM-3 expression characterises regulatory T cells in tumour tissues and is associated with lung cancer progression. *PLoS ONE*. 2012; 7(2): e30676.

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Kaynakça

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 17163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gorris MAJ, Hallilovic A, Rabold K et al. Eight color multiplex immunohistochemistry for simultaneous detection of multiple immune checkpoint molecules within the tumor microenvironment. The Journal of Immunology. 2017; 200, 1701262.
5. Wu J, Lin G, Zhu Y et al. Low TIM-3 expression indicated poor prognosis of metastatic prostate cancer and acts as an independent predictor of castration resistant status. Scientific Reports 2017; 7:8869.
6. Anderson AC. TIM-3: an emerging target in the cancer immunotherapy landscape. Cancer Immunology Research. 2014; 2:393- 398
7. Gao X, Zhu Y, Gang L et al. TIM-3 expression characterises regulatory T cells in tumour tissues and is associated with lung cancer progression. PLoS ONE. 2012; 7(2): e30676.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.

4. Gorris MAJ, Hallilovic A, Rabold K et al. Eight color multiplex immunohistochemistry for simultaneous detection of multiple immune checkpoint molecules within the tumor microenvironment. *The Journal of Immunology*. 2017; 200, 1701262.
5. Wu J, Lin G, Zhu Y et al. Low TIM-3 expression indicated poor prognosis of metastatic prostate cancer and acts as an independent predictor of castration resistant status. *Scientific Reports* 2017; 7:8869.
6. Anderson AC. TIM-3: an emerging target in the cancer immunotherapy landscape. *Cancer Immunology Research*. 2014; 2:393- 398
7. Gao X, Zhu Y, Gang L et al. TIM-3 expression characterises regulatory T cells in tumour tissues and is associated with lung cancer progression. *PLoS ONE*. 2012; 7(2): e30676.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 1763 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Gorris MAJ, Hallilovic A, Rabold K et al. Eight color multiplex immunohistochemistry for simultaneous detection of multiple immune checkpoint molecules within the tumor microenvironment. The Journal of Immunology. 2017; 200, 1701262.
5. Wu J, Lin G, Zhu Y et al. Low TIM-3 expression indicated poor prognosis of metastatic prostate cancer and acts as an independent predictor of castration resistant status. Scientific Reports 2017; 7:8869.
6. Anderson AC. TIM-3: an emerging target in the cancer immunotherapy landscape. Cancer Immunology Research. 2014; 2:393- 398
7. Gao X, Zhu Y, Gang L et al. TIM-3 expression characterises regulatory T cells in tumour tissues and is associated with lung cancer progression. PLoS ONE. 2012; 7(2): e30676.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd