Congo Red, Amyloid Special Stain Kit

English	2
العربية (Arabic)	6
简体中文 (Chinese Simplified)	9
中國傳統的 (Chinese Traditional)	13
Dansk (Danish)	17
Nederlands (Dutch)	21
Français (French – Canada)	25
Français (French – France)	29
Deutsch (German)	33
Italiano (Italian)	37
日本語 (Japanese)	41
한국어 (Korean)	45
Norsk (Norwegian)	49
Polski (Polish)	53
Português (Portuguese – Brazil)	57
Português (Portuguese – Portugal)	61
Română (Romanian)	65
Русский (Russian)	69
Slovenski (Slovenian)	73
Español (Spanish – Central America)	77
Español (Spanish – Spain)	81
Svensk (Swedish)	85
ไทย (Thai)	89
Türk (Turkish)	93
Tiếng Việt (Vietnamese)	97

Congo Red, Amyloid Special Stain Kit

REF 38016SS9

Product Name

Congo Red, Amyloid Stain, Special Stain Kit

Intended Use

Detection/Measurement

The Leica Biosystems Congo Red, Amyloid Stain, Special Stain Kit does not detect or measure an analyte or marker. The Congo Red, Amyloid Stain Kit when used with appropriate histological protocols may be used to provide visualization of amyloid deposits in tissue sections.

Product Function

The results obtained through use of Congo Red, Amyloid Stain, Special Stain Kit do not provide objective medical evidence. The coloration and contrast the Leica Biosystems Stain Kit provide to histologic specimens allows visualization of microscopic anatomy. This visualization, when interpreted by a trained professional, is utilized alongside other information such as the patient's medical history, physical condition, as well as results from other medical testing to render a medical diagnosis.

Specific Information Provided

Congo Red, Amyloid Stain Special Stain Kit is not intended for the detection, definition or differentiation of a specific disorder, condition, or risk factor. The staining demonstrated with use of these products, when used as intended, provides trained professionals information which may define the physiological or pathological state of the tissue specimen.

Automation

Congo Red, Amyloid Stain, Special Stain Kit is not automated but can be used on automated staining platforms. Use on an automated platform should be validated at the point of use.

Qualitative/Quantitative

The Leica Biosystems Congo Red, Amyloid Stain, Special Stain kit is qualitative stain.

Specimen Type

Congo Red, Amyloid Stain, Special Stain Kit can be used with any paraffin embedded human or animal specimen.

Testing Population

The Leica Biosystems Congo Red, Amyloid Stain, Special Stain kit is intended for use with any patient requiring evaluation of biopsy or resection tissue for the assessment of a suspected pathology or disease.

Intended User

Congo Red, Amyloid Stain, Special Stain Kit is intended for use by qualified laboratory personnel and/or designee of the laboratory.

In Vitro Diagnostic

Congo Red, Amyloid Stain, Special Stain Kit is intended for *in vitro* diagnostics use only.

Test Principle

Amyloid is an amorphous eosinophilic material that is composed primarily of fibrillary proteins arranged in a linear β pleated sheet orientation. Deposits of amyloid occur in various organs including spleen, liver and kidneys as a consequence of chronic inflammatory diseases. Other diseases associated with amyloid deposits include Multiple Myeloma and Alzheimer's dementia.

Although the exact mechanism by which Congo Red stains amyloid is unknown, it is generally believed that the Congo Red molecule binds the linearly arranged amyloid fibrils through hydrogen bonds. When viewed with polarized light microscopy, red to green birefringence occurs due to the parallel alignment of the dye molecules on these linearly arranged amyloid fibrils.

The specificity of this technique is enhanced by using a staining solution with an elevated pH and high salt concentration. These conditions reduce the likelihood of non-specific electrostatic interactions of Congo Red with molecules other than amyloid.

Positive staining with Congo Red and the resultant apple green birefringence with cross polarization, is generally considered the most specific method available to the light microscopist for the detection of amyloid.

Calibrators & Controls

Congo Red, Amyloid Stain, Special Stain Kit does not require the use of any calibrators or controls.

Reagent Limitations

No reagent limitations are applicable to this product.

Applicable Products

Product Code	Material Description
38016SS9	Congo Red, Amyloid Stain, Special Stain Kit
38016SS9A	Congo Red Solution, 500 mL
38016SS9B	Potassium Hydroxide Solution, 100 mL
38016SS4C	Gill II Hematoxylin, 500 mL

Congo Red, Amyloid Special Stain Kit

REF 38016SS9

Materials Not Included

Congo Red, Amyloid Stain, Special Stain Kit protocol requires the use of graded alcohols, xylene, or xylene substitutes, deionized or distilled water. Positive amyloid control slide(s), not included in this kit, should be included in each run.

Devices Required

Congo Red, Amyloid Stain, Special Stain Kit can be used on any automated staining platform or with a manual staining method.

Storage and Stability

Store at room temperature (15-30 °C).

CAUTION: Do not use after the expiration date.

In Use Stability

User discretion should be utilized when determining in-use stability.

Sterility

Congo Red, Amyloid Stain, Special Stain Kit components are not sterile products.

Warnings/Precautions

Normal precautions exercised in handling laboratory reagents should be followed. Dispose of waste observing all local, state, provincial or national regulations. Refer to Material Safety Data Sheet and product labeling for any updated risk, hazard, or safety information.

Infectious Material Status

Congo Red, Amyloid Stain, Special Stain Kit does not include any infectious material. However, specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

Special Facilities

Congo Red, Amyloid Stain, Special Stain Kit should be used per facility guidelines.

Specimen Handling

- Fixation
Neutral buffered formalin is recommended for fixation. Prolonged storage in formalin-based fixatives may result in a gradual decrease in staining intensity of amyloid deposits. Other fixatives recommended in literature citations include absolute alcohol, Carnoy's, and Bouin's fixative.
- Paraffin Sections
Following processing and paraffin embedding, section tissue specimens at 6–10 microns. Amyloid deposits in sections less than 6 microns in thickness may fail to demonstrate birefringence.

Preparation for Use

- Add 0.5 mL of the Potassium Hydroxide Solution to 50 mL of the Congo Red Solution.
- Mix well, filter and store in tightly capped bottle.
- The Congo Red working solution is stable for up to 7 days.

Note: Exercise caution when handling Potassium Hydroxide Solution.

Direction for Use

Conventional Staining Protocol

1. Deparaffinize tissue sections with xylene and hydrate through graded alcohols to deionized or distilled water.
2. Rinse in 3 changes of deionized water.
3. Place in 95% ethanol for 5-10 seconds.
4. Stain in the Congo Red working solution for 20 minutes. (A)
5. Rinse in 5-8 changes of deionized water.
6. Stain in Gill II Hematoxylin for 1-3 minutes.
7. Wash in running tap water for 2 minutes.
8. Rinse in 3 changes of deionized water.
9. Dehydrate in two changes of 95% alcohol and three changes of absolute alcohol for one minute each. (B)
10. Clear in two changes of xylene and mount in a xylene miscible medium.

Notes: (A). Background staining may increase if the Congo Red working solution is applied for more than 25 minutes.

(B). Dehydrating alcohols must be eosin free, as eosin may result in nonspecific birefringence.

Congo Red, Amyloid Special Stain Kit

REF 38016SS9

Table 1. Example of Conventional Congo Red, Amyloid Stain Protocol.

Steps	Action	Chemical	Time (mm: ss)
1-3	Deparaffinize	Xylene	3:00
4-5	Hydration	100% Alcohol	2:00
6	Hydration	80% or 95% Alcohol	1:00
7-9	Rinse	Deionized Water	0:30
10	Hydration	95% Ethanol	0:05 to 0:10
11	Stain	Congo Red working solution	20:00
12-16	Rinse	Deionized Water	0:30
17	Counterstain	Gill II Hematoxylin	1:00 to 3:00
18	Wash	Water Wash	2:00
19-21	Rinse	Deionized Water	0:30
22-23	Dehydration	95% Alcohol	1:00
24-25	Dehydration	100% Alcohol	1:00
26-27	Clearing	Xylene	2:00

Microwave Staining Protocol

Exercise caution when using the microwave to heat any solution or reagent. The microwave must be properly ventilated to prevent the accumulation of fumes in the laboratory. Microwave transparent coplin jars and caps should be used during the staining process. The caps should be loosely applied to prevent spills. Caps with ventilation holes also may be used. All microwaves should be used in accordance with the manufacturer's instructions.

1. Deparaffinize tissue sections with xylene and hydrate through graded alcohols to deionized or distilled water.
2. Rinse in 3 changes of deionized water.
3. Place in 95% ethanol for 5-10 seconds.
4. Place sections in a plastic Coplin jar containing the Congo Red working solution and microwave at 500 watts for 30 seconds.
5. Rinse in 5-8 changes of deionized water.
6. Stain in Gill II Hematoxylin for 1-3 minutes.
7. Wash in running tap water for 2 minutes
8. Rinse in 3 changes of deionized water.
9. Dehydrate in two changes of 95% alcohol and three changes of absolute alcohol for one minute each. (A)
10. Clear in two changes of xylene and mount in a xylene miscible medium.

Note: (A). Dehydrating alcohols must be eosin free, as eosin may result in nonspecific birefringence.

Readiness for Use

Once appropriate staining protocol is chosen and bath-layout is created, pour all the reagent into the reagent vessel. Place the reagent vessel back into the respective station.

Quality Control

Tissue sections containing amyloid, fixed and processed in a similar manner to the test specimens, should be included in each staining assay to ensure Special Staining Kit is performing as intended.

A loss of staining intensity has been reported to occur in tissue sections stored for an extended period of time.

Expected Results

Conventional light microscopy

- Amyloid, eosinophils, elastin — dull to brick red
- Nuclei — purple/blue

Polarized microscopy

- Amyloid — apple green birefringence

Apple green birefringence detected by polarized microscopy is considered the most specific indicator of amyloid.

Massive, long standing amyloid deposits may display diminished birefringence as compared to newly formed deposits.

Congo Red, Amyloid Special Stain Kit

REF 38016SS9

Analytical Performance

The Leica Biosystems Congo Red Stain Kits are not used to detect a specific analyte or marker. These products are used to provide visualization of amyloid deposits in tissue sections. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

Clinical Performance

The Leica Biosystems Congo Red are not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the Leica Biosystems Bluing Agents in a clinical setting.

Disposal

Spent or expired components of the Congo Red, Amyloid Stain, Special Stain Kit should be discarded in accordance with organizational, local, state, and federal regulations.

Leica Biosystems Richmond, Inc.
5205 Route 12
Richmond, IL 60071
USA
(1-844-534-2262)

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