

Kreatech™ FISH probes

Product Information Sheet

10C002D590

SE 10 (D10Z1)-002


Dark Red

25 µl

DANGER



FORMAMIDE




Kreatech Biotechnology B.V.
Vlierweg 20
1032 LG Amsterdam
The Netherlands
www.LeicaBiosystems.com

Analyte Specific Reagent

Analytical and performance characteristics are not established

PI-10C002D590_D2.1

Published Sept 2015



Not to scale

Kreatech™ SE 10 (D10Z1)-002 FISH probe

Introduction: The **Satellite Enumeration (SE) 10 (D10Z1)-002** FISH probe is optimized to identify chromosomes 10.

Critical region (dark red): The **SE 10 (D10Z1)-002** FISH probe is direct-labeled with PlatinumBright™590.

Reagent: The **SE 10 (D10Z1)-002** FISH probe is a direct-labeled DNA probe provided at four times the concentration than the amount of DNA probe used in the QC test procedure.

Kreatech FISH probes are REPEAT-FREE™ and therefore do not contain Cot-1 DNA. Hybridization efficiency is increased and background, due to unspecific binding, is highly reduced.

Patterns: Two dark red (2DR) signals will identify the non-aberrant SE 10 loci..

	Signal Pattern
Expected Signals	2DR

Warning and precautions: In case of emergencies check SDS sheets for medical advice. SDS sheets may be obtained by either contacting Leica Technical Support or visiting www.LeicaBiosystems.com. DNA probes contain formamide which is a teratogen; do not inhale or allow skin contact. Wear gloves and a lab coat when handling DNA probes. All materials should be disposed of according to your institution's guidelines for hospital waste disposal.

Reagent Storage and Handling: Store at 2-8 °C. Reagents should not be used after the expiration date on the vial label.

TECHNICAL SUPPORT Technical support is available at www.LeicaBiosystems.com/service-support/technical-support/ or toll free at 800-248-0123 or via e-mail: kreatech-support@leicabiosystems.com.

CUSTOMER SERVICE Kreatech probes may be ordered through Leica Customer Service toll free at 800-248-0123 or order via e-mail: purchase.orders@leica-microsystems.com.