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ChromoPlex™ 1 Dual Detection for BOND - 50 Test

Catalog No: DS9665

Intended Use

This detection system is for *in vitro* diagnostic use.

ChromoPlex™ 1 Dual Detection is a biotin-free, polymeric horseradish peroxidase (HRP)-linker and polymeric alkaline phosphatase (AP)-linker antibody conjugate system for the detection of tissue-bound mouse and rabbit IgG primary antibodies. It is intended for staining sections of formalin-fixed, paraffin-embedded tissue on the BOND automated system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies. Proper controls should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

The ChromoPlex 1 Dual Detection must be used with laboratory best practice in the use of tissue controls. For assurance, laboratories should stain each patient sample in conjunction with positive, negative and other tissue specific controls as needed.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation).

ChromoPlex 1 Dual Detection utilizes a novel controlled polymerization technology to prepare polymeric HRP-linker and AP-linker antibody conjugates. The detection system avoids the use of streptavidin and biotin, and therefore eliminates non-specific staining as a result of endogenous biotin.

ChromoPlex 1 Dual Detection works as follows:

- The specimen is incubated with hydrogen peroxide to quench endogenous peroxidase activity.
- A user-supplied specific primary mouse and rabbit antibody cocktail is applied.
- Poly-HRP IgG reagent localizes mouse antibodies.
- Poly-AP IgG reagent localizes rabbit antibodies.
- The first substrate chromogen, 3,3'-Diaminobenzidine tetrahydrochloride hydrate (DAB), visualizes mouse antibodies via a brown precipitate.
- The second substrate chromogen (Fast Red) visualizes rabbit antibodies via a red precipitate.
- Hematoxylin (blue) counterstaining allows the visualization of cell nuclei.

Using ChromoPlex 1 Dual Detection in combination with the BOND automated system (includes Leica BOND-MAX system and Leica BOND-III system) reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

The reagents provided are sufficient for 25 individual BOND staining runs, or a maximum of 50 slides.

To achieve a maximum of 50 slides from this detection system, slides must be batched in quantities of 2 or greater, per Slide Staining Assembly. Batching in quantities of less than 2 will result in fewer stained slides.

1. Peroxide Block (7.5 mL) 3-4% Hydrogen peroxide.
2. Polymer mHRP (7.5 mL) Poly-HRP anti-mouse containing 10% (v/v) animal serum in Tris-buffered saline and 0.09% ProClin™ 950.
3. Polymer rAP (7.5 mL) Poly-AP anti-rabbit containing 10% (v/v) animal serum in Tris-buffered saline and 0.09% ProClin 950.
4. DAB Part 1 (1 mL) 66 mM 3,3'-Diaminobenzidine tetrahydrochloride hydrate in a stabilizer solution.
5. DAB Part B (20 mL) ≤0.1% (v/v) Hydrogen peroxide in a stabilizer solution.
6. Red Part A (5.5 mL) Activator containing 0.5% ProClin 950.
7. Red Part B (1.0 mL) Substrate.
8. Red Part C (1.0 mL) Substrate.
9. Red Part D (35 mL) Buffer solution containing 0.5% ProClin 950.
10. Hematoxylin DS9665 (7.5 mL) <0.1% Hematoxylin.

Dilution and Mixing

ChromoPlex 1 Dual Detection is optimized for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

Reconstitution, mixing, dilution, or titration of these reagents is not required.

Materials Required But Not Provided

Mounting media compatible with ChromoPlex 1 Dual Detection (see **Precautions**).

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system.

Storage and Stability

Store at 2–8 °C. Do not freeze. Do not use after the expiration date indicated on the tray handle label. Return to 2–8 °C immediately after use.

There are no obvious signs to indicate instability of this product; therefore positive and negative controls should be run simultaneously with unknown specimens (refer to "Quality Control" in the "Using BOND Reagents" section of your BOND user documentation).

If unexpected staining is observed that cannot be explained by variations in laboratory procedures, and a problem with the detection system is suspected, contact your local distributor or the regional office of Leica Biosystems immediately.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- WARNING - some mounting media may cause the red chromogen to crystallize or fade.
- DO NOT USE DPX, Entellan[®], Leica CV Mount, Eukitt[®] quick-hardening mounting medium, Surgipath MM 24[®], Surgipath Sub-X[®] and Pertex mounting media as they are not suitable for use with this product.
- Fading or crystallization has not been observed with: Leica CV Ultra Mounting Media (Cat 14070937891), Poly-Mount[®] (Polysciences Inc. Cat 08381), VectaMount[™] (Vector Laboratories Cat H-5000) and Limonene Mount (Electron Microscopy Sciences Cat 17987-01).
- The compatibility of other mounting media should be fully validated by the user prior to use. Failure to do so may result in a reduction of staining due to fading or crystallization.

DAB Part 1

Contains Ethylene Glycol (>90%) and 66 mM (<10%) 3,3'-diaminobenzidine tetrahydrochloride hydrate.

GHS07: Exclamation mark.

GHS08: Health hazard.

Signal words: Danger.

H302: Harmful if swallowed.
H341: Suspected of causing genetic defects.

H350: May cause cancer.

P201: Obtain special instructions before use.

P202: Do not handle until all safety precautions have been read and understood.

P280: Wear protective gloves/protective clothing/eye protection/face protection.

P264: Wash hands thoroughly after handling.

P270: Do not eat, drink or smoke when using this product.

P301+312: IF SWALLOWED: Call a POISON CENTER or doctor if you feel unwell.

P330: Rinse mouth.

P308+313: IF exposed or concerned: Get medical advice/attention.

Red Part C

Contains 1.0 M (<10%) Hydrochloric Acid.

GHS05: Corrosion.

Signal words: Warning.

H290: May be corrosive to metals.

P234: Keep only in original container.

P390: Absorb spillage to prevent material damage.

- The activator and substrate solutions may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable nitrile gloves when handling reagents. Dispose of according to local regulations.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems. Alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com
- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user¹.
- Do not mix reagents from different detection systems.

Instructions for Use

The ChromoPlex 1 Dual Detection has been developed for use on the BOND automated system (includes Leica BOND-MAX system and Leica BOND-III system) using *IHC Protocol K - 50 Test. Operating parameters for application of the detection system reagents on the BOND Processing Module have been optimized by Leica Biosystems. These can be displayed by following the instructions in your BOND user documentation.

Use rapid alcohol/xylene dehydration.

After completing or aborting *IHC Protocol K - 50 Test, there is a delay of approximately 15 minutes before the BOND system will start another run on the same Slide Staining Assembly.

Product Specific Limitations

ChromoPlex 1 Dual Detection has been optimized at Leica Biosystems for use with BOND ancillary reagents. Laboratories may use their own primary antibody cocktails provided they have been diluted to an appropriate concentration with BOND Primary Antibody Diluent (Catalog No. AR9352). Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances.

The appropriate concentration of user's own primary antibodies may vary due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and primary antibody concentrations.

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls. They should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Some mounting media are not compatible with ChromoPlex 1 Dual Detection (see **Precautions**).

Contact your local distributor or regional office of Leica Biosystems for further information.

Troubleshooting

Refer to reference 4 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

Date of Issue

24 September 2018

ChromoPlex™ 1 Dual Detection for BOND - 50 Test

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Bibliographie

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
 3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
- ProClin 950 est une marque commerciale de Supelco, membre du groupe Sigma-Aldrich Corporation.

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Bibliografia

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
 3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
- ProClin 950 è un marchio di fabbrica di Supelco, società del gruppo Sigma-Aldrich.

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Bibliografie

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
 3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
- ProClin 950 ist eine Marke von Supelco, einem Unternehmen der Sigma-Aldrich Corporation.

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1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 es una marca comercial de Supelco, parte de Sigma-Aldrich Corporation.

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Bibliografia

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 é uma marca comercial de Supelco, parte da Sigma-Aldrich Corporation.

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1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 är ett varumärke för Supelco som ingår i Sigma-Aldrich Corporation.

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1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. Το ProClin 950 είναι εμπορικό σήμα της Supelco, τμήμα της Sigma-Aldrich Corporation.

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Bibliografi

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 er et varemærke tilhørende Supelco, en del af Sigma-Aldrich Corporation.

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Bibliografie

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 is een handelsmerk van Supelco, een onderdeel van Sigma-Aldrich Corporation.

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Litteratur

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 er et varemerke for Supelco, en del av Sigma-Aldrich Corporation.

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Kaynakça

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950, Sigma-Aldrich Şirketinin bir parçası olan Supelco'nun bir ticari markasıdır.

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1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 е търговска марка на Supelco, част от Sigma-Aldrich Corporation.

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1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
 3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
- A ProClin 950 a Sigma-Aldrich Corporation részét képező Supelco védjegye.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 este o marcă comercială a Supelco, care face parte din Sigma-Aldrich Corporation.

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1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 — торговая марка компании Supelco, входящей в состав корпорации Sigma-Aldrich.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 jest znakiem handlowym firmy Supelco, będącej częścią Sigma-Aldrich Corporation.

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Literatura

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 je blagovna znamka družbe Supelco, ki je del družbe Sigma-Aldrich Corporation.

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Literatura

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996. ProClin 950 je ochranná známka společnosti Supelco, která je součástí Sigma-Aldrich Corporation.

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Literatúra

1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
 3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
- ProClin 950 je ochranná známka spoločnosti Supelco, súčasťou spoločnosti Sigma-Aldrich Corporation.

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1. Clinical laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. List of substances which may be candidates for further scientific review and possible identification, classification and regulation as potential occupational carcinogens. Fed Reg 1980; 45:157.
3. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
4. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
يمثل ProClin 950 علامة تجارية لشركة Supelco، وهي جزء من مؤسسة Sigma-Aldrich.

Leica Biosystems Newcastle Ltd
Balliol Business Park
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

Leica Biosystems Canada
71 Four Valley Drive
Concord, Ontario L4K 4V8
Canada
☎ +1 800 248 0123

Leica Biosystems Inc
1700 Leider Lane
Buffalo Grove IL 60089
USA
☎ +1 800 248 0123

Leica Biosystems Melbourne
Pty Ltd
495 Blackburn Road
Mt Waverley VIC 3149
Australia
☎ +61 2 8870 3500