Novocastra™ Lyophilized Mouse Monoclonal Antibody Dystrophin

Product Code: NCL-DYSB

Intended Use	FOR RESEARCH USE ONLY.
Specificity	Human dystrophin molecule
Clone	34C5
Ig Class	IgG1, kappa
Antigen Used for Immunizations	Prokaryotic recombinant protein corresponding to amino acids 321 to 494 of the dystrophin molecule.
Hybridoma Partner	Mouse myeloma (p3-NS1-Ag4-1).
Preparation	Lyophilized tissue culture supernatant containing 15 mM sodium azide. Reconstitute with the volume of sterile distilled water indicated on the vial label.
Effective on Frozen Tissue	Not fully evaluated.
Effective on Paraffin Wax Embedded Tissue	Yes (using 1 mM EDTA (pH 8.0) unmasking solution combined with the high temperature antigen unmasking technique: see overleaf).
Recommendations on Use	Immunohistochemistry: Typical working dilution 1:40–1:80. High temperature antigen unmasking technique using 1 mM EDTA (pH 8.0). 60 minutes primary antibody incubation at 25 °C. Standard ABC technique. Western Blotting: Not recommended.
Positive Controls	Immunohistochemistry: Skeletal muscle. Western Blotting: Not recommended.
Staining Pattern	Membrane.
Storage and Stability	Store unopened lyophilized antibody at 4 °C. Under these conditions, there is no significant loss in product performance up to the expiry date indicated on the vial label. The reconstituted antibody is stable for at least two months when stored at 4 °C. For long term storage, it is recommended that aliquots of the antibody are frozen at -20 °C (frost-free freezers are not recommended). Repeated freezing and thawing must be avoided. Prepare working dilutions on the day of use.

Warnings and Precautions

Novocastra Lyophilised Antibodies

Contains a mixture of:
Sodium
Azide (<10%),
Benzylpenicillin Sodium
(<10%), Streptomycin
Sulphate (<10%).
Signal words: Danger.

H302: Harmful if swallowed.
H317: May cause an allergic
skin reaction.
H334: May cause allergy
or asthma symptoms or
breathing difficulties if inhaled.
H411: Toxic to aquatic life with
long lasting effects.
EUH032: Contact with acids
liberates very toxic gas.

P261: Avoid breathing dust/fumes/gas/mist/vapours/spray.
P264: Wash hands thoroughly after handling.
P270: Do not eat, drink or smoke when using this product.
P272: Contaminated work clothing should not be allowed out of the workplace.
P273: Avoid release to the environment.
P280: Wear protective gloves/protective clothing/eye protection/face protection.
P284: In case of inadequate ventilation wear respiratory protection.
P301+312: IF SWALLOWED: Call a POISON CENTER/doctor if you feel unwell.
P302+352: IF ON SKIN: Wash with plenty of soap and water.
P304+340: IF INHALED: Remove person to fresh air and keep comfortable for
breathing.
P330: Rinse mouth.
P333+313: If skin irritation or rash occurs: Get medical advice/attention.
P342+311: If experiencing respiratory symptoms: Call a POISON CENTER/doctor.
P362+364: Take off contaminated clothing and wash it before reuse.
P391: Collect spillage.
P501: Dispose of contents/container to hazardous or special waste collection
point.

Novocastra Lyophilised Antibodies

Il contiene un mélange
de: Azoture De
Sodium (<10%),
Benzylpenicillin Sodium
(<10%), Streptomycin
Sulphate (<10%).
Mentions
d'avertissement:
Danger.

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**Novocastra
Lyophilised
Antibodies**

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General Overview

Dystrophin is a 427 kD cytoskeletal protein that is a member of the β spectrin/ α actinin family. It has a molecular structure that is typical of this family. Four separate domains can be identified on the basis of sequence homologs and protein-binding capabilities. These are the N-terminal actin binding domain, the central rod domain, the cysteine-rich domain and the COOH-terminal coiled-coil domain. Dystrophin is the product of the Duchenne muscular dystrophy and Becker muscular dystrophy (DMD/BMD) gene which is located on the X-chromosome at position Xp21. Several different isoforms and splice variants of dystrophin exist including protein products of 260 kD, 140 kD, 116 kD and 71 kD.

Important: For reliable interpretation of dystrophin labelling patterns using tissue sections, the use of a SPECTRIN control is essential.

General References

Blake D J, Weir A, Newey S E, et al.. *Physiological Reviews*. 82 (2): 291–329 (2002).
Oliveira A S, Gabbai A A, Schmidt B, et al.. *Arq Neuropsiquiatr*. 50 (4): 478–485 (1992).
Haginoya K, Yamamoto K, Iinuma K, et al.. *Journal of Neurology*. 238 (7): 375–378 (1991).

Instructions for Use

High Temperature Antigen Unmasking Technique for Immunohistochemical Demonstration on Paraffin Sections

1. Cut and mount sections on slides coated with a suitable tissue adhesive.
2. Deparaffinize sections and rehydrate to distilled water.
3. Place sections in 0.5% hydrogen peroxide/methanol for 10 minutes (or use other appropriate endogenous peroxidase blocking procedure). Wash sections in tap water.
4. Heat 1500 mL of the recommended unmasking solution (0.01 M citrate buffer, pH 6.0 (or Epitope Retrieval Solution, RE7113) unless otherwise indicated overleaf) until boiling in a stainless steel pressure cooker. Cover but do not lock lid.
5. Position slides into metal staining racks (do not place slides close together as uneven staining may occur) and lower into pressure cooker ensuring slides are completely immersed in unmasking solution. Lock lid.
6. When the pressure cooker reaches operating temperature and pressure (after about 5 minutes) start a timer for 1 minute (unless otherwise indicated on the data sheet).
7. When the timer rings, remove pressure cooker from heat source and run under cold water with lid on. DO NOT OPEN LID UNTIL THE INDICATORS SHOW THAT PRESSURE HAS BEEN RELEASED. Open lid, remove slides and place immediately into a bath of tap water.
8. Wash sections in TBS* buffer (pH 7.6) for 1 x 5 minutes.
9. Place sections in diluted normal serum (or RTU Normal Horse Serum) for 10 minutes.
10. Incubate sections with primary antibody. Use Antibody Diluent RE7133 (where available).
11. Wash in TBS buffer for 2 x 5 minutes.
12. Incubate sections in an appropriate biotinylated secondary antibody.
13. Wash in TBS buffer for 2 x 5 minutes.
14. Incubate slides in ABC reagent (or RTU streptavidin/peroxidase complex).
15. Wash in TBS buffer for 2 x 5 minutes.
16. Incubate slides in DAB or other suitable peroxidase substrate.
17. Wash thoroughly in running tap water.
18. Counterstain with hematoxylin (if required), dehydrate and mount.

Solutions

0.01 M CITRATE BUFFER (pH 6.0) or RE7113 (where available).

Add 3.84 g of citric acid (anhydrous) to 1.8 L of distilled water. Adjust to pH 6.0 using concentrated NaOH. Make up to 2 L with distilled water.

1 mM EDTA (pH 8.0) or RE7116 (where available).

Add 0.37 g of EDTA (SIGMA product code E-5134) to 1 litre of distilled water. Adjust pH to 8.0 using 1.0 M NaOH.

20 mM TRIS/ 0.65 mM EDTA/ 0.005% TWEEN (pH 9.0) or RE7119 (where available).

Dissolve 14.4 g Tris (BDH product code 271197K) and 1.44 g EDTA (SIGMA product code E-5134) to 0.55 L of distilled water. Adjust pH to 9.0 with 1 M HCl and add 0.3 mL Tween 20 (SIGMA product code P-1379). Make up to 0.6 L with distilled water. This is a 10x concentrate which should be diluted with distilled water as required (eg 150 mL diluted with 1350 mL of distilled water).

* In most applications, 10 mM phosphate, 0.15 M NaCl, pH 7.6 (PBS) can be used instead of 50 mM Tris, 0.15 M NaCl, pH 7.6 (TBS).

Safety Note

To ensure the correct and safe use of your pressure cooker, PLEASE READ MANUFACTURER'S INSTRUCTIONS.

protocol/HTAUT/07/08 DYSB 03/05/2016

© Leica Biosystems Newcastle Ltd · 95.8147 Rev A