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Catalog No: PA0552

Intended Use

This reagent is for in vitro diagnostic use.

Pax-5 (1EW) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human Pax-5 in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Pax-5 (1EW) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human Pax-5 is achieved by first, allowing the binding of Pax-5 (1EW) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

Pax-5 (1EW) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin® 950 as a preservative.

Total volume = 7 mL.

Clone

1EW.

Immunogen

Recombinant protein corresponding to part of the C-terminus of the human Pax-5 molecule.

Specificity

Human Pax-5.

Subclass

IgG1.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 0.34 mg/L as determined by ELISA.

Dilution and Mixing

Pax-5 (1EW) primary antibody is optimally diluted for use on the BOND system. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system.

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Pax-5 (1EW) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin® 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Pax-5 (1EW) primary antibody was developed for use on the automated BOND system in combination with BOND Polymer Refine Detection. The use of Peroxide Block as the first step in IHC Staining Protocol F will affect the staining achieved with this antibody, for this reason the recommended staining protocol for Pax-5 (1EW) is an edited version of IHC Staining Protocol F, where the Peroxide Block step before the MARKER has been deleted and a peroxide block step has been inserted after the MARKER. This modified staining protocol needs to be created by the user. For instructions on how to edit protocols Refer to "Adding and removing protocol steps" in your BOND user documentation. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 2 for 20 minutes.

Results Expected

Normal Tissues

Pax-5 is a B cell specific nuclear transcription factor that is expressed in pro B, pre B and mature B cells. Clone 1EW detected the Pax-5 protein in follicular and non-follicular B cells of tonsil 3/3, spleen 3/3, colon 3/3 and small bowel 3/3. Staining of occasional B lymphocytes was also seen in a variety other tissues 33/69, including adrenal, kidney, esophagus, stomach, liver, pancreas, lung, skin, prostate, cervix, breast, thyroid, thymus, bone marrow, and pituitary.

Tumor Tissues

Clone 1EW stained 14/14 diffuse large B cell lymphomas, 6/6 follicle center lymphomas, 3/3 MALTomas, 1/1 mantle cell lymphomas, 1/1 Burkitt's lymphoma, 1/1 T cell rich B cell lymphoma, 5/5 Hodgkin's disease and 0/11 T cell lymphomas. No staining was observed in a variety of non-hematological malignancies (n=41).

Pax-5 (1EW) is recommended for the detection of human Pax-5 protein in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.

Product Specific Limitations

Pax-5 (1EW) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Hansson M, Jerkeman M and Dictor M. Biphentotypic bigenotypic lymphoma with simultaneous expression of PAX5/BSAP and B- and T-cell markers. European Journal of Haematology. 2007; 79:159–165.

ProClin: 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Hansson M, Jerkeman M and Dictor M. Biphentypic bigenotypic lymphoma with simultaneous expression of PAX5/BSAP and B- and T-cell markers. European Journal of Haematology. 2007; 79:159-165.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Hansson M, Jerkeman M and Dictor M. Biphentotypic bigenotypic lymphoma with simultaneous expression of PAX5/BSAP and B- and T-cell markers. European Journal of Haematology. 2007; 79:159-165.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD och Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Hansson M, Jerkeman M and Dictor M. Biphenotypic bigenotypic lymphoma with simultaneous expression of PAX5/BSAP and B- and T-cell markers. European Journal of Haematology. 2007; 79:159–165.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Hansson M, Jerkeman M and Dictor M. Biphenotypic bigenotypic lymphoma with simultaneous expression of PAX5/BSAP and B- and T-cell markers. European Journal of Haematology. 2007; 79:159–165.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Hansson M, Jerkeman M and Dictor M. Biphenotypic bigenotypic lymphoma with simultaneous expression of PAX5/BSAP and B- and T-cell markers. European Journal of Haematology. 2007; 79: 159–165.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Hansson M, Jerkeman M and Dictor M. Biphentotypic biogenotypic lymphoma with simultaneous expression of PAX5/BSAP and B- and T-cell markers. European Journal of Haematology. 2007; 79:159–165.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd