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Catalog No: PA0799

Intended Use

This reagent is for in vitro diagnostic use.

Cytokeratin 19 (b170) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of human cytokeratin 19 intermediate filament protein in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Cytokeratin 19 (b170) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of human cytokeratin 19 intermediate filament protein is achieved by first, allowing the binding of Cytokeratin 19 (b170) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system, reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

Reagents Provided

Cytokeratin 19 (b170) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35% ProClin™ 950 as a preservative.

Total volume = 7 mL.

Clone

b170.

Immunogen

Cytoskeletal preparation from human MCF-7 cells.

Specificity

Human cytokeratin 19 intermediate filament protein.

Subclass

IgG1.

Total Protein Concentration

Approx 10 mg/mL.

Antibody Concentration

Greater than or equal to 0.07 mg/L as determined by ELISA.

Dilution and Mixing

Cytokeratin 19 (b170) primary antibody is optimally diluted for use on the BOND system. Reconstitution, mixing, dilution or titration of this reagent is not required.

Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system.

Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Cytokeratin 19 (b170) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user¹.

Precautions

- This product is intended for in vitro diagnostic use.
- The concentration of ProClin™ 950 is 0.35%. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com.

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions². Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Instructions for Use

Cytokeratin 19 (b170) primary antibody was developed for use on the automated BOND system in combination with BOND Polymer Refine Detection. The recommended staining protocol for Cytokeratin 19 (b170) primary antibody is IHC Protocol F. Enzyme pretreatment is recommended using BOND Enzyme 1 for 10 minutes.

Results Expected

Normal Tissues

Clone b170 stained epithelium of the gastrointestinal tract, alveolar lining cells in lung, mesothelium, ductal and acinar cells of salivary gland, bile duct epithelium in liver, ductal and acinar epithelium of pancreas, epithelium of Bowman's capsule and tubules in kidney, basal and secretory epithelium of prostate, follicular epithelium of the thyroid gland, thymic epithelium and Hassall's corpuscles, eccrine sweat glands and ducts in skin, endometrial glandular epithelium, ductal epithelia in breast, basal cells of non-keratinizing stratified epithelia of tonsil and cervix, and epithelial cells in parathyroid, pituitary and ovary (Total number of cases stained = 62).

Tumor Tissues

Clone b170 stained 37/39 ductal carcinomas of breast, squamous cell carcinoma of esophagus 2/2, larynx 2/2, lung 2/2 and uterus 2/2, adenocarcinomas of stomach 2/2, pancreas 2/2, prostate 2/2, colon 2/2, rectum 1/1 and lung 1/1. Also stained were 2/3 thyroid carcinomas, 1/1 non small cell carcinoma of lung, 1/1 lymph node metastatic carcinoma, 2/2 cholangiocarcinomas, 0/1 hepatocellular carcinoma (HCC), 1/1 HCC + cholangiocarcinoma, 1/2 renal cell carcinoma, 1/1 clear cell carcinoma of ovary, 1/1 mucinous cystadenocarcinoma of ovary and 2/2 transitional cell carcinoma of ovary. No staining was seen on lymphoma 0/2, maltoma 0/1, glioblastoma 0/2 or squamous cell carcinoma of skin 0/2 (Total number of cases stained = 79).

Cytokeratin 19 (b170) is recommended for use as part of an antibody panel for the characterization of normal and malignant cells of epithelial origin.

Product Specific Limitations

Cytokeratin 19 (b170) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kurokawa I, Nishimura K, Hakamada A, et al. Rippled-pattern sebaceoma with an immunohistochemical study of cytokeratins. Letters to the Editor. Journal of the European Academy of Dermatology and Venereology. 2007; 21:133–134.
5. Saggiorato E, De Pompa R, Volante M, et al. Characterization of thyroid 'follicular neoplasms' in one-needle aspiration cytological specimens using a panel of immunohistochemical markers: a proposal for clinical application. Endocrine-Related Cancer. 2005; 12:305–317.
6. Lewindon PJ, Pereira TN, Hoskins AC, et al. The role of hepatic stellate cells and transforming growth factor-β1 in cystic fibrosis liver disease. American Journal of Pathology. 2002;160(5):1705–1715.
7. Cheung CC, Ezzat S, Freeman JL, et al. Immunohistochemical diagnosis of papillary thyroid carcinoma. Modern Pathology. 2001;14(4):338–342.
8. Fujita M, Furukawa H, Hattori M, et al. Sequential observation of liver cell regeneration after massive hepatic necrosis in auxiliary partial orthotopic liver transplantation. Modern Pathology 2000;13(2):152–157.

ProClin[®] 950 is a trademark of Supelco, a part of Sigma-Aldrich Corporation.

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Bibliographie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kurokawa I, Nishimura K, Hakamada A, et al. Rippled-pattern sebaceoma with an immunohistochemical study of cytokeratins. Letters to the Editor. Journal of the European Academy of Dermatology and Venereology. 2007; 21:133-134.
5. Saggiorato E, De Pompa R, Volante M, et al. Characterization of thyroid 'follicular neoplasms' in one-needle aspiration cytological specimens using a panel of immunohistochemical markers: a proposal for clinical application. Endocrine-Related Cancer. 2005; 12:305-317.
6. Lewindon PJ, Pereira TN, Hoskins AC, et al. The role of hepatic stellate cells and transforming growth factor-β1 in cystic fibrosis liver disease. American Journal of Pathology. 2002;160(5):1705-1715.
7. Cheung CC, Ezzat S, Freeman JL, et al. Immunohistochemical diagnosis of papillary thyroid carcinoma. Modern Pathology. 2001;14(4):338-342.
8. Fujita M, Furukawa H, Hattori M, et al. Sequential observation of liver cell regeneration after massive hepatic necrosis in auxiliary partial orthotopic liver transplantation. Modern Pathology 2000;13(2):152-157.

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Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kurokawa I, Nishimura K, Hakamada A, et al. Rippled-pattern sebaceoma with an immunohistochemical study of cytokeratins. Letters to the Editor. Journal of the European Academy of Dermatology and Venereology. 2007; 21:133–134.
5. Saggiolato E, De Pomba R, Volante M, et al. Characterization of thyroid 'follicular neoplasms' in one-needle aspiration cytological specimens using a panel of immunohistochemical markers: a proposal for clinical application. Endocrine-Related Cancer. 2005; 12:305–317.
6. Lewindon PJ, Pereira TN , Hoskins AC, et al. The role of hepatic stellate cells and transforming growth factor-β1 in cystic fibrosis liver disease. American Journal of Pathology. 2002;160(5):1705–1715.
7. Cheung CC, Ezzat S, Freeman JL, et al. Immunohistochemical diagnosis of papillary thyroid carcinoma. Modern Pathology. 2001;14(4):338–342.
8. Fujita M, Furukawa H, Hattori M, et al. Sequential observation of liver cell regeneration after massive hepatic necrosis in auxiliary partial orthotopic liver transplantation. Modern Pathology 2000;13(2):152–157.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kurokawa I, Nishimura K, Hakamada A, et al. Rippled-pattern sebaceoma with an immunohistochemical study of cytokeratins. Letters to the Editor. Journal of the European Academy of Dermatology and Venereology. 2007; 21:133–134.
5. Saggiolato E, De Pompa R, Volante M, et al. Characterization of thyroid "follicular neoplasms" in one-needle aspiration cytological specimens using a panel of immunohistochemical markers: a proposal for clinical application. Endocrine-Related Cancer. 2005; 12:305–317.
6. Lewindon PJ, Pereira TN, Hoskins AC, et al. The role of hepatic stellate cells and transforming growth factor- β 1 in cystic fibrosis liver disease. American Journal of Pathology. 2002;160(5):1705–1715.
7. Cheung CC, Ezzat S, Freeman JL, et al. Immunohistochemical diagnosis of papillary thyroid carcinoma. Modern Pathology. 2001;14(4):338–342.
8. Fujita M, Furukawa H, Hattori M, et al. Sequential observation of liver cell regeneration after massive hepatic necrosis in auxiliary partial orthotopic liver transplantation. Modern Pathology 2000;13(2):152–157.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kurokawa I, Nishimura K, Hakamada A, et al. Rippled-pattern sebaceoma with an immunohistochemical study of cytokeratins. Letters to the Editor. Journal of the European Academy of Dermatology and Venereology. 2007; 21:133–134.
5. Saggiolato E, De Pompa R, Volante M, et al. Characterization of thyroid 'follicular neoplasms' in one-needle aspiration cytological specimens using a panel of immunohistochemical markers: a proposal for clinical application. Endocrine-Related Cancer. 2005; 12:305–317.
6. Lewindon PJ, Pereira TN, Hoskins AC, et al. The role of hepatic stellate cells and transforming growth factor-β1 in cystic fibrosis liver disease. American Journal of Pathology. 2002;160(5):1705–1715.
7. Cheung CC, Ezzat S, Freeman JL, et al. Immunohistochemical diagnosis of papillary thyroid carcinoma. Modern Pathology. 2001;14(4):338–342.
8. Fujita M, Furukawa H, Hattori M, et al. Sequential observation of liver cell regeneration after massive hepatic necrosis in auxiliary partial orthotopic liver transplantation. Modern Pathology 2000;13(2):152–157.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kurokawa I, Nishimura K, Hakamada A, et al. Rippled-pattern sebaceoma with an immunohistochemical study of cytokeratins. Letters to the Editor. Journal of the European Academy of Dermatology and Venereology. 2007; 21:133–134.
5. Saggiarato E, De Pompa R, Volante M, et al. Characterization of thyroid 'follicular neoplasms' in one-needle aspiration cytological specimens using a panel of immunohistochemical markers: a proposal for clinical application. Endocrine-Related Cancer. 2005; 12:305–317.

6. Lewindon PJ, Pereira TN, Hoskins AC, et al. The role of hepatic stellate cells and transforming growth factor- β 1 in cystic fibrosis liver disease. *American Journal of Pathology*. 2002;160(5):1705–1715.
 7. Cheung CC, Ezzat S, Freeman JL, et al. Immunohistochemical diagnosis of papillary thyroid carcinoma. *Modern Pathology*. 2001;14(4):338–342.
 8. Fujita M, Furukawa H, Hattori M, et al. Sequential observation of liver cell regeneration after massive hepatic necrosis in auxiliary partial orthotopic liver transplantation. *Modern Pathology* 2000;13(2):152–157.
- To ProClin[®] 950 είναι εμπορικό σήμα της Supelco, τμήμα της Sigma-Aldrich Corporation.

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Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kurokawa I, Nishimura K, Hakamada A, et al. Rippled-pattern sebaceoma with an immunohistochemical study of cytokeratins. Letters to the Editor. Journal of the European Academy of Dermatology and Venerology. 2007; 21:133-134.
5. Saggiorato E, De Pomba R, Volante M, et al. Characterization of thyroid 'follicular neoplasms' in one-needle aspiration cytological specimens using a panel of immunohistochemical markers: a proposal for clinical application. Endocrine-Related Cancer. 2005; 12:305-317.
6. Lewindon PJ, Pereira TN, Hoskins AC, et al. The role of hepatic stellate cells and transforming growth factor-β1 in cystic fibrosis liver disease. American Journal of Pathology. 2002;160(5):1705-1715.
7. Cheung CC, Ezzat S, Freeman JL, et al. Immunohistochemical diagnosis of papillary thyroid carcinoma. Modern Pathology. 2001;14(4):338-342.
8. Fujita M, Furukawa H, Hattori M, et al. Sequential observation of liver cell regeneration after massive hepatic necrosis in auxiliary partial orthotopic liver transplantation. Modern Pathology 2000;13(2):152-157.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
 2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
 3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
 4. Kurokawa I, Nishimura K, Hakamada A, et al. Rippled-pattern sebaceoma with an immunohistochemical study of cytokeratins. Letters to the Editor. Journal of the European Academy of Dermatology and Venereology. 2007; 21:133–134.
 5. Saggiorato E, De Pompa R, Volante M, et al. Characterization of thyroid 'follicular neoplasms' in one-needle aspiration cytological specimens using a panel of immunohistochemical markers: a proposal for clinical application. Endocrine-Related Cancer. 2005; 12:305–317.
 6. Lewindon PJ, Pereira TN, Hoskins AC, et al. The role of hepatic stellate cells and transforming growth factor-β1 in cystic fibrosis liver disease. American Journal of Pathology. 2002;160(5):1705–1715.
 7. Cheung CC, Ezzat S, Freeman JL, et al. Immunohistochemical diagnosis of papillary thyroid carcinoma. Modern Pathology. 2001;14(4):338–342.
 8. Fujita M, Furukawa H, Hattori M, et al. Sequential observation of liver cell regeneration after massive hepatic necrosis in auxiliary partial orthotopic liver transplantation. Modern Pathology 2000;13(2):152–157.
- A ProClin™ 950 a Sigma-Aldrich Corporation részét képező Supelco védjegye.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kurokawa I, Nishimura K, Hakamada A, et al. Rippled-pattern sebaceous with an immunohistochemical study of cytokeratins. Letters to the Editor. Journal of the European Academy of Dermatology and Venereology. 2007; 21:133-134.
5. Saggiolato E, De Pompa R, Volante M, et al. Characterization of thyroid 'follicular neoplasms' in one-needle aspiration cytological specimens using a panel of immunohistochemical markers: a proposal for clinical application. Endocrine-Related Cancer. 2005; 12:305-317.

6. Lewindon PJ, Pereira TN, Hoskins AC, et al. The role of hepatic stellate cells and transforming growth factor- β 1 in cystic fibrosis liver disease. *American Journal of Pathology*. 2002;160(5):1705–1715.
 7. Cheung CC, Ezzat S, Freeman JL, et al. Immunohistochemical diagnosis of papillary thyroid carcinoma. *Modern Pathology*. 2001;14(4):338–342.
 8. Fujita M, Furukawa H, Hattori M, et al. Sequential observation of liver cell regeneration after massive hepatic necrosis in auxiliary partial orthotopic liver transplantation. *Modern Pathology* 2000;13(2):152–157.
- ProClin® 950 — торговая марка компании Supelco, входящей в состав корпорации Sigma-Aldrich.

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Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kurokawa I, Nishimura K, Hakamada A, et al. Rippled-pattern sebaceoma with an immunohistochemical study of cytokeratins. Letters to the Editor. Journal of the European Academy of Dermatology and Venereology. 2007; 21:133–134.
5. Saggiatoro E, De Pompa R, Volante M, et al. Characterization of thyroid 'follicular neoplasms' in one-needle aspiration cytological specimens using a panel of immunohistochemical markers: a proposal for clinical application. Endocrine-Related Cancer. 2005; 12:305–317.
6. Lewindon PJ, Pereira TN, Hoskins AC, et al. The role of hepatic stellate cells and transforming growth factor-β1 in cystic fibrosis liver disease. American Journal of Pathology. 2002;160(5):1705–1715.
7. Cheung CC, Ezzat S, Freeman JL, et al. Immunohistochemical diagnosis of papillary thyroid carcinoma. Modern Pathology. 2001;14(4):338–342.
8. Fujita M, Furukawa H, Hattori M, et al. Sequential observation of liver cell regeneration after massive hepatic necrosis in auxiliary partial orthotopic liver transplantation. Modern Pathology 2000;13(2):152–157.

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Kurokawa I, Nishimura K, Hakamada A, et al. Rippled-pattern sebaceoma with an immunohistochemical study of cytokeratins. Letters to the Editor. Journal of the European Academy of Dermatology and Venereology. 2007; 21:133-134.
5. Saggiorato E, De Pompa R, Volante M, et al. Characterization of thyroid "follicular neoplasms" in one-needle aspiration cytological specimens using a panel of immunohistochemical markers: a proposal for clinical application. Endocrine-Related Cancer. 2005; 12:305-317.
6. Lewindon PJ, Pereira TN, Hoskins AC, et al. The role of hepatic stellate cells and transforming growth factor-β1 in cystic fibrosis liver disease. American Journal of Pathology. 2002;160(5):1705-1715.
7. Cheung CC, Ezzat S, Freeman JL, et al. Immunohistochemical diagnosis of papillary thyroid carcinoma. Modern Pathology. 2001;14(4):338-342.
8. Fujita M, Furukawa H, Hattori M, et al. Sequential observation of liver cell regeneration after massive hepatic necrosis in auxiliary partial orthotopic liver transplantation. Modern Pathology 2000;13(2):152-157.

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Leica Biosystems Newcastle Ltd
Leica Biosystems Inc
Pty Ltd