BIOSYSTEMS

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### Catalog No: PA0831

#### Intended Use

This reagent is for *in vitro* diagnostic use.

Anaplastic Lymphoma Kinase (5A4) monoclonal antibody is intended to be used for the qualitative identification by light microscopy of endogenous ALK and ALK fusion proteins in formalin-fixed, paraffin-embedded tissue by immunohistochemical staining using the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

The clinical interpretation of any staining or its absence should be complemented by morphological studies and proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

#### Summary and Explanation

Immunohistochemical techniques can be used to demonstrate the presence of antigens in tissue and cells (see "Using BOND Reagents" in your BOND user documentation). Anaplastic Lymphoma Kinase (5A4) primary antibody is a ready to use product that has been specifically optimized for use with BOND Polymer Refine Detection. The demonstration of endogenous ALK and ALK fusion proteins is achieved by first allowing the binding of Anaplastic Lymphoma Kinase (5A4) to the section, and then visualizing this binding using the reagents provided in the detection system. The use of these products, in combination with the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system), reduces the possibility of human error and inherent variability resulting from individual reagent dilution, manual pipetting and reagent application.

#### Reagents Provided

Anaplastic Lymphoma Kinase (5A4) is a mouse anti-human monoclonal antibody produced as a tissue culture supernatant, and supplied in Tris buffered saline with carrier protein, containing 0.35 % ProClin™ 950 as a preservative.

Total volume = 7 mL.

#### Clone

5A4

#### Immunogen

Prokaryotic recombinant protein corresponding to a region which spans the tyrosine kinase catalytic domain and part of the C-terminus of the NPM-ALK transcript (419-520 amino acids).

#### Specificity

Human ALK protein and C-terminal ALK fusion proteins containing the immunized target sequence.

#### Ig Class

IgG1

#### Total Protein Concentration

Approx 10 mg/mL.

#### Antibody Concentration

Greater than or equal to 2.6 mg/L as determined by ELISA.

#### Dilution and Mixing

Anaplastic Lymphoma Kinase (5A4) primary antibody is optimally diluted for use on the BOND system (includes Leica BOND-MAX system and Leica BOND-III system). Reconstitution, mixing, dilution or titration of this reagent is not required.

#### Materials Required But Not Provided

Refer to "Using BOND Reagents" in your BOND user documentation for a complete list of materials required for specimen treatment and immunohistochemical staining using the BOND system (includes Leica BOND-MAX system and Leica BOND-III system).

#### Storage and Stability

Store at 2–8 °C. Do not use after the expiration date indicated on the container label.

The signs indicating contamination and/or instability of Anaplastic Lymphoma Kinase (5A4) are: turbidity of the solution, odor development, and presence of precipitate.

Return to 2–8 °C immediately after use.

Storage conditions other than those specified above must be verified by the user<sup>1</sup>.

#### Precautions

- This product is intended for *in vitro* diagnostic use.
- The concentration of ProClin™ 950 is 0.35 %. It contains the active ingredient 2-methyl-4-isothiazolin-3-one, and may cause irritation to the skin, eyes, mucous membranes and upper respiratory tract. Wear disposable gloves when handling reagents.
- To obtain a copy of the Material Safety Data Sheet contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, [www.LeicaBiosystems.com](http://www.LeicaBiosystems.com)

- Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions<sup>2</sup>. Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents or specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.
- Consult Federal, State or local regulations for disposal of any potentially toxic components.
- Minimize microbial contamination of reagents or an increase in non-specific staining may occur.
- Retrieval, incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

## Instructions for Use

Anaplastic Lymphoma Kinase (5A4) primary antibody was developed for use on the automated BOND system (includes Leica BOND-MAX system and Leica BOND-III system) in combination with BOND Polymer Refine Detection. The recommended staining protocol for Anaplastic Lymphoma Kinase (5A4) primary antibody is IHC Protocol F. Heat induced epitope retrieval is recommended using BOND Epitope Retrieval Solution 2 for 30 minutes.

## Results Expected

### Normal Tissues

Clone 5A4 detected Anaplastic Lymphoma Kinase protein in some neuroendocrine cells of the pituitary gland and small intestine, alveolar macrophages in lung and a proportion of prostatic stromal cells. All other tissues tested were negative. (Total number of normal cases evaluated = 130).

### Tumor Tissues

Clone 5A4 stained 5/27 anaplastic large cell lymphomas and 2/5 large cell carcinomas of the lung. No staining was observed in a variety of additional abnormal tissues evaluated, including lung tumours (0/68), breast tumours (0/5), thyroid tumours (0/5), tumours of the colon (0/4), liver tumours (0/4), metastatic tumours (0/4), tumours of the oesophagus (0/3), stomach (0/3), and head and neck (0/3), ovarian tumours (0/3), tumours of the adrenal gland (0/2), bladder tumours (0/2), brain tumours (0/2), tumours of the small intestine (0/2), rectal tumours (0/2), kidney tumours (0/2), tumours of the lymph node (0/2), salivary gland (0/2), testis (0/2), cervix (0/2) and endometrium (0/2), an osteosarcoma (0/1), a pancreatic adenocarcinoma (0/1), an adenocarcinoma of the prostate (0/1), a hyperplastic prostate, (0/1), a squamous cell carcinoma of the skin (0/1) and a melanoma (0/1). (Total number of abnormal tissues evaluated = 162).

**Anaplastic Lymphoma Kinase (5A4) is recommended to aid in the evaluation and expression of endogenous ALK and ALK fusion proteins, in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.**

## Product Specific Limitations

Anaplastic Lymphoma Kinase (5A4) has been optimized at Leica Biosystems for use with BOND Polymer Refine Detection and BOND ancillary reagents. Users who deviate from recommended test procedures must accept responsibility for interpretation of patient results under these circumstances. The protocol times may vary, due to variation in tissue fixation and the effectiveness of antigen enhancement, and must be determined empirically. Negative reagent controls should be used when optimizing retrieval conditions and protocol times.

## Troubleshooting

Refer to reference 3 for remedial action.

Contact your local distributor or the regional office of Leica Biosystems to report unusual staining.

## Further Information

Further information on immunostaining with BOND reagents, under the headings Principle of the Procedure, Materials Required, Specimen Preparation, Quality Control, Assay Verification, Interpretation of Staining, Key to Symbols on Labels, and General Limitations can be found in "Using BOND Reagents" in your BOND user documentation.

## Bibliography

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Li XQ, Hisaoka M, Shi DR et al. Expression of anaplastic lymphoma kinase in soft tissue tumors: an immunohistochemical and molecular study of 249 cases. Human Pathology. 2004; 35:711-721
5. Savic S, Bode B, Diebold J et al. Detection of ALK-positive non-small-cell lung cancers on cytological specimens: high accuracy of immunocytochemistry with the 5A4 clone. Journal of Thoracic Oncology. 2013; 8(8):1004-11
6. Savic S, Diebold J, Zimmermann AK et al. Screening for ALK in non-small cell lung carcinomas: 5A4 and D5F3 antibodies perform equally well, but combined use with FISH is recommended. Lung Cancer. 2015; 89 (2):104-9
7. Williams AS, Greer W, Bethune D et al. ALK+ lung adenocarcinoma in never smokers and long-term ex-smokers: prevalence and detection by immunohistochemistry and fluorescence in situ hybridisation. Virchows Archiv: an international journal of pathology. 2016; 469(5):533-540
8. Barreca A, Lasorsa E, Riera L et al. Anaplastic lymphoma kinase in human cancer. Journal of Molecular Endocrinology. 2011; 47(1):R11-R23
9. Zhao Z, Verma V. and Zhang M. Anaplastic lymphoma kinase: Role in cancer and therapy perspective. Cancer Biology and Therapy. 2015; 16(12):1691-701

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## Bibliografie

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 1163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Li XQ, Hisaoka M, Shi DR et al. Expression of anaplastic lymphoma kinase in soft tissue tumors: an immunohistochemical and molecular study of 249 cases. Human Pathology. 2004; 35:711-721
5. Savic S, Bode B, Diebold J et al. Detection of ALK-positive non-small-cell lung cancers on cytological specimens: high accuracy of immunocytochemistry with the 5A4 clone. Journal of Thoracic Oncology. 2013; 8(8):1004-11
6. Savic S, Diebold J, Zimmermann AK et al. Screening for ALK in non-small cell lung carcinomas: 5A4 and D5F3 antibodies perform equally well, but combined use with FISH is recommended. Lung Cancer. 2015; 89 (2):104-9
7. Williams AS, Greer W, Bethune D et al. ALK+ lung adenocarcinoma in never smokers and long-term ex-smokers: prevalence and detection by immunohistochemistry and fluorescence in situ hybridisation. Virchows Archiv: an international journal of pathology. 2016; 469(5):533-540
8. Barreca A, Lasorsa E, Riera L et al. Anaplastic lymphoma kinase in human cancer. Journal of Molecular Endocrinology. 2011; 47(1):R11-R23
9. Zhao Z, Verma V. and Zhang M. Anaplastic lymphoma kinase: Role in cancer and therapy perspective. Cancer Biology and Therapy. 2015; 16(12):1691-701

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## Bibliografi

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Li XQ, Hisaoka M, Shi DR et al. Expression of anaplastic lymphoma kinase in soft tissue tumors: an immunohistochemical and molecular study of 249 cases. Human Pathology. 2004; 35:711-721
5. Savic S, Bode B, Diebold J et al. Detection of ALK-positive non-small-cell lung cancers on cytological specimens: high accuracy of immunocytochemistry with the 5A4 clone. Journal of Thoracic Oncology. 2013; 8(8):1004-11
6. Savic S, Diebold J, Zimmermann AK et al. Screening for ALK in non-small cell lung carcinomas: 5A4 and D5F3 antibodies perform equally well, but combined use with FISH is recommended. Lung Cancer. 2015; 89 (2):104-9
7. Williams AS, Greer W, Bethune D et al. ALK+ Lung adenocarcinoma in never smokers and long-term ex-smokers: prevalence and detection by immunohistochemistry and fluorescence in situ hybridisation. Virchows Archiv: an international journal of pathology. 2016; 469(5):533-540
8. Barreca A, Lasorsa E, Riera L et al. Anaplastic lymphoma kinase in human cancer. Journal of Molecular Endocrinology. 2011; 47(1):R11-R23
9. Zhao Z, Verma V. and Zhang M. Anaplastic lymphoma kinase: Role in cancer and therapy perspective. Cancer Biology and Therapy. 2015; 16(12):1691-701

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## Kaynakça

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Li XQ, Haisaka M, Shi DR et al. Expression of anaplastic lymphoma kinase in soft tissue tumors: an immunohistochemical and molecular study of 249 cases. Human Pathology. 2004; 35:711-721
5. Savic S, Bode B, Diebold J et al. Detection of ALK-positive non-small-cell lung cancers on cytological specimens: high accuracy of immunocytochemistry with the 5A4 clone. Journal of Thoracic Oncology. 2013; 8(8):1004-11
6. Savic S, Diebold J, Zimmermann AK et al. Screening for ALK in non-small cell lung carcinomas: 5A4 and D5F3 antibodies perform equally well, but combined use with FISH is recommended. Lung Cancer. 2015; 89 (2):104-9
7. Williams AS, Greer W, Bethune D et al. ALK+ lung adenocarcinoma in never smokers and long-term ex-smokers: prevalence and detection by immunohistochemistry and fluorescence in situ hybridisation. Virchows Archiv: an international journal of pathology. 2016; 469(5):533-540
8. Barreca A, Lasorsa E, Riera L et al. Anaplastic lymphoma kinase in human cancer. Journal of Molecular Endocrinology. 2011; 47(1):R11-R23
9. Zhao Z, Verma V. and Zhang M. Anaplastic lymphoma kinase: Role in cancer and therapy perspective. Cancer Biology and Therapy. 2015; 16(12):1691-701

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1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Li XQ, Hisoaka M, Shi DR et al. Expression of anaplastic lymphoma kinase in soft tissue tumors: an immunohistochemical and molecular study of 249 cases. Human Pathology. 2004; 35:711-721

5. Savic S, Bode B, Diebold J et al. Detection of ALK-positive non-small-cell lung cancers on cytological specimens: high accuracy of immunocytochemistry with the 5A4 clone. *Journal of Thoracic Oncology*. 2013; 8(8):1004-11
6. Savic S, Diebold J, Zimmermann AK et al. Screening for ALK in non-small cell lung carcinomas: 5A4 and D5F3 antibodies perform equally well, but combined use with FISH is recommended. *Lung Cancer*. 2015; 89 (2):104-9
7. Williams AS, Greer W, Bethune D et al. ALK+ lung adenocarcinoma in never smokers and long-term ex-smokers: prevalence and detection by immunohistochemistry and fluorescence in situ hybridisation. *Virchows Archiv: an international journal of pathology*. 2016; 469(5):533-540
8. Barreca A, Lasorsa E, Riera L et al. Anaplastic lymphoma kinase in human cancer. *Journal of Molecular Endocrinology*. 2011; 47(1):R11-R23
9. Zhao Z, Verma V. and Zhang M. Anaplastic lymphoma kinase: Role in cancer and therapy perspective. *Cancer Biology and Therapy*. 2015; 16(12):1691-701

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## Bibliografia

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 7163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Li XQ, Hisaoka M, Shi DR et al. Expression of anaplastic lymphoma kinase in soft tissue tumors: an immunohistochemical and molecular study of 249 cases. Human Pathology. 2004; 35:711-721
5. Savic S, Bode B, Diebold J et al. Detection of ALK-positive non-small-cell lung cancers on cytological specimens: high accuracy of immunocytochemistry with the 5A4 clone. Journal of Thoracic Oncology. 2013; 8(8):1004-11
6. Savic S, Diebold J, Zimmermann AK et al. Screening for ALK in non-small cell lung carcinomas: 5A4 and D5F3 antibodies perform equally well, but combined use with FISH is recommended. Lung Cancer. 2015; 89 (2):104-9
7. Williams AS, Greer W, Bethune D et al. ALK+ lung adenocarcinoma in never smokers and long-term ex-smokers: prevalence and detection by immunohistochemistry and fluorescence in situ hybridisation. Virchows Archiv: an international journal of pathology. 2016; 469(5):533-540

8. Barreca A, Lasorsa E, Riera L et al. Anaplastic lymphoma kinase in human cancer. *Journal of Molecular Endocrinology*. 2011; 47(1):R11-R23
9. Zhao Z, Verma V. and Zhang M. Anaplastic lymphoma kinase: Role in cancer and therapy perspective. *Cancer Biology and Therapy*. 2015; 16(12):1691-701

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## Literatúra

1. Clinical Laboratory Improvement Amendments of 1988, Final Rule 57 FR 1163 February 28, 1992.
2. Villanova PA. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. 1991; 7(9). Order code M29-P.
3. Bancroft JD and Stevens A. Theory and Practice of Histological Techniques. 4th Edition. Churchill Livingstone, New York. 1996.
4. Li XQ, Hisaoka M, Shi DR et al. Expression of anaplastic lymphoma kinase in soft tissue tumours: an immunohistochemical and molecular study of 249 cases. Human Pathology. 2004; 35:711-721
5. Savić S, Bode B, Diebold J et al. Detection of ALK-positive non-small-cell lung cancers on cytological specimens: high accuracy of immunocytochemistry with the 5A4 clone. Journal of Thoracic Oncology. 2013; 8(8):1004-11
6. Savić S, Diebold J, Zimmermann AK et al. Screening for ALK in non-small cell lung carcinomas: 5A4 and D5F3 antibodies perform equally well, but combined use with FISH is recommended. Lung Cancer. 2015; 89 (2):104-9
7. Williams AS, Greer W, Bethune D et al. ALK+ lung adenocarcinoma in never smokers and long-term ex-smokers: prevalence and detection by immunohistochemistry and fluorescence in situ hybridisation. Virchows Archiv: an international journal of pathology. 2016; 469(5):533-540
8. Barreca A, Lasorsa E, Riera L et al. Anaplastic lymphoma kinase in human cancer. Journal of Molecular Endocrinology. 2011; 47(1):R11-R23
9. Zhao Z, Verma V. and Zhang M. Anaplastic lymphoma kinase: Role in cancer and therapy perspective. Cancer Biology and Therapy. 2015; 16(12):1691-701

## Dátum vydania

30 november 2018

Leica Biosystems Newcastle Ltd  
Balliol Business Park  
Benton Lane  
Newcastle Upon Tyne NE12 8EW  
United Kingdom  
☎ +44 191 215 4242

Leica Biosystems Canada  
71 Four Valley Drive  
Concord, Ontario L4K 4V8  
Canada  
☎ +1 800 248 0123

Leica Biosystems Inc  
1700 Leider Lane  
Buffalo Grove IL 60089  
USA  
☎ +1 800 248 0123

Leica Biosystems Melbourne  
Pty Ltd  
495 Blackburn Road  
Mt Waverley VIC 3149  
Australia  
☎ +61 2 8870 3500