Novocastra™ Epitope Retrieval Solutions

Product Code: RE7113, RE7114, RE7116, RE7119

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

EN FR IT DE ES PT SV EL DA NL CS

Instructions for Use

Please read before using this product.

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Check the integrity of the packaging before use.

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Intended Use

For *in vitro* diagnostic use.

The Novocastra™ Epitope Retrieval Solutions RE7113, RE7114, RE7116 and RE7119 are intended for Heat Induced Epitope Retrieval on formalin-fixed, paraffin-embedded tissue sections as part of an immunohistochemical procedure. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Principle of Procedure

Heat Induced Epitope Retrieval of formalin-fixed, paraffin-embedded tissue sections using an appropriate pH solution improves the staining of some antibodies by exposing epitopes within tissue that have been masked during fixation. The development of Epitope Retrieval using heat began in 1991 with the report by Shi et al.¹ Since then numerous studies have been published looking at the effects of Epitope Retrieval solutions, molarity, pH, and heating methods.² A universal Heat Induced Epitope Retrieval technique suitable for all epitopes does not exist so a number of different heating methods and Epitope Retrieval solutions, including those listed below may be used.

These products are used in an immunohistochemical (IHC) procedure that allows the qualitative identification by light microscopy of antigens in sections of formalin-fixed, paraffin-embedded tissue, via sequential steps with interposed washing steps (for IHC staining Principle of Procedure see appropriate detection system instructions for Use). Heat Induced Epitope Retrieval is not recommended for all antibodies (see primary antibody Recommendations on Use). Optimum conditions for Epitope Retrieval should be validated by the user as these are dependant upon tissue, fixation and/or primary antibody.

Reagents Provided

One of the following:

RE7113 Epitope Retrieval Solution pH 6 (x10 Concentrate) 1L RE7114 Epitope Retrieval Solution pH 6 (x10 Concentrate) 500 ml	Citrate-based buffer containing surfactant
RE7116 Epitope Retrieval Solution pH 8 (x10 Concentrate) 1L	EDTA-based buffer containing surfactant
RE7119 Epitope Retrieval Solution pH 9 (x10 Concentrate) 1L	Tris/EDTA-based buffer containing surfactant

Reconstitution, Mixing, Dilution, Titration

The Epitope Retrieval Solutions (x10 Concentrate) RE7113, RE7114, RE7116 and RE7119 require dilution with de-ionized water to prepare working solutions (see Methodology). Further dilution may result in loss of antigen staining. The user must validate any such change.

Storage and Stability

Store at 2–8 °C. Do not freeze. Return to 2–8 °C immediately after use. Do not use after expiration date indicated on the product label. Storage conditions other than those specified must be verified by the user. There are no obvious signs to indicate instability of this product, therefore, positive and negative controls should be run simultaneously with patient samples.

Specimen Preparation

The recommended fixative is 10% neutral-buffered formalin for paraffin-embedded tissue sections.

Warnings and Precautions

One or more components in the product is hazardous.

A Material Safety Data Sheet is available upon request or available from www.LeicaBiosystems.com

For professional users. Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions.³

Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water.

Minimize microbial contamination of reagents or an increase in non-specific staining may occur. Incubation times or temperatures other than those specified may give erroneous results. Any such change must be validated by the user.

Procedure

A. Reagents required but not supplied

See primary antibody Instructions for Use.

B. Equipment required but not supplied

1. Stainless steel pressure cooker (it is recommended that the gaskets are changed at regular intervals to maintain optimum retrieval conditions). To ensure safe and correct use of the pressure cooker users must read the manufacturer's instructions.
2. General immunohistochemistry laboratory equipment.

C. Methodology

Prior to undertaking this methodology, users must be trained in immunohistochemical techniques.

The combination of the primary antibody, its dilution and optimum conditions for Epitope Retrieval, together with the detection system should be validated by the user on a series of known positive and negative controls.

1. If using Novocastra™ Epitope Retrieval Solutions (x10 Concentrate) RE7113, RE7114, RE7116 or RE7119 prepare a working solution by diluting 1 part concentrate with 9 parts de-ionized water.
2. Heat 1.5L of the working solution until boiling in a pressure cooker. Cover but do not lock lid. Position slides into metal staining racks (do not place slides close together as uneven staining may occur) and lower into pressure cooker ensuring slides are completely immersed in retrieval solution. Lock lid.
3. When the pressure cooker reaches operating temperature and pressure, time for 1 minute (optimum time should be validated by the user, as this is dependant upon tissue, fixation and/or primary antibody).
4. Remove pressure cooker from heat source and run under cold water with lid on. DO NOT OPEN LID UNTIL THE INDICATORS SHOW THAT PRESSURE HAS BEEN RELEASED. Open lid, remove slides and place immediately in cool tap water.
5. Proceed with IHC protocol according to manufacturers' Instructions for Use for the primary antibody and detection system.

Quality Control

Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls in addition to the following procedures. Controls should be fresh autopsy/biopsy/ surgical specimens, formalin-fixed, processed and paraffin wax-embedded as soon as possible in the same manner as the patient sample(s).

Positive Tissue Control

Used to indicate correctly prepared tissues and proper staining techniques. One positive tissue control should be included for each set of test conditions/primary antibody in each staining run. A tissue with weak positive staining is more suitable than a tissue with strong positive staining for optimal quality control and to detect minor levels of reagent degradation.⁴ For recommended positive control tissue see primary antibody Instructions for Use. If the positive tissue control fails to demonstrate positive staining, results with the test specimens should be considered invalid.

Negative Tissue Control

Should be examined after the positive tissue control to verify the specificity of the labeling of the target antigen by the primary antibody. For recommended negative control tissue see primary antibody Instructions for Use. Alternatively, the variety of different cell types present in most tissue sections frequently offers negative control sites, but this should be verified by the user. Non-specific staining, if present, usually has a diffuse appearance. Sporadic staining of connective tissue may also be observed in sections from excessively formalin-fixed tissues. Use intact cells for interpretation of staining results. Necrotic or degenerated cells often stain non-specifically.⁵ False-positive results may be seen due to non-immunological binding of proteins or substrate reaction products. They may also be caused by endogenous enzymes such as pseudoperoxidase (erythrocytes), endogenous peroxidase (cytochrome C), or endogenous biotin⁶ (eg. liver, breast, brain, kidney). To differentiate endogenous enzyme activity or non-specific binding of enzymes from specific immunoreactivity, additional patient tissues may be stained exclusively with substrate-chromogen, Streptavidin-HRP or labeled polymer and substrate-chromogen, respectively. If specific staining occurs in the negative tissue control, results with the patient specimens should be considered invalid.

Negative Reagent Control

Use a non-specific negative reagent control in place of the primary antibody with a section of each patient specimen to evaluate non-specific staining and allow better interpretation of specific staining at the antigen site.

Patient Tissue

Examine stained patient specimens last. Positive staining intensity should be assessed within the context of any non-specific background staining of the negative reagent control. As with any immunohistochemical test, a negative result means that the antigen was not detected, not that the antigen was absent in the cells/tissue assayed. If necessary, use a panel of antibodies to identify false-negative reactions.

Limitations

Immunohistochemistry is a multistep diagnostic process that consists of specialized training in the selection of the appropriate reagents; tissue selection, fixation, and processing; preparation of the IHC slide; and interpretation of the staining results.

Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning or contamination with other tissues or fluids may produce artifacts, antibody trapping, or false negative results. Inconsistent results may be due to variations in fixation and embedding methods, or to inherent irregularities within the tissue.⁷

Excessive or incomplete counterstaining may compromise proper interpretation of results.

The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Novocastra™ Epitope Retrieval Solutions are for use on paraffin-embedded sections with specific fixation requirements. Unexpected antigen expression may occur, especially in neoplasms. The clinical interpretation of any stained tissue section must include morphological analysis and the evaluation of appropriate controls.

Performance Characteristics

The performance of Novocastra™ Epitope Retrieval Solutions RE7113, RE7114, RE7116 and RE7119 have been validated using a range of Novocastra™ mouse IgG, mouse IgM and rabbit IgG primary antibodies.

These products are stable up to the expiry date(s) indicated on the product labels.

Bibliography

1. Shi S–R, Key ME and Kalra KL. Antigen retrieval in formalin-fixed, paraffin-embedded tissues: an enhancement method for immunohistochemical staining based upon microwave oven heating of tissue sections. *Journal of Histochemistry and Cytochemistry*. 1991; 39:741–748.
2. Shi S–R, Cote RJ, Taylor CR. Antigen Retrieval Immunohistochemistry: Past, Present, and Future. *Journal of Histochemistry and Cytochemistry*. 1997; 45(3):327–343.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29–P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio–Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL, Cooper K. Beware of biotin: a source of false–positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161–67.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

Amendments to Previous Issue

New information has been added to the Warnings and Precautions section.

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Bibliographie

1. Shi S-R, Key ME and Kalra KL. Antigen retrieval in formalin-fixed,paraffin-embedded tissues: an enhancement method for immunohistochemical staining based upon microwave oven heating of tissue sections. *Journal of Histochemistry and Cytochemistry*. 1991; 39:741-748.
2. Shi S-R, Cote RJ, Taylor CR. Antigen Retrieval Immunohistochemistry: Past, Present, and Future. *Journal of Histochemistry and Cytochemistry*. 1997; 45(3):327-343.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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1. Shi S-R, Key ME and Kalra KL. Antigen retrieval in formalin-fixed, paraffin-embedded tissues: an enhancement method for immunohistochemical staining based upon microwave oven heating of tissue sections. *Journal of Histochemistry and Cytochemistry*. 1991; 39:741-748.
2. Shi S-R, Cote RJ, Taylor CR. Antigen Retrieval Immunohistochemistry: Past, Present, and Future. *Journal of Histochemistry and Cytochemistry*. 1997; 45(3):327-343.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Literatur

1. Shi S–R, Key RE and Kalra KL. Antigen retrieval in formalin–fixed, paraffin–embedded tissues: an enhancement method for immunohistochemical staining based upon microwave oven heating of tissue sections. *Journal of Histochemistry and Cytochemistry*. 1991; 39:741–748.
2. Shi S–R, Cote RJ, Taylor CR. Antigen Retrieval Immunohistochemistry: Past, Present, and Future. *Journal of Histochemistry and Cytochemistry*. 1997; 45(3):327–343.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29–P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio–Preiser C, Wolff CM, Rilke F. Field & Wood, Inc. Philadelphia.
5. P5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL, Cooper K. Beware of biotin: a source of false–positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161–167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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1. Shi S-R, Key RE and Kalra KL. Antigen retrieval in formalin-fixed, paraffin-embedded tissues: an enhancement method for immunohistochemical staining based upon microwave oven heating of tissue sections. *Journal of Histochemistry and Cytochemistry*. 1991; 39:741-748.
2. Shi S-R, Cote RJ, Taylor CR. Antigen Retrieval Immunohistochemistry: Past, Present, and Future. *Journal of Histochemistry and Cytochemistry*. 1997; 45(3):327-343.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc. Philadelphia.
5. P5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Bibliografia

1. Shi S-R, Key RE and Kalra KL. Antigen retrieval in formalin-fixed, paraffin-embedded tissues: an enhancement method for immunohistochemical staining based upon microwave oven heating of tissue sections. *Journal of Histochemistry and Cytochemistry*. 1991; 39:741-748.
2. Shi S-R, Cote RJ, Taylor CR. Antigen Retrieval Immunohistochemistry: Past, Present, and Future. *Journal of Histochemistry and Cytochemistry*. 1997; 45(3):327-343.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1-15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL, Cooper K. Beware of biotin: a source of false-positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161-167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Bibliografi

1. Shi S–R, Key ME and Kalra KL. Antigen retrieval in formalin-fixed,paraffin-embedded tissues: an enhancement method for immunohistochemical staining based upon microwave oven heating of tissue sections. *Journal of Histochemistry and Cytochemistry*. 1991; 39:741–748.
2. Shi S–R, Cote RJ, Taylor CR. Antigen Retrieval Immunohistochemistry: Past, Present, and Future. *Journal of Histochemistry and Cytochemistry*. 1997; 45(3):327–343.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29–P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio–Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL, Cooper K. Beware of biotin: a source of false–positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161–167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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1. Shi S–R, Key RE and Kalra KL. Antigen retrieval in formalin–fixed,paraffin–embedded tissues: an enhancement method for immunohistochemical staining based upon microwave oven heating of tissue sections. *Journal of Histochemistry and Cytochemistry*. 1991; 39:741–748.
2. Shi S–R, Cote RJ, Taylor CR. Antigen Retrieval Immunohistochemistry: Past, Present, and Future. *Journal of Histochemistry and Cytochemistry*. 1997; 45(3):327–343.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29–P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio–Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL, Cooper K. Beware of biotin: a source of false–positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161–167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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1. Shi S–R, Key ME and Kalra KL. Antigen retrieval in formalin–fixed,paraffin–embedded tissues: an enhancement method for immunohistochemical staining based upon microwave oven heating of tissue sections. *Journal of Histochemistry and Cytochemistry*. 1991; 39:741–748.
2. Shi S–R, Cote RJ, Taylor CR. Antigen Retrieval Immunohistochemistry: Past, Present, and Future. *Journal of Histochemistry and Cytochemistry*. 1997; 45(3):327–343.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29–P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio–Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL, Cooper K. Beware of biotin: a source of false–positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161–167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Literatura

1. Shi S–R, Key ME and Kalra KL. Antigen retrieval in formalin-fixed,paraffin-embedded tissues: an enhancement method for immunohistochemical staining based upon microwave oven heating of tissue sections. *Journal of Histochemistry and Cytochemistry*. 1991; 39:741–748.
2. Shi S–R, Cote RJ, Taylor CR. Antigen Retrieval Immunohistochemistry: Past, Present, and Future. *Journal of Histochemistry and Cytochemistry*. 1997; 45(3):327–343.
3. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29–P.
4. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio–Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
5. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
6. Mount SL, Cooper K. Beware of biotin: a source of false–positive immunohistochemistry. *Current Diagnostic Pathology*. 2001; 7:161–167.
7. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.

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Nepoužity.

Datum vydání

05 Listopadu 2010

Leica Biosystems Newcastle Ltd
Balliol Business Park West
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242