Novocastra™ Liquid Mouse Monoclonal Antibody ZAP-70

Product Code: NCL-L-ZAP-70

Leica Biosystems Newcastle Ltd
Balliol Business Park
Benton Lane
Newcastle Upon Tyne NE12 8EW
United Kingdom
☎ +44 191 215 4242

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Intended Use

For in vitro diagnostic use.

NCL-L-ZAP-70 is intended for the qualitative identification by light microscopy of ZAP-70 molecules in paraffin sections. The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Principle of Procedure

Immunohistochemical (IHC) staining techniques allow for the visualization of antigens via the sequential application of a specific antibody to the antigen (primary antibody), a secondary antibody to the primary antibody and an enzyme complex with a chromogenic substrate with interposed washing steps. The enzymatic activation of the chromogen results in a visible reaction product at the antigen site. The specimen may then be counterstained and coverslipped. Results are interpreted using a light microscope and aid in the differential diagnosis of pathophysiological processes, which may or may not be associated with a particular antigen.

Clone

L453R

Immunogen

A prokaryotic recombinant protein corresponding to a portion of the N-terminus.

Specificity

Human ZAP-70.

Reagent Composition

NCL-L-ZAP-70 is a liquid tissue culture supernatant, dialyzed and diluted in PBS containing sodium azide as a preservative.

Ig Class

IgG2b (kappa)

Total Protein Concentration

Total Protein

Refer to vial label for lot specific total protein concentration.

Antibody Concentration

Greater than or equal to 449 mg/L. Refer to vial label for lot specific Ig concentration.

Recommendations On Use

Immunohistochemistry on paraffin sections.

Heat Induced Epitope Retrieval (HIER): Please follow the instructions for use in Novocastra Epitope Retrieval Solution pH 6.

Suggested dilution: 1:100 for 30 minutes at 25 °C. This is provided as a guide and users should determine their own optimal working dilutions.

Visualization: Please follow the instructions for use in the Novolink™ Polymer Detection Systems. For further product information or support, contact your local distributor or regional office of Leica Biosystems, or alternatively, visit the Leica Biosystems' Web site, www.LeicaBiosystems.com

Storage and Stability

Store at 2–8 °C. Do not freeze. Return to 2–8 °C immediately after use. Do not use after expiration date indicated on the vial label. Storage conditions other than those specified above must be verified by the user.

Specimen Preparation

The recommended fixative is 10% neutral-buffered formalin for paraffin-embedded tissue sections.

Warnings and Precautions

This reagent has been prepared from the supernatant of cell culture. As it is a biological product, reasonable care should be taken when handling it.

This reagent contains sodium azide. A Material Safety Data Sheet is available upon request or available from www.LeicaBiosystems.com

Consult federal, state or local regulations for disposal of any potentially toxic components.

Specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions.† Never pipette reagents by mouth and avoid contacting the skin and mucous membranes with reagents and specimens. If reagents or specimens come in contact with sensitive areas, wash with copious amounts of water. Seek medical advice.

Minimize microbial contamination of reagents or an increase in non-specific staining may occur.

Incubation times or temperatures, other than those specified, may give erroneous results. Any such changes must be validated by the user.

Quality Control

Differences in tissue processing and technical procedures in the user's laboratory may produce significant variability in results, necessitating regular performance of in-house controls in addition to the following procedures.

Controls should be fresh autopsies/biopsy/surgical specimens, formalin-fixed, processed and paraffin wax-embedded as soon as possible in the same manner as the patient sample(s).

Positive Tissue Control

Used to indicate correctly prepared tissues and proper staining techniques.

One positive tissue control should be included for each set of test conditions in each staining run.

A tissue with weak positive staining is more suitable than a tissue with strong positive staining for optimal quality control and to detect minor levels of reagent degradation.²

Recommended positive control tissue is tonsil.

If the positive tissue control fails to demonstrate positive staining, results with the test specimens should be considered invalid.

Negative Tissue Control

Should be examined after the positive tissue control to verify the specificity of the labeling of the target antigen by the primary antibody.

Recommended negative control tissue is cerebellum.

Alternatively, the variety of different cell types present in most tissue sections frequently offers negative control sites, but this should be verified by the user.

Non-specific staining, if present, usually has a diffuse appearance. Sporadic staining of connective tissue may also be observed in sections from excessively formalin-fixed tissues. Use intact cells for interpretation of staining results. Necrotic or degenerated cells often stain non-specifically.³ False-positive results may be seen due to non-immunological binding of proteins or substrate reaction products. They may also be caused by endogenous enzymes such as pseudoperoxidase (erythrocytes), endogenous peroxidase (cytochrome C), or endogenous biotin (eg. liver, breast, brain, kidney) depending on the type of immunostain used. To differentiate endogenous enzyme activity or non-specific binding of enzymes from specific immunoreactivity, additional patient tissues may be stained exclusively with substrate chromogen or enzyme complexes (avidin-biotin, streptavidin, labeled polymer) and substrate-chromogen, respectively. If specific staining occurs in the negative tissue control, results with the patient specimens should be considered invalid.

Negative Reagent Control

Use a non-specific negative reagent control in place of the primary antibody with a section of each patient specimen to evaluate non-specific staining and allow better interpretation of specific staining at the antigen site.

Patient Tissue

Examine patient specimens stained with NCL-L-ZAP-70 last. Positive staining intensity should be assessed within the context of any non-specific background staining of the negative reagent control. As with any immunohistochemical test, a negative result means that the antigen was not detected, not that the antigen was absent in the cells/tissue assayed. If necessary, use a panel of antibodies to identify false-negative reactions.

Results Expected

Normal Tissues

Clone L453R detected the ZAP-70 protein in the cytoplasm and nucleus of T cells and NK cells. In benign lymphoid tissue paracortical T-lymphocytes and rare scattered small lymphocytes in the mantle zones and follicle centers were positive. In the normal tissue ZAP-70 cytoplasmic and/or nuclear staining was predominantly observed in infiltrating lymphocytes. Occasionally NK cells were stained and there was no staining of normal mature B-cells. Other than infiltrating lymphocytes, IHC studies have shown no unexpected reactivity, suggesting the clone is specific to ZAP-70 and not cross-reactive with Syk. (n=39).

Abnormal Tissues

Clone L453R stained 6/7 peripheral T-cell lymphomas, 1/1 T-lymphoblastic lymphoma, 1/1 angioimmunoblastic T-cell lymphoma, 1/1 NK/T cell lymphoma and 0/2 myeloid leukemias. In B-cell neoplasms, clone L453R stained 13/27 chronic lymphocytic leukemias, 2/2 mantle cell lymphomas, 1/4 anaplastic large cell lymphomas, 0/18 diffuse large cell B-cell lymphomas, 0/7 follicular lymphomas, 0/2 extranodal marginal zone B-cell lymphomas of mucosa-associated lymphoid tissue and 0/1 Burkitt's lymphoma. In Hodgkin's lymphoma, the neoplastic cells were negative in all cases studied (0/5). Except for infiltrating lymphocytes, no staining was seen in non-hematological malignancies (n=68).

NCL-L-ZAP-70 is recommended for the detection of human ZAP-70 protein in normal and neoplastic tissues, as an adjunct to conventional histopathology using non-immunologic histochemical stains.

General Limitations

Immunohistochemistry is a multistep diagnostic process that consists of specialized training in the selection of the appropriate reagents; tissue selection, fixation, and processing; preparation of the IHC slide; and interpretation of the staining results.

Tissue staining is dependent on the handling and processing of the tissue prior to staining. Improper fixation, freezing, thawing, washing, drying, heating, sectioning or contamination with other tissues or fluids may produce artifacts, antibody trapping, or false negative results. Inconsistent results may be due to variations in fixation and embedding methods, or to inherent irregularities within the tissue.⁴ Excessive or incomplete counterstaining may compromise proper interpretation of results.

The clinical interpretation of any staining or its absence should be complemented by morphological studies using proper controls and should be evaluated within the context of the patient's clinical history and other diagnostic tests by a qualified pathologist.

Antibodies from Leica Biosystems Newcastle Ltd are for use, as indicated, on either frozen or paraffin-embedded sections with specific fixation requirements. Unexpected antigen expression may occur, especially in neoplasms. The clinical interpretation of any stained tissue section must include morphological analysis and the evaluation of appropriate controls.

Bibliography - General

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.
5. Orchard J, Ibbotson R, Best G, et al. ZAP-70 in B cell malignancies. *Leukaemia and Lymphoma* 2005; 46(12):1689–1698.
6. Wang J, Young L, Win W, et al. Distribution and ZAP-70 expression of WHO categories in Shanxi, China. *Applied Immunohistochemistry and Molecular Morphology* 2005; 13(4):323–332.
7. Mustelin T and Tasken K. Positive and negative regulation of T-cell activation through kinases and phosphatases. *Biochemical Journal* 2003; 371:15–27.
8. Van Oers N and Weiss A. The Syk/ZAP-70 protein tyrosine kinase connection to antigen receptor signalling processes. *Seminars in Immunology* 1995; 7:227–236.

Amendments to Previous Issue

Reagent Composition, Total Protein Concentration, Recommendations On Use, Warnings and Precautions.

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1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.
5. Orchard J, Ibbotson R, Best G, et al. ZAP-70 in B cell malignancies. Leukaemia and Lymphoma 2005; 46(12):1689–1698.
6. Wang J, Young L, Win W, et al. Distribution and ZAP-70 expression of WHO categories in Shanxi, China. Applied Immunohistochemistry and Molecular Morphology 2005; 13(4):323–332.
7. Mustelin T and Tasken K. Positive and negative regulation of T-cell activation through kinases and phosphatases. Biochemical Journal 2003; 371:15–27.
8. Van Oers N and Weiss A. The Syk/ZAP-70 protein tyrosine kinase connection to antigen receptor signalling processes. Seminars in Immunology 1995; 7:227–236.

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Literatur - Allgemein

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991;7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.
5. Orchard J, Ibbotson R, Best G, et al. ZAP-70 in B cell malignancies. Leukaemia and Lymphoma 2005; 46(12):1689–1698.
6. Wang J, Young L, Win W, et al. Distribution and ZAP-70 expression of WHO categories in Shanxi, China. Applied Immunohistochemistry and Molecular Morphology 2005; 13(4):323–332.
7. Mustelin T and Tasken K. Positive and negative regulation of T-cell activation through kinases and phosphatases. Biochemical Journal 2003; 371:15–27.
8. Van Oers N and Weiss A. The Syk/ZAP-70 protein tyrosine kinase connection to antigen receptor signalling processes. Seminars in Immunology 1995; 7:227–236.

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Bibliografi - Allmän

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F, Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and its pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.
5. Orchard J, Ibbotson R, Best G, et al. ZAP-70 in B cell malignancies. *Leukaemia and Lymphoma* 2005; 46(12):1689–1698.
6. Wang J, Young L, Win W, et al. Distribution and ZAP-70 expression of WHO categories in Shanxi, China. *Applied Immunohistochemistry and Molecular Morphology* 2005; 13(4):323–332.
7. Mustelin T and Tasken K. Positive and negative regulation of T-cell activation through kinases and phosphatases. *Biochemical Journal* 2003; 371:15–27.
8. Van Oers N and Weiss A. The Syk/ZAP-70 protein tyrosine kinase connection to antigen receptor signalling processes. *Seminars in Immunology* 1995; 7:227–236.

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1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavaï M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.
5. Orchard J, Ibbotson R, Best G, et al. ZAP-70 in B cell malignancies. *Leukaemia and Lymphoma* 2005; 46(12):1689–1698.
6. Wang J, Young L, Win W, et al. Distribution and ZAP-70 expression of WHO categories in Shanxi, China. *Applied Immunohistochemistry and Molecular Morphology* 2005; 13(4):323–332.
7. Mustelin T and Tasken K. Positive and negative regulation of T-cell activation through kinases and phosphatases. *Biochemical Journal* 2003; 371:15–27.
8. Van Oers N and Weiss A. The Syk/ZAP-70 protein tyrosine kinase connection to antigen receptor signalling processes. *Seminars in Immunology* 1995; 7:227–236.

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Bibliografi – generell

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.

2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.
5. Orchard J, Ibbotson R, Best G, et al. ZAP-70 in B cell malignancies. *Leukaemia and Lymphoma* 2005; 46(12):1689–1698.
6. Wang J, Young L, Win W, et al. Distribution and ZAP-70 expression of WHO categories in Shanxi, China. *Applied Immunohistochemistry and Molecular Morphology* 2005; 13(4):323–332.
7. Mustelin T and Tasken K. Positive and negative regulation of T-cell activation through kinases and phosphatases. *Biochemical Journal* 2003; 371:15–27.
8. Van Oers N and Weiss A. The Syk/ZAP-70 protein tyrosine kinase connection to antigen receptor signalling processes. *Seminars in Immunology* 1995; 7:227–236.

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Bibliografie - General

1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. *Progress in Surgical Pathology*. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F, Field & Wood, Inc., Philadelphia.
3. Nadji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. *Laboratory Medicine*. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. *American Journal of Clinical Pathology*. 1980; 73:626.
5. Orchard J, Ibbotson R, Best G, et al. ZAP-70 in B cell malignancies. *Leukaemia and Lymphoma* 2005; 46(12):1689–1698.
6. Wang J, Young L, Win W, et al. Distribution and ZAP-70 expression of WHO categories in Shanxi, China. *Applied Immunohistochemistry and Molecular Morphology* 2005; 13(4):323–332.
7. Mustelin T and Tasken K. Positive and negative regulation of T-cell activation through kinases and phosphatases. *Biochemical Journal* 2003; 371:15–27.
8. Van Oers N and Weiss A. The Syk/ZAP-70 protein tyrosine kinase connection to antigen receptor signalling processes. *Seminars in Immunology* 1995; 7:227–236.

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1. National Committee for Clinical Laboratory Standards (NCCLS). Protection of laboratory workers from infectious diseases transmitted by blood and tissue; proposed guideline. Villanova, P.A. 1991; 7(9). Order code M29-P.
2. Battifora H. Diagnostic uses of antibodies to keratins: a review and immunohistochemical comparison of seven monoclonal and three polyclonal antibodies. Progress in Surgical Pathology. 6:1–15. eds. Fenoglio-Preiser C, Wolff CM, Rilke F. Field & Wood, Inc., Philadelphia.
3. Nadjji M, Morales AR. Immunoperoxidase, part I: the techniques and pitfalls. Laboratory Medicine. 1983; 14:767.
4. Omata M, Liew CT, Ashcavai M, Peters RL. Nonimmunologic binding of horseradish peroxidase to hepatitis B surface antigen: a possible source of error in immunohistochemistry. American Journal of Clinical Pathology. 1980; 73:626.
5. Orchard J, Ibbotson R, Best G, et al. ZAP-70 in B cell malignancies. Leukaemia and Lymphoma 2005; 46(12):1689–1698.
6. Wang J, Young L, Win W, et al. Distribution and ZAP-70 expression of WHO categories in Shanxi, China. Applied Immunohistochemistry and Molecular Morphology 2005; 13(4):323–332.
7. Mustelin T and Tasken K. Positive and negative regulation of T-cell activation through kinases and phosphatases. Biochemical Journal 2003; 371:15–27.
8. Van Oers N and Weiss A. The Syk/ZAP-70 protein tyrosine kinase connection to antigen receptor signalling processes. Seminars in Immunology 1995; 7:227–236.

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Leica Biosystems Inc
Pty Ltd