Leica Biosystems Cytology Stains.

English.....	2
العربية (Arabic).....	5
简体中文 (Chinese Simplified).....	8
繁體中文 (Chinese Traditional).....	11
Dansk (Danish).....	14
Nederlands (Dutch).....	17
Français (French – Canada).....	21
Français (French – France).....	25
Deutsch (German).....	29
Italiano (Italian).....	33
日本語 (Japanese).....	36
한국어 (Korean).....	39
Norsk (Norwegian).....	42
Polski (Polish).....	45
Português do Brasil (Portuguese – Brazil).....	49
Português (Portuguese – Portugal).....	52
Română (Romanian).....	55
Русский (Russian).....	58
Slovenščina (Slovenian).....	62
Español (Spanish – Central America).....	65
Español (Spanish – Spain).....	69
Svenska (Swedish).....	73
ภาษาไทย (Thai).....	76
Türkçe (Turkish).....	80
Tiếng Việt (Vietnamese).....	83

Cytology Stains

- REF**
- **Orange G-6: 3801660, 3801660E, 3801660BBE, 3801661, 3801661E, 3801662, 3801662E**
 - **EA-50: 3801620, 3801620E, 3801620BBE, 3801621, 3801621E, 3801622, 3801622E**
 - **EA-50 Modified: 3801680, 3801680E, 3801681, 3801681E, 3801682, 3801682E**
 - **EA-65: 3801640, 3801640E, 3801640BBE, 3801641, 3801641E, 3801642, 3801642E**

Product Name

Leica Biosystems Cytology Stains.

Intended Use

Detection/Measurement

The Leica Biosystems Cytology Stains do not detect or measure an analyte or marker. The Leica Biosystems Cytology Stains are used to produce various colorations of the cytoplasm within the cellular elements of cytology specimens.

Product Function

The results obtained through use of the Leica Biosystems Cytology Stains do not provide objective medical evidence. The Leica Biosystems Cytology Stains are used in a cytology staining protocol to provide contrast and coloration the cytoplasm of the cellular elements. The cytology stained specimen, when interpreted by a trained professional, is utilized alongside other information such as the patient’s medical history, physical condition, as well as results from other medical testing to render a medical diagnosis.

Specific Information Provided

Leica Biosystems Cytology Stains are not intended for the detection, definition or differentiation of a specific disorder, condition, or risk factor. The staining demonstrated with use of these products, when used as intended, provides trained professionals information which may define the physiological or pathological state of the tissue specimen.

Automation

The Leica Biosystems Cytology Stains are not automated but can be used on automated staining platforms. Use on an automated platform should be validated at the point of use.

Qualitative/Quantitative

The Leica Biosystems Cytology Stains are qualitative stains.

Specimen Type

The Leica Biosystems Cytology Stains may be used with cytologic specimens: smears, tissue imprints, cell cultures, monolayers.

Testing Population

The Leica Cytology Stains are intended for use with any patient requiring evaluation of a cytology or biopsy specimen for the assessment of a suspected pathology or disease.

Intended User

The Leica Biosystems Cytology Stains are intended for use by qualified laboratory personnel and/or designee of the laboratory.

In Vitro Diagnostic

The Leica Biosystems Cytology Stains are intended for *in vitro* diagnostics use only.

Intended User

The Leica Biosystems Cytology Stains are intended for use by qualified laboratory personnel and/or designee.

Test Principle

Orange G-6 stain will stain keratinized cells. EA stains will stain the cytoplasm of cytology specimens’ preparations.

Calibrators & Controls

The Leica Cytology Stains do not require the use of any calibrators or controls.

Reagent Limitations

Not Applicable.

Applicable Products

Product Code	Material Description
3801660	Orange G-6 (1 gal)
3801660E	Orange G-6 (5 L)
3801660BBE	Orange G-6 (5 L bag in box)
3801661	Orange G-6 (1/2 gal)
3801661E	Orange G-6 (2.5 L)
3801662	Orange G-6 (1 QT)
3801662E	Orange G-6 (1 L)
3801620	EA-50 (1 gal)
3801620E	EA-50 (5 L)
3801620BBE	EA-50 (5L bag in a box)
3801621	EA-50 (1/2 gal)

Cytology Stains

- REF**
- **Orange G-6: 3801660, 3801660E, 3801660BBE, 3801661, 3801661E, 3801662, 3801662E**
 - **EA-50: 3801620, 3801620E, 3801620BBE, 3801621, 3801621E, 3801622, 3801622E**
 - **EA-50 Modified: 3801680, 3801680E, 3801681, 3801681E, 3801682, 3801682E**
 - **EA-65: 3801640, 3801640E, 3801640BBE, 3801641, 3801641E, 3801642, 3801642E**

3801621E	EA-50 (2.5 L)
3801622	EA-50 (1 QT)
3801622E	EA-50 (1 L)
3801680	EA-50 Modified (1 gal)
3801680E	EA-50 Modified (5 L)
3801681	EA-50 Modified (1/2 gal)
3801681E	EA-50 Modified (2.5 L)
3801682	EA-50 Modified (1 QT)
3801682E	EA-50 Modified (1 L)
3801640	EA-65 (1 gal)
3801640E	EA-65 (5 L)
3801640BBE	EA-65 (5 L bag in a box)
3801641	EA-65 (1/2 gal)
3801641E	EA-65 (2.5 L)
3801642	EA-65 (1 QT)
3801642E	EA-65 (1 L)

NOTE: Products listed here may not be available in all regions.

Materials Not Included

The Leica Biosystems Cytology stains are commonly used with conjunction with hematoxylin and ancillary reagents as well as additional staining reagents, which are not included.

Devices Required

Leica Biosystems Cytology Stains can be used on any automated staining platform or with a manual staining method.

Storage and Stability

The product shall be stable for 24 months post-production when stored at ambient temperature.

In cold weather, chilling may cause some separation. Warm to room temperature, agitate and filter before using.

Store reagents at room temperature (15-30°C) in a well-ventilated place.

CAUTION: Do not use after the expiration date.

In Use Stability

User discretion should be utilized when determining in-use stability.

Sterility

The Leica Biosystems Cytology Stains are not sterile products.

Warnings/Precautions

1. Use the following personal protective equipment when working with this product.

- Eyes - Wear safety glasses or goggles.
- Body - Wear laboratory coat.
- Respiratory - Use NIOSH/MSHA approved respirator when ventilation is not adequate.
- Hands - Wear latex or nitrile gloves.

2. Emergency Overview: Highly flammable liquids and vapors. Vapor may cause flash fire. Cannot be made non-poisonous. May be fatal or cause blindness if swallowed. Contains material that may cause blood, nervous system, reproductive system, liver, gastrointestinal tract, respiratory tract, skin, and eye damage. Keep away from heat, sparks, and flame. Keep container closed. Use only with adequate ventilation.

Infectious Material Status

The Leica Biosystems Cytology Stains do not include any infectious material. However, specimens, before and after fixation, and all materials exposed to them, should be handled as if capable of transmitting infection and disposed of with proper precautions per facility guidelines.

Special Facilities

The Leica Biosystems Cytology Stains should be used per facility guidelines.

Cytology Stains

- REF**
- **Orange G-6: 3801660, 3801660E, 3801660BBE, 3801661, 3801661E, 3801662, 3801662E**
 - **EA-50: 3801620, 3801620E, 3801620BBE, 3801621, 3801621E, 3801622, 3801622E**
 - **EA-50 Modified: 3801680, 3801680E, 3801681, 3801681E, 3801682, 3801682E**
 - **EA-65: 3801640, 3801640E, 3801640BBE, 3801641, 3801641E, 3801642, 3801642E**

Specimen Handling

Suggested fixatives include ethanol or routine cytology fixative. Ensure that the material is fixed before air drying occurs. Ensure that any polyethylene glycol has been fully removed in ethanol before rehydration.

Preparation for Use

The Cytology Stains products are available in a ready-to-use format.

1. Rehydrate as necessary using routine procedure. Immerse in Hematoxylin until desired intensity is reached.
2. Wash in tap water.
3. Differentiate as necessary when utilizing a regressive stain method. The stain should now only remain in the nuclei. This step should be controlled microscopically for optimum results. Under or over staining can be corrected at this stage by repeating the above steps as appropriate.
4. Wash well in water (the pH of the tap water should be above 7 to ensure correct color of nuclei).
5. Rinse in 95% alcohol.
6. Blue if necessary in alkaline tap water or commercial bluing agent.
7. Immerse in Orange G-6 until desired intensity is reached.
8. Rinse twice in 95% alcohol.
9. Stain in EA-50, Modified EA-50, or EA-65 for 3 minutes or until desired intensity is reached.
10. Rinse well in 95% alcohol.
11. Dehydrate, clear and mount in synthetic resin.
12. To avoid contamination, filter solutions daily and change at least weekly.

Results: The nuclei and their components should be stained in varying shades of blue/purple.

- Keratinized cells - orange.
- Cytoplasm of superficial cells - pink.
- Cytoplasm of intermediate cells - blue/green.
- Cytoplasm of parabasal cells - blue/green.

Readiness for Use

Once appropriate staining protocol is chosen and bath layout is created, pour all the reagent into the reagent vessel. Place the reagent vessel back into the respective station.

Quality Control

A routine quality control slide(s) containing tissue fixed and processed in a similar manner to the test specimens should be performed prior to routine use to ensure Cytology Stains are performing as intended.

Expected Results

The nuclei and their components should be stained in varying shades of blue/purple. Keratinized cells – orange, cytoplasm of superficial cells – pink, cytoplasm of intermediate cells - blue/green, cytoplasm of parabasal cells - blue/green.

Analytical Performance

The Leica Biosystems Cytology Stains are not used to detect a specific analyte or marker. These products are used in conjunction with other products in a Papanicolaou stain protocol to stain cytology specimens cytoplasm various colors. Analytical parameters such as analytical sensitivity, analytical specificity, trueness (bias), precision (repeatability and reproducibility), accuracy (resulting from trueness and precision), limits of detection and quantitation, measuring range, linearity, cut-off, including determination of appropriate criteria for specimen collection and handling and control of known relevant endogenous and exogenous interference, cross-reactions do not apply to the performance of this system.

Clinical Performance

The Leica Biosystems Cytology Stains are not intended for use as a means of detecting a specific disease or pathological process or state. Clinical performance indices such as diagnostic sensitivity, diagnostic specificity, positive predictive value, negative predictive value, likelihood ratio as well as expected values in normal and affected populations do not apply to the use of the Leica Biosystems Bluing Agents in a clinical setting


Disposal

Spent or excess of Cytology Stains should be discarded in accordance with organizational, local, state, and federal regulations.

Leica Biosystems Richmond, Inc.
5205 Route 12
Richmond, IL 60071
USA
(1-844-534-2262)

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